

Sep 08, 2022

# 🌐 S-Trap™ plate digestion protocol (Protifi) of proteins for LC-MS / proteomics

📁 In 1 collection

DOI

[dx.doi.org/10.17504/protocols.io.kxygxzd2zv8j/v1](https://dx.doi.org/10.17504/protocols.io.kxygxzd2zv8j/v1)

ronan o'cualain<sup>1</sup>, David Knight<sup>1</sup>, Stacey Warwood<sup>1</sup>, James Allsey<sup>1</sup>, Emmakeevill<sup>1</sup>

<sup>1</sup>University of Manchester

BioMS CRF, UoM



ronan o'cualain

Inoviv



## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygxzd2zv8j/v1>

**Protocol Citation:** ronan o'cualain, David Knight, Stacey Warwood, James Allsey, Emmakeevill 2022. S-Trap™ plate digestion protocol (Protifi) of proteins for LC-MS / proteomics. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kxygxzd2zv8j/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** July 20, 2022

**Last Modified:** September 08, 2022

**Protocol Integer ID:** 67156

**Keywords:** S-trap™ 96-well plate, Digestion, Elution, Digestion protocol, Eppendorf thermomixer, Mass spec analysis, proteomics, quartz, s-trapping, S-trap, Protifi, desalting, clean-up, LC-MS, off-line, offline, enzymatic digestion, automatable, proteomics this protocol, proteomic, proteins for lc, well plate protein clean, long protocol from protifi, protein, protifi, digestion, protocol detail, house bioms procedure

## Abstract

This protocol details the in-house BioMS procedure of S-Trap™ 96-well plate protein clean-up and digestion.

It is adapted from the long protocol from Protifi (as on August 2022) - <https://files.protifi.com/protocols/s-trap-96-well-plate-long-1-4.pdf>










## Attachments






[iiaebptmp.docx](#)

173KB

## Guidelines

- You may purchase S-trap plate positions (or full plates) and other consumables directly from us at the BioMS core facility on PPMS - <https://corefacilities.manchester.ac.uk/?BioMS>).
- Allow  02:30:00 for this process.
- This plate based protocol is to be used for samples with a protein load between  50 µg to  250 µg of starting material.<sup>1</sup>
- You have cell or tissue lysates in  50 µL of S-Trap lysis buffer (  5 % (v/v) SDS with  50 millimolar (mM) TEAB  7.5 ).<sup>2</sup>
- Protein lysates have been sheared of nucleic acid (e.g. from BioMS's LE220+ AFA sonication system).<sup>3</sup>
- Protein lysates have been reduced and alkylated and clarified by centrifuging at  14000 x g for  00:10:00 (see the protocol "*Reduction and alkylation of protein lysates for LC-MS (proteomics) using dithiothreitol (DTT) and iodoacetamide (IAM)*" in this collection).<sup>4</sup>
- The protein concentration of a pool of the protein lysates has been quantified and samples prepared to be 100 ug of protein in 50ul of 5% SDS (see notes).

### Notes:



1. If you have less than  50 µg of protein - use the "*S-Trap™ column digestion protocol (Protifi) of proteins for LC-MS / proteomics*" in this collection. Speak with a member of BioMS if you are unsure which option to use.
2. Alternatively, if your lysate volume is larger than  50 µL simply remove a  25 µL sample aliquot for processing.
3. This step is essential as omitting it may clog the S-Trap column. If you do not have access to ultrasonication, see protocol from **EMDmillipore** for details.
4. If it is thought that the insoluble pellet contains proteins of interest, then proceed with homogenising (with a pipette and tip) the insoluble pellet as much as possible with fresh 5% S-trap lysis buffer. A suspension of insoluble particles may be formed this way. Proceed with a separate S-trap column to process this additional sample.

## Materials

### Locate the following buffers and reagents:

A	B
Location	Buffer/reagent
Bench	-- 12% phosphoric acid -- S-Trap binding buffer (90% aqueous methanol containing a final concentration of 100 mM TEAB, pH 7.1) -- 0.1% formic acid in water -- 0.1% formic acid in 30% acetonitrile -- MTBE / Methanol solution (10/3 (v/v))
Freezer 3	Aliquots of trypsin (10 uL at 2ug.uL <sup>-1</sup> )

### Locate the following consumables:





- 2x S-Trap plates (one for the processing, and another to use as a balance in the centrifuge - the balance one will be next to the centrifuge you will be using). Please purchase the use of the S-trap plate in advance from BioMS core facility on PPMS - <https://corefacilities.manchester.ac.uk/?BioMS>).
- 2x fresh collection plates: one for sample flow-through and washes, use the other clean plate for elution and collection of peptides generated by the process.
-  1.5 mL <sup>(i)</sup>, and  2 mL <sup>(ii)</sup> Eppendorf tubes

for





<sup>(i)</sup> diluting samples (if necessary), and for the digestion step.



and

<sup>(ii)</sup> collecting unbound material from the S-trapping process.

- You may need additional S-Trap binding buffer, there are aliquots of  5 mL of 1M 100 millimolar (mM) TEAB,  7.1 stored in freezer 3, thaw and add  45 mL of methanol to make a final volume of  50 mL to use.

### Identify the following equipment that you will use:

- one each of  10 µL or  20 µL , a  200 µL pipette, and a  1 mL pipette.
- A plate centrifuge such as a Thermo Megafuge 16 with plate rotor fitted. )
- A Vortex mixer (if diluting samples).
- Eppendorf Thermomixer.
- An Eppendorf Thermomixer with 1.5 mL thermoblock and set it to 47 °C, 01:00:00, and a speed of 0 rpm (i.e. no shaking).

Set the Eppendorf Thermomixer to  47 °C ,  01:00:00 , with the PCR 96 thermoblock, and a speed of 0

**RPM (i.e. no shaking).**



## Catalogue numbers:

- ⊗ ortho-Phosphoric acid 85% **Merck MilliporeSigma (Sigma-Aldrich) Catalog #1005732500**
- ⊗ Methanol Optima™ LC/MS Grade Fisher Chemical **Fisher Scientific Catalog #A456-4**
- ⊗ Pierce&trade; 0.1% Formic Acid (v/v) in Acetonitrile, LC-MS Grade **Thermo Fisher Catalog #85174**
- ⊗ Pierce&trade; 0.1% Formic Acid (v/v) in Water, LC-MS Grade **Thermo Fisher Catalog #85170**
- ⊗ tert-Butyl methyl ether 99% **Alfa Aesar Catalog #L14030**
- ⊗ Trypsin TPCK Treated **Worthington Biochemical Corporation Catalog #LS003740**
- ⊗ SDS sodium dodecyl sulfate **Catalog #75746-250G**
- ⊗ Triethylammonium bicarbonate (TEAB) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T7408**
- ⊗ Protifi S-trap 96-well plate **Catalog #C02-96well-1**

## Troubleshooting





## Safety warnings




- ⚠ Please refer to the copies of Risk Assessment Forms held in both B2071 and B2075 for hazards to health, and other identified hazards and risks, associated with the use of this protocol before starting.

## Before start

Prepare your protein samples using the other protocols in this collection.

The following steps are optimised for volumes of  50  $\mu\text{L}$  and  100  $\mu\text{L}$  of protein.

For other **volumes** and **amounts** of protein, adjust accordingly, by dilution into S-trap lysis buffer. It is recommended that the final concentration of SDS be at least greater than  3 % (v/v) , and up to  15 % (v/v) , and a protein load between  50  $\mu\text{g}$  to  250  $\mu\text{g}$  for the process to work successfully.

Use the 1.5 mL adaptor for the Eppendorf Thermomixer, and set the thermomixer to  47 °C ,  01:15:00 , and a speed of  0 rpm (i.e. no shaking).

## Sample preparation

10m

- 1 To the 50  $\mu\text{L}$  volume of sample in S-trap lysis buffer, add 5  $\mu\text{L}$  of 12 % (v/v) aqueous phosphoric acid at 1:10 for a final concentration of 1.2 % (v/v) phosphoric acid and vortex mix.



### Note

#### Notes:

1.To create a 50  $\mu\text{L}$  sample with a concentration of 50  $\mu\text{g}/\mu\text{L}$  protein, You can estimate the amount of lysate required using the following calculation:

amount lysate ( $\mu\text{L}$ ) =  $50\mu\text{L}/\text{calculated pooled lysate concentration}(\text{in } \mu\text{g}/\mu\text{L})$  and make up the volume to 50  $\mu\text{L}$  with 1x SDS solubilization buffer,

**e.g.** if the pooled lysate was determined to have a concentration of 1.6  $\text{mg}/\text{mL}$  then take  $50/1.6 =$  31  $\mu\text{L}$  of each sample and add 19  $\mu\text{L}$  of 1x SDS solubilization buffer.

If your samples are dilute, i.e. less than 0.5  $\mu\text{g}/\mu\text{L}$  , it is be a good idea to concentrate your samples before proceeding with the S-trap plate process. Methods to do this include the use of a speed-vac or lyophilisation.

Speak with a member of the BioMS team before doing so.

### Note

This step is essential to completely denature proteins and trap them efficiently.

The pH will be  $\leq$  1.0 . If the sample pH is not  $\leq$  1.0 , add additional phosphoric acid to reach  $\text{pH} \leq$  1.0 .

A quick way to check the pH is to spot 2  $\mu\text{L}$  of the acidified lysate on a strip of filter paper.

The final phosphoric acid concentration is different between S-Trap micros, and minis/midis.

- 2 Add 350  $\mu\text{L}$  of S-Trap binding buffer to the acidified lysis buffer and mix.



- Put the S-Trap plate on top of a clean 96 well plate, add the acidified methanolic SDS lysate into the plate.

#### Note

No plate pre-equilibration is necessary. Solution typically begins to drip through immediately.

## Sample Trapping


30m

- Locate an S-trap balance plate, with a receiver 96 well plate beneath. Centrifuge the plate at 1000 x g for 00:02:00 in the Thermo megafuge 16 centrifuge.

2m



#### Note



- Repeat the previous two steps until all sample has been applied to the S-Trap plate.

#### Note





Protein should be trapped within the protein-trapping matrix of the plate.





## Sample Washing

30m





- 6 Wash captured protein with one wash of  200  $\mu$ L of MTBE solution, simply add  200  $\mu$ L of the MTBE solution to the column, and spin at  1000 x g for  00:02:00 .

2m



### Note

This will remove methanol insoluble biomolecules from the quartz filter.

- 7 Following this, perform three washes of  200  $\mu$ L of S-Trap binding buffer, again, add  200  $\mu$ L of the S-trap binding buffer, and centrifuge at  1000 x g for  00:02:00 .








2m



### Note

**Note:** If you wish, you may transfer the flow through and washes back into an eppendorf sample tube after each centrifugation step, otherwise empty the collection tube so that the washes do not come in contact with the binding matrix. If discarding the washes then collect in a beaker and put in acetonitrile/solvent waste when finished.

### Note

- **Note:** Depending on the number of protein samples you need to process, you may find that you need additional S-Trap binding buffer.
- If so, there are aliquots of  5 mL of  100 millimolar (mM) TEAB at  7.1 stored in  50 mL Falcon tubes in freezer 3.
- - take one out, thaw at  Room temperature , and add  45 mL of methanol (located in fume hood) to make a final volume of  50 mL , mix, and use.

## Digest proteins

1h 15m

- 8 Move S-Trap digestion plate on top of a clean receiver plate.





9 Locate the trypsin aliquots. They are in the top shelf of freezer 3.



Note

Trypsin must be added to the protein at a ratio of 1:10 wt:wt (enzyme:protein).

10

The frozen aliquots are at a volume of  10  $\mu\text{L}$  containing  20  $\mu\text{g}$  of trypsin (concentration of  $[M] 2 \mu\text{g}/\mu\text{L}$  1).

11

Add  250  $\mu\text{L}$  of digestion buffer to the aliquot. This gives a total volume of  260  $\mu\text{L}$ , enough for 2 S-Trap digestions.



12

Add  125  $\mu\text{L}$  of digestion buffer containing protease into the top of the wells.



13

Place cover over the stacked plates.

14

Incubate in the thermomixer for  01:15:00 at  47  $^{\circ}\text{C}$  for trypsin.

1h 15m

Note

**Some dripping may occur during incubation; this is not of concern. REMEMBER - DO NOT SHAKE.**



## Note



- 15 **OPTIONALSTEP:** If you wish, you may also set up this digestion step **overnight**, with no impact on the S-trap process.
- To do this, set the Thermomixer to  $37^{\circ}\text{C}$  and incubate overnight, again with no shaking.



## Elute peptides:

20m

- 16 Add  $80\ \mu\text{L}$  of digestion buffer to all wells of the S-Trap digestion plate.
- 17 Centrifuge the plate at  $1000\ \text{x g}$  for  $00:02:00$  or until all solution has passed through.






2m



## Note

Do not centrifuge the plate prior to addition of  $80\ \mu\text{L}$  of digestion buffer used in this first elution.



18 Add  80  $\mu\text{L}$  of 0.1% aqueous formic acid to all wells of the S-Trap digestion plate and spin through at  1000 x g for  00:02:00 .

2m



19 Further elute peptides with  55  $\mu\text{L}$  of 30% aqueous acetonitrile containing 0.1% formic acid and spin through at  1000 x g for  00:02:00 .


2m



20 This elution assists in recovery of hydrophobic peptides.

#### Note


The final acetonitrile concentration will be around  5 % (v/v) .

21 Proceed to **R3 plate desalting** or store in a refrigerator  Overnight .

2m



## Note



When you are ready, please dispose of any solvent waste in the non-chlorinated waste drum.