

Aug 02, 2023

Version 1

sample_prep_serum.nan V.1

Forked from a private protocol

 In 1 collection



DOI

<https://dx.doi.org/10.17504/protocols.io.6qpvr34epvmk/v1>

NAN KB¹, Mario Uchimiya², John Glushka², Christopher Esselman², Leandro Ponce², Laura Morris², Arthur Edison²

¹Network for Advanced NMR (NAN); ²University of Georgia

Christopher Esselman: Protocol review;

Leandro Ponce: Protocol review



NAN-KB UGA

Network for Advanced NMR

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvr34epvmk/v1>

Protocol Citation: NAN KB, Mario Uchimiya, John Glushka, Christopher Esselman, Leandro Ponce, Laura Morris, Arthur Edison 2023. sample_prep_serum.nan. protocols.io <https://dx.doi.org/10.17504/protocols.io.6qpvr34epvmk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 02, 2023

Last Modified: August 02, 2023

Protocol Integer ID: 85858

Keywords: serum sample, protein precipitation method, serum, protein, sample

Funders Acknowledgements:

National Science Foundation

Grant ID: 194670

Disclaimer

This protocol is developed and maintained by Network for Advanced NMR (NAN). The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to this protocol is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with this protocol, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.



Abstract

This is a modified protocol for a protein precipitation method for plasma/serum samples. This protocol was originally proposed by:

Citation

Nagana Gowda GA, Raftery D (2014)

. Quantitating metabolites in protein precipitated serum using NMR spectroscopy. Analytical chemistry.

<https://doi.org/10.1021/ac5005103>

LINK

See also:

Citation

Beckonert O, Keun HC, Ebbels TM, Bundy J, Holmes E, Lindon JC, Nicholson JK (2007)

. Metabolic profiling, metabolomic and metabonomic procedures for NMR spectroscopy of urine, plasma, serum and tissue extracts.

Nature protocols.

Citation

Nagana Gowda GA, Raftery D (2019)

. Analysis of Plasma, Serum, and Whole Blood Metabolites Using ¹H NMR Spectroscopy. Methods in molecular biology (Clifton, N.J.).

https://doi.org/10.1007/978-1-4939-9690-2_2

LINK



Materials

1. Chemicals and reagents

- Phosphate buffer in D₂O (100 mM, pH 7.4, 1/3 mM DSS-d₆)
- Methanol

2. Equipment

- Calibrated micropipettes (100 μ L, 200 μ L, and 1000 μ L)
- Pippette tips
- 1.5-mL Eppendorf tubes
- 5-mm SampleJet NMR tubes (Bruker)
- Centrifuge
- Vortex mixer
- Speed-vac concentrator



1 Day-1/2

1.1 Thaw samples On ice or at 4 °C

1.2 Add 600 µL of 100% cold methanol to 300 µL of samples On ice

- Use 1.5-mL Eppendorf tubes
- Keep methanol cold On ice



1.3 Vortex the samples for 00:00:10

10s

1.4 Incubate the samples at -20 °C for 00:20:00

20m

1.5 Centrifuge the samples at 4 °C at 16000 rcf for 00:30:00

30m



1.6 Transfer the supernatants to new 1.5-mL Eppendorf tubes

1.7 Dry the samples in a speed-vac concentrator



Overnight (time varies depending on the sample condition)

Note

Store dried samples at -80 °C until the following step if needed

2 Day-2/2

2.1 Thaw the samples On ice or in 4 °C



2.2 Add  600 µL of the phosphate buffer to each sample



2.3 Vortex the samples at  4 °C for  00:10:00

10m

2.4 Centrifuge the samples at  4 °C for  00:00:10

10s



2.5 Transfer  580 µL to 5-mm NMR tubes

Note

No stickers/labels on the caps and tubes allowed

Citations

Nagana Gowda GA, Raftery D. Quantitating metabolites in protein precipitated serum using NMR spectroscopy.
<https://doi.org/10.1021/ac5005103>

Beckonert O, Keun HC, Ebbels TM, Bundy J, Holmes E, Lindon JC, Nicholson JK. Metabolic profiling, metabolomic and metabonomic procedures for NMR spectroscopy of urine, plasma, serum and tissue extracts.

Nagana Gowda GA, Raftery D. Analysis of Plasma, Serum, and Whole Blood Metabolites Using 1H NMR Spectroscopy.
https://doi.org/10.1007/978-1-4939-9690-2_2