

Oct 05, 2022

Sample preparation and lysis of homogenized malaise trap samples

 In 4 collections

DOI

<https://dx.doi.org/10.17504/protocols.io.dm6gpjrmjgzp/v1>

Dominik Buchner¹

¹University of Duisburg-Essen, Aquatic Ecosystem Research



Dominik Buchner

University of Duisburg-Essen, Aquatic Ecosystem Research



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gpjrmjgzp/v1>

Protocol Citation: Dominik Buchner 2022. Sample preparation and lysis of homogenized malaise trap samples. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.dm6gpjrmjgzp/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 29, 2022

Last Modified: October 05, 2022

Protocol Integer ID: 70661

Keywords: lysis of homogenized malaise trap sample, lysis before dna extraction, homogenized malaise trap sample, dna extraction, steps of sample preparation, protocol of the leeselab, sample preparation, malaise trap, leeselab, extraction, lysi, sample, dna

Abstract

This protocol describes the steps of sample preparation and lysis before DNA extraction for the Malaise trap metabarcoding protocol of the LeeseLab.

Guidelines

Follow general lab etiquette. Wear gloves to prevent contaminating the samples. Clean the workspace before starting with 80% EtOH.



Materials

Materials required:

Below all materials needed for the protocol are listed. Vendors and part numbers are listed but interchangeable depending on the supply situation.

Chemicals:

Tris ultrapure 99.9%  Tris ultrapure 99.9% **Diagonal Catalog #A1086.1000**

Sodium chloride  Sodium chloride **Fisher Scientific Catalog #10112640**

EDTA disodium salt  EDTA disodium salt **Merck MilliporeSigma (Sigma-Aldrich) Catalog #E5134-50G**

SDS ultrapure  Sodium dodecyl sulfate **Diagonal Catalog #A1112.0500**

Proteinase K  Proteinase K **7BioScience Catalog #RP100B**

Hydrochloric acid fuming 37%

 Hydrochloric acid fuming 37% **Merck MilliporeSigma (Sigma-Aldrich) Catalog #1003171011**

Sodium hydroxide  Sodium hydroxide - pellets **Fisher Scientific Catalog #S/4920/60**

Calcium chloride  Calcium chloride 94% **Carl Roth Catalog #A119.1**

Glycerol 87%  Glycerol 87 % for molecular biology **Panreac AppliChem Catalog #A3739,1000**

Labware:

2 mm dia zirconia beads  Zirconia Beads 2 mm dia **BioSpec Products Catalog #11079124zx**

2 mL screwcap tubes T  2 mL screwcap tube **Sarstedt Catalog #72.693**

Wide-bore tips  ART Wide Bore Tip 1000 uL **Thermo Fisher Scientific Catalog #2079G**

Disposable PES Bottle Top Filters

 Fisherbrand Disposable PES Bottle Top Filters **Fisher Scientific Catalog #15973307**

Stock solutions:

 1 L Tris stock solution **[M] 1 Mass Percent**  7.5

- Add  121.14 g Tris ultrapure 99.9% to a beaker
- Adjust volume to  800 mL with ddH₂O
- Adjust pH to  7.5 with HCl
- Adjust volume to  1 L
- Sterilize by filtering and store at  Room temperature



🧪 1 L NaCl stock solution [M] 5 Mass Percent

- Add 🧪 292.2 g Sodium chloride to a beaker
- Adjust volume to 🧪 1 L with ddH₂O
- Sterilize by filtering and store at 🌡 Room temperature

🧪 1 L EDTA stock solution [M] 0.5 Mass Percent 📏 pH 8

- Add 🧪 186.12 g EDTA disodium salt to a beaker
- Adjust volume to 🧪 1 L with ddH₂O
- Adjust pH to 📏 pH 8 with sodium hydroxide
- Sterilize by filtering and store at 🌡 Room temperature

🧪 1 L SDS stock solution [M] 10 Mass Percent

- Add 🧪 100 g SDS ultrapure to a beaker
- Adjust volume to 🧪 1 L with ddH₂O
- Sterilize by filtering and store at 🌡 Room temperature

🧪 1 L Proteinase K storage buffer ([M] 50 millimolar (mM) Tris , [M] 3 millimolar (mM) CaCl₂ [M] 50 % (v/v) glycerol) 📏 pH 7.8

- Add 🧪 50 mL of [M] 1 Mass Percent Tris stock solution 📏 pH 7.5
- Add 🧪 333 mg calcium chloride
- Add 🧪 500 mL of glycerol 87%
- Adjust volume to 🧪 900 mL with ddH₂O
- Adjust pH to 📏 pH 7.8 with sodium hydroxide
- Adjust volume to 🧪 1 L with ddH₂O
- Sterilize by filtering and store at 🌡 Room temperature

Working solutions:

🧪 1 L **TNES buffer** ([M] 50 millimolar (mM) Tris , [M] 400 millimolar (mM) NaCl , [M] 20 millimolar (mM) EDTA , [M] 0.5 Mass / % volume SDS , 📏 pH 7.5)

- Add 🧪 50 mL of [M] 1 Mass Percent Tris stock solution 📏 pH 7.5
- Add 🧪 80 mL of [M] 5 Mass Percent NaCl stock solution



- Add  40 mL of  0.5 Mass Percent EDTA stock solution  8
- Add  50 mL of  10 Mass / % volume SDS stock solution
- Adjust volume to  800 mL with ddH₂O
- Adjust pH to  7.5 with HCl
- Sterilize by filtering and store at  Room temperature

 200 mL Proteinase K working solution ( 10 mg/mL Proteinase K)

- Dissolve  2 g of Proteinase K in  200 mL Proteinase K storage buffer
- Store at  -20 °C

Safety warnings

 Reagents are potentially damaging to the environment. Dispose waste responsibly.

Before start

Make sure all buffers are prepared before starting.



- 1 For each sample prepare a screwcap tube pre-filled with a few 2 mm zirconia beads.
- 2 Shake the sample well.
Transfer  800 μL of the small size fraction and  200 μL of the large size fraction to a 2 mL screwcap tube. It might be beneficial to use wide-bore tips or cut off the tip when using regular pipette tips.

- 3  11.000 x g, 00:03:00

3m

- 4 Remove as much ethanol as possible with a  1000 μL pipette.

- 5 Add  900 μL of **TNES buffer** and  100 μL of **Proteinase K working solution**.
Vortex shortly.

Note

Depending on the amount of samples this can be prepared as a mastermix. We usually prepare TNES + Proteinase K in batches for 24 samples. Proteinase K tends to self-digest if the time for samples preparation takes too long.

- 6 Bead-beat for  00:02:00 at  2400 rpm

2m

- 7 Incubate  1400 rpm, 56°C, 00:20:00

- 8 Store at  -20 °C until DNA extraction.