

Feb 19, 2026

Version 2

Sample Preparation for Automated Extraction in Periphyton Metabolomic Analysis V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.81wgbwjoogpk/v2>



Lin Zi¹, Mélissa EON¹, Nicolas Creusot¹

¹INRAE EABX

MetaboHUB-Bordeaux



Lin Zi

INRAE

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.81wgbwjoogpk/v2>

Protocol Citation: Lin Zi, Mélissa EON, Nicolas Creusot 2026. Sample Preparation for Automated Extraction in Periphyton Metabolomic Analysis. [protocols.io](https://dx.doi.org/10.17504/protocols.io.81wgbwjoogpk/v2) <https://dx.doi.org/10.17504/protocols.io.81wgbwjoogpk/v2> Version created by [Lin Zi](#)

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 17, 2026

Last Modified: February 19, 2026

Protocol Integer ID: 243444

Keywords: chemical risk assessment, periphyton, metabolomics, automated extraction, robotic extraction in periphyton metabolomic, periphyton metabolomic analysis this protocol, periphyton metabolomic analysis, periphyton metabolomic, automated extraction, periphyton, sample preparation process, automated robotic extraction, sample preparation, using micronic tube, procedure, sample preparation for automated extraction

Funders Acknowledgements:

PARC: Partnership for the Assessment of Risks from Chemicals

Grant ID: 101057014

Bordeaux Metabolome Platform

Grant ID: <https://doi.org/10.15454/1.5572412770331912E12>

MetaboHUB

Grant ID: ANR-11-INBS-0010

MetaboHUB

Grant ID: ANR-21-ESRE-0035

Abstract

This protocol details the sample preparation process using micronic tubes for automated robotic extraction in periphyton metabolomics. The procedure utilizes periphyton previously established on 9-mm glass discs and exposed to chemical stressors in a 48-well plate format to evaluate ecotoxicological risk.

Guidelines

- Since the procedure of measuring dry biomass without caps, it's better to do the measurement fast and close the lid of the micronic plate whenever it's possible to avoid contamination of the samples.
- Try to limit freeze and de-freeze samples before the samples are dry.

Materials



micronic 96-well plate



micronic 1.1 mL tubes



micronic screw caps



metal bead




GENEVAC EZ2

Safety warnings

- ! Safety: Wear a lab coat, gloves, and safety glasses; utilize the fume hood and follow chemical safety rules. Follow safety rules for using liquid nitrogen.
- Waste: Recycle all waste according to quality document EABX_GEN_706.
- Traceability: Optimize labeling and storage for all samples, reagents, and stock solutions.
- Monitoring: Perform gravimetric monitoring for all stock solutions.
- Documentation: Record the weight, volume, and chemical nature of all internal standards in the laboratory notebook.

Before start

Periphyton samples were colonized on 9-mm glass discs and subsequently exposed to chemicals within a 48-well plate format to evaluate chemical risk. The colonization procedure has been described in "Periphyton growing and sampling in situ or in mesocosm for metabolomic analysis.**protocols.io** <https://protocols.io/view/periphyton-growing-and-sampling-in-situ-or-in-meso-dmxt47nn>". The contamination procedure has been described in "Contamination procedure for chemical risk assessment in periphyton.**protocols.io** <https://protocols.io/view/contamination-procedure-for-chemical-risk-assessme-hgxyb3xpx>".



Periphyton grown on 9-mm glass discs in 48-well plate



Quenching

5m

- 1 After chemical exposure in the 48-well plate, discard  850 μL medium from each well, leaving approximately  150 μL inside the well to prevent periphyton desiccation. If plan to perform exometabolome analysis, transfer  850 μL medium from each well to 1.1 mL Micronic tube. To ensure sufficient space for capping, discard  150 μL of the transferred volume (leaving  700 μL medium).
- 2 Seal the microplate with parafilm and use a tape to stick the microplate to the holder. Then quench the microplate in liquid N_2 for  00:05:00 .
- 3 After taking out from liquid N_2 , seal the microplate with parafilm and store the microplate in freezer  -80 $^{\circ}\text{C}$ before further use.

5m

Transfer samples from microplate into micronic tubes

4h 5m

- 4 Design the micronic plate plan in excel. Each 96 micronic plate contains at least two blank tubes. One blank tube can be replaced by QC during injection. Label the plate on the right side of the rack: "Project_date, number of the plate" e.g. "Parc_12/2025, 1"



Perc. 01/12/2025, plate 1											
A	20A1	20A2	20A3	20A4	20A5	20A6	20A7	20A8	20A9	20B1	20B2
B	20B6	20B5	20B6	20B7	20B8	20B9	20C1	20C2	20C3	20C4	20C5
C	20C7	20C8	20C9	20D1	20D2	20D3	20D4	20D5	20D6	20D7	20D8
D	20E1	20E2	20E3	20E4	20E5	20E6	20E7	20E8	20F1	20F2	20F3
E	20F4	20F5	20F6	20F7	20F8	20F9	21A1	21A2	21A3	21A4	21A5
F	21A7	21A8	21A9	21B1	21B2	21B3	21B4	21B5	21B6	21B7	21B8
G	21C1	21C2	21C3	21C4	21C5	21C6	21C7	21C8	21C9	21D1	21D2
H	21D4	21D5	21D6	21D7	21D8	21E1	21E2	21E3	21E4	21E5	21E6
Perc. 01/12/2025, plate 2											
A	21F7	21F8	21F9	21F1	21F2	21F3	21F4	21F5	21F6	21F7	21F8
B	22A1	22A2	22A3	22A4	22A5	22A6	22A7	22A8	22A9	22B1	22B2
C	22B4	22B5	22B6	22B7	22B8	22C1	22C2	22C3	22C4	22C5	22C6
D	22C7	22C8	22C9	22D1	22D2	22D3	22D4	22D5	22D6	22D7	22D8
E	22E1	22E2	22E3	22E4	22E5	22E6	22E7	22E8	22F1	22F2	22F3
F	22F4	22F5	22F6	22F7	22F8	22F9	23A1	23A2	23A3	23A4	23A5
G	23A7	23A8	23A9	23B1	23B2	23B3	23B4	23B5	23B6	23B7	23B8
H	23C1	23C2	23C3	23C4	23C5	23C6	23C7	23C8	23C9	23D1	23D2
Perc. 01/12/2025, plate 3											
A	23D4	23D5	23D6	23D7	23D8	23E1	23E2	23E3	23E4	23E5	23E6
B	23F1	23F2	23F3	23F4	23F5	23F6	23F7	23F8	23F9	23G1	23G2
C	24A1	24A2	24A3	24A4	24A5	24A6	24A7	24A8	24A9	24B1	24B2
D	24B4	24B5	24B6	24B7	24B8	24C1	24C2	24C3	24C4	24C5	24C6
E	24C7	24C8	24C9	24D1	24D2	24D3	24D4	24D5	24D6	24D7	24D8
F	24E1	24E2	24E3	24E4	24E5	24E6	24E7	24E8	24E9	24F1	24F2
G	24F4	24F5	24F6	24F7	24F8	24F9	25A1	25A2	25A3	25A4	25A5
H	25A7	25A8	25A9	25B1	25B2	25B3	25B4	25B5	25B6	25B7	25B8
Perc. 01/12/2025, plate 4											
A	25C1	25C2	25C3	25C4	25C5	25C6	25C7	25C8	25C9	25D1	25D2
B	25D4	25D5	25D6	25D7	25D8	25E1	25E2	25E3	25E4	25E5	25E6
C	25E7	25E8	25E9	25F1	25F2	25F3	25F4	25F5	25F6	25F7	25F8
D	26A1	26A2	26A3	26A4	26A5	26A6	26A7	26A8	26A9	26B1	26B2
E	26B4	26B5	26B6	26B7	26B8	26C1	26C2	26C3	26C4	26C5	26C6
F	26C7	26C8	26C9	26D1	26D2	26D3	26D4	26D5	26D6	26D7	26D8
G	26E1	26E2	26E3	26E4	26E5	26E6	26E7	26E8	26E9	26F1	26F2
H	26F4	26F5	26F6	26F7	26F8	26F9	27A1	27A2	27A3	27A4	27A5
Perc. 01/12/2025, plate 5											
A	27A7	27A8	27B1	27B2	27B3	27B4	27B5	27B6	27B7	27B8	27B9
B	27C1	27C2	27C3	27C4	27C5	27C6	27C7	27C8	27C9	27D1	27D2
C	27D4	27D5	27D6	27D7	27D8	27E1	27E2	27E3	27E4	27E5	27E6
D	27E7	27E8	27F1	27F2	27F3	27F4	27F5	27F6	27F7	27F8	27F9
E	28A1	28A2	28A3	28A4	28A5	28A6	28A7	28A8	28A9	28B1	28B2
F	28B4	28B5	28B6	28B7	28B8	28C1	28C2	28C3	28C4	28C5	28C6
G	28C7	28C8	28C9	28D1	28D2	28D3	28D4	28D5	28D6	28D7	28D8
H	28E1	28E2	28E3	28E4	28E5	28E6	28E7	28E8	28F1	28F2	28F3
Perc. 01/12/2025, plate 6											
A	28F4	28F5	28F6	28F7	28F8	28F9	29A1	29A2	29A3	29A4	29A5
B	29A7	29A8	29B1	29B2	29B3	29B4	29B5	29B6	29B7	29B8	29B9
C	29C1	29C2	29C3	29C4	29C5	29C6	29C7	29C8	29C9	29D1	29D2
D	29D4	29D5	29D6	29D7	29D8	29E1	29E2	29E3	29E4	29E5	29E6
E	29E7	29E8	29F1	29F2	29F3	29F4	29F5	29F6	29F7	29F8	29F9
F	30A1	30A2	30A3	30A4	30A5	30A6	30A7	30A8	30A9	30B1	30B2
G	30B4	30B5	30B6	30B7	30B8	30C1	30C2	30C3	30C4	30C5	30C6
H	30C7	30C8	30C9	30D1	30D2	30D3	30D4	30D5	30D6	30D7	30D8
Perc. 01/12/2025, plate 7											
A	30E1	30E2	30E3	30E4	30E5	30E6	30E7	30E8	30E9	30F1	30F2
B	30F4	30F5	30F6	30F7	30F8	30F9	31A1	31A2	31A3	31A4	31A5
C	31A7	31A8	31B1	31B2	31B3	31B4	31B5	31B6	31B7	31B8	31B9
D	31C1	31C2	31C3	31C4	31C5	31C6	31C7	31C8	31C9	31D1	31D2
E	31D4	31D5	31D6	31D7	31D8	31E1	31E2	31E3	31E4	31E5	31E6
F	31F1	31F2	31F3	31F4	31F5	31F6	31F7	31F8	31F9	31G1	31G2
G	31G4	31G5	31G6	31G7	31G8	31G9	32A1	32A2	32A3	32A4	32A5
H	32A7	32A8	32B1	32B2	32B3	32B4	32B5	32B6	32B7	32B8	32B9
Perc. 01/12/2025, plate 8											
A	32C7	32C8	32D1	32D2	32D3	32D4	32D5	32D6	32D7	32D8	32D9
B	32E1	32E2	32E3	32E4	32E5	32E6	32E7	32E8	32E9	32F1	32F2
C	32F4	32F5	32F6	32F7	32F8	32F9	33A1	33A2	33A3	33A4	33A5
D	33A7	33A8	33B1	33B2	33B3	33B4	33B5	33B6	33B7	33B8	33B9
E	33C1	33C2	33C3	33C4	33C5	33C6	33C7	33C8	33C9	33D1	33D2
F	33D4	33D5	33D6	33D7	33D8	33D9	33E1	33E2	33E3	33E4	33E5
G	33E7	33E8	33E9	33F1	33F2	33F3	33F4	33F5	33F6	33F7	33F8
H	33F9	33G1	33G2	33G3	33G4	33G5	33G6	33G7	33G8	33G9	33H1
Perc. 01/12/2025, plate 9											
A	33H4	33H5	33H6	33H7	33H8	33H9	34A1	34A2	34A3	34A4	34A5
B	34A7	34A8	34B1	34B2	34B3	34B4	34B5	34B6	34B7	34B8	34B9
C	34C1	34C2	34C3	34C4	34C5	34C6	34C7	34C8	34C9	34D1	34D2
D	34D4	34D5	34D6	34D7	34D8	34D9	34E1	34E2	34E3	34E4	34E5
E	34E7	34E8	34E9	34F1	34F2	34F3	34F4	34F5	34F6	34F7	34F8
F	34F9	34G1	34G2	34G3	34G4	34G5	34G6	34G7	34G8	34G9	34H1
G	34H4	34H5	34H6	34H7	34H8	34H9	34I1	34I2	34I3	34I4	34I5
H	34I7	34I8	34I9	34J1	34J2	34J3	34J4	34J5	34J6	34J7	34J8
Perc. 01/12/2025, plate 10											
A	34J9	34K1	34K2	34K3	34K4	34K5	34K6	34K7	34K8	34K9	34L1
B	34L4	34L5	34L6	34L7	34L8	34L9	34M1	34M2	34M3	34M4	34M5
C	34M7	34M8	34M9	34N1	34N2	34N3	34N4	34N5	34N6	34N7	34N8
D	34N9	34O1	34O2	34O3	34O4	34O5	34O6	34O7	34O8	34O9	34P1
E	34P4	34P5	34P6	34P7	34P8	34P9	34Q1	34Q2	34Q3	34Q4	34Q5
F	34Q7	34Q8	34Q9	34R1	34R2	34R3	34R4	34R5	34R6	34R7	34R8
G	34R9	34S1	34S2	34S3	34S4	34S5	34S6	34S7	34S8	34S9	34T1
H	34T4	34T5	34T6	34T7	34T8	34T9	34U1	34U2	34U3	34U4	34U5

dry ground BM
ground in solvent

Exemple of micronic 96-well plate plan

- 5 Weight empty micronic 1.1 mL tubes (NBS MP32033L, without caps) and put the tubes at 96 well micronic rack. Close the lid of the rack before further use.

 **MICRONIC**
INNOVATION IN SAMPLE STORAGE

1.1 mL micronic tube

 **MICRONIC**
INNOVATION IN SAMPLE STORAGE

Micronic 96-well rack

- 6 Add  200 μL cold EtOH (100%, Grade UPLC) in each well of 48 well plate.
- 7 Place the plate in a beaker which contained ice water. Then, put the beaker in sonication for  00:05:00 . Don't place the plate directly in the sonication bath, as the power of sonication is not even which can cause well to well contamination.

5m



Ultra-sonication bath: microplate in a beaker contained ice cold water

- 8 After sonication, using pipette to wash the discs and transfer the EtOH with periphyton into pre-weighted micronic tubes.
- 9 Dry EtOH by using GENEVAC (Hydro mode), until no more liquid and the sample is completely dry. Probably it will take more than  04:00:00 . If you will using GENEVAC to dry samples in the next day after adding EtOH, put the samples in  4 °C . If you will dry the samples a few days later, put the samples in  -20 °C .
- 10 After GENEVAC, measure the weight of tubes plus dry samples (without caps). Avoid putting the plate on ice during the measurement of dry weight as ice can increase humidity then affect the dry weight measurement.
- 11 Add one metal bead into each micronic tube for sample grinding before close the cap.
- 12 Store the micronic plates in  -80 °C before metabolites extraction.

4h



Protocol references

Lin Zi, Zahrasadat Alavikakhki, Mélissa Eon, Chloé Bonnineau, Nicolas Creusot 2026. Periphyton growing and sampling in situ or in mesocosm for metabolomic analysis. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.4r3l2qjd4l1y/v2>

Lin Zi, Mélissa EON, Zahrasadat Alavikakhki, Chloé Bonnineau, Nicolas Creusot 2026. Contamination procedure for chemical risk assessment in periphyton. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bp2l6zjpdgge/v2>