

Apr 21, 2025

Version 1

Sample Preparation for Quantitative Whole Cell Proteome Analysis by Mass Spectrometry V.1

 [Nature Cell Biology](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.n92ld5918v5b/v1>

Markus Hartl¹

¹Max Perutz Labs – MS Facility



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ld5918v5b/v1>

External link: <https://doi.org/10.1038/s41556-025-01712-y>

Protocol Citation: Markus Hartl 2025. Sample Preparation for Quantitative Whole Cell Proteome Analysis by Mass Spectrometry. protocols.io <https://dx.doi.org/10.17504/protocols.io.n92ld5918v5b/v1>

**Manuscript citation:**

Adriaenssens, E., Schaar, S., Cook, A.S.I. *et al.* Reconstitution of BNIP3/NIX-mitophagy initiation reveals hierarchical flexibility of the autophagy machinery. *Nat Cell Biol* **27**, 1272–1287 (2025). <https://doi.org/10.1038/s41556-025-01712-y>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: April 09, 2025

Last Modified: April 21, 2025

Protocol Integer ID: 126559

Keywords: ASAPCRN, sample preparation for quantitative whole cell proteome analysis, quantitative whole cell proteome analysis, mass spectrometry, sample preparation, whole cell

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-000350

Marie Skłodowska-Curie MSCA Postdoctoral fellowship

Grant ID: 101062916

Abstract

-

Guidelines

GENERAL NOTE: The protocols are copied from the internal electronic lab notebook and contain specific notes, for example, sample-specific volumes and incubation times. They should be omitted from a general protocol and rather end concentrations and total amounts should be provided where necessary.



Materials

Buffer: 2% SDC [M] 100 millimolar (mM) Tris/HCl pH=  8.8 ( 20 mg SDC in  1 mL)

- PBS
- sodium deoxycholate
- Tris/HCl pH 8.8
- benzonase
- micro BCA assay (Thermo Fisher)
- iodoacetamide
- DTT
- ammonium bicarbonate
- LysC (mass spectrometry grade, FUJIFILM Wako chemicals)
- trypsin (Trypsin Gold, Promega)
- trifluoroacetic acid



Cell Lysis

33m

- 1 Pre-heat 2% sodium deoxycholate (SDC) [1M] 100 millimolar (mM) Tris/HCl pH 8.8 at 95 °C .
- 2 Add 50 µL hot 2% SDC [1M] 100 millimolar (mM) Tris to cell pellet (add more if not soluble).
- 3 Incubate samples for 00:10:00 at 95 °C . 10m
- 4 Resuspend thoroughly with pipet, put to 95 °C for 00:03:00 . 3m
- 5 Cool down, add 1 µL benzonase and incubate On ice for 10-30 min.
- 6 Sonicate 5 cycles 30/30 sec level H in Bioruptor.
- 7 Let sit for 00:10:00 for extra benzonase treatment. 10m
- 8 Centrifuge at 15000 x g, 10°C, 00:10:00 to pellet cell debris. 10m
- 9 Transfer supernatant to fresh 0.6mL tube.
- 10 Determine protein concentration with micro BCA assay.

Note

NOTE: we use the Micro BCA kit according to manufacturer's instructions
<https://www.thermofisher.com/order/catalog/product/23235>.



In Solution Digest

4h 39m

- 11 Protein concentration: ~ 5 $\mu\text{g}/\mu\text{L}$.
- 12 Transfer 50 μg protein to 0.2mL tube, fill up to 15 μL .
- 13 Reduction: 10 millimolar (mM) DTT (250 millimolar (mM) stock) 00:30:00 50 $^{\circ}\text{C}$. Here we used 0.6 μL of DTT. 30m
- 14 Alkylation: 20 millimolar (mM) iodoacetamide (IAA) (500 millimolar (mM) stock = 0.092 g/mL) 00:30:00 Room temperature in the dark. Here we used 0.6 μL of IAA. 30m
- 15 Quench the remaining iodoacetamide (IAA) with half amount of DTT used in step 14 and incubate for 00:10:00 at Room temperature . Here we used 0.3 μL DTT to quench. 10m
- 16 Dilute to 1% final SDC concentration with 100 millimolar (mM) Tris pH 8.5 . Here we used 15 μL of SDC.
- 17 Add 1:100 LysC (100 ng/ μL aliquot) – incubate at Room temperature for up to 03:00:00 . Here we used 0.5 μL from 1 $\mu\text{g}/\mu\text{L}$. 3h
- 18 Add 1:50 Trypsin (100 ng/ μL aliquot) – incubate at 37 $^{\circ}\text{C}$ Overnight . Here we used 1 μL from 1 $\mu\text{g}/\mu\text{L}$. Here, the total volume was 33 μL at this point of the protocol.
- 19 Add 10% TFA to reach final 1% TFA, vortex, spin down at 14000 rpm, 00:10:00 . For us here, we added 3.3 μL to reach 10% TFA. 10m
- 20 Transfer supernatant to new tube.
- 21 Add 10% TFA to reach final 2% TFA, vortex, spin down at 14000 rpm, 00:10:00 . Here we added 3.3 μL , to reach a total volume of 39.6 μL . 10m



22 Load directly onto prepared C18 STAGETip or transfer to new tube and load 20%. For us here, 20% equals  10 μg $\rightarrow$  8 μL .

23 Desalt: Load the peptides on a C18 STAGETip.

Note

3-5 plugs depending on protein amount, 3 plugs for max.  20 μg .

24 Equilibrate:



24.1  100 μL 100% MeOH, spin  2000 rpm, 00:03:00 .

3m

24.2  100 μL 80% ACN, 0.1% TFA, spin  2000 rpm, 00:03:00 .

3m

24.3  100 μL 0.1% TFA, spin  2000 rpm, 00:03:00 .

3m

25 Load sample: spin  1700 rpm .



26 Wash:  100 μL 0.1% TFA $\rightarrow$ keep flow through (FT) and wash together.



27 Elute: 2x  30 μL 60% ACN, 0.1% TFA in PCR tube.

28 Speedvac to remove ACN and take up sample in  20 μL 0.1% TFA 2% ACN.