



Apr 23, 2025

Version 2

Sample Preparation for Quantitative Whole Cell Proteome Analysis by Mass Spectrometry V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.n92ld5918v5b/v2>

Elias Adriaenssens¹

¹Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ld5918v5b/v2>

Protocol Citation: Elias Adriaenssens 2025. Sample Preparation for Quantitative Whole Cell Proteome Analysis by Mass Spectrometry. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n92ld5918v5b/v2> Version created by [Elias Adriaenssens](#)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working



Created: April 23, 2025

Last Modified: April 23, 2025

Protocol Integer ID: 141123

Keywords: ASAPCRN, sample preparation for quantitative whole cell proteome analysis, quantitative whole cell proteome analysis, mass spectrometry, mass spectrometry this protocol, sample preparation, samplespecific volume, experiment specific note, internal electronic lab notebook, sample, preparation, automatic generation process, cell

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-000350

Marie Skłodowska-Curie MSCA Postdoctoral fellowship

Grant ID: 101062916

Abstract

This protocol was generated from an internal electronic lab notebook and contains experiment specific notes, for example, samplespecific volumes and incubation times. Due to the automatic generation process it might also contain errors and should be interpreted cautiously. An updated protocols will still be uploaded. This protocol was generated from an internal electronic lab notebook and contains experiment specific notes, for example, sample-specific volumes and incubation times. Due to the automatic generation process it might also contain errors and should be interpreted cautiously. An updated protocols will still be uploaded.

Guidelines

GENERAL NOTE: This protocol was generated from an internal electronic lab notebook and contains experiment specific notes, for example, sample-specific volumes and incubation times. Due to the automatic generation process it might also contain errors and should be interpreted cautiously. An updated protocols will still be uploaded.



Materials

Buffer: 2% SDC [M] 100 millimolar (mM) Tris/HCl pH=  8.8 ( 20 mg SDC in  1 mL)

- PBS
- sodium deoxycholate
- Tris/HCl pH 8.8
- benzonase
- micro BCA assay (Thermo Fisher)
- iodoacetamide
- DTT
- ammonium bicarbonate
- LysC (mass spectrometry grade, FUJIFILM Wako chemicals)
- trypsin (Trypsin Gold, Promega)
- trifluoroacetic acid

Safety warnings

 This protocol was generated from an internal electronic lab notebook and contains experiment specific notes, for example, sample-specific volumes and incubation times. Due to the automatic generation process it might also contain errors and should be interpreted cautiously. An updated protocols will still be uploaded. This protocol was generated from an internal electronic lab notebook and contains experiment specific notes, for example, sample-specific volumes and incubation times. Due to the automatic generation process it might also contain errors and should be interpreted cautiously. An updated protocols will still be uploaded.



Cell Lysis

33m

- 1 Pre-heat 2% sodium deoxycholate (SDC) [M] 100 millimolar (mM) Tris/HCl pH 8.8 at 95 °C .
- 2 Add 50 µL hot 2% SDC [M] 100 millimolar (mM) Tris to cell pellet (add more if not soluble).
- 3 Incubate samples for 00:10:00 at 95 °C . 10m
- 4 Resuspend thoroughly with pipet, put to 95 °C for 00:03:00 . 3m
- 5 Cool down, add 1 µL benzonase and incubate On ice for 10-30 min.
- 6 Sonicate 5 cycles 30/30 sec level H in Bioruptor.
- 7 Let sit for 00:10:00 for extra benzonase treatment. 10m
- 8 Centrifuge at 15000 x g, 10°C, 00:10:00 to pellet cell debris. 10m
- 9 Transfer supernatant to fresh 0.6mL tube.
- 10 Determine protein concentration with micro BCA assay.

Note

NOTE: we use the Micro BCA kit according to manufacturer's instructions
<https://www.thermofisher.com/order/catalog/product/23235>.



In Solution Digest

4h 39m

- 11 Protein concentration: ~ 5 μL .
- 12 Transfer 50 μg protein to 0.2mL tube, fill up to 15 μL .
- 13 Reduction: 10 millimolar (mM) DTT (250 millimolar (mM) stock) 00:30:00 50 $^{\circ}\text{C}$. Here we used 0.6 μL of DTT. 30m
- 14 Alkylation: 20 millimolar (mM) iodoacetamide (IAA) (500 millimolar (mM) stock = 0.092 μL) 00:30:00 Room temperature in the dark. Here we used 0.6 μL of IAA. 30m
- 15 Quench the remaining iodoacetamide (IAA) with half amount of DTT used in step 14 and incubate for 00:10:00 at Room temperature . Here we used 0.3 μL DTT to quench. 10m
- 16 Dilute to 1% final SDC concentration with 100 millimolar (mM) Tris pH 8.5 . Here we used 15 μL of SDC.
- 17 Add 1:100 LysC (100 μL aliquot) – incubate at Room temperature for up to 03:00:00 . Here we used 0.5 μL from 1 μL . 3h
- 18 Add 1:50 Trypsin (100 μL aliquot) – incubate at 37 $^{\circ}\text{C}$ Overnight . Here we used 1 μL from 1 μL . Here, the total volume was 33 μL at this point of the protocol.
- 19 Add 10% TFA to reach final 1% TFA, vortex, spin down at 14000 rpm, 00:10:00 . For us here, we added 3.3 μL to reach 10% TFA. 10m
- 20 Transfer supernatant to new tube.
- 21 Add 10% TFA to reach final 2% TFA, vortex, spin down at 14000 rpm, 00:10:00 . Here we added 3.3 μL , to reach a total volume of 39.6 μL . 10m



22 Load directly onto prepared C18 STAGETip or transfer to new tube and load 20%. For us here, 20% equals  10 μg $\rightarrow$  8 μL .

23 Desalt: Load the peptides on a C18 STAGETip.

Note

3-5 plugs depending on protein amount, 3 plugs for max.  20 μg .

24 Equilibrate:



24.1  100 μL 100% MeOH, spin  2000 rpm, 00:03:00 .

3m

24.2  100 μL 80% ACN, 0.1% TFA, spin  2000 rpm, 00:03:00 .

3m

24.3  100 μL 0.1% TFA, spin  2000 rpm, 00:03:00 .

3m

25 Load sample: spin  1700 rpm .



26 Wash:  100 μL 0.1% TFA $\rightarrow$ keep flow through (FT) and wash together.



27 Elute: 2x  30 μL 60% ACN, 0.1% TFA in PCR tube.

28 Speedvac to remove ACN and take up sample in  20 μL 0.1% TFA 2% ACN.