

Jan 20, 2026

Sanger sequencing of barcoded beads

 [Nature Communications](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.5jyl8xpzdv2w/v1>

Koen Theunis^{1,2,3,4}

¹Laboratory of Computational Biology, VIB Center for AI & Computational Biology, Leuven, Belgium;

²VIB-KU Leuven Center for Brain & Disease Research, Leuven, Belgium;

³Department of Human Genetics, KU Leuven, Leuven, Belgium;

⁴Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, United States

Aertslab



Koen Theunis

Laboratory of Computational Biology, VIB Center for AI & Com...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl8xpzdv2w/v1>



External link: <https://doi.org/10.1038/s41467-026-68742-4>

Protocol Citation: Koen Theunis 2026. Sanger sequencing of barcoded beads. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.5jyl8xpzdv2w/v1>

Manuscript citation:

Dickmanken H, Wojno M, Mahieu L, Theunis K, Ekşi EC, Christiaens V, Kempynck N, Rop FVD, Roels N, Spanier KI, Vandepoel R, Hulselmans G, Poovathingal S, Aerts S (2026) Evaluating single-cell ATAC-seq atlasing technologies using sequence-to-function modeling. Nature Communications 17(). doi: [10.1038/s41467-026-68742-4](https://doi.org/10.1038/s41467-026-68742-4)

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 19, 2026

Last Modified: January 20, 2026

Protocol Integer ID: 238928

Keywords: ASAPCRN, hydrogel bead barcode qc, barcoded hydrogel bead, barcoded bead, hydrogel bead, beads this protocol

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000430

Aligning Science Across Parkinson's

Grant ID: ASAP-025179

Abstract

This protocol is used to test the purity of barcoded hydrogel beads like: <https://www.protocols.io/view/hydrop-v2-bead-generation-amp-ligation-barcoding-8epv5n11dv1b/v2>

Materials

Oligo's for HyDropATAC v2 Sanger QC

PCR	Name	Sequence (5'>3')	Purification
first	ATAC-QC oligo	GTCTCGTGGGCTCGGCTG TCTCTTATACACATCTGAC GCTGCCGACGA	Standard desalted
second	P7 indexing primer	CAAGCAGAAGACGGCATA CGAGATXXXXXXXXXXGTC TCGTGGGCTCGGAGATGT G	Standard desalted
second	P5-Partial	AATGATACGGCGACCACCGA	Standard desalted

Oligo's for HyDropATAC v1 Sanger QC

PCR	Name	Sequence (5'>3')	Purification
first	ATAC-QC oligo	GTGACTGGAGTTCAGACGT GTGCTCTTCCGATCTGCC GAGCCCACGAGAC	Standard desalted
second	P7 indexing primer	CAAGCAGAAGACGGCATA CGAGATXXXXXXXXXXGTG ACTGGAGTTCAGACGTGTG	Standard desalted
second	Dropseq-P5-oligo	AATGATACGGCGACCACCGAGATCTACACGCCTGTCC GCGGAAGCAGTGGTATCA ACGCAGAGT*A*C	Standard desalted

Wash finalized barcoded beads

- 1 Take finalized barcoded beads from storage (-80°C)
- 2 Add 5 μL beads to a 1.5mL microcentrifuge tube containing 500 μL of 0.04% BSA in PBS.
- 3 Spin 1 minute at 500 x g and remove supernatant but leave about 20 μL .
From this supernatant a negative control can be taken in order to investigate oligo bleeding from the beads.
- 4 Repeat these washes for a total of 3 washes

Pick individual beads

- 5 Prepare the **lid** of a petri dish (Falcon 351008) with plenty of 5 μL PBS/BSA drops as washing medium for picked beads later. Use 20-30 μL drops in the petri dish for washing and priming of the 75 μm capillaries (MXL3-STR and MXL3-75, CooperSurgical). Adjust the Stripper[®] pipette to 0.5 μL and prime with PBS/BSA.
- 6 Dilute beads by pipetting in 10-20 μL drops of PBS/BSA on a petri dish **lid** under a stereo microscope until individual beads are distinguishable from each other.
- 7 Now pick beads with Stripper[®] pipette and 75 μm capillaries (MXL3-STR and MXL3-75, CooperSurgical) and transfer them over three wash drops in order to verify that 1 bead was picked. The tip can be rinsed and primed from large drops of PBS/BSA in between the bead transfers.
- 8 Transfer a single bead to a PCR tube containing 1.5 μL PBS/BSA.
- 9 Wash the capillary 3 times in a large volume of PBS/BSA. Now additional beads can be picked.
- 10 A blank can be taken from an empty droplet that was used for washing.
- 11 At the end of the experiment, a few multi-droplet controls ($n = 2$, $n = 5$, $n = 10$) can be picked. Don't do this before picking single beads because beads can stay in the dead



volume of the capillary.

First PCR amplification of the barcoded bead oligo

- 12 Add 18 μL PCR master mix 1 to the PCR tubes containing beads and/or blanks.
Final concentrations: 1x KAPA HiFi hot start ready mix (Roche Cat. No. KK2602), 30mM DTT, 1 μM ATAC-QC oligo.

13

	A	B	C
	72 °C	3 m	
	98 °C	30 s	
	98 °C	10 s	
	59 °C	30 s	22 cycles
	72 °C	30 s	
	4-15 °C	hold	

- 14 Perform the purification according to the manufacturers recommendation with 1.5X of Ampure beads (add 30 μL). Elute in 10 μL EB buffer and transfer 8 μL to the next PCR reaction.

Second PCR reaction

- 15 Add 12 μL PCR master mix 2 to 8 μL sample.
Final concentrations: 1X KAPA HiFi hot start ready mix, 1 μM P7 indexing primer, 1 μM P5-Partial primer.

16

	A	B	C
	95 °C	3 m	
	98 °C	20 s	
	67 °C	30 s	22 cycles
	72 °C	20 s	
	72 °C	1 m	



A	B	C
4-15 °C	hold	

- 17 Perform the purification according to the manufacturers recommendation with 1.5X of Ampure beads (add 30 μ L). Elute beads in 17 μ L EB buffer and transfer 15 μ L to a fresh tube.

Sanger sequencing

- 18 Send for Sanger sequencing (eg LGC Genomics) using a P7 primer (for HyDropV2) or P5 primer (for HyDropV1).

Expected result

19

