

Jan 18, 2019

## Sanger Sequencing

DOI

[dx.doi.org/10.17504/protocols.io.pxddpi6](https://dx.doi.org/10.17504/protocols.io.pxddpi6)

Livia Sacchetto, Izabela M Rezende<sup>1</sup>

<sup>1</sup>Stanford University School of Medicine

MRCA



Livia Sacchetto

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.pxddpi6](https://dx.doi.org/10.17504/protocols.io.pxddpi6)

**Protocol Citation:** Livia Sacchetto, Izabela M Rezende 2019. Sanger Sequencing. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.pxddpi6>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** May 05, 2018

**Last Modified:** January 18, 2019

**Protocol Integer ID:** 11973

## Abstract


For sequence by dideoxy-method using the ABI3130 platform (Applied Biosystems, USA) and specific primers.




## Materials

### MATERIALS


 Microcentrifuge Tubes

 Filter Tips


 Nuclease-Free Water, 1000ml **Promega Catalog #P1199**


 Primer

 Micropipettors

 Sequencing plate

 PCR purified template

 Sequencing buffer

 Sanger sequencer

 adhesive plate seal



- 1 In a 0.2 mL microtube, place the items below:

PCR product	20 to 30 ng of DNA per 100 bp
Primer Forward or Reverse	10 pmol
BigDye™ Terminator v3.1 5X Sequencing Buffer (Applied Biosystems)	2 µL
Nuclease-free water	to a final volume total of 10,0 µL

- 2 Pipette each reaction from the tube into the corresponding well of the sequencing plate (10 µL).
- 3 Seal the plate; centrifuge briefly to collect the contents of the wells at the bottom. Protect from extended light exposure or elevated temperatures before cycling.
- 4 Use Sanger sequencing with the Applied Biosystems® 3730 DNA Analyser, according to the manufacturer's instructions.