

Sep 29, 2023

🌍 Sanger Tree of Life HMW DNA Extraction: Manual Nucleated Blood Nanobind®

📁 In 3 collections

DOI

dx.doi.org/10.17504/protocols.io.5jyl8p2w8g2w/v1

Amy Denton¹, Graeme Oatley¹, Pacific Biosciences², Caroline Howard¹

¹Tree of Life, Wellcome Sanger Institute, Hinxton, Cambridgeshire, CB10 1SA; ²Pacific Biosciences

Tree of Life at the Wellc...

Earth BioGenome Project



Tree of Life Genome Note Editor

Wellcome Sanger Institute



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl8p2w8g2w/v1>

Protocol Citation: Amy Denton, Graeme Oatley, Pacific Biosciences, Caroline Howard 2023. Sanger Tree of Life HMW DNA Extraction: Manual Nucleated Blood Nanobind®. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5jyl8p2w8g2w/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 04, 2023

Last Modified: September 29, 2023

Protocol Integer ID: 87304

Keywords: HMW DNA extraction, nucleated blood extraction, Nanobind, manual DNA extraction, reference genome, long read sequencing, hmw dna from nucleated blood, sanger tree of life hmw dna extraction, extracting hmw dna, life hmw dna extraction, manual extraction of hmw dna, hmw dna fragmentation, hmw dna, quantity hmw dna, dna extraction, nucleated blood sample, nucleated blood within the chordate group, nucleated blood, acronyms hmw, li pacbio protocol, procedure from pacbio, dna, read sequencing, molecular weight li, high molecular weight li, manual extraction, extraction, sanger tree

Funders Acknowledgements:

Wellcome Trust

Grant ID: 218328

Wellcome Trust

Grant ID: 206194

Gordon and Betty Moore Foundation

Grant ID: GBMF8897

Abstract

This protocol describes the manual extraction of HMW DNA from nucleated blood samples intended for long-read sequencing using the Nanobind® tissue kit and following the 'Extracting HMW DNA from nucleated blood using Nanobind® kits' procedure from PacBio. This process is effective for any species with nucleated blood within the Chordate group, covered by the Tree of Life Programme, however difficulties can arise with samples where preservation has not been optimal. The output of this protocol is high quality and quantity HMW DNA, which can be directed towards the HMW DNA Fragmentation: Diagenode Megaruptor®3 for LI PacBio protocol.

Acronyms:

HMW: high-molecular weight

LI: low input



Guidelines

- This protocol uses the 'Extracting HMW DNA from nucleated blood using Nanobind® kits' procedure, with the inclusion of the standard sample inputs and the standard elution volumes used by Sanger Tree of Life.
- This protocol is suitable for nucleated blood from birds, amphibians or fish that has been either flash frozen or stored in ethanol.
- Before starting the protocol, keep samples on dry ice to maintain temperature and prevent nucleic acid degradation.
- An experienced operator can expect to comfortably process 4 samples with a start to finish period of 1 hour. This estimation excludes overnight incubation at room temperature and subsequent QC checks.

Additional Notes:


- FluidX tubes are used throughout the Tree of Life programme in order to track samples, therefore rather than the microcentrifuge tubes which have been mentioned in this protocol for DNA storage, all routine DNA extracts are stored in FluidX tubes.

Materials

- 1.5 mL Protein LoBind microcentrifuge tubes (Eppendorf Cat. no. 0030 108.116)
- Nanobind® tissue kit (Cat. no. 102-302-100)
- 100% absolute ethanol
- 100% absolute isopropanol
- 1x phosphate-buffered saline (PBS)
- 15 mL or 50 mL centrifuge tubes

Equipment:

- Pipettes for 0.5 to 1000 µL and filtered tips
- Wide-bore pipette tips (200 µL, filtered if available)
- DynaMag™-2 magnetic rack (Cat. no. 12321D) or similar
- Eppendorf ThermoMixer C (Cat. no. 5382000031)
- Eppendorf SmartBlock 2.0 ml (Cat. no. 5362000035)
- HulaMixer Sample Mixer (Cat. no. 15920D)
- Mini centrifuge (Cat. no. SS-6050)
- Eppendorf Centrifuge 5425/5425 R (Cat. no. 5405000263)
- Vortexer (Vortex Genie™ 2 SI-0266)
- Timer

Protocol PDF:  Sanger Tree of Life HMW DNA Extr... 79KB

Troubleshooting



Safety warnings

- The operator must wear a lab coat, powder-free nitrile gloves and safety specs to perform the laboratory procedures in this protocol. Cotton glove liners are strongly recommended when handling the samples on dry ice.
- Waste needs to be collected in a suitable container (e.g. plastic screw-top jar or Biobin) and disposed of in accordance with local regulations.
- Liquid waste needs to be collected in a suitable container (e.g. glass screw-top jar) and disposed of in accordance with local regulations.

Before start

- Add 100% ethanol to the Buffers CW1 and CW2 as per manufacturer's instructions.
- Set a heat block to 55°C.

Laboratory protocol

- 1 Add 10–20 μL of nucleated blood to a 1.5 mL Protein LoBind tube.
- 2 Add 180–190 μL of 1 X PBS for a total volume of 200 μL .
- 3 Add 20 μL of Proteinase K.
- 4 Add 20 μL of RNase A.
- 5 Pulse vortex the blood sample 10 times for 1 second at the maximum setting, and then spin the blood sample tube on a mini-centrifuge for 2 seconds.
- 6 Incubate the blood sample at room temperature for 3 minutes.
- 7 Add 200 μL of Buffer BL3 and pulse vortex 10 times for 1 second at the maximum setting, then spin the blood sample on a mini-centrifuge for 2 seconds.
- 8 Incubate the blood sample on a heat block at 55 °C at 900 rpm for 10 minutes.
- 9 Pulse vortex the blood sample 3 times for 1 second at the maximum setting, and then spin the blood sample on a mini-centrifuge for 2 seconds.
- 10 Add a Nanobind disk to the blood sample tube.
- 11 Add 350 μL of 100% isopropanol to the blood sample and Nanobind disk.
- 12 Inversion mix the blood sample 5 times, then place the sample on a rotating mixer at 9 rpm for 15 minutes at room temperature.

- 13 Place the blood sample tube on the magnetic rack to allow for disk (and therefore HMW DNA) capture.
- 14 Discard the supernatant with a pipette, trying not to disturb the Nanobind disk or pipette the DNA. To minimise carryover contamination, remove any excess liquid in the tube cap.
- 15 Add 700 μL of Buffer CW1, then remove the blood sample tube from the magnetic rack, inversion mix 4 times, before placing the blood sample tube back on the magnetic rack and discarding the supernatant. To ensure that the Nanobind disk is captured towards the top of the tube, first place the tube in the rack separated from the magnet, then flip the rack so that the lid of the tube is on the benchtop. Add the magnet into the rack, then flip the now-assembled magnetic rack so that the sample is the right side up.
- 16 Add 500 μL of Buffer CW2, then remove the blood sample tube from the magnetic rack, inversion mix 4 times, before placing the blood sample tube back on the magnetic rack and discarding the supernatant. To ensure that the Nanobind disk is captured towards the top of the tube, first place the tube in the rack separated from the magnet, then flip the rack so that the lid of the tube is on the benchtop. Add the magnet into the rack, then flip the now-assembled magnetic rack so that the sample is the right side up.
- 17 Repeat step 16 for a total of two Buffer CW2 washes.
- 18 Spin the blood sample tube on the mini-centrifuge for 2 seconds, then place the blood sample tube back on the magnetic rack.
- 19 Remove any residual buffer that may be present in the tube; if the Nanobind disk is blocking access to the bottom of the tube to remove this liquid, it can be gently pushed towards the magnet with the tip of the pipette.
- 20 Repeat steps 17 and 18 to ensure that there is no residual buffer in the tube.
- 21 Remove the sample from the magnetic rack, then add 200 μL of Buffer EB directly onto the Nanobind disk and incubate the sample for 10 minutes at room temperature.
- 22 Transfer the eluate containing the HMW DNA to a new 1.5 mL microcentrifuge tube with a standard p200 pipette and 200 μL pipette tip. Repeat until all the eluate is transferred.
- 23 Spin the sample tube containing the Nanobind disk on a centrifuge at $10,000 \times g$ for 15 seconds, then combine any additional liquid that comes off the disk with the previously transferred eluate.



- 24 Repeat step 23 for a total of 2 spins on the centrifuge if any visible DNA remains on the disk.
- 25 Pipette mix the eluate 10 times with a standard p200 pipette and 200 μ L pipette tip to homogenise the DNA within the sample and disrupt any viscous regions.
- 26 Allow the eluate to rest overnight at room temperature to allow the DNA to solubilise.
- 27 Following the overnight rest, pipette mix 10 times with a standard p200 pipette and 200 μ L pipette tip, then perform the required QC.
- 28 Store the DNA at 4 °C.

Protocol references

[Procedure-checklist-Extracting-HMW-DNA-from-nucleated-blood-using-Nanobind-kits.pdf \(pacb.com\)](#)