

Aug 17, 2023

sCD40L ELISA Assay

DOI

dx.doi.org/10.17504/protocols.io.5jyl8pp2dg2w/v1

Fang Zhou^{1,2}

¹Massachusetts General Hospital; ²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network



Daniel Choo

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl8pp2dg2w/v1>

Protocol Citation: Fang Zhou 2023. sCD40L ELISA Assay. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.5jyl8pp2dg2w/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 17, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 86576

Keywords: ASAPCRN, protocol for the scd40l elisa kit, scd40l elisa kit, concentration of scd40l, scd40l, csf samples in mice, invitrogen

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-000312

Abstract

This protocol for the sCD40L ELISA kit (Cat#BMS6010, Invitrogen) is used to measure the concentration of sCD40L in serum and CSF samples in mice.

Troubleshooting

- 1 Prepare Wash Buffer.
- 2 Wash the microwell strips twice with 300 μ L Wash Buffer per well. Allow the Wash buffer to sit for 10–15s before harshly drop it in the strain. Wash again.
- 3 Empty wells on the paper towel, use the strip immediately after washing.
- 4 Prepare the standard dilution as the protocol. Using Sample Diluent (provided) to dilute samples.
- 5 Add 100 μ L of Sample Diluent in duplicate to the blank wells.
- 6 For serum samples, add 50 μ L of Sample Diluent to the sample wells and then add 50 μ L of serum samples to the sample well. For CSF samples, add 5 μ L of Sample Diluent to the sample wells and then add 95 μ L of serum samples to the sample well.
- 7 Prepare Bio-Conjugate (fresh) and add 50 μ L of Bio-Conjugate to all wells.
- 8 Apply an adhesive film to the plate and let it incubate at room temperature (RT) for 2 hrs on a microplate shaker.
- 9 Remove adhesive film and empty wells. Wash 4 times using Wash Buffer.
- 10 Prepare Streptavidin-HRP and add 100 μ L of diluted Streptavidin-HRP to all wells.
- 11 Apply an adhesive film to the plate and let it incubate at RT for 1 hr on a microplate shaker.
- 12 Remove adhesive film and empty wells. Wash 4 times using Wash Buffer.
- 13 Add 100 μ L TMB substrate solution to all wells. Incubate the strips at RT until the color changes. Avoid light.



- 14 Add 100 μ l the stop solution to all wells.

- 15 Measure the absorbance of each microwell on a spectro-photometer at a wavelength of 450 nm.