


Aug 11, 2025

Version 1

🌐 scDNA-scRNA-seq (SDR-seq) V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.6qpvr9q43vmk/v1>



Dominik Lindenhofer^{1,2}, Julia Bauman³, John A. Hawkins^{1,4}, Donnacha Fitzgerald^{1,5,6}, Umut Yildiz¹, Haeyeon Jung¹, Anastasiia Korosteleva^{1,7}, Mikael Marttinen^{8,9}, Moritz Kueblbeck¹, Judith B. Zaugg^{10,11}, Kyung-Min Noh^{1,12}, Sascha Dietrich^{5,13,14,6}, Wolfgang Huber^{1,6}, Oliver Stegle^{1,4}, Lars M. Steinmetz^{1,2,3,15}

¹Genome Biology Unit, European Molecular Biology Laboratory (EMBL), Heidelberg, Germany;

²DZHK (German Centre for Cardiovascular Research), partner site Heidelberg/Mannheim, Heidelberg, Germany;

³Department of Genetics, Stanford University School of Medicine, Stanford, CA, USA;

⁴Division of Computational Genomics and Systems Genetics, German Cancer Research Center (DKFZ), Heidelberg, Germany;

⁵Department of Hematology, Oncology and Rheumatology, Heidelberg University Hospital, Heidelberg, Germany;

⁶Molecular Medicine Partnership Unit (MMPU), European Molecular Biology Laboratory (EMBL), Heidelberg, Germany;

⁷Faculty of Biosciences, Heidelberg University, Heidelberg, Germany;

⁸Molecular Systems Biology Unit, European Molecular Biology Laboratory (EMBL), Heidelberg, Germany;

⁹Faculty of Medicine and Health Technology, Tampere University, Tampere, Finland;

¹⁰Structural and Computational Biology Unit, European Molecular Biology Laboratory, Heidelberg, Germany;

¹¹Department of Biomedicine, University of Basel | University Hospital Basel, Basel, Switzerland;

¹²Department of Biomedicine, Aarhus University, Aarhus, Denmark;

¹³Department of Oncology, Hematology and Clinical Immunology, Medical Faculty of Heinrich-Heine-Universität, Düsseldorf University Hospital, Düsseldorf, Germany;

¹⁴Center for Integrated Oncology Aachen-Bonn-Cologne-Düsseldorf (CIO ABCD), Aachen Bonn Cologne Düsseldorf, Germany;

¹⁵Stanford Genome Technology Center, Palo Alto, CA, USA

Protocols



Dominik Lindenhofer

EMBL Heidelberg

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvr9q43vmk/v1>

Protocol Citation: Dominik Lindenhofer, Julia Bauman, John A. Hawkins, Donnacha Fitzgerald, Umut Yildiz, Haeyeon Jung, Anastasiia Korosteleva, Mikael Marttinen, Moritz Kueblbeck, Judith B. Zaugg, Kyung-Min Noh, Sascha Dietrich, Wolfgang Huber, Oliver Stegle, Lars M. Steinmetz 2025. scDNA-scRNA-seq (SDR-seq). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.6qpvr9q43vmk/v1>

Manuscript citation:

Lindenhofer, Dominik, Julia R. Bauman, John A. Hawkins, Donnacha Fitzgerald, Umut Yildiz, Jan M. Marttinen, Moritz Kueblbeck, et al. 2024. "Functional Phenotyping of Genomic Variants Using Multiomic scDNA-scRNA-Seq." bioRxiv. <https://doi.org/10.1101/2024.05.31.596895>.

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 30, 2025

Last Modified: August 11, 2025

Protocol Integer ID: 119358

Keywords: SDR-seq, single-cell, scDNA-scRNA-seq, scRNA-seq, scDNA-seq, Multiomics, capture sequence to the cdna, rna target, reverse transcription, situ reverse transcription, rna, transcriptomic target, reverse primers for each gdna, multiplexed pcr amplifies gdna, transcriptomic targets in thousand, r2 truseq for rna, containing cdna, scrna, scdna, distinct cell barcode oligo, nextera for gdna, cell barcode oligo, microfluidic tapestri device from missionbio, unique molecular identifier, sdr, microfluidic tapestri device, adding unique molecular identifier, cell barcoding, gdna, single cells per experiment, pcr reagent, cell suspension, sample barcode, reverse primer, single cell, multiplexed pcr, capture sequence, cdna, sequencing depth, cell

Funders Acknowledgements:

HFSP postdoctoral fellowship

Grant ID: LT0023/2022-L

NSF graduate research fellowship

Grant ID: DGE-1656518

Stanford Graduate Fellowship

Grant ID: -

Artificial Intelligence in Health Innovation Cluster

Grant ID: -

German Federal Ministry of Education and Research

Grant ID: SIMONA, 031L0263A

EMBL-GSK collaboration Fund

Grant ID: -

EMBL predoctoral fund

Grant ID: -

Dieter Schwarz Foundation Professorship

Grant ID: -

Abstract

scDNA-scRNA-seq (SDR-seq) is a method that enables a targeted readout of both genomic and transcriptomic targets in thousands of single cells per experiment. It combines an *in-situ* reverse transcription (RT) step of fixed cells with a multiplexed PCR in an emulsion-based droplets using the microfluidic Tapestri device from MissionBio. A single-cell suspension is fixed with glyoxal, permeabilized and used to perform *in-situ* reverse transcription. This converts RNA to cDNA while adding unique molecular identifiers (UMIs), a sample barcode (sample BC) and a capture sequence to the cDNA. These cells containing cDNA and gDNA are used as an input for the Tapestri machine where cells are lysed, treated with proteinase K and mixed with reverse primers for each gDNA or RNA target in a first droplet. This droplet is then fused into a second droplet that contains PCR reagents, forward primers containing a capture sequence (CS) overhang and a cell barcoding bead with distinct cell barcode oligos and a matching CS. The multiplexed PCR amplifies gDNA and RNA targets simultaneously in each droplet while cell barcoding is ensured by the complementary CS in the forward primers and the cell barcode oligos. After breaking emulsions, separate NGS libraries for both gDNA and RNA can be generated by having distinct overhangs on the reverse primers (R2N Nextera for gDNA; R2 TruSeq for RNA). This enables to sequence each library at different sequencing depth and length.

Materials

Glyoxal – Sigma #128465

Acetic acid (Glacial) – Sigma #695092

RNasin® Ribonuclease Inhibitor – Promega #N2615

Enzymatics Rnase Inhibitor - Enzymatics #280520

IGEPAL CA-630 – Sigma #I8896

Digitonin – ThermoFisher Scientific #BN2006

dNTP-Mix - ThermoFisher Scientific #R0192

Maxima H Minus Reverse Transcriptase - ThermoFisher Scientific #EP0751

Tapestri Single-Cell DNA CoreAmbient Kit v2 - MB51-0007

Tapestri Single-Cell DNA Core -20 Kit v2 - MB51-0010

Tapestri Single-Cell DNA Bead Kit - MB51-0009

Troubleshooting

Design and order oligos

1 In-situ RT primers

The following oligos need to be ordered to perform the *in-situ* RT reactions. We recommend ordering at least 8 RT primers for optimal performance of SDR-seq. More RT primers can be ordered in case more samples are required or to improve removal of contaminating ambient RNA reads during analysis.

A	B
Name	Sequence
RT_primer_1	GTACTCGCAGTAGTCGACACGTCTCGCCTTANNNNNNNNTTTTTTTT TTTTTTTTVN
RT_primer_2	GTACTCGCAGTAGTCGACACGTCCTAGTACGNNNNNNNNNTTTTTTTT TTTTTTTTVN
RT_primer_3	GTACTCGCAGTAGTCGACACGTCTTCTGCCTNNNNNNNNNTTTTTTTT TTTTTTTTVN
RT_primer_4	GTACTCGCAGTAGTCGACACGTCGCTCAGGANNNNNNNNTTTTTTTT TTTTTTTTVN
RT_primer_5	GTACTCGCAGTAGTCGACACGTCAGGAGTCCNNNNNNNNNTTTTTTTT TTTTTTTTVN
RT_primer_6	GTACTCGCAGTAGTCGACACGTCCATGCCTANNNNNNNNTTTTTTTT TTTTTTTTVN
RT_primer_7	GTACTCGCAGTAGTCGACACGTCGTAGAGAGNNNNNNNNNTTTTTTTT TTTTTTTTVN
RT_primer_8	GTACTCGCAGTAGTCGACACGTCCCTCTCTGNNNNNNNNNTTTTTTTT TTTTTTTTVN

Overview of 8 RT primers used for in-situ RT reaction. CS: GTACTCGCAGTAGTC; CS extension for universal primer: GACACGTC. Sample barcodes: N701-706

2 Targeted primer design for gDNA and RNA primers

Targeted gDNA primers were designed using the Tapestry Designer online tool (<https://designer.missionbio.com>). Targeted RNA primers were selected by taking the “inner primer” sequences designed using the TAP-seq primer prediction tool with a targeted Tm of 60 ° C (min 58 ° C and max 62 ° C) and a product size range from 150-300 bp (<https://github.com/argschwind/TAPseq>). Two target-specific primers are needed for each gDNA target. Only one target-specific primer is needed for RNA targets, because all RNA targets are amplified by a universal primer binding to the CS + extension introduced via the *in-situ* RT. Primers designed by the respective design tools require overhangs as described in the table below.

A	B
Name	Sequence
Droplet_gDNA_fwd	GTACTCGCAGTAGTC-gDNA_specific_sequence
Droplet_gDNA_rev	GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAG-gDNA_specific_sequence
Droplet_RNA_fwd	GTACTCGCAGTAGTCGACACGTC
Droplet_RNA_rev	GTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT-cDNA_specific_sequence

3 Primers for separate library generation

After breaking emulsions following the multiplexed PCR in emulsion-based droplets library primers with the respective overhangs are required to give specificity for gDNA (R2N) and RNA (R2) targets. gDNA targets will be amplified with library primers containing R1N and R2N binding sequences. RNA targets will be amplified with library primers containing R1N and R2 binding sequences.

A	B
Name	Sequence
R1N_i5_S502	AATGATACGGCGACCACCGAGATCTACACCTCTCTATTCGTCGGCAGCGTC

A	B
R1N_i5_S503	AATGATACGGCGACCACCGAGATCTACACTATCCTCTTCGTCCGGCA GCGTC
R1N_i5_S505	AATGATACGGCGACCACCGAGATCTACACGTAAGGAGTCGTCCGGCA GCGTC
R1N_i5_S506	AATGATACGGCGACCACCGAGATCTACACACTGCATATCGTCCGGCA GCGTC
R1N_i5_S507	AATGATACGGCGACCACCGAGATCTACACAAGGAGTATCGTCCGGCA GCGTC
R1N_i5_S508	AATGATACGGCGACCACCGAGATCTACACCTAAGCCTTCGTCCGGCA GCGTC
R2N_i7_N701	CAAGCAGAAGACGGCATAACGAGATTTCGCCTTAGTCTCGTGGGCTCG G
R2N_i7_N702	CAAGCAGAAGACGGCATAACGAGATCTAGTACGGTCTCGTGGGCTCG G
R2N_i7_N703	CAAGCAGAAGACGGCATAACGAGATTTCTGCCTGTCTCGTGGGCTCG G
R2N_i7_N704	CAAGCAGAAGACGGCATAACGAGATGCTCAGGAGTCTCGTGGGCTCG G
R2N_i7_N705	CAAGCAGAAGACGGCATAACGAGATAGGAGTCCGTCTCGTGGGCTCG G
R2N_i7_N706	CAAGCAGAAGACGGCATAACGAGATCATGCCTAGTCTCGTGGGCTCG G
R2_i7_N701	CAAGCAGAAGACGGCATAACGAGATTTCGCCTTAGTGACTGGAGTTCA



A	B
	GACGTGT
R2_i7_N702	CAAGCAGAAGACGGCATAACGAGATCTAGTACGGTGACTGGAGTTCA GACGTGT
R2_i7_N703	CAAGCAGAAGACGGCATAACGAGATTTCTGCCTGTGACTGGAGTTCA GACGTGT
R2_i7_N704	CAAGCAGAAGACGGCATAACGAGATGCTCAGGAGTGACTGGAGTTCA GACGTGT
R2_i7_N705	CAAGCAGAAGACGGCATAACGAGATAGGAGTCCGTGACTGGAGTTCA GACGTGT
R2_i7_N706	CAAGCAGAAGACGGCATAACGAGATCATGCCTAGTGACTGGAGTTCA GACGTGT

Prepare primer pools for targeted gDNA and RNA amplification

4 Prepare oligo pools

- 4.1 Forward and reverse oligo pools containing the targeted gDNA and RNA primers were ordered as oPools from IDT (50 pmol/oligo).
- 4.2 Reverse gDNA and RNA oligo pools were diluted to a concentration of 120 μM in H₂O each and then mixed 1:1 for usage in Step 2.6 of the **Tapestri V2 User Guide**.
- 4.3 Forward gDNA oligo pools were diluted to 20 μM . The universal Droplet_RNA_fwd oligo that binds at all RNA targets was also diluted to 20 μM . The forward gDNA oligo pools and the universal Droplet_RNA_fwd oligo were mixed 1:1 for usage in Step 4.8 of the **Tapestri V2 User Guide**.

Prepare Buffers

- 5 All Buffer reactions are calculated for fixation of two sampels. Adjust volumes according to sample number.

Glyoxal fixation buffer pH 0

	A	B	C	D	E
	Reagent	Stock	Final (1x)	Volume (µl)	Comment
	EtOH (%)	100	20	400	
	Glyoxal (%)	40	3	150	
	Acetic acid (Glacial) (%)	100	0.75	15	
	NaCl (mM)	5000	120	48	
	H2O			1362	
	NaOH (1M)			25	to pH 4 (adjust volume)
	Final Volume			2000	

Wash Buffer 1

	A	B	C	D
	Reagent	Stock	Final (1x)	Volume (µl)
	DTT (mM)	100	1	45
	BSA (%)	10	2	900
	RNasin Ribonuclease Inhibitor (U/µl)	20	0.04	9
	PBS (1x)			3546

A	B	C	D
Final Volume			4500

Wash Buffer 2

A	B	C	D
Reagent	Stock	Final (1x)	Volume (µl)
TRIS pH7.5 (mM)	1000	10	40
NaCl (mM)	5000	10	8
MgCl ₂ (mM)	1000	3	12
Tween20 (%)	10	0.1	40
RNasin Ribonuclease Inhibitor (U/µl)	20	0.1	20
DTT (mM)	100	1	40
BSA (%)	10	2	800
H ₂ O			3040
Final Volume			4000

Permeabilization Buffer

A	B	C	D	E
Reagent	Stock	Final (1x)	Volume (µl)	Comment
IGEPAL CA-630 (%)	10	0.1	6	
Digitonin (%)	5	0.01	1.2	preheat to 65 C for 10 min

	A	B	C	D	E
	Wash Buffer 2			592.8	
	Final Volume			600	


Resuspension buffer

	A	B	C	D
	Reagent	Stock	Final (1x)	Volume (µl)
	DTT(mM)	100	1	10
	BSA (%)	10	2	200
	RNasin Ribonuclease Inhibitor (U/µl)	20	0.3	15
	1x PBS			775
	Final Volume			1000

Dissociation - Fixation - Permeablization

6 Dissociation

- 6.1 Prepare single cell suspension according to cell line or tissue requirements as for other single cell technologies. Filter cell suspension through a 40 µm cell strainer.
- 6.2 We recommend using 1.5×10^6 cells as input for fixation. We tested performance of SDR-seq with a minimum of 0.35×10^6 cells as input for fixation. Resuspend cells in 1 ml DPBS -/-. To proceed with samples viability should be >85 %.



Example of cells before fixation using Countess II

7 Fixation

- 7.1 We recommend using 15 ml Polypropylene centrifuge tubes and swinging bucket rotors in the appropriate centrifuges to minimize cell loss during this process.
- 7.2 Spin cells at 300g for 5 min and remove supernatant.
- 7.3 Resuspend cells in 200 μ l of glyoxal fixation buffer (pH adjusted) and incubate at room temperature for 7 min.
- 7.4 Add 1 ml of ice-cold wash buffer 1 and gently invert the tube 2-3 times to mix.
- 7.5 Spin at 400 g for 5 min at 4 °C and carefully take off supernatant.
- 7.6 Repeat step 7.4 and 7.5 for a total of 2x washes.

- 7.7 Resuspend cells in 175 μ l ice-cold permeabilization buffer and incubate on ice for 4 min.
- 7.8 Add 1 ml of ice-cold wash buffer 2 and gently invert the tube 2-3 times to mix.
- 7.9 Spin at 400 g for 5 min at 4 °C and carefully take off supernatant.
- 7.10 Resuspend cells in 200 μ l resuspension buffer.
- 7.11 Filter cell suspension through a 40 μ m cell strainer.
- 7.12 Count cells.



Example of cells after fixation using Countess II

- 7.13 Dilute cell suspension to 1.4×10⁶ cells/ml

In-situ reverse transcription

8 We recommend to use a total of 48 wells of a 96-well PCR plate with 10,000 cells/well per experiment during in-situ reverse transcription (RT). This is calculated to have enough surplus for the Tapestri microfluidic device (optimal input ranges from 1.05-2×10⁵ cells). Multiple samples can be multiplexed per Tapestri lane by using distinct sample barcodes. We recommend at using at least 8 RT sample barcode primers per SDR-seq experiment to improve removal of ambient RNA during analysis.

8.1 Prepare the following RT MM on ice. Mix by pipetting up and down for 10x times.

A	B	C	D	E
Reagent	Stock	Final	1x (μl)	52x (μl)
5x Maxima H Minus RT Buffer	5	1	4	208
Enzymatics Rnase Inhibitor (U/μl)	40	0.25	0.125	6.5
RNasin Ribonuclease Inhibitor (U/μl)	20	0.25	0.25	13
dNTPs (μM per base)	10000	500	1	52
Maxima H Minus RT (U/μl)	200	20	2	104
H2O	NA	NA	0.625	32.5
Total volume per reaction			8	416
Final reaction volume			20	1040

8.2 If using 8 sample barcodes as we recommend prepare the following master mix per sample barcode on ice. Otherwise adjust accordingly.

A	B	C
	1x (µl)	6.25x (µl)
MM RT	8	50
RT Primer (12.5 µM)	4	25
Fixed cells (1.4 M/ml)	8	50

8.3 Add cells last, mix by pipetting up and down 10x.

8.4 Aliquot 20 µl of the master mix containing RT primers and fixed cells into the appropriate number of wells into a 96-well plate. In the above example this would be 6 wells per sample barcode.

8.5 Put the RT-plate into a Thermocycler that is preheated to 50 °C and run the following program to perform the *in-situ* RT:

A	B	C
Temperature °C	Time	Cycle
50	10 min	
8	12 sec	3x
15	45 sec	
20	45 sec	
30	30 sec	
42	2 min	
50	3 min	
50	5 min	

	A	B	C
	4	-	

- 8.6 Put the RT plate on ice.
- 8.7 Pool all RT reactions into a 15 ml Polypropylene centrifuge tube containing 10 ml DPBS +/- with 1 % BSA. Make sure to pipet up and down in the RT plate 2-3 times without introducing bubbles to resuspend cells in the RT plate before pooling. Wash each well with 200 µl DPBS +/- with 1 % BSA to reduce cell loss during this step.
- 8.8 Spin cells at 500 g for 10 min at 4 °C and carefully take off supernatant.
- 8.9 Resuspend in 35-50 µl ice-cold Cell Buffer provided by MissionBio.
- 8.10 Count and dilute cell suspension to 4,500 cells/µl.

Encapsulation, Lysis, Digestion, Barcoding, Targeted PCR and CleanUp

- 9 We then followed the [Tapestri V2 User Guide](#) Step2 (Encapsulate Cells), Step 3 (Lyse and Digest Cells), Step 4 (Barcode Cells), Step 5 (UV Treatment and Targeted PCR Amplification) and Step 6 (Cleanup PCR Products). gDNA and RNA targets combined were treated as the total amplicon panel size and the specifications of the multiplexed PCR was adjusted according to the recommendations of the [Tapestri V2 User Guide](#).

Library PCR

- 10 Library PCR primers were used to amplify separate gDNA (R1N + R2N) and RNA (R1N +R2) NGS libraries.
- 10.1 Set up the following PCR reactions:

PCR – gDNA library

A	B
Reagent	μl
MisssionBio Library PCR MM (2x)	25
Cleaned up sample	15
Primer R1N – i5 (2.5 μM)	5
Primer R2N – i7 (2.5 μM)	5

PCR – RNA library


A	B
Reagent	μl
MisssionBio Library PCR MM (2x)	25
Cleaned up sample	15
Primer R1N – i5 (2.5 μM)	5
Primer R2N – i7 (2.5 μM)	5

10.2 Run the following PCR cycling program for both libraries on a thermocycler.


A	B	C
Tempearture °C	Time	Cycle
95	3 min	
98	20 sec	10x
62	20 sec	

A	B	C
72	45 sec	
72	2 min	
4	-	

10.3 Follow the library cleanup described in Step 7.7 and the subsequent recommendations in Step 8 (Quantify and Normalize Sequencing Library) in the [Tapestri V2 User Guide](#).



Example Bioanalyzer track of gDNA library.



Example Bioanalyzer track of RNA library.

Sequencing

- 11 The separate gDNA and RNA libraries are sequenced with the following read length specifications on Illumina sequencers. Read depth can be adjusted according to experimental needs. We recommend a read depth of at least 10000 reads/cell per 100 amplicons for each library.

gDNA library

	A	B	C
	Read	Cycles	Type
	R1N	161	Cell barcode + gDNA amplicon
	R2N	161	gDNA amplicon
	i5	8	Library PCR sample BC
	i7	8	Library PCR sample BC

RNA library

	A	B	C
	Read	Cycles	Type
	R1N	76	Cell barcode + <i>in-situ</i> RT sample BC +UMI
	R2	46	RNA amplicon
	i5	8	LibraryPCR sample BC
	i7	8	Library PCR sample BC

Analysis

- 12 To generate count matrices for gDNA and RNA libraries please use the package **SDRranger**. We recommend selecting for high-quality cells first followed by variant calling using the GATK HaplotypeCaller (v4.2.3.0) in each cell separately.

Acknowledgements

We thank all members of the Steinmetz lab for feedback and technical expertise. We thank the Genomics Core Facility (GeneCore) at EMBL for consultation and sequencing. We thank Jan-Philipp Mallm and the Single-cell Open Lab (scOpenLab) at DKFZ for expertise and access to the Tapestry instrument from MissionBio. Support by the DKFZ scOpenLab is gratefully acknowledged. We thank the Flow Cytometry Core Facility at EMBL for support and flow cytometry services. **Funding:** D.L. received funding from a HFSP postdoctoral fellowship (LT0023/2022-L). J.R.B. is supported by an NSF graduate research fellowship (DGE-1656518) and a Stanford Graduate Fellowship (Smith Fellowship). J.A.H. is supported by the European Research Council (Synergy Grant DECODE under grant agreement no. 810296) and a fellowship through state funds approved by the State Parliament of Baden-Württemberg for the Innovation Campus Health + Life Science Alliance Heidelberg Mannheim. D.F. received funding from a grant from the German Federal Ministry of Education and Research (SIMONA, 031L0263A). M.K. received funding from an EMBL-GSK collaboration Fund. U.Y. is supported by the EMBL predoctoral fund. L.S.M. is supported by the Dieter Schwarz Foundation Professorship. Work in the Steinmetz lab is supported by a grant from the Else Kröner Fresenius Stiftung (2021_EKS.:61).