

Sep 14, 2022

scNMT-seq v2

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g7yrz3gwz/v1>

Liyun Zhao¹, Thomas Lefevre², Thierry Voet², Bernard Thienpont¹

¹Laboratory for Functional Epigenetics, Department of Human Genetics, KU Leuven;

²Laboratory of Reproductive Genomics, Department of Human Genetics, KU Leuven



Liyun Zhao

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7yrz3gwz/v1>

Protocol Citation: Liyun Zhao, Thomas Lefevre, Thierry Voet, Bernard Thienpont 2022. scNMT-seq v2. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.q26g7yrz3gwz/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: September 08, 2022

Last Modified: September 14, 2022

Protocol Integer ID: 69707

Keywords: Single-cell, scRNA-seq, G&T-seq, BS-seq, chromatin accessibility, DNA methylation, gene expression, rna, single cell, single, cell, ASAPCRN, scNMTseq, open chromatin, seq v2 scnmt, testing modifications of the scnmt, original scnmtseq, single cell chromatin status, scnmt, transcriptome sequencing, single cell nucleosome, methylation profile, seq pipeline, methylome, strand synthesis, problem of unmapped read, unmapped read, seq,

Abstract

scNMT-seq (single cell Nucleosome, Methylome, and Transcriptome sequencing) allows the parallel study of a single cell chromatin status, methylation profile, and transcriptome.

Here, we are developing and testing modifications of the scNMT-seq pipeline. The protocol is carried out in 96w plates and typically takes 4-5 days to complete.

The number of pre-amplification cycles is adjusted to tackle the problem of poor recovery after BS conversion. Primers are optimized for first-strand and second-strand synthesis to solve the problem of unmapped reads and poor amplification. Both are testified as compatible with the original the original scNMTseq.

Materials

- ✕ GpC Methyltransferase (M.CviPI) - 1,000 units **New England Biolabs Catalog #M0227L**
- ✕ IGEPAL-CA630 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #I3021 SIGMA-ALDRICH**
- ✕ RNase Inhibitor **Thermo Fisher Catalog #N8080119**
- ✕ RLT Plus Buffer **Qiagen**
- ✕ Dynabeads MyOne Streptavidin T1 **Thermo Fisher Scientific Catalog #65601**
- ✕ Superscript II **Invitrogen - Thermo Fisher Catalog #18064014**
- ✕ ERCC RNA Spike-In Mix or order custom-made synthetic sequences **Thermo Fisher Scientific Catalog #4456740**
- ✕ dNTP Mix (10 mM ea) **Thermo Fisher Catalog #18427013**
- ✕ Nextera XT DNA Sample Preparation Kit, 96 samples **Illumina, Inc. Catalog #FC-131-1096**
- ✕ Nextera XT Index Kit v2 Set A (96 indexes 384 samples) **Illumina, Inc. Catalog #FC-131-2001**
- ✕ Betaine 5M **Merck MilliporeSigma (Sigma-Aldrich) Catalog #B0300**
- ✕ Magnesium Chloride (MgCl₂) Solution - 6.0 ml **New England Biolabs Catalog #B9021S**
- ✕ DTT **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D0632**
- ✕ Kapa HiFi Hotstart ReadyMix (2x) **Kapa Biosystems Catalog #KK2612**
- ✕ Ampure XP beads **Beckman Coulter Catalog #A63881**
- ✕ Lambda DNA **Thermo Fisher Catalog #SD0011**
- ✕ EZ-96 DNA Methylation-Direct MagPrep **Zymo Research Catalog #D5044**
- ✕ Klenow (3' → 5' exo-) (High Concentration) **Enzymatics Catalog #P7010-HC-L**
- ✕ QuantiFluor® dsDNA System **Promega Catalog #E2670**

	A	B
	Name	Sequence (5' to 3')
	Oligo-dT primer	Biotin-TEG-AAGCAGTGGTATCAACGCAGAGTACTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTVN
	TSO	AAGCAGTGGTATCAACGCAGAGTACATrGrG+G
	IS PCR	AAGCAGTGGTATCAACGCAGAGT
	Pre-amplification primer	Biotin-TGACTGGAGTTCAGACGTGTGCTCTTCCGATCTHHHHHHH*H
	Adapter2 primer	ACACTCTTCCCTACACGACGCTCTTCCGATCTDDDDDDD*D



	A	B
	PE 1.0	AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATC* T
	iPCRTag primer	CAAGCAGAAGACGGCATACGAGATXXXXXXXXXXGTGACTGGAGTTCAGACGTGTGCTC TTCCGATC*T

All oligos should be ordered with HPLC purification



Cell isolation and GpC methylation

1

Prepare plates containing **2.5µl** GpC methylase reaction mixture in each well:

	A	B	C	D
	Component (initial)	Component(final)	Volume (µl)	Master mix (110 samples)
	M.CviPI reaction buffer (10x)	1x	0.25	27.5
	M.CviPI (4U/µl)	2U	0.5	55
	SAM (320µM)	160µM	1.25	137.5
	IGEPAL (10%)	0.1%	0.025	2.75
	RNasein (20U/µl)	1U/µl	0.125	13.75
	Nuclease-free water		0.35	38.5

2

Isolate cells manually or using FACS in  2.5 µL of GpC methylase reaction buffer

3

After cell isolation, spin down plates at **≥1000g** for **≥10s**  4 °C

4

Incubate the samples at  37 °C for  00:15:00

15m



5

Stop reaction by adding  5 μL **RLT plus buffer**

6

Store the plates at  $-80\text{ }^{\circ}\text{C}$ until processed

Oligo-dT₃₀VN bead preparation

30m

7

Add  55 μL Dynabeads into a new Eppendorf tube. Place the tube on a magnet for  00:02:00 and discard supernatant

2m

8

Resuspend beads in  200 μL **Dynabead solution A** ( 0.1 Mass Percent **NaOH**,  0.05 Mass Percent **NaCl**). Place the tube on a magnet for  00:02:00 and discard supernatant

2m

9

Repeat step 8 once

10

Resuspend beads in  200 μL **Dynabead solution B** ( 0.1 Mass Percent **NaCl**). Place on a magnet for  00:02:00 and discard supernatant

2m

11

Resuspend the beads in  55 μL of **2x B&W** ( 2 Mass Percent **NaCl**,  10 millimolar (mM) **Tris-HCl**,  1 millimolar (mM) **EDTA**) and  55 μL Biotinylated Oligo-dT30VN ( 100 micromolar (μM)). Incubate  00:20:00 on a thermomixer while shaking at **2000rpm** at  Room temperature

20m

*In the meantime, prepare the **bead resuspension buffer***



A	B
Superscript FS buffer (5x)	220µl
Nuclease-free Water	825µl
RNase inhibitor (20U/µl)	55µl

After adding RNase inhibitor, use beads within 30min

In the meantime, prepare **1x B&W** buffer by mixing  440 µL Nuclease-free water with

 440 µL **2x B&W** buffer

12

Place beads on a magnet for  00:02:00 and discard supernatant

2m

13

Resuspend the beads in  200 µL **1x B&W** buffer. Place beads on a magnet for

 00:02:00 and discard supernatant

2m

14

Repeat step 13 **three** more times

15

Resuspend the beads in the **bead resuspension buffer**

Physical separation of mRNA and gDNA

16

Thaw the 96-well plate containing the single cell lysates on ice

17

Add  1 µL **ERCC spike-ins** at 1:1Million - 1:128Million dilution to each sample using a multi-dispensing pipette. Run the pulse centrifugation program to spin ERCCs down to the bottom

18

Take 4 tubes( 1069 µL per tube) of **G&T wash buffer**([M] 50 millimolar (mM) **Tris-HCl**  8.3 , [M] 75 millimolar (mM) **KCl**, [M] 3 millimolar (mM) **MgCl₂**, **0.5% Tween**

20 Solution) and add to each tube  137.5 µL **DTT** and  25 µL **RNaseIn**

19

Add  50 µL **G&T-Seq wash buffer** per well to the "G&T-Seq wash plate"

20

Add  10 µL **Oligo-dT beads** per well to the "bead plate"

21

Add an empty **non-skirted** 96 well plate labeled "gDNA collection"

22

Spin all plates and run the adapted G&T-separation program robotically or manually.

*While the separation program is running, prepare the **RT master mix***

	A	B	C	D
	Component (C_initial)	C_final	Volume(µl)	Mastermix (110 samples)
	dNTP (10mM)	1mM	0.5	55
	TSO (100µM)	1µM	0.05	5.5
	MgCl ₂ (1M)	6mM	0.03	3.3
	Betain (5M)	1M	1	110
	S II First strand buffer (5x)	1x	1	110
	DTT (100mM)	5mM	0.25	27.5
	Nuclease-free water		1.8	198
	RNase inh (20U/µl)	0.5U/µl	0.125	13.75
	Superscript reverse transcriptase II (200U/µl)	10U/µl	0.25	27.5

Adding enzyme within less than 30 min before running the Reverse Transcription program

**Note**

Separation is performed robotically on the Hamilton platform in this protocol. If performed manually, steps should be as follows

22.1

Manually pipette  10 μL of prepared **oligo-dT beads** to each well of the **sample plate** using a multichannel pipette

22.2

Mix at maximum speed for  00:20:00

20m

22.3

Place on magnet for  00:05:00 . Aspirate  17.5 μL and transfer to the empty low-bind plate for gDNA collection

5m

22.4

Add  15 μL of G&T-seq wash buffer off magnet.

22.5

Mix at maximum speed for  00:10:00

10m

22.6

Place on magnet for  00:02:00 . Aspirate  15 μL and transfer to the empty low-bind plate for gDNA collection

2m

22.7

Repeat steps 22.3-22.6 twice more

Note

Lysate (17.5ul) combined with 3 washes (15ul each) should now have been collected into the gDNA plate

Reverse transcription

1h 45m



23

Collect the polyA(+) mRNA plate and using the multi-dispenser dispense  5 μ L **RT master mix** into each well of the bead-containing 96-well plate

24

Seal the mRNA and gDNA plates and spin.

Store gDNA at  -80 °C until processed

25

Incubate the polyA(+) mRNA 96-well plate on a **thermomixer C** using the program below (approx. duration  01:45:00)

1h 45m

	A	B	C	D
	Cycle	Temp (°C)	Time	Mixing (rpm)
	1	42	2 min	2000
	2	42	60 min	1500
	3	50	30 min	1500
	4	60	10 min	1500

26

In the meantime prepare **PCR mastermix**

	A	B	C
	Component	Volume(μ l)	Mastermix (110 samples)
	KAPA HiFi HotStart ReadyMix (2x)	6.25	687.5
	IS PCR primer (10 μ M)	0.124	13.64
	Nuclease-free water	1.13	124.3



PCR amplification of cDNA

30s

27

Add  7.5 μ L **PCR reaction mastermix**, seal the plate and centrifuge

28

Resuspend the beads for  00:00:30 at **2000rpm** using the Thermomixer

30s

29

Perform cDNA amplification as follows

	A	B	C
	Cycles	Temperature (°C)	Time
	1	98	3 min
	18-25	98	20 s
		67	15 s
		72	6 min
	1	72	5 min
	1	4	Hold

Amplification cycles differ

PCR cleanup of amplified cDNA

22m 10s

30

Add  12.5 μ L **Agencourt AMPure beads** (1:1 ratio), mix thoroughly by pipetting up and down

31

Incubate  00:05:00 at  Room temperature

5m



32

Pellet the beads on a Low-elution magnet for  00:05:00

5m

33

Remove the supernatant without disturbing the beads

34

Wash the beads twice with  150 μL of freshly prepared 80% ethanol for

 00:00:10

10s

35

Allow the beads to dry for approximately  00:05:00 . Resuspend in  25 μL nuclease-free water. Incubate for  00:02:00  Room temperature

7m

36

Return the 96-well plate to the magnet and allow the Agencourt AMPure beads to settle for  00:05:00

5m

37

Carefully transfer the supernatant to a new 96-well plate

Note

Quality control: QUBIT+BIOANALYZER
expected cDNA concentration: $\geq 1\text{ng}/\mu\text{l}$
expected cDNA length: 500-2000bp, peaking at 1-1.5kb

Library preparation of cDNA (Nextera XT)

9m

38

Dilute the **cDNA** of each sample to **0.2ng/ μl** with nuclease-free water

39

Add  2.5 μL **Tagment DNA(TD) buffer** to a new Hard-Shell skirted 96-well plate



40

Add  1.25 μL **diluted cDNA** and  1.25 μL **amplicon tagment mix (ATM)** to TD buffer and mix

41

Centrifuge the plate at  280 x g, 20°C, 00:01:00

1m

42

Incubate on a thermal cycler

	A	B	C
	Segment	Temp(°C)	Duration(min)
	1	55	10
	2	10	Hold

43

Add  1.25 μL Neutralize Tagment Buffer (NT)

44

Vortex & spin down at  800 x g, 20°C, 00:01:00

1m

45

Incubate at  Room temperature  00:05:00

5m

46

Add  1.25 μL **Index (i7) adapter** to each column and  1.25 μL **Index 2 (i5) adapter** to each row

47

Add  3.75 μL **Nextera PCR mastermix** and mix

48

Centrifuge the samples at  280 x g, 20°C, 00:01:00 and amplify as follows:

1m



	A	B	C
	Cycle	Temp (°C)	Duration
	1	72	3min
	2	95	30s
	3-14	95	10s
		55	30s
		72	30s
	15	72	5min
	16	4	Hold

49

Centrifuge the plate at  280 x g, 20°C, 00:01:00

1m

50

Purify libraries at a **0.66:1 ratio** and elute in  12.5 µL **EB buffer****Note***Libraries can be stored for at least a year at -20°C*

51

Pool libraries and quantify using qPCR

Note*expected pool concentration: 4nM
expected pool size: 250-1500bp*



scBS-seq library preparation (gDNA)

52

Bisulfite conversion

39m

53

19m

Prepare the CT conversion reagent by mixing  7.9 mL **M-Solubilisation buffer** and  3 mL **M-Dilution buffer** and  00:15:00 vortexing at  Room temperature. Finally, add  1.6 mL **M-Reaction buffer** and vortex  00:04:00 at  Room temperature.

54

Add  32.5 μ L **AMPure XP beads** to the gDNA plate (0.65:1 ratio)

55

20m

Incubate  00:20:00  Room temperature

56

20m

Place the plate on the magnet for  00:20:00 and discard the supernatant

57

Wash the beads twice with  200 μ L 80% ethanol

58

Resuspend the beads in  10 μ L elution buffer, optionally containing 60fg unmethylated lambda DNA

**Note**

Do not transfer the samples from the beads
Do not dry the beads after the second wash, a dry step when purifying gDNA lowers recovery

59

Add  65 μL CT conversion reagent **without mixing**

Note

Watch out for bubbles, centrifuge shortly if necessary

60

Incubate the mixture as follows:

	A	B	C
	Segment	Temperature($^{\circ}\text{C}$)	Duration(min)
	1	98	8
	2	65	180
	3	4	Hold

Note

BS converted DNA is stable for 3 days at -20°C or 20h at 4°C

Purification of the bisulfite converted DNA

33m

61



Mix  300 μ L M-binding buffer and  5 μ L MagBinding beads

Note

*Tip: to minimize loss of sample due to pipetting use a thermomixer to mix instead of pipetting
Use a deep-well plate*

62

Transfer the samples to the M-binding buffer - MagBinding beads mix and incubate

 00:05:00

 Room temperature

5m

63

Pellet the beads on a magnet for  00:03:00 and discard the supernatant

3m

64

Resuspend the beads in  200 μ L M-Wash buffer

65

Pellet beads on the magnet and discard the supernatant. Resuspend the beads in

 100 μ L M-Desulphonation buffer and incubate  00:15:00

 Room temperature

15m

Note

The beads sink quite fast to the bottom, during these 15 mins you can slowly mix on regular basis with the thermomixer

66

Pellet beads on the magnet and discard the supernatant. Wash the beads twice with

 200 μ L M-Wash buffer

67

10m



Dry the beads on a heating element at 55 °C for 00:10:00

In the meantime, prepare the **pre-amplification mix** as follows

	A	B	C	D
	Component	Amount (μl)	Final concentration	Mastermix (110 samples)
	Blue buffer (10×)	4	1x	440
	dNTP mix (10mM)	1.6	0.4mM	176
	Preamp Oligo (10 μM)	1.6	0.4 μM	176
	H2O	32.8		3608
	Total volume	40		4400

Pre-amplification

8m 5s

68

Resuspend the beads in a 40 μL **pre-amplification mix**

69

Incubate the mixture at 55 °C for 00:04:00 and place it on the magnet

4m

70

After the beads are pelleted transfer 39 μL to a new plate

71

Incubate the samples 00:03:00 at 65 °C and immediately cool on a pre-cooled aluminum rack

3m 10s

Centrifuge the plate at 500 x g, Room temperature, 00:00:10

72

Add 1 μL **klenow exo- polymerase** (50U/μl)

Vortex the samples and amplify as follows:



	A	B	C	D
	Segment	Temp (°C)	Duration (min)	Ramp speed (°C/min)
	1	4	5	-
	2	4-37	8.25	4
	3	37	30	-
	4	4	Hold	

73

In the meantime, prepare 6 tubes of **pre-amplification mix**

Note

Only add klenow exo to the mix before use

	A	B	C	D
	Component	Amount (µl)	Final concentration	Mastermix (samples)
	Blue buffer (10×)	0.25	1x	
	dNTP mix (10mM)	0.1	0.4mM	
	Preamp Oligo (10 µM)	1	4 µM	
	Klenow exo- (50 U/µl)	0.5	10 U/µl	
	H2O	0.65		
	Total volume	2.5		

74

Heat the plate to 95 °C for 00:00:45 and transfer it to an aluminum rack pre-cooled on ice

45s



75

Centrifuge the plate at 500g for  00:00:10 at 15-25°C

10s

76

Add  2.5 µL of the **pre-amplification mix**

77

Repeat steps 72-76 **five more times**

78

Incubate as follows:

	A	B	C	D
	Segment	Temp (°C)	Duration (min)	Ramp speed (°C/min)
	1	4	5	-
	2	4-37	8.25	4
	3	37	90	-
	4	4	Hold	

Note

The first-strand product can be stored ON at 4°C or for at least a month at -20°C

Exonuclease I treatment

1h

79

Dilute the samples to a volume of  98 µL with **nuclease-free water**

80

1h

Add  2 μL **exonuclease I** (20U/ μL) to the pre-amplified product and incubate

 01:00:00 at  37 °C with the **heated lid set to**  50 °C

Purification

18m

81

Add  75 μL **AMPure XP beads** (0.75:1 ratio) and mix thoroughly by pipetting up and down

Note

Tip check the volume of some samples first and adjust volumes of beads to add accordingly

82

10m

Incubate  00:10:00  Room temperature

*In the meantime, prepare **Adaptor 2 mix***

A	B	C	D
Component	Amount (μL)	Final concentration	Mastermix (samples)
Blue buffer (10 $\times$)	4.7	1x	
dNTP mix (10mM)	1.9	0.4mM	
Adapter 2 Oligo (10 μM)	1.9	0.4 μM	
H ₂ O	38		
Total volume	46.5		

83

3m

Place on the magnet for  00:03:00 and discard the supernatant

84

Add  200 μL of 80% (vol/vol) ethanol while keeping the plate on the magnet then discard ethanol after $\pm 10\text{sec}$

85

Repeat 84 once. Dry the AMPure XP beads for  00:05:00  Room temperature

5m

Adapter 2 tagging

12m

86

Resuspend the beads in  46.5 μL **Adapter 2 mix**

87

Incubate for  00:10:00  Room temperature

10m

88

Transfer samples to a new plate

89

Heat mixture to  95 $^{\circ}\text{C}$ for  00:00:45 then immediately cool on ice using an aluminum rack

45s

90

Spin down at 500g for  00:00:10 at 15–25 $^{\circ}\text{C}$

10s

91

Add  1 μL **Klenow exo-** (50 U/ μL), vortex gently, and spin down at 500g for  00:00:10 at 15–25 $^{\circ}\text{C}$

10s

92

incubate as follows:

	A	B	C	D
	Segment	Temp ($^{\circ}\text{C}$)	Duration (min)	Ramp speed ($^{\circ}\text{C}/\text{min}$)
	1	4	5	-



	A	B	C	D
	2	4-37	8.25	4
	3	37	30	-
	4	4	Hold	

93

In the meantime, prepare 1 tube of **Adapter 2 mix**

	A	B	C	D
	Component	Amount (μl)	Final concentration	Mastermix (samples)
	Blue buffer (10×)	0.25	1x	
	dNTP mix (10mM)	0.1	0.4mM	
	Preamp Oligo (10 μM)	1	4 μM	
	Klenow exo- (50 U/μl)	0.5	10 U/μl	
	H2O	0.65		
	Total volume	2.5		

94

Heat the plate to 95 °C for 00:00:45 and transfer it to an aluminum rack pre-cooled on ice

45s

95

Centrifuge the plate at 500g for 00:00:10 at 15-25°C

10s

96

Add 2.5 μL of the **Adapter 2 mix**

97

Incubate as follows:



	A	B	C	D
	Segment	Temp (°C)	Duration (min)	Ramp speed (°C/min)
	1	4	5	-
	2	4-37	8.25	4
	3	37	90	-
	4	4	Hold	

Purification

21m 10s

98

Add  37.5 µL AMPure XP beads (**0.75:1 ratio**)

99

Incubate  00:10:00 at room temperature

10m

*In the meantime, prepare the **library amplification mix***

A	B	C	D
Component	Amount (µl)	Final Concentration	Mastermix (samples)
KAPA HIFI HotStart ReadyMix (2x)	25	1x	
PE1.0 (10µM)	1	0.2µM	
Nuclease-free water	23		
Total volume	49		

100

Place on a magnet for  00:03:00

3m

101

Place on a magnet for  00:03:00 and discard the supernatant

3m

102

Add  200 μ L ethanol (70%) without disturbing the beads. After  00:00:10 remove ethanol

10s

103

Repeat step 102 once then dry beads  00:05:00 at room temperature

5m

Library amplification

10m

104

Resuspend the beads in  49 μ L **library amplification mix**

105

Incubate the mixture  00:10:00  Room temperature

10m

106

Place on a magnet and transfer supernatant to a new plate

107

Add  1 μ L **10 μ M reverse iPCRTag primer**(containing a sample-specific index)

Amplify as follows:

	A	B	C
	Cycles	Temperature (°C)	Time
	1	95	3 min
	17-20	98	80 s
		65	30 s
		72	30 s
	1	72	3 min
	1	4	Hold

**Note**

The PCR product can be stored ON at 4°C or for at least a month at -20°C

Purification of amplified libraries**28m**

108

Add 37.5µl AMPure XP beads (0.75:1 ratio) and mix well

109

Incubate  00:10:00  Room temperature

10m

110

Place on the magnet for  00:03:00 and discard supernatant

3m

111

Add  200 µL ethanol (70%) without removing the plate from the magnet then discard the ethanol

112

Repeat step 111 once then dry beads  00:05:00

5m

113

Resuspend the beads in  15 µL **EB**

114

Incubate  00:10:00  Room temperature

10m

115

Place on a magnet then transfer supernatant to a new plate

Note

Library quantity and quality can be checked using Qubit HS Assay and Bioanalyzer
expected gDNA concentration: $\geq 1\text{ ng}/\mu\text{l}$
expected fragment length: $>200\text{bp}$ and on average $400\text{-}600\text{bp}$

Libraries can be stored for at least a year at -20°C

