

Jun 20, 2025

Screening compounds that inhibit SARS-CoV-2 nsp3 macrodomain: ADP-Ribose quantification by LC/MS-MS

DOI

<https://dx.doi.org/10.17504/protocols.io.8epv5rpydg1b/v1>

Noa Lahav^{1,2}, Haim Barr^{3,2}

¹Weizmann Institute; ²ASAP Discovery Consortium; ³Weizmann Institute of Science

ASAP Discovery



Noa Lahav

Weizmann Institute

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5rpydg1b/v1>

Protocol Citation: Noa Lahav, Haim Barr 2025. Screening compounds that inhibit SARS-CoV-2 nsp3 macrodomain: ADP-Ribose quantification by LC/MS-MS . protocols.io <https://dx.doi.org/10.17504/protocols.io.8epv5rpydg1b/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: September 16, 2024

Last Modified: June 20, 2025

Protocol Integer ID: 107677

Keywords: SARS-CoV-2, screening, inhibition, Mac1, Coronaviridae, ADPr, LC/MS-MS, nsp3, macrodomain, treatments for viral infectious disease, viral infectious disease, presence of potential inhibitor, potential inhibitor, inhibitor, mass spectrometry, relative potency of inhibitor, functional biochemical assay, ribose quantification by lc, mac1 ability, sar, screening compound, liquid chromatography, utilizing liquid chromatography

Funders Acknowledgements:

National Institutes of Health/National Institute Of Allergy and Infectious Diseases (NIH/NIAID)

Grant ID: U19AI171399

Disclaimer

The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institutes of Health.

Abstract

Functional biochemical assay to identify treatments for viral infectious disease that targets SARS-CoV-2 nsp3 macrodomain (MAC1). Utilizing Liquid Chromatography-Mass Spectrometry/Mass Spectrometry (LC-MS/MS) analysis, MAC1 ability to hydrolyze α NAD⁺ in the presence of potential inhibitors was measured by quantifying ADPR and α NAD⁺ levels. This assay is characterized by the relative potency of inhibitors as IC₅₀, it is modified from <https://doi.org/10.1126/sciadv.abf8711>. In this version, we have added the Addgene id.

Materials

Construct used

Addgene: Sars Cov 2 Mac domain  Sars Cov 2 Mac domain Plasmid **addgene Catalog #204787**

Assay Reagents (Concentration listed are from Stock Solutions)

1. **[M] 250 millimolar (mM)**  HEPES 0.5M buffer soln. pH 7.0 **Fisher Scientific Catalog #AAJ60064AE** (or similar)
2. **[M] 200 millimolar (mM)**  Sodium Chloride **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S9888**
3. **[M] 0.5 mg/mL**  Bovine Serum Albumin **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A7030**
4.  α -Nicotinamide adenine dinucleotide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #N6754**
5.  Adenosine 5'-diphosphoribose sodium salt **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A0752**

Protocol materials

 α -Nicotinamide adenine dinucleotide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #N6754**

 Sodium Chloride **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S9888**

 Bovine Serum Albumin **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A7030**

 Adenosine 5'-diphosphoribose sodium salt **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A0752**

 Sars Cov 2 Mac domain Plasmid **addgene Catalog #204787**

 HEPES 0.5M buffer soln. pH 7.0 **Fisher Scientific Catalog #AAJ60064AE**

Safety warnings

 Please be sure to wear proper Personal Protective Equipment (PPE) while performing this experiment.

Before start

Note: Inhibitor compounds stock concentration is 20 mM. Compounds are pre-dispensed into assay plates and stored at -20°C until use.



SARS-CoV-2 nsp3 expression and purification

1 The protein used in this assay was expressed and purified according to this protocol:

We used the following plasmid in the protein expression and purification protocol.

 Sars Cov 2 Mac domain Plasmid **addgene Catalog #204787**

Protocol

NAME

SARS-CoV-2 nsp3 macrodomain His-tagged expression and purification protocol for assay

CREATED BY

korvus.wang

Preview

Assay reagents preparation

2 Assay buffer:

	A	B	D	E
		Stock	Final in assay plate	Units
	HEPES pH=7.0	1000	25	mM
	NaCl	5000	100	mM
	BSA	10	1	mg/ml

	A	B	C	D	E
	Reagent	Stock	dispensed to plate	Final in assay plate	Units
	His-SARS-CoV-2 MAC1	183000	800	200	nM



	A	B	C	D	E
	Substrate (NAD ⁺)	Dry (dissolve fresh for each assay)	32	8	μM

Prepare 96 well-plate

1m

- 3 **PRIME** the dispenser with **Assay Buffer** by selecting the **PRIME** button until the tubes are filled completely.
- 3.1 **DISPENSE**  75 μL Mac1 assay Buffer to control wells of assay plate and  50 μL Mac1 Buffer to the rest of the plate.
 - **Note:** Control wells are for **NAD⁺ alone** and **MAC1+NAD⁺**, without compounds
- 3.2 **EMPTY** the dispenser tubes.
 - Discard the liquid discharged from the tubes.
- 4 **PRIME** the dispenser with **4x MAC1 enzyme** (800 nM) by selecting the **PRIME** button until the tubes are filled completely.
- 4.1 **DISPENSE**  25 μL 4xMac1 to all wells of assay plate besides **NAD⁺ alone** control well.
- 4.2 **EMPTY** the dispenser tubes.
 - Discard the liquid discharged from the tubes.
- 5 **CENTRIFUGE** plate  1500 rpm, Room temperature, 00:01:00 to remove bubbles
- 6 **INCUBATE** plate for  00:15:00 at  Room temperature
- 7 **PRIME** the dispenser with **4x NAD⁺** (32 μM) by selecting the **PRIME** button until the tubes are filled completely.
- 7.1 **DISPENSE**  25 μL 4xNAD⁺ to all wells of assay plate.

1m

15m



7.2 **EMPTY** the dispenser tubes.

- Discard the liquid discharged from the tubes.

8 **CENTRIFUGE** plate  1500 rpm, Room temperature, 00:01:00 to remove bubbles

1m

9 **INCUBATE** plate on an Orbi-plate shaker at  300 rpm for  02:45:00 at
 Room temperature

2h 45m

10 **STOP** the reaction by freezing the plate and store in -80 until LC/MS analysis.

LC/MS-MS analysis

- 11 Quantitative analysis of ADPR and α NAD⁺ was performed as was recently published at Lee, Alpha A et al. "Discovery of potent SARS-CoV-2 nsp3 macrodomain inhibitors uncovers lack of translation to cellular antiviral response." *bioRxiv : the preprint server for biology* 2024.08.19.608619. 21 Aug. 2024, doi:10.1101/2024.08.19.608619. Preprint. <https://pmc.ncbi.nlm.nih.gov/articles/PMC11370477/>

Protocol references

<https://doi.org/10.1126/sciadv.abf8711>