



Aug 21, 2025

Screening for L1014F mutation in *B. rufimanus*

DOI

<https://dx.doi.org/10.17504/protocols.io.kxygwx9kwv8j/v1>

Antoine Pichon¹, Louise McNamara¹, Sheila Alves¹, Stephen Byrne¹

¹Teagasc



Stephen Byrne

Teagasc

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygwx9kwv8j/v1>

Protocol Citation: Antoine Pichon, Louise McNamara, Sheila Alves, Stephen Byrne 2025. Screening for L1014F mutation in *B. rufimanus*. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kxygwx9kwv8j/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 03, 2024



Last Modified: August 21, 2025

Protocol Integer ID: 113480

Keywords: L1014F, bruchid beetle, non-destructive DNA isolation, absence of the L1014f mutation, L1014f mutation, insecticide resistance in other coleoptera species, broad bean bruchid beetle, insecticide resistance, coleoptera, gated sodium channel gene, other coleoptera species, sodium channel gene, beetle, insect species, chrysomelidae

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

This protocol describes the process of screening the broad bean bruchid beetle, *Bruchus rufimanus* (Coleoptera; Chrysomelidae) for the L1014F mutation that has been associated with insecticide resistance in other Coleoptera species (and more widely across insect species). The starting point is a quick and high-throughput non-destructive DNA isolation from beetles, followed by amplification of a region of the *Voltage-gated sodium channel* gene (VGSC) that harbours the L1014F mutation. The resulting PCR products can be sent for Sanger sequencing to confirm presence/absence of the L1014F mutation.

Materials

QuickExtract™ DNA Extraction Solution (Biosearch Technologies, Germany)

Classic++ Taq DNA Polymerase (Tonobo)

3-bruchid-L1014F (forward and reverse primers)

Safety warnings

⚠ Please read SDS associated with various consumables and kits used in this protocol and wear appropriate PPE. A site specific procedural risk assessment should be carried out prior to introducing this protocol to the lab.



Before start

This protocol assumes that you are working with *B. rufimanus* adults. The primers were designed with *B. rufimanus* sequence information and may not result in successful amplification in other *Bruchus* species. The beetles used in this protocol had been stored in 99 percent ethanol in the freezer prior to DNA extraction.



DNA Isolation

- 1 Individual beetles were placed in 1.5ml micro-centrifuge tubes with 100µl of QuickExtract DNA Extraction Solution.
- 2 These samples were vortexed for  00:00:15 using a benchtop vortex.
- 3 Samples were then placed in a ThermoMixer (Eppendorf) at  65 °C and  1000 rpm for  00:15:00 .
- 4 Samples were transferred to a heating block at  98 °C for  00:02:00 .
- 5 The samples were then placed on ice and the beetles removed and placed back into the 99 percent ethanol for longer term storage.
- 6 The 1.5ml micro-centrifuge tubes minus the beetles were briefly centrifuged using a bench-top centrifuge and the extraction solution (with DNA) was transferred to a PCR plate (semi-skirt).

PCR Amplification

- 7 A number of primers were tested and the following pair produced a single band at the correct size of ca. 500 bp.

	A	B
	Primer Name	Primer Sequence
	3-bruchid-L1014F-for	5'- CGGACCATGAACTTCCTAGGTGG-3'
	3-bruchid-L1014F-rev	5'-GTGTGTAGCAAACCGAGGAGCA-3'

- 8 Prepare a master-mix for the number of samples to be run so that each well will contain:
 -  12.5 µL of Classic++ Taq DNA Polymerase

-  1 μ L of 3-bruchid-L1014F-for
-  1 μ L of 3-bruchid-L1014F-rev
-  9.5 μ L of molecular biology grade water

9 Pipette  1 μ L of DNA (from step 6) into each well of a PCR plate and keep on ice

10 Pipette  24 μ L of master-mix (from step 8) into each well, seal the plate carefully with foil lid, and centrifuge the PCR plate briefly to ensure no droplets adhere to sides of sample wells.

11 Place the PCR plate in a PCR machine (e.g MiniAmp Plus Thermal Cycler) and run the following PCR cycle:

-  95 °C  00:01:00

followed by 35 cycles of

-  95 °C  00:00:15
-  65 °C  00:00:15
-  72 °C  00:00:15

and a final elongation step of

-  72 °C  00:03:00

12 Check a selection of samples by running on a 1.2 percent TBE agarose gel alongside a 100 bp ladder. If using GelRed carry out post-run staining to avoid incorrect sizing of PCR products. A band should be observed at ca. 500 bp.

Sanger Sequencing and Analysis

13 Ship plate of PCR products to Sanger Sequencing core facility or service provider along with an aliquot of the **3-bruchid-L1014F-for** primer for sequencing. Sequencing from the forward primer is sufficient to extract clean data covering the region of the L1014F mutation.



- 14 The sequence files can be visualised in Geneious Prime or FinchTV. Poor quality regions at start and end of sequences can be trimmed off and each sequence can be annotated with the longest open reading frame and sequences aligned (with traces) to identify the L1014F mutation.