

Feb 25, 2019 Version 1

## scTHS-seq V.1

DOI

[dx.doi.org/10.17504/protocols.io.xf8fjrw](https://dx.doi.org/10.17504/protocols.io.xf8fjrw)

Kun Zhang<sup>1</sup>, Brandon Sos<sup>1</sup>, Dinh H Diep<sup>1</sup>, Thu E Duong<sup>1</sup>

<sup>1</sup>University of California, San Diego



Dinh H Diep

University of California, San Diego

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.xf8fjrw](https://dx.doi.org/10.17504/protocols.io.xf8fjrw)

**Protocol Citation:** Kun Zhang, Brandon Sos, Dinh H Diep, Thu E Duong 2019. scTHS-seq . **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.xf8fjrw>

**Manuscript citation:**

<https://www.nature.com/articles/nbt.4038>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development

**We are still developing and optimizing this protocol**

**Created:** January 24, 2019

**Last Modified:** February 25, 2019

**Protocol Integer ID:** 19680



## Abstract

Single cell transposome hypersensitive site sequencing (scTHS-seq) combined transposome hypersensitive site sequencing (THS-seq) with combinatorial cellular indexing using customized barcoded transposomes. Compared to a similar method, ATAC-seq, THS-seq have improved coverage of distal enhancer sites which is due to the utilization of an engineered super-mutant Tn5 transposase as well as in vitro transcription to perform linear amplification from a single insertion site of the Tn5 transposase instead of two insertion sites required in ATAC-seq. scTHS-seq was used to generate single nucleus accessibility profile of 36,869 cells with an average of 10,168 unique reads per cell that are also associated with a predicted accessible region. The doublet rate of scTHS-seq is estimated to be between 20-25% which is comparable to similar combinatorial indexing methods.

## Materials

### MATERIALS

- ⊗ HiScribe T7 High Yield RNA Synthesis Kit - 50 rxns **New England Biolabs Catalog #E2040S**
- ⊗ Taq 5X Master Mix - 500 rxns **New England Biolabs Catalog #M0285L**
- ⊗ Pierce™ Dimethylformamide (DMF), Sequencing grade **Thermo Fisher Scientific Catalog #20672**
- ⊗ 2x Kapa HiFi Hotstart Readymix **Kapa Biosystems Catalog #KK2602**
- ⊗ Tango Buffer **Thermo Fisher Scientific Catalog #BY5**
- ⊗ Ultrapure SMART MMLV Reverse Transcriptase for RT-PCR **Takara Bio Inc. Catalog #639524**
- ⊗ Advantage® UltraPure dNTP Combination Kit (100 mM each dNTP) **Takara Bio Inc. Catalog #639125**
- ⊗ RNase H **Enzymatics Catalog #Y9220L**
- ⊗ KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**
- ⊗ twin.tec PCR Plate 96 LoBind semi-shirted clear 25 pcs. **Eppendorf Catalog #30129504**

### STEP MATERIALS

- ⊗ KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**
- ⊗ KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**
- ⊗ KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**
- ⊗ KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**
- ⊗ KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**
- ⊗ Taq 5X Master Mix - 500 rxns **New England Biolabs Catalog #M0285L**
- ⊗ KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**
- ⊗ KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**
- ⊗ KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**
- ⊗ KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**

For purification of engineered super-mutant Tn5-059, see previous publications (<https://genomebiology.biomedcentral.com/articles/10.1186/s13059-016-0882-7>) for details.



## Protocol materials

☒ KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**

☒ KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**

☒ Pierce™ Dimethylformamide (DMF), Sequencing grade **Thermo Fisher Scientific Catalog #20672**

☒ Tango Buffer **Thermo Fisher Scientific Catalog #BY5**

☒ RNase H **Enzymatics Catalog #Y9220L**

☒ twin.tec PCR Plate 96 LoBind semi-shirted clear 25 pcs. **Eppendorf Catalog #30129504**

☒ HiScribe T7 High Yield RNA Synthesis Kit - 50 rxns **New England Biolabs Catalog #E2040S**

☒ KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**

☒ KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**

☒ KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**

☒ Advantage® UltraPure dNTP Combination Kit (100 mM each dNTP) **Takara Bio Inc. Catalog #639125**

☒ KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**

☒ Taq 5X Master Mix - 500 rxns **New England Biolabs Catalog #M0285L**

☒ Ultrapure SMART MMLV Reverse Transcriptase for RT-PCR **Takara Bio Inc. Catalog #639524**

☒ KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**

☒ KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**

☒ 2x Kapa HiFi Hotstart Readymix **Kapa Biosystems Catalog #KK2602**

☒ KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**

☒ Taq 5X Master Mix - 500 rxns **New England Biolabs Catalog #M0285L**

☒ KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**

☒ KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**

☒ Taq 5X Master Mix - 500 rxns **New England Biolabs Catalog #M0285L**

☒ KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**

☒ KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**

☒ KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**

☒ KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**

## Tagmentation

- 1 Make 10 mL 75% glycerol with 0.3% tritonX-100

🧴 9.45 g Glycerol (100%)

🧴 300  $\mu$ L TritonX-100, 10%

🧴 2200  $\mu$ L nuclease free water

- 2 Make 3.045 mL transposons dilution buffer

🧴 2345  $\mu$ L 75% glycerol/0.3% TritonX-100

🧴 700  $\mu$ L Qiagen EB

- 3 Make 1000  $\mu$ L transposase dilution buffer

🧴 667  $\mu$ L 75% glycerol/0.3% TritonX-100

🧴 1  $\mu$ L DTT, 1M

🧴 6.25  $\mu$ L Tris-HCl, 4M, pH 7.5

🧴 325.75  $\mu$ L nuclease free water

- 4 Make 10 mL of 1.5X D6 buffer

🧴 1500  $\mu$ L Tango Buffer, 10X

🧴 2400  $\mu$ L DMF, 100%

🧴 6100  $\mu$ L nuclease free water

- 5 Make 10 mL of 1X lysis buffer

🧴 25  $\mu$ L Tris-HCl, 4M, pH 7.5

🧴 20  $\mu$ L Sodium chloride, 5M

🧴 15  $\mu$ L Magnesium chloride, 2M

🧴 100  $\mu$ L IGEPAL-630, 10%



🧪 2000  $\mu$ L BSA, 10%

🧪 1000  $\mu$ L PBS, 10X

🧪 6840  $\mu$ L nuclease free water

Finally, add 1 protease inhibitor tablet (EDTA-free).

## 6 Make 10 mL of 2X FACS buffer

🧪 1000  $\mu$ L PBS, 10X

🧪 2000  $\mu$ L BSA, 10%

🧪 80  $\mu$ L EDTA, 0.5M

🧪 6920  $\mu$ L nuclease free water

## 7 Generate 384 annealed transposons

🧪 7  $\mu$ L Indexed Tn5 top strand (100uM)

🧪 7  $\mu$ L 5' phosphorylated mosaic end sequence bottom strand (100uM)

Combine top and bottom strand oligos in a PCR tube. In a thermocycler, incubate for 2 minutes at 🌡️ 95 °C then ramp to 🌡️ 14 °C at a rate of 🌡️ -0.1 °C per second.

## 8 Dilute annealed transposons

🧪 1  $\mu$ L Annealed transposons

🧪 6.7  $\mu$ L transposon dilution buffer

## 9 Make 420 uL of 4.8 uM transposase Tn5-059

🧪 252  $\mu$ L transposase, Tn5-059, 8 uM



🧪 168  $\mu$ L transposase dilution buffer

## 10 Load Tn5-059

🧪 1  $\mu$ L 4.8 uM transposase Tn5-059

🧪 1  $\mu$ L diluted annealed transposons (6.5 uM)

11 Incubate

 22 °C Thermocycler  00:30:00

12 Wash nuclei or cells suspension with ice-cold PBS

13 Count cells

14 Pellet 400,000 cells / nuclei by centrifuging at 330 RCF (cells) or 900 RCF (nuclei)

15 Resuspend 400,000 cells/nuclei in lysis buffer. Use 100 uL per 100,000 cells/nuclei

16 Incubate



 4 °C On ice  00:05:00

17 Pellet by centrifuging at 900 RCF

18 Resuspend nuclei in 100 uL 1.5X D6 buffer initially. Count and add more 1.5X D6 buffer to create ~250 nuclei/uL

19 Aliquot 4.0 uL of cells in 1.5X D6 buffer into wells containing 2.0 uL of prepared transposomes

20 Incubation without shaking

 37 °C Thermocycler  00:30:00

21 Terminate the tagmentation reaction with EDTA

 4 µL EDTA, 50 mM



22 Incubate

37 °C Thermocycler 00:15:00

23 Store cells at -20C overnight or until ready to sort

## FACS (BD Influx)

24 Coat 12 × 200 uL, 1 × 15 mL conical using 1X FACS buffer

01:00:00 Let tubes sit in FACS buffer

25 Add 10 uL PBS to each well of 10 to 15 × 96-wells plate for sorting

### Note

Cover and store all plates on ice.

26 Add 1 volume (10 uL) ice cold 2X FACS buffer to each sample well with cells in tagmentation mix

27 When 15 mL conical have been coated for 1 hr, pool all samples together using multichannel pipette and adding to 12 × 200 uL tubes

28 Centrifuge at 500 RCF 4 °C

29 Resuspend nuclei in 1.5 mL 1X FACS buffer

30 Add stain for sorting

75 uL Propidium Iodide


31 Sort 100 cells into each well using BD Influx, remove doublets using forward and side scatter plots and selecting for PI-staining events





## In vitro transcription (10 plates of sorted nuclei)

- 32 Add 11 uL of 8M guanidine chloride to each well using multi-channel to lyse nuclei and release DNA


 14 mL Guanidine hydrochloride, 8M

- 33 Add 30 uL of KAPA Pure Beads to each well to bind DNA


 42.24 mL KAPA Pure Beads

 KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**

- 34 Mix well and bind for at least 8 minutes

 00:08:00 Room temperature

- 35 Place plate on magnet for at least 3 minutes then discard the supernatant

 00:03:00 Room temperature


- 36 Prepare fresh 80% ethanol solution. Wash beads (while still on magnet) twice with 170 uL of 80% ethanol

 384 mL Freshly prepared 80% ethanol

- 37 After last wash, completely remove all ethanol and transfer to a heated thermomixer to dry

 37 °C  00:01:00

- 38 Add 8 uL of nuclease free water to each well using a multichannel pipette

 8640 µL nuclease-free water

- 39 Add 2 uL of 5X Taq Master Mix to perform gap filling



🧪 2160  $\mu$ L 5X Taq Master Mix

🔗 Taq 5X Master Mix - 500 rxns **New England Biolabs Catalog #M0285L**

40 Incubate plates for gap filling

🌡️ 72 °C ⌚ 00:08:00

41 Add 2  $\mu$ L of HiScribe T7 IVT reaction components to each well

🧪 2112  $\mu$ L T7 buffer, 10X

🧪 2112  $\mu$ L UTP

🧪 2112  $\mu$ L GTP

🧪 2112  $\mu$ L ATP

🧪 2112  $\mu$ L CTP

🧪 2112  $\mu$ L T7 RNA Polymerase

42 Incubate

🌡️ 37 °C Thermocycler ⌚ 16:00:00

43 Add 40  $\mu$ L of KAPA Pure Beads to each well to bind DNA

🧪 42.24 mL KAPA Pure Beads

🔗 KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**

44 Mix well and bind for at least 8 minutes

⌚ 00:08:00 Room temperature

45 Place plate on magnet for at least 3 minutes then discard the supernatant

⌚ 00:03:00 Room temperature



- 46 Prepare fresh 80% ethanol solution. Wash beads (while still on magnet) twice with 170 uL of 80% ethanol

🧪 384 mL Freshly prepared 80% ethanol

- 47 After last wash, completely remove all ethanol and transfer to a heated thermomixer to dry

🌡️ 37 °C ⌚ 00:01:00

- 48 Add 9 uL of nuclease free water to each well using a multichannel pipette

🧪 9720 µL nuclease-free water

## cDNA synthesis (10 plates of sorted nuclei)

- 49 Add 2.5 uL of random hexamers

🧪 3000 µL random hexamers, 20 uM

- 50 Incubate

🌡️ 70 °C Thermocycler ⌚ 00:03:00

### Note

Cool immediately on ice after.

- 51 Add 9.5 uL of MMLV RT mixture

🧪 4224 µL First strand buffer, 5X

🧪 2112 µL Advantage UltraPure dNTPs, 100 uM each

🧪 2112 µL DTT, 100mM

🧪 528 µL MMLV RT

**52** Incubate in a thermocycler

22 °C 00:10:00

42 °C 01:00:00

70 °C 00:10:00

**53** Add 1 uL RNase H

1127.52 µL nuclease-free water

125.28 µL RNase H, 5 Units/uL

**54** Incubate in thermocycler

37 °C 00:20:00

**55** Add 2.5 uL second strand synthesis primer

3000 µL sss\_scnXTv2, 20 uM

**56** Incubate in thermocycler

65 °C 00:02:00

**Note**

Cool immediately on ice after incubation

**57** Add 5.9 uL 5X Taq Master Mix

7363 µL Taq Master Mix, 5X



58 Incubate in thermocycler

72 °C 00:08:00

## Tagmentation on beads (10 plates of sorted nuclei)

59 Add 40 uL of KAPA Pure Beads to each well to bind DNA

42.24 mL KAPA Pure Beads

KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**

60 Mix well and bind for at least 8 minutes

00:08:00 Room temperature

61 Place plate on magnet for at least 3 minutes then discard the supernatant

00:03:00 Room temperature

62 Prepare fresh 80% ethanol solution. Wash beads (while still on magnet) twice with 170 uL of 80% ethanol

384 mL Freshly prepared 80% ethanol

63 After last wash, completely remove all ethanol and transfer to a heated thermomixer to dry

37 °C 00:01:00

64 Add 4 uL of 1.5X tagmentation buffer to each well using a multichannel pipette

3172 µL nuclease-free water

780 µL Tango Buffer, 10X

1248 µL DMF, 100%

**65** Generate annealed transposons 40  $\mu\text{L}$  i7 top strand (100uM) 40  $\mu\text{L}$  5' phosphorylated mosaic end sequence bottom strand (100uM)

Combine top and bottom strand oligos in a PCR tube. In a thermocycler, incubate for 2 minutes at 95 °C then ramp to 14 °C at a rate of -0.1 °C per second.

**66** Dilute annealed transposons 74  $\mu\text{L}$  Annealed transposons 1026  $\mu\text{L}$  transposon dilution buffer**67** Make 420  $\mu\text{L}$  of 2.4  $\mu\text{M}$  transposase Tn5-059 330  $\mu\text{L}$  transposase, Tn5-059, 8  $\mu\text{M}$  770  $\mu\text{L}$  transposase dilution buffer**68** Load Tn5-059 1100  $\mu\text{L}$  2.4  $\mu\text{M}$  transposase Tn5-059 1100  $\mu\text{L}$  diluted annealed transposons (6.5  $\mu\text{M}$ )**69** Incubate

22 °C Thermocycler 00:30:00

**70** Add 2  $\mu\text{L}$  of loaded transposomes to each well**71** Incubate

55 °C Thermocycler 00:06:00

**Note**


Cool immediately on ice after.

72 Add 6 uL guanidine hydrochloride

 7488 µL Guanidine hydrochloride, 8M

73 Add 30 uL of KAPA Pure Beads to each well to bind DNA


 31.68 mL KAPA Pure Beads

 KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**

74 Mix well and bind for at least 8 minutes

 00:08:00 Room temperature

75 Place plate on magnet for at least 3 minutes then discard the supernatant

 00:03:00 Room temperature

76 Prepare fresh 80% ethanol solution. Wash beads (while still on magnet) twice with 170 uL of 80% ethanol

 384 mL Freshly prepared 80% ethanol

77 After last wash, completely remove all ethanol and transfer to a heated thermomixer to dry

 37 °C  00:01:00

78 Add 16 uL nuclease free water



🧪 21.6 mL nuclease-free water

## Test amplification (10 plates of sorted nuclei)

79 Preheat hotstart polymerase mix

🧪 250  $\mu$ L KAPA HiFi Hotstart ReadyMix, 2X

🌡️ 98 °C Thermocycler ⌚ 00:00:45

### Note

Cool immediately on ice

80 Add primers to make a master mix to 417 nM final each

🧪 42  $\mu$ L i5XX, 5 uM

🧪 42  $\mu$ L i7XX, 5 uM

81 Add SYBR Green for qPCR to 0.8X final

🧪 4  $\mu$ L SYBR Green, 100X

82 Aliquot 8 uL of amplification mix to each well to qPCR plate

83 Sample 5 wells from each plate, adding 4 uL of each sample. Include an no nuclei control well.

84 Perform qPCR. Monitor the sample progress, terminate when all samples appears to be in exponential phase (usually 6-8 cycles). No nuclei control may appear to amplify but majority of those products are primer dimers.

1. 72C x 5 min








2. 98C x 30 s
3. 59C x 30 s
4. 72C x 1 min
- Plate read
5. Go to 2 for 20 cycles
6. 72C for 5 min

- 85 Assess PCR products on a 6% TBE PAGE gel. If lack of products, repeat test amplification on a different sampling of wells but adding 3-4 more cycles of amplification. If samples appear to be over-amplified, reduce the number of cycles by 1.

## Library amplification (10 plates of sorted nuclei)

- 86 Preheat hotstart polymerase mix

 12000 µL KAPA HiFi Hotstart ReadyMix, 2X


 98 °C Thermocycler  00:00:45

### Note


Cool immediately on ice

- 87 Aliquot 12 uL of preheated KAPA HiFi ReadyMix to each well of a new PCR plate








- 88 Aliquot 2 uL of 16 different s5XX primers to each set of 2 plates (one per row)

 150 µL each s5XX primer, 5 uM

- 89 Aliquot 2 uL of 12 different s7XX primers to each column of every plate (one per column)

 100 µL s7XX primer, 5 uM

- 90 Transfer 8 uL of sample to each well of PCR plate

- 91 Perform PCR using the number of cycles determined previously with test PCR
1. 72C x 5 min
  2. 98C x 30 s
  3. 59C x 30 s
  4. 72C x 1 min
  5. Go to 2 for XX cycles
  6. 72C for 5 min
- 92 Pool 4 uL of each well from 2 plates with non-overlapping s5XX primers. Final volume is ~800 uL per pool. 5 pools total
- 93 Add 800 uL of KAPA Pure Beads to each pool to bind DNA
-  3 mL KAPA Pure Beads
-  KAPA Pure Beads (60 mL) **Roche Catalog #07983298001**
- 94 Mix well and bind for at least 8 minutes
-  00:08:00 Room temperature
- 95 Place plate on magnet for at least 3 minutes then discard the supernatant
-  00:03:00 Room temperature
- 96 Prepare fresh 80% ethanol solution. Wash beads (while still on magnet) twice with 1000 uL of 80% ethanol
-  12 mL Freshly prepared 80% ethanol
- 97 After last wash, completely remove all ethanol and transfer to a heated thermomixer to dry
-  37 °C  00:01:00
- 98 Add 20 uL of nuclease free water to each pool



🧪 100 µL nuclease-free water

99 Quantify each pool using dsDNA HS Qubit

100 Perform gel size selection for 220-1000 bp using 6% TBE PAGE. Use 5-wells gel, and loading ~1 ug of library per well.