

Sep 23, 2023

SDS-PAGE and western blot analysis

DOI

<https://dx.doi.org/10.17504/protocols.io.eq2lyj33plx9/v1>

Elias Adriaenssens¹

¹Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2lyj33plx9/v1>

Protocol Citation: Elias Adriaenssens 2023. SDS-PAGE and western blot analysis. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.eq2lyj33plx9/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 28, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 84147

Keywords: SDS-PAGE, western blot analysis, ASAPCRN, western blot analysis this protocol, western blot analysis, sd, page, analysis, protocol

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-000350

Marie Skłodowska-Curie MSCA Postdoctoral fellowship

Grant ID: 101062916

Abstract

This protocol describes SDS-PAGE and western blot analysis.

Attachments



[758-1926.pdf](#)

47KB

Materials

Materials

- 100 mM DTT
- PageRuler Prestained protein marker (Thermo Fisher)
- nitrocellulose membranes (RPN132D, GE Healthcare)
- Mini Trans-Blot Cell (Bio-Rad)
- 0.1% Tween 20
- horseradish peroxidase (HRP)-coupled antibodies

RIPA buffer

	A	B
	Tris-HCl pH 8.0	50 mM
	NaCl	150 mM
	sodium deoxycholate	0.50%
	SDS	0.10%
	NP-40	1%

 cOmplete™, Mini, EDTA-free (Protease Inhibitor) **Roche Catalog ##11836170001**

 Roche PhosSTOP™ **Merck MilliporeSigma (Sigma-Aldrich) Catalog #4906837001**

 Pierce™ Detergent Compatible Bradford Assay Kit **Thermo Fisher Catalog #23246**

 SuperSignal™ West Femto Maximum Sensitivity Substrate **Thermo Fisher Catalog #34096**



SDS-PAGE and western blot analysis

10m

- 1 For SDS-PAGE and western blot analysis, collect cells by trypsinization and subsequent centrifugation at 300 x g, 4°C, 00:05:00 . 5m
- 2 Wash the cell pellets in PBS and centrifuge once more at 300 x g, 4°C, 00:05:00 . 5m
- 3 Remove the supernatant and lyse the cell pellets in RIPA buffer supplemented by cOmplete EDTA-free protease inhibitors (11836170001, Roche) and phosphatase inhibitors (Phospho-STOP, 4906837001, Roche).
- 4 After incubating in RIPA buffer for 00:20:00 On ice , clear the samples by centrifugation at 20000 x g, 4°C, 00:10:00 . 30m
- 5 Collect the soluble supernatant fraction and measure the protein concentrations using the Pierce Detergent Compatible Bradford Assay Kit (23246, Thermo Fisher).
- 6 Then, adjust the samples for equal loading and mix it with 6x protein loading dye, supplemented with 100 mM DTT. Then, boil it for 00:05:00 at 95 °C . 5m
- 7 Load the samples on 4-12% SDS-PAGE gels (NP0321BOX, NP0322BOX, or NP0323BOX, Thermo Fisher) with PageRuler Prestained protein marker (Thermo Fisher).
- 8 Transfer the proteins onto nitrocellulose membranes (RPN132D, GE Healthcare) for 01:00:00 at 4 °C using the Mini Trans-Blot Cell (Bio-Rad). 1h
- 9 After the transfer, block the membranes with 5% milk powder dissolved in PBS-Tween (0.1% Tween 20) for 01:00:00 at Room temperature . 1h
- 10 Incubate the membranes Overnight at 4 °C with primary antibodies dissolved in the blocking buffer and wash. 1h
- 10.1 Wash for 00:05:00 .(1/3) 5m



10.2 Wash for  00:05:00 .(2/3)

5m

10.3 Wash for  00:05:00 .(3/3)

5m

10.4 Then, incubate with species-matched secondary horseradish peroxidase (HRP)-coupled antibodies diluted 1:10,000 in blocking buffer for  01:00:00 at  Room temperature .

1h

11 Afterwards wash the membranes three times with PBS-T and further process for western blot detection.



12 Incubate the membranes with SuperSignal West Femto Maximum Sensitivity Substrate (34096, Thermo Fisher) and image it with a ChemiDoc MP system (Bio-Rad).



13 Analyze the images with ImageJ [57].

