

Oct 09, 2022

SDS-PAGE


DOI

dx.doi.org/10.17504/protocols.io.kqdg3pm11l25/v1



Ana Belem García González¹, Georgina Diego¹, Irán Alessandra Chaparro Rodríguez¹, Jair Alexis Gardea Sáenz¹

¹Tecnologico de Monterrey Campus Chihuahua

 Denisse Chacón Ramírez

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg3pm11l25/v1>

Protocol Citation: Ana Belem García González, Georgina Diego, Irán Alessandra Chaparro Rodríguez, Jair Alexis Gardea Sáenz 2022. SDS-PAGE. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kqdg3pm11l25/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: August 09, 2022

Last Modified: October 09, 2022

Protocol Integer ID: 68368

Keywords: page this protocol, sd, team tec, protocol, page, chihuahua, step

Abstract

This protocol shows the steps carried out by team Tec-Chihuahua to perform SDS-PAGE

Materials

SOLUTIONS:

- **Acrylamide/Bis 30%**

- 29.2 g of Acrylamide.
- 0.8 g of Bis Acrylamide
- Dissolve in 50 mL of distilled water in constant agitation, gauge to a volume of 100 mL.
- Store at 4 °C degrees while protected from light.

- **Tris-HCl/ SDS 4X pH 8.8 (1.5 M Tris-HCl, 0.4% SDS)**

- 18.17 g of TRIZMA® base
- 0.4 g of SDS
- Dissolve in 80 mL of distilled water
- Adjust to a pH de 8.8 with HCl and gauge to a 100 mL
- Filtrate the solution
- Store at 4 °C degrees

- **Tris-HCl/ SDS 4X pH 6.8 (0.5 M Tris-HCl, 0.4% SDS)**

- 6.06 g de TRIZMA® base
- 0.4 g de SDS
- Dissolve en 80 mL of distilled water.
- Adjust to a pH de 6.8 with HCl and gauge to a 100 mL
- Filtrate the solution
- Store at 4 °C degrees

- **10% Ammonium Persulfate**

- 500 mg of APS (Ammonium Persulfate)
- Dissolve en 5 mL of distilled ultra pure water
- Store at -20 °C for up to two weeks
- Store at 4 °C degrees for one use only

- **2X Loading Buffer**

- 2 mL of glycerol
- 400 uL of Mercapto
- 0.02 g of bromophenol blue
- 0.4 g of SDS
- 2.5 mL of 8.6 pH buffer

- **1X Running Buffer**

- 14 g of glycine.
- 3 g of Trizma base.
- 1 g of SDS.



- Gauge to a 1 L.

- **Staining Solution**

- 1.25 g of brilliant blue R in 250 mL of methanol
- 200 mL of distilled water.
- 50 mL of Glacial Acetic Acid
- Store in room temperature while protected from light

- **Destaining Solution**

- 250 mL of methanol
- 62.5 mL of Glacial Acetic Acid
- 312.5 mL of distilled water.
- Store in room temperature

Usage: Pour solution into a container and allow to stir until the SDS gel is clear or electrophoresis bands are visible.

Troubleshooting


Protocol

3h 45m

- 1 Clean the components of the electrophoresis camera with 70% ethanol and gauze 's.
- 2 Assemble the chamber and check that there are no leaks by pouring distilled water between the glasses.
- 3 Prepare polyacrylamide gels:

	A	B
	Reagents	2 minigels
	Distilled water	3.4 mL
	Acrylamide/ Bis 30%	4 mL
	Tris- HCl/SDS 4X pH 8.8 (1.5 M Tris-HCl, 0.4% SDS)	2.5 mL
	ASP 10%	100 µL
	TEMED	4 µL

Separation gel (12%)

- 4 Pour the solution between the glasses with a 1 mL micropipette, leaving a space of 1.5 cm for the concentrating gel. To level, distilled water is added and allowed to settle for  00:30:00 or until a line is seen between the gel and the water.

30m




	A	B
	Reagents	2 minigels
	Distilled water	2.7 mL



A	B
Acrylamide/ Bis 30%	1 mL
Tris-HCl/ SDS 4X pH 8.8 (1.5 M Tris-HCl, 0.4% SDS)	1.3 mL
ASP 10%	50 µL
TEMED	4 µL

- Concentrating gel (6%)


5 Pour solution onto separating gel using a 1 mL micropipette.

6  00:30:00  00:30:00 Insert the comb (carefully avoiding the formation of bubbles) and leave to solidify for  00:30:00

1h 30m


7 When the gels are polymerized, prepare an electrophoresis chamber with 1X running buffer until it covers the gels and it reaches the line of two gels.



8 Sample PTake the pellets contained in Eppendorf tubes

- Take the pellets contained in Eppendorf tubes
- Add  300 µL of 1X loading buffer and resuspend the pellet.reparation:

9 To denature proteins, heat samples in boiling water for 5 min.



10 Load gels with the hot sample.

11 In the first well add  7 µL of the molecular weight marker.

12 Once the samples are loaded, run the gel at 80 V for  00:20:00 and then at 180 V for  00:45:00

1h 5m



- 13 Turn off the camera and disarm it.
- 14 Remove the gels from the glasses and place in a container with staining solution. Leave stirring for one hour.
- 15 Remove staining solution after one hour, add destaining solution and leave stirring for  00:20:00 Change the destaining solution and leave again for  00:20:00
- 16 Leave stirring until the gel is transparent.
- 17 Analyze the gels.

40m