



Jun 13, 2025

Version 1

🌐 Seasonal and zoonotic influenza WGS using ONT rapid barcoding technology V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl41onjlo5/v1>

Richard Chalvignac¹, Quentin Semanas¹, Hadrien Regue¹, Antonin Bal¹, Laurence Josset¹

¹Hospices Civils de Lyon - Genepii Sequencing Platform



Richard Chalvignac

Hospices Civils de Lyon

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl41onjlo5/v1>

Protocol Citation: Richard Chalvignac, Quentin Semanas, Hadrien Regue, Antonin Bal, Laurence Josset 2025. Seasonal and zoonotic influenza WGS using ONT rapid barcoding technology. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.3byl41onjlo5/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: June 11, 2025

Last Modified: June 13, 2025

Protocol Integer ID: 219941

Keywords: zoonotic rapid whole genome, zoonotic influenza wg, ont rapid barcoding technology this protocol, clade resolution, sequencing, tierce bioinformatic tool, oxford nanopore technology, sequencing result, genome, using ont, rapid barcoding technology, whole genome sequencing

Abstract

This protocol is suitable for seasonal and zoonotic rapid whole genome sequencing using universal primers for influenza A and Oxford Nanopore Technologies (ONT) rapid chemistry (SQK-RBK114.24 or SQK-RBK114.96). Using ONT and tierce bioinformatic tools, it enables sub-typing and clade resolution sequencing result in one working-day.



Nucleic acid extraction

1h

- 1 RNA extraction is performed using either Maxwell[®] CSC Pathogen Total Nucleic Acid Kit **AS1860** or EMAG[®] Nucleic Acid Extraction System Generic program, allowing the management of up to 16 or 24 samples simultaneously depending on the extraction system.

1h

Universal Influenza A Primers

- 2 Univ primers are used for the 8 segments amplification. Primers sequences are listed here :

Primers	5'-3' Sequences
Uni12 Forward (Fw)	GGG-GGG-AGC-AAA-AGC-AGG
Uni12b Forward (Fw) (Uni12_inf3)	GGG-GGG-AGC-GAA-AGC-AGG
Uni13 Reverse (Rv)	CGG-GTT-ATT-AGT-AGA-AAC-AAG-G

Table 1 : Universal Influenza A Primers

Pre-PCR and amplification

4h

- 3 Pre-PCR is prepared as follows :
 - ezDNase **11766051** (Thermofisher Scientific) is used for nucleic-acid extracts DNase treatment

Reagent	Volume / Sample (μL)
nuclease-free water	0.625
ezDNase buffer (10X)	1.25
ezDNase enzyme	0.625
nucleic acid extract	10
Total	12.5

Table 2 : Nucleic-acid extract DNase treatment

DNase treatment is performed by incubating the prepared mix at  37 °C for

 00:10:00

- 3.1 SuperScript III One Step RT-PCR System **12574035** (Thermofisher Scientific) combined with RNasin **N2511** (Promega) is then used as follows :

	Reagent	Volume / Sample (μL)
	2X Reaction Mix	12.5
	ezDNase - nucleic acid extract mix	5
	nuclease-free water	4.75
	Primer Uni12 forward (10μM)	0.5
	Primer Uni12b forward (10μM)	0.5
	Primer Uni13 reverse (10μM)	1
	RNasin	0.25
	SuperScript III RT/ Platinum Taq High Fidelity Enzyme Mix	0.5
	Total	20

Table 3 : Mix preparation for one-step RT-PCR

One-step RT-PCR is performed on a thermocycler as follows  03:30:00 :

	Temp (C°)	Time	Cycle(s)
	45	01:00:00	1
	94	00:02:00	
	94	00:00:30	5
	45	00:00:30	

	Temp (C°)	Time	Cycle(s)
	68	00:03:00	
	94	00:00:30	31
	57	00:00:30	
	68	00:03:00	

Table 4 : Thermocycler one-step RT-PCR program

Library preparation

- 4 The library construction is performed according to manufacturer's instruction for the Rapid Barcoding Kit (SQK-RBK114.24/96).
Before starting this step it's possible to perform a flow cell check on the ONT device to gain some time on the whole experiment.

Briefly, each sample is quantified using QuBit fluorometer and normalised to obtain

 50 ng of amplified DNA in a  10 µL volume.

 1.5 µL of rapid barcode are added to the reaction volume and the mix is then incubated as follows :

 30 °C	 00:02:00
 80 °C	 00:02:00

- 4.1 Once the incubation time is over, each sample is pooled in a single 1.5mL LoBind Eppendorf tube and a 1:1 AMPure XP beads purification is performed. Two 80% ethanol washes are performed on a magnetic rack. Elution is done in  11 µL of Elution Buffer (EB)

- 4.2 Rapid Adapter mix is performed as follows :

	Reagent	Volume (µL)
	Adapter Buffer (ADB)	3.5



	Reagent	Volume (µL)
	Rapid Adapter (RA)	1.5
	Total	5

Table 5 : Rapid Adapter Mix (unused mix can be stored at -20°C for further use)

 1 µL of Rapid Adapter mix is added to the  11 µL of previously barcoded and purified samples.

The resulting mix is incubated  00:05:00 at room temperature and then stocked at 4°C prior to loading on the flow cell

Sequencing and analysis

1h 30m

- Sequencing is performed on any ONT sequencing device. Preferentially use one with computing capacity allowing for real-time High Accuracy basecalling mode (HAC).

1h

R10.4.1 flow cell is loaded according to manufacturer's instruction and sequencing run is launched for at least  01:00:00 for small batches of samples (<8). The time required for successful sequencing is depending on samples multiplexing but also individual samples overall quality and viral load thus should be tested to ensure good results. Once the short run is over, it can be re-launched for a longer time to use as supplementary data if needed.

- 5.1 Analysis can be performed immediately using Epi2Me wf-flu (<https://github.com/epi2me-labs/wf-flu>).

30m

The RBK kit usage must be specified while setting-up the analysis run.

The wf-flu is not suited for zoonotic influenza type so the samples will always be Undetermined in the subtyping results of the report. If there is a suspicion for zoonotic influenza, an alternative solution can be used as described after.

Using intermediate files generated in the analysis directory, it's still possible to obtain the consensus for the samples of interest (wf-flu_ID\output\SampleID\consensus\draft.consensus.fasta).

The consensus can then be uploaded on Nextclade with the automatically suggested reference dataset, allowing for clade determination based on the HA sequence