

Jun 04, 2020

🌐 Selecting the correct transfer membrane for Western blot enhances detection of alpha-synuclein and tau proteins

📁 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.bg6kjjcw>

Gabriele Kaminski Schierle¹

¹University of Cambridge


Molecular Neuroscience Group

Tech. support email: gsk20@cam.ac.uk



Gabriele Kaminski Schierle

University of Cambridge



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bg6kjjcw>

Protocol Citation: Gabriele Kaminski Schierle 2020. Selecting the correct transfer membrane for Western blot enhances detection of alpha-synuclein and tau proteins. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bg6kjjcw>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 04, 2020

Last Modified: June 04, 2020

Protocol Integer ID: 37804

Keywords: Western blot, gel transfer, alpha-synuclein, tau, PVDF, nitrocellulose, fixation, synuclein fixation of the membrane, synuclein from tissue, synuclein fixation, synuclein, tau protein, tau proteins the ability, correct transfer membrane for western blot, correct transfer membrane, transfer membrane, recombinant human tau, proteins by western blot, polyvinylidene difluoride, protein, membrane, choice of membrane, detection of alpha,

Abstract

The ability to detect proteins by Western blot, particularly from tissue samples, can be hampered by the choice of membrane and whether a fixation step is used. Here we compare three transfer membranes, polyvinylidene difluoride (PVDF) 0.2 μm , PVDF 0.45 μm and nitrocellulose 0.45 μm for the detection of rat α -synuclein from tissue and recombinant human tau. For α -synuclein fixation of the membrane directly after transfer was imperative for detection of the protein and use of the PVDF 0.45 μm membrane gave the highest detection signal. For tau, the signal was highest on nitrocellulose 0.45 μm and fixation did not enhance the signal.



Materials

STEP MATERIALS

- ☒ NuPAGE® LDS Sample Buffer (4X) **Thermo Fisher Catalog #NP0008**
- ☒ NuPAGE® 4-12% Bis-Tris Protein Gels, 1.5 mm, 15-well **Thermo Fisher Catalog #NP0336PK2**
- ☒ NuPAGE® MES SDS Running Buffer (20X) **Thermo Fisher Catalog #NP0002**
- ☒ Immobilon-PSQ PVDF Membrane 0.2um roll **Merck MilliporeSigma (Sigma-Aldrich) Catalog #ISEQ00005**
- ☒ Immobilon-P PVDF Membrane, 0.45um, roll **Merck MilliporeSigma (Sigma-Aldrich) Catalog #IPVH00010**
- ☒ Nitrocellulose Transfer Membrane **Amersham plc Catalog #10600002**
- ☒ Methanol
- ☒ NuPAGE™ transfer buffer **Thermo Fisher Scientific Catalog #NP0006**
- ☒ PBS
- ☒ Tween-20 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P9416**
- ☒ Rabbit IgG horse radish peroxidase (HRP) linked Whole Ab **Catalog #GENA934-1ML**
- ☒ SuperSignal® West Pico PLUS Chemiluminescent Substrate **Thermo Fisher Catalog #34579**
- ☒ Polyclonal Rabbit Anti-Human Tau Unconjugated Ig fraction **Agilent Technologies Catalog #A002401-2**
- ☒ α -Synuclein (D37A6) XP® Rabbit mAb **Cell Signaling Technology Catalog #4179**
- ☒ Rabbit IgG horse radish peroxidase (HRP) linked Whole Ab **Catalog #GENA934-1ML**



Protocol materials

- ⊗ Polyclonal Rabbit Anti-Human Tau Unconjugated Ig fraction **Agilent Technologies Catalog #A002401-2**
- ⊗ NuPAGE[®] LDS Sample Buffer (4X) **Thermo Fisher Catalog #NP0008**
- ⊗ NuPAGE[®]; 4-12% Bis-Tris Protein Gels, 1.5 mm, 15-well **Thermo Fisher Catalog #NP0336PK2**
- ⊗ Rabbit IgG horse radish peroxidase (HRP) linked Whole Ab **Catalog #GENA934-1ML**
- ⊗ Immobilon-PSQ PVDF Membrane 0.2um roll **Merck MilliporeSigma (Sigma-Aldrich) Catalog #ISEQ00005**
- ⊗ Immobilon-P PVDF Membrane, 0.45um, roll **Merck MilliporeSigma (Sigma-Aldrich) Catalog #IPVH00010**
- ⊗ Methanol
- ⊗ SuperSignal[®]; West Pico PLUS Chemiluminescent Substrate **Thermo Fisher Catalog #34579**
- ⊗ NuPAGE[®]; MES SDS Running Buffer (20X) **Thermo Fisher Catalog #NP0002**
- ⊗ NuPAGE[™] transfer buffer **Thermo Fisher Scientific Catalog #NP0006**
- ⊗ α -Synuclein (D37A6) XP[®] Rabbit mAb **Cell Signaling Technology Catalog #4179**
- ⊗ Nitrocellulose Transfer Membrane **Amersham plc Catalog #10600002**
- ⊗ PBS
- ⊗ Tween-20 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P9416**

Safety warnings

- ⚠ Follow safety guidelines and wear correct PPE when handling fixation solution containing 4% paraformaldehyde and 0.1% glutaraldehyde

1 Protein samples were boiled for 00:03:00 in 1x

NuPAGE[®] LDS Sample Buffer (4X) Thermo Fisher Catalog #NP0008

and ran on

NuPAGE[®] 4-12% Bis-Tris Protein Gels, 1.5 mm, 15-well Thermo Fisher Catalog #NP0336PK2

in

NuPAGE[®] MES SDS Running Buffer (20X) Thermo Fisher Catalog #NP0002

using the

Equipment

XCell SureLock Mini-Cell Electrophoresis System NAME

Electrophoresis System TYPE

Invitrogen BRAND

EI0001 SKU

at 200 V for 00:30:00 .

2 Three membranes were selected to test Western blot transfer and detection, PVDF 0.2 µm

Immobilon-PSQ PVDF Membrane 0.2um roll Merck MilliporeSigma (Sigma-Aldrich) Catalog #ISEQ00005

PVDF 0.45 µm

Immobilon-P PVDF Membrane, 0.45um, roll Merck MilliporeSigma (Sigma-Aldrich) Catalog #IPVH00010

and nitrocellulose 0.45 µm

Nitrocellulose Transfer Membrane Amersham plc Catalog #10600002

.

3 PVDF membranes were first activated with methanol by incubating for 00:03:00 . Nitrocellulose membranes do not need to be activated first.



4 The gel was transferred onto each membrane using the

Equipment

Invitrogen XCell II™ Blot Module

NAME

Western blot transfer module

TYPE

Invitrogen

BRAND

EI9051

SKU

with

 NuPAGE™ transfer buffer Thermo Fisher Scientific Catalog #NP0006

+ 20%

 Methanol

at 30 V for  01:30:00 . The blot module was kept  On ice and surrounded by ice to keep cool and prevent protein damage.

5

Safety information

PFA and glutaraldehyde are toxic and should be handled in a chemical hood with appropriate PPE. The fixation solution should never be disposed of down the sink, but disposed of in a separate container for disposal through correct waste systems.

For detecting α -synuclein in tissue samples the membranes were fixed immediately after transfer with 4% paraformaldehyde (PFA), 0.1% glutaraldehyde in PBS for  00:30:00 . Lack of fixation leads to poor/no α -synuclein detection and fixation with only 4% PFA also leads to lower detection than with the addition of 0.1% glutaraldehyde.

6 Following disposal of the fixation solution into properly designated containers for removal, the membranes were blocked for  01:00:00 in 5% BSA in

 PBS

+ 0.05%

 Tween-20 Merck MilliporeSigma (Sigma-Aldrich) Catalog #P9416



(PBS-T) at Room temperature .

- 7 The membranes were then incubated either overnight or for 01:00:00 at Room temperature in primary antibody probing for rat α -synuclein 1:1000 dilution

α -Synuclein (D37A6) XP® Rabbit mAb **Cell Signaling Technology Catalog #4179**

- 8 The membranes were washed three times for 00:05:00 in PBS-T.

- 9 The membranes were incubated for 01:00:00 at Room temperature in secondary antibody 1:4000

Rabbit IgG horse radish peroxidase (HRP) linked Whole Ab **Catalog #GENA934-1ML**

- 10 The membranes were washed five times for 00:05:00 in PBS-T.

11 The membranes were developed in


 SuperSignal[®] West Pico PLUS Chemiluminescent Substrate **ThermoFisher** **Catalog #34579**

and imaged in a

Equipment	
G:Box Chemi	NAME
Imaging system	TYPE
Syngene	BRAND
XX6	SKU

12 .

Expected result



The best signal for detection of alpha-synuclein from rat tissue was obtained on PVDF 0.45 μ m with fixation in 4% PFA and 0.1% glutaraldehyde.

- 13 For recombinant human tau, fixation was not required and the membranes were instead immediately blocked for 01:00:00 in 5% BSA in PBS-T at Room temperature .
- 14 The membranes were then incubated either overnight or for 01:00:00 at Room temperature in primary antibody probing for human tau 1:200 dilution
 Polyclonal Rabbit Anti-Human Tau Unconjugated Ig fraction **Agilent Technologies Catalog #A002401-2**
- 15 The membranes were washed three times for 00:05:00 in PBS-T.

16 The membranes were incubated for ⌚ 01:00:00 at 🌡 Room temperature in secondary antibody, 1:4000 dilution

🔬 Rabbit IgG horse radish peroxidase (HRP) linked Whole Ab **Catalog #GENA934-1ML**

17 The membranes were washed five times for ⌚ 00:05:00 in PBS-T.


18 The membranes were developed in Supersignal West Pico (#34580, ThermoFisher) and imaged in a

Equipment

G:Box Chemi	NAME
Imaging system	TYPE
Syngene	BRAND
XX6	SKU

19

Expected result



The best signal for recombinant Tau was obtained on nitrocellulose 0.45 µm.

- 20 This optimisation was originally performed by Dr Colin Hockings, a former postdoc in the Molecular Neuroscience Group.