


Aug 18, 2025

Selective Isolation of Lipid Classes from cold-water corals via SPE for Quantitative Fatty Acid Profiling

DOI

<https://dx.doi.org/10.17504/protocols.io.5qpvdx37g4o/v1>



Audrey M. Pruski¹, Gilles Vétion¹

¹Laboratoire d'Ecogéochimie des environnements benthiques, Sorbonne Université

Hussain's den



Audrey M. PRUSKI

Sorbonne University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpvdx37g4o/v1>

Protocol Citation: Audrey M. Pruski, Gilles Vétion 2025. Selective Isolation of Lipid Classes from cold-water corals via SPE for Quantitative Fatty Acid Profiling. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5qpvdx37g4o/v1>

Manuscript citation:

In situ diet patterns and health status of cold-water corals in the Lacaze-Duthiers canyon (NW Mediterranean Sea): insights from fatty acid biomarkers on lipid classes. Audrey M. Pruski, Gilles Vétion, Franck Lartaud, Erwan Peru, and Nadine Le Bris. Deep Sea Research Part I, accepted with minor modifications

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 30, 2025


Last Modified: August 18, 2025

Protocol Integer ID: 221340

Keywords: Lipids, Fatty acids, Solid Phase Extraction, Corals, extraction of total lipid, water corals via spe, coral polyp sample, selective isolation of lipid class, separation of lipid class, quantitative fatty acid profiling, quantitative fatty acid profiling this protocol, fatty acid methyl ester, water coral, total lipid, lipid class, extraction, solid phase extraction, step acidic transesterification

Abstract

This protocol describes the extraction of total lipids from coral polyp samples using the method of Bligh and Dyer (1959), followed by the separation of lipid classes via solid phase extraction (SPE) on aminopropyl (NH₂) silica cartridges as described by Kaluzny et al. (1985). A one-step acidic transesterification, adapted from Indarti et al. (2005) and Christie (2003), is then applied to produce fatty acid methyl esters (FAMES) for GC-MS analysis.



Protocol flowchart

Image Attribution

Pruski Audrey M.



Guidelines

Sample Handling: Handle sample quickly after collection to prevent lipolysis

Equipment: Ensure all glassware, tubes, and Pasteur pipettes are pre-burnt ( 450 °C) to avoid contamination

Chemicals and solvents: Use GC-MS or HPLC reagent grade products

Internal Standard: Nonadecanoic acid (C19:0) is used to correct for extraction losses during transesterification

Temperature Control: Maintain precise heating and cooling conditions to ensure consistent results

Storage: Keep eluates and FAME extracts at  -20 °C if not analyzed immediately

Materials

Equipments:

- Freeze dryer (Christ Alpha 1-4)
- TissueLyser II (QIAGEN)
- Solid phase extraction manifold (Phenomenex®)
- Hand-operated vacuum pump
- Rotary shaker (Bioblock)
- Ultracentrifuge (Sigma 6K-15)
- Speed-Vac concentrator (Thermo Scientific™ Savant™)
- GC-MS system (Agilent)

Materials:

- Strata NH₂ cartridges (Phenomenex®)
- 20 and 30 ml Borosilicate glass tubes round bottom with PBT red screw cap & PTFE liner
- 20 ml Borosilicate glass tubes conical bottom with PBT red screw cap & PTFE liner
- 2ml GC vials with 200µl insert
- Polar capillary column (e.g. VF-Wax-MS or ZB-WAX-MS)

Chemicals and reagents (HPLC or GC grade):

- Acetic acid (CH₃CO₂H)
- Butylated hydroxytoluene (BHT)
- Chloroform (CH₃Cl)
- Dichloromethane (CH₂Cl₂)
- Diethyl ether ((C₂H₅)₂O)
- Ethyl acetate (C₄H₈O₂)
- Hexane (C₆H₁₄)
- Hydrochloric acid (HCl)
- Methanol (CH₃OH)
- Nonadecanoic acid (C19:0= internal standard) (C₁₉H₃₈O₂)
- Phosphoric acid (H₃PO₄)
- Potassium bicarbonate (KHCO₃)
- 2-Propanol (C₃H₈O)
- Sodium chloride (NaCl)
- Sodium hydroxide (NaOH)
- Sulphuric acid (H₂SO₄)
- Water (ultrapure)

- Commercial FAME standards (Qualmix Fish Synthetic, Ladoran Fine Chemicals, INTERCHIM, France; Supelco 37, PUFA No. 1 and No. 3, SUPELCO France) for GC-MS calibration

Note: use GC/HPLC grade reagents

Biological samples:

- Coral nubbins

Safety warnings



Safety information

- Always wear safety goggles and gloves
- Work in a fume hood when handling solvents
- Dispose of waste in appropriate containers



Lipid Extraction from the Biological Matrix (Cold-water coral nubbins)

1d 9h 6m

1 Sample Preparation

1.1 Flash-freeze  coral nubbins on board to prevent degradation

10m

1.2 Freeze-dry samples and store at  -80 °C until analysis

1d

1.3 Grind freeze-dried  coral nubbins into a fine powder using a bead mill homogenizer while keeping the chamber filled with liquid nitrogen to prevent heating

10m

2 Preparation of Working Solutions

4h

2.1 Chloroform/methanol/water (1:2:0.8) reagent: 5.265 mL CHCl_3 , 10.525 mL CH_3OH and 4.21 mL H_2O . Keep at  4 °C in a flammable storage refrigerator

2.2 1M NaCl solution in ultrapure water ($58\text{g}\cdot\text{L}^{-1}$). Keep at  4 °C in a refrigerator

2.3 0.3N NaOH solution ($12\text{g}\cdot\text{L}^{-1}$). Keep at  4 °C in a refrigerator

2.4 0.2M H_3PO_4 solution ($19.6\text{g}\cdot\text{L}^{-1}$). Keep at  4 °C in a refrigerator

2.5 BHT (butylated hydroxytoluene) at $50\text{mg}\cdot\text{mL}^{-1}$ in Chloroform. Keep at  4 °C in a flammable storage refrigerator

3 Total lipid extraction

4h

**Note**

Total lipids are extracted using a modified version of the Bligh and Dyer method (1959), optimized to improve extraction efficiency.

3.1 Place the freeze-dried  coral nubbins (~  3 g) in a 30 mL glass tube (round bottom, screw tap with PTFE liner)

3.2 Add  20 mL of $\text{CHCl}_3/\text{CH}_3\text{OH}/\text{H}_2\text{O}$ (1:2:0.8) with 1M NaCl + 0.2M H_3PO_4

Note

The solution of NaCl 1M- H_3PO_4 0.2M increases protein unwinding and lipid recovery

3.3 Sonicate  00:05:00 ( On ice if needed), then rotate  00:20:00 at  Room temperature

25m

3.4 Add  2.5 mL CHCl_3 , then  2.5 mL H_2O . Shake  00:05:00

5m

3.5  2500 rpm, 4°C, 00:05:00 . Recover the organic phase (bottom) using a Pasteur pipette

5m

3.6 Repeat extraction 3 times with  5 mL CHCl_3

3.7 Pool all organic phases (final volume ≈ 22.765 mL) and vortex.  Keep organic phase cold  On ice

3.8 Wash organic phase with  2.5 mL cold 0.3N NaOH. Shake  00:01:00

1m

3.9  2500 rpm, 4°C, 00:05:00 and recover the organic phase

5m



- 3.10 Take  20 mL or more into a new glass tube (conical bottom, screw tap with PTFE liner) and evaporate at  40 °C under vacuum (rotary evaporator)

Note

⚠ Lipid extracts should never be allowed to dry out!

- 3.11 Cool at  -20 °C and redissolve in  400 µL CHCl₃ + BHT (50 mg.L⁻¹)

Note

BHT is added optionally to chloroform to prevent lipid oxidation during extraction!

- 3.12 Centrifuge  2000 rpm, 4°C, 00:05:00

5m

Store at  -20 °C

Lipid Class Separation by Solid Phase Extraction (SPE)

1d

4 Solvent Preparation

- 4.1 Solvent A: Chloroform:2-propanol (2:1 v/v)
Solvent B: 2% Acetic acid in diethyl ether
Solvent C: Methanol
Solvent D: Hexane
Solvent E: 1% Diethyl ether + 10% Dichloromethane in hexane
Solvent F: 5% Ethyl acetate in hexane
Solvent G: 15% Ethyl acetate in hexane
Solvent H: Chloroform:Methanol (2:1 v/v)

**Note**

Keep all solvents at  4 °C in a flammable storage refrigerator!

5 Separation of neutral lipids, fatty acids and phospholipids

1d

Note

3 Strata NH₂ (aminopropyl silica) cartridges (Phenomenex) are needed per sample. 7 lipid fractions are efficiently recovered from the initial chloroform extract with minimal impurities.

5.1 Condition the first cartridge (C1) with 2 ×  2 mL hexane, vacuum dry 2–5 min

5.2 Load  390 µL lipid extract. Vacuum to dryness

5.3 Elute neutral lipids (E1): 4 ×  2 mL solvent A

5.4 Elute free fatty acids (E2): 4 ×  2 mL solvent B

5.5 Elute phospholipids (E3): 4 ×  2 mL solvent C

5.6 Store eluates at  -20 °C .

6 Neutral Lipid Fractionation**Note**

In the original method by Kaluzny et al. (1985), the elution of steryl esters was demonstrated using commercially available lipid standards, which did not include wax esters. However, based on our experimental results, wax esters appear to co-elute with steryl esters under the same solvent conditions used for fraction E4.



6.1 Condition a second cartridge (C2) as above with hexane

6.2 Load  390 μ L from E1. No need to dry

6.3 Elute wax esters and sterol esters (E4): 4 $\times$  2 mL solvent D

Note

Note: Mono-, di-, and triglycerides and sterols are retained on C2.

6.4 Condition a third cartridge (C3) as above

6.5 Stack C2 onto C3 using an adapter

6.6 Elute triglycerides (E5): 6 $\times$  2 mL solvent E

6.7 Elute sterols (E6): 9 $\times$  2 mL solvent F

6.8 Remove and discard C3, place C2 back on manifold

6.9 Elute diglycerides (E7): 4 $\times$  2 mL solvent G

6.10 Elute monoglycerides (E8): 4 $\times$  2 mL solvent H

7 Evaporation of lipid fractions

7.1 Evaporate under vacuum all eluates ( 45 $^{\circ}$ C)



7.2 Redissolve in  100 μL CHCl_3 + BHT. Store at  $-20\text{ }^\circ\text{C}$

Note

Eluates for GC-MS analysis of fatty acids are E2 (Free Fatty Acids), E3 (Phospholipids), E4 (Wax esters and Sterol esters), E5 (Triglycerides), E7 (Diglycerides), and E8 (Monoglycerides). E6 may be used for sterol analysis!

One-Step Acidic Transesterification

1d

8 Lipid extraction and Transesterification

8.1 Preparation of Working Solutions

- Transmethylation reagent:  3.4 mL Methanol +  0.6 mL H_2SO_4 (98%), and  4.0 mL CHCl_3 with BHT (50 mg/L). Keep at  $4\text{ }^\circ\text{C}$ in a flammable storage refrigerator

Note

Reagent should be cold and freshly prepared

- C19:0 (internal standard) stock solution $\sim 1\text{ mg}\cdot\text{mL}^{-1}$ in hexane. The solution is diluted before use (1/9). Keep at  $-20\text{ }^\circ\text{C}$

Note

Adjust the concentration of the internal standard so that its peak on the chromatogram falls within the range of peak intensities, ideally between the most abundant and the least abundant fatty acids. This ensures accurate quantification across the entire dynamic range

8.2 To the eluates: add 20 μL internal standard C19:0 + 8 mL of transmethylation reagent. Cap and shake

8.3 Incubate  01:30:00 at  $90\text{ }^\circ\text{C}$. Shake occasionally (3-4 times)

1h 30m



9 Phase Separation

- 9.1 Cool (On ice or -20 °C). Open carefully
- 9.2 Add 2 mL cold H₂O. Vortex. 300 rpm, 4°C, 00:05:00 Centrifuge 5 min, 3000 rpm, 4°C
- 9.3 Recover lower organic phase in a new glass tube

10 Organic Phase Wash

- 10.1 Re-extract aqueous phase 2 times with 2 mL cold hexane:CHCl₃ (4:1)
- 10.2 Combine all organic phases (~8–9 mL)

11 Final Purification

- 11.1 Wash with 4 mL cold 2% K₂CO₃. 3000 rpm, 4°C, 00:05:00
- 11.2 Discard the upper polar layer and transfer 6–9 mL of the lower organic phase to a clean screw-cap tube

12 Drying and Dissolving FAMES

- 12.1 Evaporate solvent in a SpeedVac concentrator at Room temperature . Cool at -20 °C
- 12.2 Redissolve in 30–50 µL pure hexane. Vortex, 3000 rpm, 4°C, 00:05:00

10m

12.3 Transfer to GC vials with a 200 µl insert. Store at  -80 °C

12.4 Fatty acid methyl extracts are ready for GC-MS analysis

Note

In this study, the analyses were carried out using the instrumental conditions described in a companion protocol developed by our group (Pruski & Vétion, 2024).

Protocol references

Bligh, E.G., Dyer, W.J., 1959. A rapid method of total lipid extraction and purification. *Canadian Journal of Biochemistry and Physiology*, 37, 911–917.

Christie, W.W., 2003. Lipid analysis, volume 15 : isolation, separation, identification and structural analysis of lipids. Third edition. The Oily Press, Bridgawater, England.

Indarti, E., Majid, M.I.A., Hashim, R., Chong, A., 2005. Direct FAME synthesis for rapid total lipid analysis from fish oil and cod liver oil. *Journal of Food Composition and Analysis*, 18, 161–170. <https://doi.org/10.1016/j.jfca.2003.12.007>

Kaluzny, M. A., Duncan, L. A., Merritt, M. V., & Epps, D. E., 1985. Rapid separation of lipid classes in high yield and purity using bonded phase columns. *Journal of Lipid Research*, 26, 135–140

Pruski, A.M., Vétion, G. 2024. Comparison of direct transesterification and conventional multi-step methods for fatty acid analysis in marine sediments. *protocols.io*. <https://dx.doi.org/10.17504/protocols.io.81wgbz3ylpk/v1>

Simpson, N.J.K. 2000. Solid-Phase Extraction: Principles, Techniques, and Applications. CRC Press.