


May 05, 2020 Version 1

Separation and purification of human PBMC from FRESH BLOOD V.1

DOI

dx.doi.org/10.17504/protocols.io.bfxmjpk6



Marco Cosentino¹, Elisa Storelli², Alessandra Luini¹, Massimiliano LM Legnaro³, Emanuela Rasini³, Marco Ferrari³, Franca Marino³

¹Center for Research in Medical Pharmacology, University of Insubria (Varese, Italy);

²Center for Research in Medical Pharmacology, University of Insubria; ³University of Insubria



Elisa Storelli

Center for Research in Medical Pharmacology, University of I...

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bfxmjpk6

Protocol Citation: Marco Cosentino, Elisa Storelli, Alessandra Luini, Massimiliano LM Legnaro, Emanuela Rasini, Marco Ferrari, Franca Marino 2020. Separation and purification of human PBMC from FRESH BLOOD. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bfxmjpk6>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 04, 2020

Last Modified: May 05, 2020

Protocol Integer ID: 36557

Keywords: PBMC, Fresh Blood, Neuroimmune-Pharmacology, Parkinson's Disease, Cell isolation, Primary cell culture,

Abstract

Separation and purification of PBMC from FRESH BLOOD: list of published work using this protocol

Kustrimovic, N., Comi, C., Magistrelli, L., Rasini, E., Legnaro, M., Bombelli, R., Aleksic, I., Blandini, F., Minafra, B., Riboldazzi, G., Sturchio, A., Mauri, M., Bono, G., Marino, F., & Cosentino, M. (2018). Parkinson's disease patients have a complex phenotypic and functional Th1 bias: cross-sectional studies of CD4+ Th1/Th2/T17 and Treg in drug-naïve and drug-treated patients. *Journal of neuroinflammation*, 15(1), 205. <https://doi.org/10.1186/s12974-018-1248-8>

Kustrimovic, N., Rasini, E., Legnaro, M., Bombelli, R., Aleksic, I., Blandini, F., Comi, C., Mauri, M., Minafra, B., Riboldazzi, G., Sanchez-Guajardo, V., Marino, F., & Cosentino, M. (2016). Dopaminergic Receptors on CD4+ T Naive and Memory Lymphocytes Correlate with Motor Impairment in Patients with Parkinson's Disease. *Scientific reports*, 6, 33738. <https://doi.org/10.1038/srep33738>

Cosentino M., Ferrari M., Kustrimovic N., Rasini E., Marino F. (2015). Influence of dopamine receptor gene polymorphisms on circulating T lymphocytes: A pilot study in healthy subjects. *Human immunology*, 76, 10, 747-752. <https://doi.org/10.1016/j.humimm.2015.09.032>

Materials

MATERIALS

 Ficoll Paque PLUS **GE Healthcare Catalog #17144003-500 ml**

 Fetal Bovine Serum (FBS) **EuroClone Catalog #ECS0180L-500 ml**

 RPMI 1640 **EuroClone Catalog #ECM 0495L- 500 ml**

 Trypan Blue solution 0.4% **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T8154- 100 ml**

Instrumentation required:

- Laminar flow hood
- Centrifuge
- Cellometer (automated cell counter) or Optical Microscope (manual cell count)
- Flow Cytometer
- Autoclave




Before start

If you need to obtain **PBMC for cell culture**, make sure you are using **sterile PBS, culture medium, filtered Lysis Buffer and sterile plastic disposables as well**. Moreover, work under laminar flow hood when you are processing samples. Otherwise, use non-sterile solutions and plastic disposables, and process samples in cell isolation laboratory.

ALL REAGENTS USED IN THIS PROTOCOL MUST BE AT ROOM TEMPERATURE!






- 1 Put the needed amount of blood sampl into a  50 mL conical tube.
- 2 Add an equal volume of **PBS 1X** and mix well.

Document


 NAME
PBS 1X


CREATED BY
Marco Ferrari

PREVIEW

- 3 Place  3 mL of FICOLL in a  15 mL conical tube.
- 4 Carefully layer  12 mL of diluted blood on FICOLL with a glass Pasteur Pipette to a final volume of 15 ml as shown in the figure below.





- 5 Centrifuge samples  400 x g, 00:40:00 at room temperature (RT) without break.

Equipment

Allegra AVANTI 30

NAME

Centrifuge

TYPE

Beckman Coulter


BRAND




Beckman Italy

SKU

- 6 After centrifugation, take out the tubes carefully to not disturb the mononuclear cell layer that appears as a white, cloudy band between the plasma and FICOLL as shown in the

figure below.



- 7 Carefully with a glass Pasteur pipette transfer the mononuclear lymphocyte cell layer to another 15 ml conical tube. 
- 8 Wash the isolated PBMC with **PBS/FBS 2%** to a final volume of  10 mL and centrifuge at  600 x g, 00:10:00 at RT.



Document






NAME

RPMI - FBS and PBS - FBS

CREATED BY

Elisa Storelli

PREVIEW

- 9 Remove supernatants, resuspend pellet in  1 mL of **Lysis Buffer** and add another  9 mL of **Lysis Buffer**. Immediately centrifuge tubes at  300 x g, 00:10:00 at RT.

Document





NAME

Lysis Buffer

CREATED BY

Elisa Storelli

PREVIEW

- 10 Remove supernatant and resuspend pellet in  10 mL of **PBS/FBS 2%** and centrifuge at  600 x g, 00:10:00 at RT.



Document




NAME

RPMI - FBS and PBS - FBS

CREATED BY

Elisa Storelli

PREVIEW

- 11 Remove supernatant and resuspend the obtained pellet in  10 mL of **RPMI/FBS 10%** for cell counting.

Document



NAME


RPMI - FBS and PBS - FBS


CREATED BY

Elisa Storelli

PREVIEW

- 12 **For manual cell count use Türk solution for checking purity.**

Mix  10 μ L of cell suspension with an equal amount of Türk solution (dilution factor = 2), allow mixture 3 min at room temperature.

Take  10 μ L of the mixture and place it inside a Bürker chamber and view under an optical microscope using 40X magnification.

Count the cells in each square found in the four corners and in the central square (see figure 1 below), including those that lie on the bottom and left-hand perimeters, but not those that lie on the top and right hand perimeters (see figure 2 below).

Total number of cells per ml = mean number of cells x dilution factor x 10^4 (hemacytometer volume).




Figure 1





Figure 2

Document



NAME
Trypan blue and Turk solution

CREATED BY
Marco Ferrari

PREVIEW


13 OPTIONAL STEP

For automatic cell count with Cellometer machine use Trypan Blue. The machine will calculate the n° of cells/ml and the % of viability.

Take  10 μ L of cell suspension and add an equal amount of Trypan Blue. Use all the volume to place it in a counting chamber. Place the chamber inside Cellometer and

count.

Equipment	
Cellometer Auto T4	NAME
Automated Cell Counter	TYPE
Nexcelom Bioscience	BRAND
Euroclone	SKU

Document	
	NAME
Trypan blue and Turk solution	
CREATED BY	
Marco Ferrari	PREVIEW

- 14 If needed, check the purity of PBMC suspension by using morphological parameter of the flow cytometer.

For this test 0.5×10^6 PBMC in 500 μ l of PBS are enough.





Equipment

BD FACS Celesta

NAME

Flow Cytometer

TYPE

Becton Dickinson

BRAND

Milan Italy BD

SKU

15 Expected results

Expected result

VIABILITY - The expected viability by Trypan Blue should be $\geq 90\%$.

PURITY - The PBMC suspension obtained should contain at least 80% of lymphocytes, 10-15% of monocytes and few contaminant PMN cells ($\leq 5\%$) as confirmed by flow cytometry.

YIELD - The expected amount of PBMCs should be $\pm 28,5 \times 10^6$ starting from 25 ml of fresh blood.