

Sep 18, 2019

Version 3

seqFISH Tissue Preservation V.3

DOI

dx.doi.org/10.17504/protocols.io.7e7hjhn

Long Cai¹, Nina Dar¹, Yiing Lin², Shin Lin³

¹California Institute of Technology; ²Washington University, Saint Louis; ³University of Washington

Human BioMolecular Atlas Program (HuBMAP) Method Development Community

Tech. support email: Jeff.spraggins@vanderbilt.edu



Shin Lin

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.7e7hjhn>

Protocol Citation: Long Cai, Nina Dar, Yiing Lin, Shin Lin 2019. seqFISH Tissue Preservation. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.7e7hjhn>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: In development

We are still developing and optimizing this protocol

Created: September 18, 2019

Last Modified: September 18, 2019

Protocol Integer ID: 27839

Keywords: seqFISH, HUBMAP, seqfish tissue preservation tissue preparation, tissue preservation tissue preparation, tissue, preparation

Abstract

Tissue Preparation for seqFISH



Materials

	Reagent	company	Cat#
	16% PFA RNAse free ampules	Electron Microscopy	100503-917
	10x PBS, RNAse free	invitrogen	AM9625
	Ultra pure water, DNAse and RNAse free	invitrogen	10977-015
	RNA zap	Invitrogen	AM9780
	Ultra pure Sucrose, RNAse free	VWR	0335

	Material	Company	Cat#
	50 ml tube, DNAse and	Corning	430290



	RNA se free		
	15 ml tube s, DNA se and RNA se free	Corn ing	430 052
	1.7 ml micr ocen trifu ge tube s, DNA se and RNA se free certi fied	VWR	870 03- 294

Troubleshooting



Fixation

- 1 Clean surface with RNAzap. Make the following solutions.

1L PBS 1x, RNase-free

1. take 1 L RNase-free water bottle
2. use serologic pipet and draw up 100 mL water into another container
3. add 100 mL of 10x RNase-free PBS
4. mix thoroughly

50% sucrose

1. tare 1 L RNase-free disposable container
2. take 500 g RNase-free sucrose container
3. add RNase-free water to sucrose and transfer to 1 L container
4. continue to add RNase-free water until mass to 1 kg
5. mix thoroughly, may need agitation on platform at RT for some time
6. store at 4°C

NOTE: for fixation, minimum PFA to tissue volume is 20:1, ideal closer to 50:1

- 2 Prepare 4% PFA in 1x PBS, 120 mL total for 4 tissues 30 mL each in 50 mL conicals.

	prep are 3x conical s for 4 tissues	
	16% PFA	10 ml
	10x PBS	4 ml
	RNA se free water	26 ml
	total	40 ml



- 3 Cut freshly dissected tissue to small pieces for heart, try cutting slab, 2 cm x 1 cm x 0.4 cm thickness: 2 cm dimension from epicardium to endocardium.
- 4 Place each piece in a 50 ml tube (or smaller if you can cut to smaller pieces) in 4% PFA leave at RT for 3 hrs.
** the sample should be completely covered.
- 5 Wash three times with equal volume (30mL) of 1x PBS to remove PFA leave a bit of extra liquid in the bottom between each wash.

Sucrose gradient

- 6 Place the tissue in 10% Sucrose at RT and wait for it to sink in 50 mL conical, add 10 mL of sucrose wait for sample to sink to bottom – minimum incubation time of 30 min remove supernatant for next bath.
- 7 Place the tissue in 20% Sucrose at RT and wait for it to sink.
- 8 Place the tissue in 30% sucrose in 4°C – ready for shipping. Make sure there is enough 30% sucrose solution in container e.g. if in 50 mL conicals, have 50 mL of 30% sucrose.