


Mar 13, 2019

Version 3

Sequential smFISH V.3

DOI

dx.doi.org/10.17504/protocols.io.y6nfzde



ZengU19 BICCN Grant¹

¹Allen Institute

BICCN / BICAN



Stephanie Mok

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.y6nfzde>

Protocol Citation: ZengU19 BICCN Grant 2019. Sequential smFISH. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.y6nfzde>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 13, 2019

Last Modified: March 13, 2019

Protocol Integer ID: 21422

Keywords: smFISH, multiplexed single molecule fish protocol for use, multiplexed single molecule fish protocol, sequential smfish, mouse tissue, human tissue, tissue, protocol, nature protocol

Abstract

We have developed a multiplexed single molecule FISH protocol for use at the Institute. This protocol was optimized on human tissue, but will work on mouse tissue as well. It was adapted from Lyubimova et. al., Nature Protocols, 2013.

Attachments



[smFISH.docx](#)

16KB

Guidelines


Ensure that all reagents are in recombinant and RNase-free format, as we have noticed RNA degradation in solutions that contain enzymes derived from whole organisms.

We filter every solution with a 0.2um syringe filter prior to use. This reduces background spots and dust that interfere with imaging of diffraction limited spots.

For the SDS treatment after fixation and permeabilization, be gentle when dropping SDS onto the section, as well as during washes. This treatment is relatively harsh and the tissue must be treated somewhat delicately.

Materials

STEP MATERIALS

 4% PFA


 PBS

 PBS


 PBS

 2X SSC

 2X SSC

 65% formamide/2X SSC

Protocol materials

 4% PFA


 PBS


 PBS

 PBS

 2X SSC

 2X SSC


 65% formamide/2X SSC

 4% PFA


 PBS

 PBS

 PBS

 2X SSC

 2X SSC

 65% formamide/2X SSC

Troubleshooting



Safety warnings


- ⚠ Please refer to the SDS (Safety Data Sheet) for hazard information and safety warnings.
Avoid exposure to formamide, DAPI

Before start

Ensure all incubators and ovens are at the appropriate temperature prior to experiment.


Tissue and Sectioning

- 1 10-14 um cryosections are taken from fresh-frozen tissue, which are collected on poly-lysine-treated #1 coverslips at room temperature (RT). After 5-10 min at RT, sections are placed at 4°C until sectioning is complete. At that point, proceed immediately to fixation and permeabilization.

 00:05:00 RT

Fixation/Permeabilization

- 2 Post-fix sections for 15 min with 4% PFA @ 4 °C.

 4% PFA

 00:15:00 Post-fixing

 4 °C Post-fixing

- 3 Wash with PBS (1/3)

 PBS

- 4 Wash with PBS (2/3)

 PBS

- 5 Wash with PBS (3/3)

 PBS

- 6 Permeabilize with cold methanol at -20 °C for 10 min.


 -20 °C Permeabilizing

 00:10:00 Permeabilizing


- 7 Air dry for 30 min in fume hood (Stopping point: store coverslips at -80°C)

 00:30:00 Air drying

- 8 Optional: Treat sections with 8% SDS/PBS for 10 minutes, followed by 3 – 5 rinses with PBS or 2XSSC

 00:10:00

- 9 Add 2ml 2X SSC

 2X SSC

 2 mL 2X SSC

Hybridization


10 Pre-heat hyb oven to 37 °C

 37 °C oven

11 Place sections in hyb buffer without probes.

12 Add 4 ul probe 400ul hyb buffer.


 4 µL probe

 400 µL hyb buffer

Note

Specific to 6-well plate format – if using perfusion chamber, this volume can be reduced.


13 Incubate at 37 C for 2H.

 37 °C Incubation


 02:00:00 Incubation


Wash

14 Add 2 ml wash buffer to each well.

 2 mL wash buffer


15 Incubate at 37 C for 15 min.


 37 °C Incubation

 00:15:00 Incubation

16 Remove wash buffer.

17 Add 2 ml fresh wash buffer and incubate at 37 C for 15 min.

 2 mL wash buffer

 37 °C Incubation

 00:15:00 Incubation

18 Replace wash buffer with fresh wash buffer + DAPI (final 5ug/mL) and incubate at 37 C for 15 min.



37 °C Incubation

00:15:00 Incubation

19 GLOX buffer step if performing antibody stain

20 Mount and image or store at 4 C in 2XSSC until imaging session

4 °C

2X SSC

Stripping

21 65% formamide/2X SSC, 10 min X 3, 30 C

65% formamide/2X SSC

00:10:00

30 °C

22 Wash in 2XSSC (1/3)

23 Wash in 2XSSC (2/3)

24 Wash in 2XSSC (3/3)

25 Following stripping, proceed to hybridization step.

go to step #10 Hybridization step