

Feb 27, 2019

Serial Dilution of Nucleofected iPSC Pools

 In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.x76frre

Celeste M M. Karch¹, Rita Martinez¹, Jacob Marsh¹


¹Washington University in St Louis

Neurodegeneration Method Development Community
Tech. support email: ndcn-help@chanzuckerberg.com



Celeste M M. Karch

Washington University in St Louis



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.x76frre>

Protocol Citation: Celeste M M. Karch, Rita Martinez, Jacob Marsh 2019. Serial Dilution of Nucleofected iPSC Pools.
protocols.io <https://dx.doi.org/10.17504/protocols.io.x76frre>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 17, 2019

Last Modified: February 27, 2019

Protocol Integer ID: 20446

Keywords: serial dilution of nucleofected, ipsc pool, serial dilution, nucleofected, ipsc

Attachments



Comprehensive

Genomi...


31KB

Guidelines







This protocol is part of the **Genomic Editing: iPSC collection**.

Troubleshooting

Safety warnings


 Please refer to the SDS (Safety Data Sheet) for information about hazards, and to obtain advice on safety precautions.



- 1 Coat 3 wells of a 6 well plate with  1 mL Matrigel (supplemented with RGD fragment)
- 2 Aspirate media from cells in culture.
- 3 Wash with  1 mL -  2 mL PBS per well.
- 4 Add  1 mL Accutase per well.
- 5 Incubate at  37 °C for  00:10:00 to achieve single cells.

Note

Individual donor lines exhibit variable sensitivity to accutase-mediated dissociation. Monitor cells regularly to identify when cells achieve single cells.

- 6 Collect cells in 5 mL DMEM/F12 and transfer to a 15mL conical tube.
- 7 Spin at 750-800 rpm for  00:03:00 .
- 8 Aspirate media.
- 9 Resuspend cells in mTesR1 supplemented with 5 uM Rock Inhibitor.
- 10 Plate several dilutions of cells over the three wells (typically 25 µl, 50 µl and 75 µl in 2 mL of mTesR1).
- 11 Freeze down the remaining cells by adding an equal volume of 2x iPSC Freezing Media (20% DMSO in FBS).



- 12 Incubate at 🌡️ 37 °C overnight.
- 13 Change mTesR1 the following day.