

Oct 16, 2019

Setting up a liquid culture and harvesting C.elegans

DOI

<https://dx.doi.org/10.17504/protocols.io.xcpfivn>

Vidur Sabharwal¹

¹Tata Institute of Fundamental Research



Vidur Sabharwal

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.xcpfivn>

Protocol Citation: Vidur Sabharwal 2019. Setting up a liquid culture and harvesting C.elegans. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.xcpfivn>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 22, 2019



Last Modified: October 16, 2019

Protocol Integer ID: 19567

Keywords: liquid culture

Materials

S Basal [5.85 g NaCl, 1 g K₂ HPO₄, 6 g KH₂PO₄, 1 ml cholesterol (5 mg/ml in ethanol), H₂O to 1 litre. Sterilize by autoclaving.]

1 M Potassium citrate pH 6.0 [20 g citric acid monohydrate, 293.5 g tri-potassium citrate monohydrate, H₂O to 1 litre. Sterilize by autoclaving.]

Trace metals solution [1.86 g disodium EDTA, 0.69 g FeSO₄ • 7 H₂O, 0.2 g MnCl₂ • 4 H₂O, 0.29 g ZnSO₄ • 7 H₂O, 0.025 g CuSO₄ • 5 H₂O, H₂O to 1 litre. Sterilize by autoclaving. Store in the dark.]

1 M CaCl₂ [55.5 g CaCl₂ in 1 litre H₂O. Sterilize by autoclaving.]

S Medium [1 litre S Basal, 10 ml 1 M potassium citrate pH 6, 10 ml trace metals solution, 3 ml 1 M CaCl₂, 3 ml 1 M MgSO₄. Add components just before use, using sterile technique; do not autoclave.]

IP buffer [15mM HEPES pH 7.6, 10mM KCl, 1.5mM MgCl₂, 0.1mM EDTA, 0.5mM EGTA, 44mM Sucrose, 100mM NaCl.

Sterilize by autoclaving

Add the following components before use, 1x Protein Inhibitor cocktail from 25x stock (Sigma cat no. #11873580001), 30mM NEM from 300mM(0.1g in 2.67ml) stock in ethanol(Sigma cat no. #E3876)]



- 1 Transfer 20 L4 hermaphrodites each to 10 90cm plates. They will starve out in approximately 6 days. Start the following steps on the day the plates are almost starved.

20 °C 144:00:00
- 2 Add 30ml S Medium to a sterilized 250 ml flask.
- 3 Inoculate the S Medium with a concentrated E. coli OP50 pellet to a final concentration of 30g/L
- 4 Wash each of 5 large plates (90cm) of C. elegans (just cleared of bacteria) with 5 ml S Medium and add to the 30 ml flask.
- 5

84:00:00 20 °C

Put the flask on a shaker at 200rpm
- 6 Put the flask on ice to allow the worms to settle.

00:15:00
- 7 Aspirate most of the liquid from the flask.
- 8 Transfer the remaining liquid to a 15 ml borosil conical tube using a glass pasteur pipette and spin for at least 2 min at 500 × g to pellet the worms. Young larvae may take longer than 2 min to pellet.

00:02:00
- 9 Aspirate the remaining liquid and wash with 10x volume of M9 twice. Spin for at least 2 min at 500 × g to pellet the worms

STEP CASE

No Contamination found 2 steps

In case of no contamination the worms will be separated from OP50 since OP 50 does not settle at 500xg

- 10 Transfer the worms to a 1.5ml eppendorf, incubate on ice for 10 mins

00:10:00

 and remove as much supernatant as possible.



- 11 Add IP Buffer to the Pellet with 2 parts homogenization buffer : 1 part worm pellet and store at -80C for upto 3 months