


Feb 24, 2020

# Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) real-time RT-PCR E gene 2020

DOI

[dx.doi.org/10.17504/protocols.io.bcv9iw96](https://dx.doi.org/10.17504/protocols.io.bcv9iw96)



Judy A Northill<sup>1</sup>, Ian M Mackay<sup>1</sup>

<sup>1</sup>Public Health Virology, Forensic and Scientific Services

Public Health Virology, Fo...

Coronavirus Method Deve...



Judy A Northill

Public Health Virology, Forensic and Scientific Services

OPEN ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.bcv9iw96](https://dx.doi.org/10.17504/protocols.io.bcv9iw96)

**Protocol Citation:** Judy A Northill, Ian M Mackay 2020. Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) real-time RT-PCR E gene 2020. [protocols.io](https://dx.doi.org/10.17504/protocols.io.bcv9iw96) <https://dx.doi.org/10.17504/protocols.io.bcv9iw96>

**Manuscript citation:**

Corman VM, Landt O, Kaiser M, Molenkamp R, Meijer A, Chu DK, et al. Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. *Eurosurveillance*. 2020;25(3):2000045. <https://www.eurosurveillance.org/content/10.2807/1560-7917.ES.2020.25.3.2000045>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** February 24, 2020

**Last Modified:** February 24, 2020

**Protocol Integer ID:** 33441

**Keywords:** CoV, coronavirus, Wuhan, Real-time, RT-PCR, PCR, virus, China, pneumonia, seafood market, WSMPV, sarbecovirus, SARS-CoV-2, COVID-19

## Abstract

A real-time RT-PCR designed to amplify a portion of the envelope gene of sequences from the *Betacoronavirus* sub-genus *Sarbecovirus*.

The probe and primers were published by Corman *et al.*, and we have slightly modified the protocol, increasing the concentration of the reverse primer, using a different kit and different cycling conditions.

This test has identified clinical positive cases of coronavirus disease 2019 (COVID-19).

## Guidelines


- If using a different brand or model of real-time thermocycler, check the concentration of ROX is adequate.
- Method assumes the user is familiar with the thermocycler and software used to run the protocol.

## Materials

### MATERIALS

 SuperScript™ III Platinum™ One-Step qRT-PCR Kit **Life Technologies Catalog #11732088**

### STEP MATERIALS

 SuperScript™ III Platinum™ One-Step qRT-PCR Kit **Life Technologies Catalog #11732088**

## Protocol materials

 SuperScript™ III Platinum™ One-Step qRT-PCR Kit **Life Technologies Catalog #11732088**

In Materials, Materials, Step 2

## Oligonucleotides


1

Oligo name	Sequence 5'-3'	Location*
E_Sarbeco_F1	ACAGGTACGTTAATAGTTAATAGCGT	26269 -26294
E_Sarbeco_R2	ATATTGCAGCAGTACGCACACA	26381 -26357
E_Sarbeco_P1	6FAM-ACACTAGCCATCCTTACTGCGCTTCG-BHQ1	26332 -26357

\*Based on numbering for GenBank accession NC\_045512 Wuhan seafood market pneumonia virus isolate Wuhan-Hu-1

## Reagents

2

 SuperScript™ III Platinum™ One-Step qRT-PCR Kit **Life Technologies** Catalog #11732088

## Synthetic controls

3 Synthetic controls are produced using the **binary synthetic template oligonucleotide positive control for in-house diagnostic real-time RT-PCR method.**

The oligonucleotide sequences required to make controls for this assay are:

Probe control:

AAAATAATACGACTCACTATAGGGTGAAGAGAATCCACAAGGAATTGAAACACTAGCCATCCTTAC  
TGCGCTTCGACAGTGTTTCAGCAGGTCCTGTTGAAAA

Primer control:

AAAATAATACGACTCACTATAGGGACAGGTACGTTAATAGTTAATAGCGTATGATCTGGCACGGGA  
CCCTCCAATGTGTGCGTACTGCTGCAATATAAAA

## Reaction Set-up

- 4
  - Assay has been designed to be used on both a Rotor-Gene 6000 / Rotor-Gene Q 5-plex using 100-place rotor discs and a ABI 7500 Fast real-time machine.
  - Total reaction volume is 20µL.
  - Prepare sufficient for number of reaction plus a 'dead volume' usually 2 extra. Adjust as necessary if using a robotic dispenser.

Reagent	Volume (µl) x1	Final reaction concentration
Nuclease-Free water	4.39	
E_Sarbeco_F1	0.04	400nM
E_Sarbeco_R2	0.09	900nM
E_Sarbeco_P1	0.04	200nM
2X Reaction mix*	10	
Superscript III/Platinum Taq enzyme mix*	0.4	
ROX reference dye (25uM)*	0.04	50nM
<b>TOTAL VOLUME</b>	<b>15</b>	

\*Superscript®III Platinum® One-Step qRT-PCR kit



Dispense 15µl to each reaction well.

Add 5µl of template, extracted RNA, controls or NTC (nuclease-free water).

Total reaction volume is 20µl.

## Amplification

### 5 PCR Amplification

<b>1 cycle</b>	<b>40 cycles</b>
50°C 5 min utes	95°C 3 sec onds
95°C 2 min utes	60°C 30 se conds *

\*Florescence acquisition step

## Result analysis

- 6 The definition used for a satisfactory positive result from a real-time fluorogenic PCR should include each of the following:
1. A sigmoidal curve – the trace travels horizontally, curves upward, continues in an exponential rise and followed by a curve towards a horizontal plateau phase
  2. A suitable level of fluorescence intensity as measured in comparison to a positive control (y-axis)
  3. A defined threshold cycle ( $C_T$ ) value which the fluorescent curve has clearly exceeded (Fig.1 arrow) and which sits early in the log-linear phase and is <40 cycles
  4. A flat or non-sigmoidal curve or a curve that crosses the threshold with a  $C_T$  value >40 cycles is considered a negative result
  5. NTCs should not produce a curve.




Figure 1. Examples of satisfactory sigmoidal amplification curve shape when considering an assay's fluorescent signal output. The crossing point or threshold cycle ( $C_T$ ) is indicated (yellow arrow); it is the value at which fluorescence levels surpass a predefined (usually set during validation, or arbitrary) threshold level as shown in this normalized linear scale depiction. LP-log-linear phase of signal generated during the exponential part of the PCR amplification; TP-a slowing of the amplification and accompanying fluorescence signal marks the transition phase; PP-the plateau phase is reached when there is little or no increase in fluorescent signal despite continued cycling.