


May 25, 2020

Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) real-time RT-PCR CCDC-ORF1ab 2020

DOI

<https://dx.doi.org/10.17504/protocols.io.bgtnjwme>



Judy Northill¹, Ian Mackay¹

¹Public Health Virology, Forensic and Scientific Services

Public Health Virology, F...

Coronavirus Method De...



Judy A Northill

Public Health Virology, Forensic and Scientific Services

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bgtnjwme>

Protocol Citation: Judy Northill, Ian Mackay 2020. Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) real-time RT-PCR CCDC-ORF1ab 2020. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bgtnjwme>

Manuscript citation:

<https://www.who.int/who-documents-detail/molecular-assays-to-diagnose-covid-19-summary-table-of-available-protocols>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 25, 2020

Last Modified: May 25, 2020

Protocol Integer ID: 37454

Keywords: CoV, coronavirus, Wuhan, Real-time, RT-PCR, PCR, virus, China, pneumonia, seafood market, WSMPV, sarbecovirus, SARS-CoV-2, COVID-19, sequences from the severe acute respiratory syndrome coronavirus, severe acute respiratory syndrome coronavirus, coronavirus disease, clinical positive cases of coronavirus disease, pcr ccdc, pcr, orf1ab gene, china cdc, portion of the orf1ab gene, sar, cov

Abstract

A real-time RT-PCR designed to amplify a portion of the ORF1ab gene of sequences from the Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2)

The probe and primers were published by the China CDC and are now included in the WHO in-house assay document.

This test has been modified for use with our standard reagents and instrumentation.

This test has identified clinical positive cases of coronavirus disease 2019 (COVID-19).

Guidelines

- If using a different brand or model of real-time thermocycler, check the concentration of ROX is adequate.
- Method assumes the user is familiar with the thermocycler and software used to run the protocol.

Materials

MATERIALS

 SuperScript™ III Platinum™ One-Step qRT-PCR Kit [Life Technologies Catalog #11732088](#)

STEP MATERIALS

 SuperScript™ III Platinum™ One-Step qRT-PCR Kit [Life Technologies Catalog #11732088](#)

Protocol materials

 SuperScript™ III Platinum™ One-Step qRT-PCR Kit [Life Technologies Catalog #11732088](#)

Oligonucleotides

1

	Oligo name	Sequence 5'-3'	Location*
	CCDC-ORF1ab-F	CCCTGTGGGTTTACACTTAA	13342-13362
	CCDC-ORF1ab-R	ACGATTGTGCATCAGCTGA	13460-13442
	CCDC-ORF1ab-P	6FAM- CCGTCTGCGGTATGTGGAAAGGTTATGG-BHQ1	13377-13404

*Based on numbering for GenBank accession NC_045512.2 Wuhan seafood market pneumonia virus isolate Wuhan-Hu-1

Reagents

2



SuperScript™ III Platinum™ One-Step qRT-PCR Kit Life Technologies Catalog #11732088

Synthetic controls

3

Synthetic controls are produced using the binary synthetic template oligonucleotide positive control for in-house diagnostic real-time RT-PCR method.

The oligonucleotide sequences required to make controls for this assay are:

Probe control:

AAAATAATACGACTCACTATAGGGTGAAGAGAATCCACAAGGAATTGAACCGTCTGCGGTATGTGGAAAGGTTATGGACAGTGTTTCAGCAGGTCCTGTTGAAAA

Primer control:

AAAATAATACGACTCACTATAGGGCCCTGTGGGTTTTACACTTAAATGATCTGGCACGGGAC
CCTCCAATCAGCTGATGCACAATCGTAAAA

Reaction Set-up

- 4
- Assay has been designed to be used on both a Rotor-Gene 6000 / Rotor-Gene Q 5-plex using 100-place rotor discs and a ABI 7500 Fast real-time machine.
 - Total reaction volume is 20 μ L.
 - Prepare sufficient for number of reaction plus a 'dead volume' usually 2 extra. Adjust as necessary if using a robotic dispenser.

Reagent	Volume (μ L) x1	Final reaction concentration
Nuclease-Free water	4.43	
CCDC-ORF1ab-F (200nM)	0.05	500 nM
CCDC-ORF1ab-R (200nM)	0.05	500 nM
CCDC-ORF1ab-P (100nM)	0.03	150nM
2X Reaction mix*	10	
Superscript III/Platinum Taq enzyme mix*	0.4	
ROX reference dye (25uM)*	0.04	50nM
TOTAL VOLUME	15	

*Superscript®III Platinum® One-Step qRT-PCR kit

Dispense 15 μ L to each reaction well.

Add 5 μ L of template, extracted RNA, controls or NTC (nuclease-free water).

Total reaction volume is 20 μ L.

Amplification

5 PCR Amplification

1 cycle	50 cycles
50°C 5 minutes	95°C 30 seconds
95°C 2 minutes	60°C 30 seconds*

*Fluorescence acquisition step

Result analysis

- 6 The definition used for a satisfactory positive result from a real-time fluorogenic PCR should include each of the following:
1. A sigmoidal curve – the trace travels horizontally, curves upward, continues in an exponential rise and followed by a curve towards a horizontal plateau phase
 2. A suitable level of fluorescence intensity as measured in comparison to a positive control (y-axis)
 3. A defined threshold cycle (C_T) value which the fluorescent curve has clearly exceeded (Fig.1 arrow) and which sits early in the log-linear phase and is <40 cycles
 4. A flat or non-sigmoidal curve or a curve that crosses the threshold with a C_T value >40 cycles is considered a negative result
 5. NTCs should not produce a curve.




Figure 1. Examples of satisfactory sigmoidal amplification curve shape when considering an assay's fluorescent signal output. The crossing point or threshold cycle (C_T) is indicated (yellow arrow); it is the value at which fluorescence levels surpass a predefined (usually set during validation, or arbitrary) threshold level as shown in this normalized linear scale depiction. LP-log-linear phase of signal generated during the exponential part of the PCR amplification; TP-a slowing of the amplification and accompanying fluorescence signal marks the transition phase; PP-the plateau phase is reached when there is little or no increase in fluorescent signal despite continued cycling.