




Jun 01, 2020

Version 3

Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) RdRp nested RT-PCR V.3

DOI

<https://dx.doi.org/10.17504/protocols.io.bgz2jx8e>



Judy Northill¹, Ian Mackay¹

¹Public Health Virology, Forensic and Scientific Services

Public Health Virology, F...

Coronavirus Method De...



Judy A Northill

Public Health Virology, Forensic and Scientific Services

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bgz2jx8e>

Protocol Citation: Judy Northill, Ian Mackay 2020. Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) RdRp nested RT-PCR. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bgz2jx8e>

**Manuscript citation:**

Hu D, Zhu C, Wang Y, Ai L, Yang L, Ye F, et al. Virome analysis for identification of novel mammalian viruses in bats from Southeast China. Scientific Reports. 2017;7(1):10917 . <https://www.nature.com/articles/s41598-017-11384-w>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 01, 2020

Last Modified: June 01, 2020

Protocol Integer ID: 37658

Keywords: COVID-19, SARS, coronavirus, Rd-Rp, RT-PCR, PCR, SARS-CoV-2, severe, coronavirus rt, severe acute respiratory syndrome coronavirus, other coronavirus, specific to sar, sar, pcr, nested rt, rdrp region, rdrp,

Abstract

A nested RT-PCR targeting the RdRp region of the sub-genus *Sarbecovirus*. The primers are modified from the pan-coronavirus RT-PCR published by Hu *et al.* 2017 to be more specific to SARS-CoV-2, though other coronaviruses may still be detected.

Assay may be used in resource poor settings where real-time cyclers are not available.

Sanger sequencing should be used to confirm SARS-CoV-2 where WGS is not available or where WGS fails due to poor quality sample.

Materials

STEP MATERIALS

☒ SuperScript III One-Step RT-PCR System with Platinum Taq **Invitrogen - Thermo Fisher Catalog #12574026**

☒ MyFi Mix **Bioline Catalog #BIO-25049**

Protocol materials

☒ SuperScript III One-Step RT-PCR System with Platinum Taq **Invitrogen - Thermo Fisher Catalog #12574026**

☒ MyFi Mix **Bioline Catalog #BIO-25049**

Oligonucleotides

1

	Name	Purpose	5'-3'	Location*
	SARS-CoV-2-15283-F	Round 1 forward primer	ATGGGTTGGGATTATCCTA AATG	15283-15305
	SARS-CoV-2-15887-R	Round 1 reverse primer	TGTTGAGAGCAAAATTCAT G	15887-15868
	SARS-CoV-2-15286-F	Round 2 forward primer	GGTTGGGATTATCCTAAAT GTGA	15286-15308
	SARS-CoV-2-15725-R	Round 2 reverse primer	GCATCGTCAGAGAGTATCA TCAT	15725-15703

*Based on numbering for GenBank accession NC_045512.2 Severe acute respiratory syndrome coronavirus 2 isolate Wuhan-Hu-1

Round 1 amplified products will produce a band of 605bp.

Round 2 amplified products will produce a band of 440bp.

Mix

2



SuperScript III One-Step RT-PCR System with Platinum Taq **Invitrogen - Thermo Fisher** Catalog #12574026



MyFi Mix **Bioline** Catalog #BIO-25049

Round 1 mix

Reagent	Volume x1 (μl)	Final Concentration
Nuclease-free water	3.2	
2 x Reaction mix*	10	1 X
Wuhan15283-F (20pmol/μl)	0.5	500nM



Wuhan15887-R (20pmol/μl)	0.5	500nM
Superscript/Platinum Taq mix*	0.8	
TOTAL	15	

*SuperScript III One-Step RT-PCR System with Platinum Taq #12574026

Dispense mix in 15μl amounts in 0.2ml PCR tubes suitable for thermocycler.

Mix maybe stored frozen if necessary.

Round 2 mix

Reagent	Volume x1 (μl)	Final Concentration
Nuclease-free water	4	
MyFi 2X mix*	10	1X
Wuhan15286-F (20pmol/μl)	0.5	500nM
Wuhan15725-R (20pmol/μl)	0.5	500nM
TOTAL	15	

*MyFi™ Mix Bioline #BIO-25049

Dispense mix in 15μl amounts in 0.2ml PCR tubes suitable for thermocycler.

Mix maybe stored frozen if necessary.

AMPLIFICATION ROUND 1

3

- Thaw required number of tubes with round 1 mix.
- Add 5μl of extract, control, or NTC (nuclease-free water) to round 1 mix above for a final reaction volume of 20μl.
- Label the tubes and record details.

The assay has been used with Eppendorf thermocyclers.

PCR cycling times

1 cycle	40 cycles	1 cycle	Hold
50°C 30 minutes	94°C 30 seconds	68°C 5 minute	15°C
94°C 2 minutes	55°C 30 seconds		
	68°C 1 minute		



AMPLIFICATION ROUND 2

- 4
 - Dilute each of the round 1 amplicons 1/100 in nuclease-free water (2ul + 198ul is ideal)
 - Thaw required number of tubes with round 2 mix.
 - Add 5µl of the diluted round 1 amplicon.
 - Total reaction volume is 20µl.

PCR cycling times

	1 cycle	40 cycles	1 cycle	Hold
	95°C 3 minutes	95°C 15 seconds	72°C 1 minute	15°C
		55°C 15 seconds		
		72°C 15 seconds		

GEL ELECTROPHORESIS


- 5

Amplified products are analysed by gel electrophoresis or equipment such as a QIAxcel or equivalent.

Round 1 amplified products will produce a band of 605bp.

Round 2 amplified products will produce a band of 440bp.

No template controls (NTC) should be not detected.



Example of result from a QIAxcel.
D3-Round 1 band, D8-round 2 band, D4 and D9 are NTC in round 1 and 2 respectively.