

Oct 29, 2023

Version 1

sgRNA library re-amplification in liquid culture V.1

DOI

dx.doi.org/10.17504/protocols.io.n92ldmr7xl5b/v1

Erik Haussner¹, Michael Böttcher¹

¹Medical Faculty of the Martin Luther University Halle-Wittenberg

Atlas of Variant Effects ...



Erik Haussner

PhD

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ldmr7xl5b/v1>

Protocol Citation: Erik Haussner, Michael Böttcher 2023. sgRNA library re-amplification in liquid culture. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n92ldmr7xl5b/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: September 15, 2023

Last Modified: October 29, 2023





Protocol Integer ID: 87835

Keywords: sgRNA, reamplification, Plasmid pool, Library, Electroporation, amplification, sgrna libraries in liquid culture, sgrna library, sgrna, liquid culture, order from addgene, distribution of library element

Abstract

In this protocol, we describe a stepwise procedure for the re-amplification of sgRNA libraries in liquid culture. In our hands, this protocol works reliably to amplify pre-cloned sgRNA libraries (e.g. order from Addgene) in a way that preserves the distribution of library elements.

Materials

-  100 µL  ElectroMAX[®]; Stbl4[®]; Competent Cells **Thermo Fisher Catalog #11635018**
- Up to 400 ng of  Sample
- Electroporator and  Electroporation cuvettes **Biozym Catalog #748010**



Protocol materials

- ✕ Electroporation cuvettes **Biozym Catalog #748010**
- ✕ ElectroMAX[®]; Stbl4[®]; Competent Cells **Thermo Fisher Catalog #11635018**
- ✕ NucleoBond Xtra Midi kit for transfection-grade plasmid DNA **Macherey-Nagel Catalog #REF 740410.50**
- ✕ GelRed[™] Nucleic Acid Gel Stain, 10,000X in Water **Gold Biotechnology Catalog #G-725**
- ✕ 1% Agarose gel **Catalog #/**
- ✕ TriTrack DNA Loading Dye (6X) **Thermo Fisher Catalog #R1161**
- ✕ 1 kb Plus DNA-Ladder **Thermo Fisher Scientific Catalog #10787018**
- ✕ 1 kb Plus DNA-Ladder **Thermo Fisher Scientific Catalog #10787018**
- ✕ ElectroMAX[®]; Stbl4[®]; Competent Cells **Thermo Fisher Catalog #11635018**
- ✕ Electroporation cuvettes **Biozym Catalog #748010**
- ✕ SOC Outgrowth Medium - 100 ml **New England Biolabs Catalog #B9020S**
- ✕ SOC Outgrowth Medium - 100 ml **New England Biolabs Catalog #B9020S**
- ✕ Liquid LB medium
- ✕ ElectroMAX[®]; Stbl4[®]; Competent Cells **Thermo Fisher Catalog #11635018**
- ✕ ElectroMAX[®]; Stbl4[®]; Competent Cells **Thermo Fisher Catalog #11635018**
- ✕ 1.5 mL Eppendorf tubes
- ✕ Liquid LB medium
- ✕ Liquid LB medium
- ✕ LB agar plates with the proper antibiotic (e.g. Kanamycin)
- ✕ Liquid LB medium
- ✕ LB agar plates with the proper antibiotic (e.g. Kanamycin)
- ✕ LB agar plates with the proper antibiotic (e.g. Kanamycin)
- ✕ 1x TBE buffer
- ✕ 1% Agarose gel **Catalog #/**
- ✕ Agarose Low Melt **Carl Roth Catalog #6351.4**
- ✕ Electroporation cuvettes **Biozym Catalog #748010**
- ✕ SOC Outgrowth Medium - 100 ml **New England Biolabs Catalog #B9020S**

Troubleshooting



Before start

Keep your original stock safe and aliquoted. For large plasmids with complementary sequences such as LTR sites, keep in mind that repeated reamplification from an already reamplified stock may lead to an accumulation of recombined plasmids and a poorer distribution of library elements.




Library transformation

25m

1 Prepare Sample

5m


 ElectroMAX[®]; Stbl4[®]; Competent Cells **Thermo Fisher Catalog #11635018**


 Electroporation cuvettes **Biozym Catalog #748010** and

 SOC Outgrowth Medium - 100 ml **New England Biolabs Catalog #B9020S** for electroporation.

1.1 Thaw

5m


 ElectroMAX[®]; Stbl4[®]; Competent Cells **Thermo Fisher Catalog #11635018**

 00:05:00 on ice.

1.2 Pre-cool Electroporation cuvettes **Biozym Catalog #748010** by placing it on ice.


1.3 Pre-warm

 SOC Outgrowth Medium - 100 ml **New England Biolabs Catalog #B9020S** at

 37 °C .

2 Add 100 ng Sample into 25 µL

10m

 ElectroMAX[®]; Stbl4[®]; Competent Cells **Thermo Fisher Catalog #11635018**

, carefully mix by pipetting up and down.



3 Add 25 µL of the plasmid/cell mix into a cuvette, electroporate at 1.2 kV, 25 uF and 200 ohm or alternative setting (see note below). Directly after electroporation, add 1 mL of pre-warmed

10m

 SOC Outgrowth Medium - 100 ml **New England Biolabs Catalog #B9020S** .





Note



The electroporator setting may vary from model to model and should be checked along with the test plasmids provided in the kit of the STBL4 cells.

Safety information

Make sure that any water or ice residue is removed from the cuvettes before inserting them into the electroporator to avoid arcing.

Library recovery

1h 15m

- 4 After electroporation, add the  1 mL resuspended cells in a 14 ml culture tube and incubate the cells in a thermoshaker  600 rpm, 37°C, 01:00:00 .

1h 15m



Note

In general, an incubation temperature of 37°C is optimal for cell recovery. Since *E. coli* tend to recombine plasmids with complementary sequences (e.g. LTRs), recovery temperature can be reduced to 30°C. This may however, result in a lower total number of recovered cells.

Determination of transformation efficiency

16h 30m

- 5 Use a small fraction of your cells to determine the electroporation efficiency of the reamplification.

**Note**

In this step, much depends on the size of the particular plasmid and the number of elements in the library. Therefore, the dilution factor must be chosen based on properties of the library and the scale of the electroporation. Smaller plasmids yield significantly more colonies than large ones, and an upscaled plasmid input at the electroporation step may result in higher dilutions being required to achieve a countable number on the respective agar plates after plating.

5.1 For 1:10,000 dilution:

5m

Prepare 1.5 mL Eppendorf tubes . Take 10 μL of recovery culture and dilute in 990 μL of Liquid LB medium (1:100 dilution). Take 100 μL of 1:100 dilution and dilute in 900 μL of Liquid LB medium (1:1,000) and plate 100 μL on LB agar plates with the proper antibiotic (e.g. Kanamycin) (1:10,000 dilution).

5.2 For 1:1,000,000 dilution:

5m

Take 10 μL of the 1:1,000 dilution and dilute in 990 μL of Liquid LB medium and plate 100 μL on a pre-warmed LB agar plates with the proper antibiotic (e.g. Kanamycin) (1:1,000,000 dilution).

Note

When preparing the dilution series, always mix stock solutions well by flicking the tube before diluting, to resuspend sedimented cells. Distribute the plated cells evenly over the plate by e.g. using glass beads.





6 Place the plates in an incubator at 37 °C Overnight .

16h

**Library extraction and quality control**

16h




- 7 Use rest of recovery to inoculate up to  500 mL of  Liquid LB medium with an added selection marker specific antibiotic like ampicillin in an Erlenmeyer flask for  Overnight culture.  600 rpm, 30°C

16h



Determination of transformation efficiency

15m

- 8 On the next day, check for overall coverage via colony counting on  LB agar plates with the proper antibiotic (e.g. Kanamycin) . The overall colony count should be 1000x the element number of your library.

15m





Note

Below we provide a simplified example for how to determine transformation coverage.

Example calculation of coverage: On the 1:10,000 dilution plate we count 100 colonies. This gives us $100 \times 10,000 = 1,000,000$ total colonies. This total colony number is divided by the number of elements (e.g. sgRNAs) in the respective library. For a library the size of 1,000 sgRNAs, the coverage would $1,000,000 / 1,000 = 1,000x$. For larger libraries, e.g. the size of the genome-wide Brunello library (80,000 sgRNAs), we would count the 1:1,000,000 dilution plate. In this case, 80 counted colonies would mean 80,000,000 total colonies which divided by the library size (80,000 sgRNAs) would again return a transformation coverage of 1,000x.

Library preparation and QC

1h

- 9 Follow the protocol instructions of the  NucleoBond Xtra Midi kit for transfection-grade plasmid DNA **Macherey-Nagel Catalog #REF 740410.50** for transfection-grade plasmid DNA for Midi Prep. Follow the protocol instructions of the for transfection-grade plasmid DNA for Midi Prep.
- 10 Determine your final  Sample concentration via NanoDrop or Qubit measurement.

25m



Equipment

new equipment

NAME

Qubit 2.0 Fluorometer instrument

BRAND

Q33226

SKU

with Qubit RNA HS Assays

SPECIFICATIONS

Equipment

NanoDrop™ One/OneC Microvolume UV-Vis Spectrophotometer

NAME

UV-Vis Spectrophotometer

TYPE

Thermo Scientific

BRAND

ND-ONE-W

SKU

- 11 Send a sample of your reamplified  Sample for Sanger sequencing.



Below we show an example chromatogram of an expected sequencing result. We recommend using sequencing primers 50-100 nt upstream of the sgRNA region. You should see clean traces up- and downstream of the SPACER region, and a noisy 20 nt signal in the SPACER region, due to the sgRNA diversity in your library.














Example of a Sanger sequencing result.

Note











! To validate the distribution of elements in your library, we strongly recommend performing next generation sequencing of your plasmid pool before proceeding with downstream experiments. To do so, follow the NGS protocol provided with your library, using the plasmid pool as template, instead of the genomic DNA (as you would in CRISPR screens).

QC: Plasmid recombination check

3h 15m

- 12 Since sgRNA library plasmids can recombine during E.coli re-amplification, it is recommended to check for recombination via linearization of  200 ng of your reamplified  Sample via a restriction digest within the backbone of your library vector.
- 13 Prepare a  1% Agarose gel **Catalog #/** by melting  1 g of  Agarose Low Melt **Carl Roth Catalog #6351.4** in  100 mL  1x TBE buffer
- 14 Let the required amount for casting cool down till it is approxametyly  50 °C and add  1 µL of  GelRed™ Nucleic Acid Gel Stain, 10,000X in Water **Gold Biotechnology Catalog #G-725** per ml of melted  1% Agarose gel **Catalog #/**
- 15 Pour the warm, still liquid gel into an electroporation chamber and wait until it has cooled down.



- 16 Mix your linearized reamplified  Sample with  5m
 TriTrack DNA Loading Dye (6X) **Thermo Fisher Catalog #R1161** and prepare
 1 kb Plus DNA-Ladder **Thermo Fisher Scientific Catalog #10787018** for gel loading. 
- 17 Add your linearized reamplified  Sample alongside with the prepared  1h
 1 kb Plus DNA-Ladder **Thermo Fisher Scientific Catalog #10787018** onto the gel and run it for  01:00:00 at 120 V.
- 18 Check the plasmid size on your gel using UV excitation.  5m

Expected result

In the best case, only one band will be present, corresponding to the size of the respective vector. In case of recombination, one or more additional bands will be present. For optimal downstream results, the band of the intact vector should be dominant.

