

Aug 23, 2018

Short insert size WGS libraries preparation for assembly of the *Lateolabrax maculatus* genome

 [GigaScience](#)

 In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.sszeeef6

Chang Li¹

¹BGI-Qingdao, BGI-Shenzhen, Qingdao, 266555, China.

GigaScience Press

BGI



Chang Li

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.sszeeef6

External link: <https://doi.org/10.1093/gigascience/gy114>

Protocol Citation: Chang Li 2018. Short insert size WGS libraries preparation for assembly of the *Lateolabrax maculatus* genome. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.sszeeef6>

Manuscript citation:

hangwei Shao, Chang Li, Na Wang, Yating Qin, Wenteng Xu, Qun Liu, Qian Zhou, Yong Zhao, Xihong Li, Shanshan Liu, Xiaowu Chen, Shahid Mahboob, Xin Liu, Songlin Chen, Chromosome-level genome assembly of the spotted sea bass, *Lateolabrax maculatus*, *GigaScience*, Volume 7, Issue 11, November 2018, gy114, <https://doi.org/10.1093/gigascience/gy114>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 22, 2018

Last Modified: August 24, 2018


Protocol Integer ID: 14905


Abstract

This protocol is used to clarify the process of the short insert size WGS libraries preparation for the *L. maculatus*.

Materials

STEP MATERIALS


 AMPure XP beads **Beckman Coulter**


 T4 DNA polymerase **Enzymatics**


 AMPure XP beads **Beckman Coulter**

 T4 DNA polymerase **Enzymatics**

Protocol materials

 AMPure XP beads **Beckman Coulter**

 T4 DNA polymerase **Enzymatics**

 AMPure XP beads **Beckman Coulter**

 T4 DNA polymerase **Enzymatics**


 AMPure XP beads **Beckman Coulter**

 T4 DNA polymerase **Enzymatics**



Genomic DNA interruption

- 1 The extracted DNA was sheared into fragments between 50 bp and ~800 bp in size using a Covaris E220 ultrasonicator (Covaris, Brighton, UK). Treat time 20s, Acoustic Duty Factor 25%, Peak Incident Power 500W, Cycles Per Burst 500, 24 cycles.

 00:00:20

Fragment selection


- 2 Fragments between 150 bp and 250 bp or 200 to 500 bp were selected using AMPure XP beads (Agencourt, Beverly, the U.S.).

 AMPure XP beads **Beckman Coulter**

End-repair

- 3 Repaired using T4 DNA polymerase, (ENZYMATICS, Beverly, the U.S.) 30 min. at 20 °C to obtain blunt ends which were then 3'-adenlyated to create sticky ends.

 T4 DNA polymerase **Enzymatics**

 00:30:00

 20 °C

Add adapter

- 4 These DNA fragments were ligated at both ends to T-tailed adapters and amplified.

PCR amplification

- 5 The temperature profile was 3 min. at 95 °C followed by 8 cycles of 20 sec. at 98 °C, 15 sec. at 60 °C, 30 sec. at 72 °C, and more 10 min. at 72 °C for further elongation.

Library purification

- 6 AMPure XP beads (Agencourt, Beverly, the U.S.) was used to purify the PCR production.