

Dec 16, 2019

# Short term effect of Aldicarb on C.elegans

DOI

dx.doi.org/10.17504/protocols.io.baivice6

Priota Islam<sup>1</sup>

<sup>1</sup>Imperial College London

**Behavioural Genomics** 



#### Priota Islam

Imperial College London

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account





DOI: https://dx.doi.org/10.17504/protocols.io.baivice6

Protocol Citation: Priota Islam 2019. Short term effect of Aldicarb on C.elegans. protocols.io

https://dx.doi.org/10.17504/protocols.io.baivice6

License: This is an open access protocol distributed under the terms of the <u>Creative Commons Attribution License</u>, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: December 16, 2019

Last Modified: December 16, 2019

Protocol Integer ID: 31029

Keywords: cholinesterase inhibitor, acetylcholine, breakdown of acetylcholine, short term effect of aldicarb, aldicarb

#### Abstract

- Aldicarb is a cholinesterase inhibitor which prevents the breakdown of acetylcholine in the synapse
- Reported to paralyze C. elegans
- Strains of interest believed to be resistant to the paralytic effect
- Strains: OW1601 & OW1603

## **Troubleshooting**



- 1 Bleach synchronize the worms on a Friday
- 2 Following Monday refeed the arrested L1s on OP50 seeded plates, make 2plates per strain keeping one plate at 20C and the other at 25C
- 3 The day before tracking (Wednesday) add 35ul of 2.5mM AK to 24 imaging plates
- 4 Also, seed those drugs treated 24 imaging plates and additional 24 non-treated imaging plates with 50ul of 1:10 OP50 and leave to dry O/N
- 5 On the day of tracking (Thursday), transfer 5 worms per strain on to the aldicarb treated plates and leave for 1.5hrs, pick onto the controlled plates too
- 6 After 1.5hrs image for 15mins on Phenix

Track 8 sets: (2 treatment sets + 2 control sets with no treatment) x 2 growth temperatures