

Dec 16, 2019

## Short term effect of Aldicarb on C.elegans

DOI

[dx.doi.org/10.17504/protocols.io.baivice6](https://dx.doi.org/10.17504/protocols.io.baivice6)

Priota Islam<sup>1</sup>

<sup>1</sup>Imperial College London

Behavioural Genomics



Priota Islam

Imperial College London

---

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.baivice6](https://dx.doi.org/10.17504/protocols.io.baivice6)

**Protocol Citation:** Priota Islam 2019. Short term effect of Aldicarb on C.elegans. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.baivice6>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** December 16, 2019

**Last Modified:** December 16, 2019

**Protocol Integer ID:** 31029



## Abstract

- Aldicarb is a cholinesterase inhibitor which prevents the breakdown of acetylcholine in the synapse
- Reported to paralyze *C. elegans*
- Strains of interest believed to be resistant to the paralytic effect
- Strains: OW1601 & OW1603



- 1 Bleach synchronize the worms on a Friday
- 2 Following Monday refeed the arrested L1s on OP50 seeded plates, make 2 plates per strain keeping one plate at 20C and the other at 25C
- 3 The day before tracking (Wednesday) add 35ul of 2.5mM AK to 24 imaging plates
- 4 Also, seed those drugs treated 24 imaging plates and additional 24 non-treated imaging plates with 50ul of 1:10 OP50 and leave to dry O/N
- 5 On the day of tracking (Thursday), transfer 5 worms per strain on to the aldicarb treated plates and leave for 1.5hrs, pick onto the controlled plates too
- 6 After 1.5hrs image for 15mins on Phenix

Track 8 sets: (2 treatment sets + 2 control sets with no treatment) x 2 growth temperatures