



Mar 06, 2023

## Sigma GeneElute total RNA extraction from fungal tissue

DOI

[dx.doi.org/10.17504/protocols.io.14egn7znyv5d/v1](https://dx.doi.org/10.17504/protocols.io.14egn7znyv5d/v1)

ALEXANDER J BRADSHAW<sup>1</sup>

<sup>1</sup>Natural History Museum of Utah



ALEXANDER J BRADSHAW

University of Utah

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn7znyv5d/v1>

**Protocol Citation:** ALEXANDER J BRADSHAW 2023. Sigma GeneElute total RNA extraction from fungal tissue. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.14egn7znyv5d/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** January 24, 2022



**Last Modified:** March 07, 2023

**Protocol Integer ID:** 57351

**Keywords:** rna extraction from fungal tissue rna extraction, fungal tissue rna extraction, sigma geneelute total rna kit, total rna extraction, rna extraction, fungal tissue, total rna, sigma geneelute, rna, extraction, sigma

## Abstract

RNA extraction from fungal tissue using the Sigma GeneElute total RNA kit.

## Troubleshooting



## Tissue and workspace preparation

15m

- 1 ALL WORK MUST BE DONE IN THE STERILE HOOD THAT HAS BEEN CLEANED FOR RNA WORK



REGULARLY CHNAGE GLOVES AND WASH HANDS WITH RNASE TO AVOID  
DEGREDEATION

MAKE SURE ALL PIPETORS, TIPS, AND ANY OTHER EQUIPMENT ARE IN THE HOOD  
AND CLEANED WITH RNASE

DO NOT MOVE ANYTHING IN OR OUT WITHOUT WASHING WITH RNASE

Reagents that need to be made:  
Fresh, day of 70% etoh

- 2 Flash Freeze tissue in Liquid Nitrogen (LN2).
- 3 Grind Fungal Tissue to a fine powder with a mortar and pestle with LN2, be careful to not allow the tissue to thaw during this Process.
- 3.1 At this point frozen ground tissue can be stored at -80C for later extraction.

## RNA extraction

15m

- 4 add 600ul of buffer RL to a 1.7ml eppendorf tube.
- 5 Place **no more than 50mg** of ground tissue into the sample tube with buffer RL.
- 6 Incubate for 15-30 minutes at room temperature Room temperature
- 7 Centrifuge at maximum speed for 2 minutes

15m



## Equipment

### Centrifuge

NAME

Benchtop Centrifuge

TYPE

Eppendorf

BRAND

5405000441

SKU

<https://online-shop.eppendorf.us/US-en/Centrifugation-44533/Centrifuges-44534/Centrifuge-5425-PF-243560.html>

LINK

Any benchtop centrifuge will suffice

SPECIFICATIONS

- 8 Transfer the lysate to a new RNase free 1.7ml tube.
- 9 Add 1:1 volume (so 600ul) of **freshly** prepared 70% etoh
- 10 Vortex for 10 seconds, or until homogenous
- 11 Assemble a column with a waste tube from the GeneElute kit.
- 12 Add 600ul of lysate/etoh solution to the assembled column
- 13 Centrifuge at 6000G for 2 minuts
- 14 If lysate is not completely passes, spin at an additional 14,000G for 1 minute
- 15 Discard waste from waste tube, and repeat step 13 and 14 if more lysate/etoh exists.



- 16 Place column into a new waste collection tube
- 17 Add 400ul of wash solution A to column
- 18 Centrifuges at max speed for 1.5 minutes
- 19 discard waste from collection tube
- 20 Repeat steps 17-19 two more times (for a total of 3x wash steps)
- 21 Centrifuge Sample for a final 2 minute dry spin at max speed
- 22 Place column in a clean eppendorf tube provided in the GeneElute kit.
- 23 add 50ul of elution solution A
- 24 Centrifuge at 200G for 2 minutes
- 25 centrifuge at 14,000G for 1 minute
- 26 If the entire 50ul has not eluted, then spin for an additional minute at 14,000G

## Quality control of RNA and storage

15m

- 27 Keep eluted samples on Ice, Prepare a bleach gel for RNA analysis
- 28 BLEACH GEL:



100ml TAE  
400ul concentrated bleach  
1g agarose

Make sure to add bleach **before** microwaving.

29 Load and run samples normally at 120V for ~20 minutes with a 100bp ladder.

30 immediately store samples at -80C after loading gel.