

Apr 18, 2020 Version 2

🌐 Simple Step to Increase Gibson Assembly Efficiency V.2

DOI

dx.doi.org/10.17504/protocols.io.be9cjh2w

Eric Danner¹

¹Max Delbrück Center

Eric D

MDC

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.be9cjh2w

Protocol Citation: Eric Danner 2020. Simple Step to Increase Gibson Assembly Efficiency . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.be9cjh2w>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 18, 2020

Last Modified: April 18, 2020

Protocol Integer ID: 35844


Abstract

The Gibson assembly solutoin is extremely toxic to some bacteria cells.

Doing a solution change after assembly and before transforming the bacteria will significantly improve bacteria survival. This increases 10-100x the amount of colonies.

Materials

MATERIALS

 Gibson Assembly Master Mix - 50 rxns **New England Biolabs Catalog #E2611L**

 NucleoSpin® Gel and PCR Clean-up **Macherey-Nagal Catalog #740609.10**



Do Normal Assembly

- 1 Do a normal Gibson assembly. Incubate in the PCR machine 60 min @ 50C
- 2 Remove the assemblies. Do a PCR quickchange/ gel extract column. This will remove the Gibson Assembly solution but keep the 100ng (or however much) of DNA.

Elute in 15u -20ull TE or water

- 3 Take the 15-20ul assembled plasmid and heat shock your bacteria. 30 second heat shock. 2 min ice. Add 200ul media and shake for 45 min. Carb can be plated immediatly.

Plate