



Sep 26, 2022

## Simultaneous detection of miRNA and mRNA at the single-cell level in plant tissues

DOI

<https://dx.doi.org/10.17504/protocols.io.x54v9deppg3e/v1>

Chi-Chih Wu<sup>1</sup>

<sup>1</sup>Biodiversity Research Center, Academia Sinica, Taipei 11529, Taiwan



Chi-Chih Wu

IPMB

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.x54v9deppg3e/v1>

**Protocol Citation:** Chi-Chih Wu 2022. Simultaneous detection of miRNA and mRNA at the single-cell level in plant tissues. [protocols.io https://dx.doi.org/10.17504/protocols.io.x54v9deppg3e/v1](https://dx.doi.org/10.17504/protocols.io.x54v9deppg3e/v1)

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** August 28, 2022

**Last Modified:** September 26, 2022

**Protocol Integer ID:** 69270

**Keywords:** mirna molecules in specific cell, mirna molecule, transgenic rice lines with different gene expression level, mirna, transgenic rice line, lna probe, cell level in plant tissue, mrna, expression levels of mrna, plant tissue, gene expression, different gene expression level, locked nucleic acid, single cell, detection signal

## Abstract

Detecting the simultaneous presence of a microRNA (miRNA) and a mRNA in a specific tissue can provide support for the prediction that the miRNA regulates the mRNA. We develop a method that uses sequence-specific miRNA-locked nucleic acid (LNA) and mRNA-LNA probes. Moreover, it augments the detection signal by rolling circle amplification, achieving a high signal-noise ratio at the single-cell level. Dot signals are counted for determining the expression levels of mRNA and miRNA molecules in specific cells. We show a high sequence specificity of our miRNA-LNA probe, revealing that it can discriminate single-base mismatches. Numerical quantification by our method is tested in transgenic rice lines with different gene expression levels.



## 1 miRNA hybridization

| A                                                                                                              |
|----------------------------------------------------------------------------------------------------------------|
| 1. The slides with sections are taken out the freezer and equilibrated to RT for 40 mins.                      |
| 2. Permeabilized in 20 ug/ml proteinase k for proper duration                                                  |
| 3. Quickly wash in DEPC-PBS                                                                                    |
| 4. Quickly dehydrate the slides in EtOH (50, 70, 99 %) and then air dry                                        |
| 5. Mount the secure seal reaction chambers onto the slides.                                                    |
| 6. Add 1x DEPC-PBS-tween 0.05 % (Wash buffer) into the chambers to keep the slides wet until RT reaction ready |

## 2 Section permeabilization

| A                     | B        | C         |
|-----------------------|----------|-----------|
|                       | Stock    | Final     |
| Formamide             | 100%     | 50%       |
| SSC                   | 20x%     | 5x        |
| tRNA                  | 10 mg/ml | 0.5 ug/ul |
| Denhardt's            | 50x      | 1x        |
| LNA probe A           | 10 uM    | 2-3 pmole |
| DEPC-H <sub>2</sub> O |          |           |

### 1. process

| A                                                                              |
|--------------------------------------------------------------------------------|
| Hybridization under the predicted melting temperature of the probe for an hour |
| Wash with 0.1X SSC three times at the temperature set in the Step 1            |



|  |                                                      |
|--|------------------------------------------------------|
|  | A                                                    |
|  | Wash with 2X SSC at RT once                          |
|  | Wahs wtih the wash buffer (PBS, 0.05% Tween-20) once |

### 3 mRNA cDNA synthesis

|  |                          |          |           |
|--|--------------------------|----------|-----------|
|  |                          |          |           |
|  | Reagent                  | Stock    | Final     |
|  | NEB Tag DNA ligase       | 40U/ul   | 0.5 U/ul  |
|  | Rnase H                  | 5 U/ul   | 0.4 U/ul  |
|  | Ribolock Rnase inhibitor | 40 U/ul  | 1 U/ul    |
|  | NEB Tag ligase buffer    | 10x      | 1x        |
|  | BSA                      | 20 ug/ul | 0.2 ug/ul |
|  | KCl                      | 1 M      | 0.05 M    |
|  | Formamide                | 100%     | 20%       |
|  | Pd_A                     | 10 uM    | 0.1 uM    |
|  | Pd_B                     | 10 uM    | 0.1 uM    |
|  | Pd_C                     | 10 uM    | 0.1 uM    |
|  |                          |          |           |
|  | DEPC-H2O                 |          |           |

#### Process

|  |                                                              |
|--|--------------------------------------------------------------|
|  |                                                              |
|  | 1. Add ligation mixture in chambers, seal with adhesive film |
|  | 2. Incubate for 30 min at 37 C followed by 45 min at 48 C    |



|                                               |
|-----------------------------------------------|
| 3. Wash 2x, 1x<br>DEPC-PBS-Tween<br>20, 0.05% |
|-----------------------------------------------|

#### 4 miRNA and mRNA padlock probe hybridization and ligation

| Reagent                  | Stock    | Final     |
|--------------------------|----------|-----------|
| NEB Tag DNA ligase       | 40U/ul   | 0.5 U/ul  |
| Rnase H                  | 5 U/ul   | 0.4 U/ul  |
| Ribolock Rnase inhibitor | 40 U/ul  | 1 U/ul    |
| NEB Tag ligase buffer    | 10x      | 1x        |
| BSA                      | 20 ug/ul | 0.2 ug/ul |
| KCl                      | 1 M      | 0.05 M    |
| Formamide                | 100%     | 20%       |
| Pd_A                     | 10 uM    | 0.1 uM    |
| Pd_B                     | 10 uM    | 0.1 uM    |
| Pd_C                     | 10 uM    | 0.1 uM    |
|                          |          |           |
| DEPC-H2O                 |          |           |

process

| A                                                            |
|--------------------------------------------------------------|
| 1. Add ligation mixture in chambers, seal with adhesive film |
| 2. Incubate for 30 min at 37 C followed by 45 min at 48 C    |
| 3. Wash 2x, 1x DEPC-PBS-Tween 20, 0.05%                      |

## 5 rolling circle amplification

|  | Reagent                  | Stock    | Final     |
|--|--------------------------|----------|-----------|
|  | Phi 29 polymerase        | 10 U/ul  | 1 U/ul    |
|  | Ribolock Rnase inhibitor | 40 U/ul  | 1 U/ul    |
|  | 10 X phi 29 buffer       | 10 x     | 1 x       |
|  | dNTP                     | 10 mM    | 0.25 mM   |
|  | BSA                      | 20 ug/ul | 0.2 ug/ul |
|  | Glycerol                 | 50%      | 5%        |
|  |                          |          |           |
|  | DEPC-H <sub>2</sub> O    |          |           |

process

|                                          |
|------------------------------------------|
| A                                        |
| 1. Add reaction mixture and seal chamber |
| 2. Incubate for over night at 37 C       |
| 3. Wash 2x, DEPC-PBS-Tween 20, 0.05 %    |

## 6 detection oligo hybridization

| A                      | B     | C      |
|------------------------|-------|--------|
| Reagent                | Stock | Final  |
| Hyb mixture            | 4 x   | 2 x    |
| Detection oligo 1-FITC | 1 uM  | 0.1 uM |



|  | A                     | B    | C      |
|--|-----------------------|------|--------|
|  | Detection oligo 2-Cy3 | 1 uM | 0.1 uM |
|  | Detection oligo 3-Cy5 | 1 uM | 0.1 uM |
|  |                       |      |        |
|  | DEPC-H2O              |      |        |

process

|  |                                                   |
|--|---------------------------------------------------|
|  |                                                   |
|  | 1. Add reaction mixture                           |
|  | 2. Incubate for 30 min at 37 C                    |
|  | 3. Wash 2x, DEPC-PBS-Tween 20, 0.05 %             |
|  | 4. Dehydrated by EtOH 50, 70, 99 %; then air dry. |
|  | 5. Mount cover slips                              |