



Mar 30, 2023

## Single cell dissociation of brain organoids

 In 2 collections

DOI

[dx.doi.org/10.17504/protocols.io.e6nvwj799lmk/v1](https://dx.doi.org/10.17504/protocols.io.e6nvwj799lmk/v1)

michela.deleidi<sup>1</sup>, Maria Jose Perez J.<sup>1</sup>, Hariam Raji<sup>1</sup>

<sup>1</sup>German Center for Neurodegenerative Diseases (DZNE), Tübingen, 72076 Germany



Federico Bertoli

ASAP

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwj799lmk/v1>

**Protocol Citation:** michela.deleidi, Maria Jose Perez J., Hariam Raji 2023. Single cell dissociation of brain organoids.  
**protocols.io** <https://dx.doi.org/10.17504/protocols.io.e6nvwj799lmk/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** March 28, 2023

**Last Modified:** May 31, 2024

**Protocol Integer ID:** 79538

**Keywords:** brain organoids, cell dissociation, Papain Dissociation System, ASAPCRN, single cell dissociation of brain organoid, brain organoid, single cell dissociation, protocol details about single cell dissociation, cell dissociation, cell, brain

## Abstract

This protocol details about single cell dissociation of brain organoids.

## Attachments



[405-879.docx](#)

15KB




## Materials

### Kit:

Papain Dissociation System.

 Papain Dissociation System **Worthington Biochemical Corporation Catalog #LK003150**

### Reconstitute powders.

- Add  5 mL Earle's medium into Papain Vial (1 Vial/2 organoids).
- Add  500 µL Earle's medium into DNase vial.
- Add  35 mL Earle's medium into Inhibitor vial (1 vial/10 organoids).

## Troubleshooting



## Single cell dissociation of brain organoids

27m

- 1 Mix  500  $\mu$ L DNase with  5 mL Papain. 


### Note

**Note:** MIX GENTLY.



- 2 Transfer single or pooled organoid to 60 mm dish.




- 3 Aspirate excess media, add  2.5 mL Papain + DNase solution. 


- 4 With a razor blade mince organoid (<1 mm).


- 5 Transfer plate to an orbital shaker  70 rpm, 00:30:00 (inside incubator).

- 6 With 1-mL pipette dissociate pieces (Mix up-down 30 times).

- 7 Put in orbital shaker  00:20:00 . 

- 8 In the meantime, add  5 mL Earle's medium +  3 mL Inhibitor to a 15-mL conical tube. 


- 9 Remove samples from the orbital shaker. With a 1-mL tip, mix up-down 30 times. 

- 10 Take  2 mL (upper part) into new tube using a 40  $\mu$ m cell strainer. Wait 1-3 min to debris to settle.





11 Transfer cell suspension to the inhibitor tube. Invert to mix 5 times.



12 Centrifuge  300 rpm, Room temperature, 00:07:00 .

7m



13 Aspirate supernatant, resuspend in  500  $\mu\text{L}$  to  1 mL 0.5% BSA-PBS (Up-down 30 times).

14 Filter the resuspended cells (  900  $\mu\text{L}$  ) with a 30  $\mu\text{m}$  cell strainer.

15 Count the cells for the final suspension and dilute. Resuspend at 1000 cells/ $\mu\text{L}$  in 0.04% BSA-PBS.