

Jun 21, 2023

Single cell dissociation of healthy paediatric skin

DOI

dx.doi.org/10.17504/protocols.io.e6nvwd6kwlmk/v1

Emily Stephenson^{1,2}, Chloe Admane^{1,2}, Keval Sidhpura^{1,2}

¹Newcastle University; ²Wellcome Sanger Institute

Human Cell Atlas Metho...



Emily Stephenson

Newcastle University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwd6kwlmk/v1>

Protocol Citation: Emily Stephenson, Chloe Admane, Keval Sidhpura 2023. Single cell dissociation of healthy paediatric skin . protocols.io <https://dx.doi.org/10.17504/protocols.io.e6nvwd6kwlmk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: June 19, 2023

Last Modified: June 21, 2023

Protocol Integer ID: 83638

Keywords: single cell dissociation of healthy paediatric skin, single cell dissociation, healthy paediatric skin, method for the enzymatic dissociation, enzymatic dissociation, single cell suspension, cell

Funders Acknowledgements:

CZI Pediatric Networks for the Human Cell Atlas

Disclaimer

This protocol has been tested on a variety of skin sites such as lip and trunk, however not all body sites has been tested. This protocol has also only been used for healthy, non-diseased skin.

Abstract

This protocol outlines the method for the enzymatic dissociation of healthy paediatric skin >3mm into a single cell suspension.

Guidelines

This protocol includes an overnight incubation step.

Materials

Petri dish

Scalpel

Forceps

100 micron fiilters

RPMI

RF-10 (RPMI plus 10%FCS, 1% Pen-strep and 1% L-glut)

48-well v bottom plate

PBS

Dispase (Roche)

Collagenase Type IV (Worthington)

50ml Falcon Tubes

Flow Buffer (PBS 2% FCS and 2mM EDTA)

Troubleshooting



Safety warnings


! This protocol uses sharp objects.

Before start

Do not forget to record the metadata for this sample. Before starting, clean MSC Class II with 70% ethanol and make up virkon. Ensure you have all materials needed to carry out the protocol.



Day 1 - Begin protocol late afternoon

- 1 Record sample meta data and assess size of sample  Sample 10m
- 1.1 If sample is <3mmx3mm freeze and embed sample in OCT (see other protocol), if sample is >3mmx3mm continue with protocol
- 2 Empty sample onto petri dish and wash sample with PBS 5m
- 3 Cut off lower dermis and fat and place in well of 48-well plate with 1ml RPMI 2m
- 4 Place epidermis and upper dermis sample in a 48 well V-bottom plate with 1ml RPMI and 20µL Dispase 5m
- 5 Add parafilm to plate and leave in 4°C fridge overnight 1m

Day 2 - Begin protocol early in the morning

4h 8m

- 6 Take plate out of fridge and empty epidermis/upper dermis onto new petri dish, remove collagenase type IV from -80 freezer 2m
- 7 Using forceps, separate epidermis and upper dermis and place each separately into new wells of the 48-well plate 10m
- 8 Place 1ml RPMI in each well of plate containing epidermis and upper dermis 5m
- 9 Add collagenase type IV 1:100 (10µL) to each of the tissues (epidermis, upper dermis and lower dermis) 5m
- 10 Place in incubator and incubate at 37°C for 3 hours 3h
- 11 Remove from incubator and, using 1ml pipette, pipette each tissue up and down to ensure the tissue has dissociated 10m



- 12 Pipette through 100micron filter into separate 50ml Falcon tubes 5m
- 13 Wash out each well an additional 3 times with 1ml RF-10 and pass through filter 2m
- 14 Wash filter with 25ml RF-10 and adjust so each Falcon tube has the same volume 2m
- 15 Centrifuge at 500g for 5 mins at 4deg (9acc/dec) 500 x g, 4°C 5m
- 16 Pour off supernatant 2m
- 17 Resuspend pellet with 1ml of Flow Buffer 5m
- 18 Take 10µL for cell count and count using Trypan Blue and C-Chip Haemocytometer, record number of cells isolated for each tissue 15m
- 18.1 If sorting for single cell RNA seq, continue with antibody staining and FACS protocol or if freezing down cells, continue with viable cell freezing protocol