


May 06, 2020

🌐 Single-Cell Dissociation of Human Trabecular Meshwork

DOI

<https://dx.doi.org/10.17504/protocols.io.bfdyji7w>



Alexi Mcadams^{1,2}, Tavé an Zyl^{1,2}

¹Department of Ophthalmology, Harvard Medical School; ²Glaucoma Service, Mass. Eye and Ear, Boston, MA

Human Cell Atlas Metho...



Alexi Mcadams

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bfdyji7w>

Protocol Citation: Alexi Mcadams, Tavé an Zyl 2020. Single-Cell Dissociation of Human Trabecular Meshwork . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bfdyji7w>

**Manuscript citation:**

van Zyl T, Yan W, McAdams A, Peng YR, Shekhar K, Regev A, Juric D, Sanes JR. Cell atlas of aqueous humor outflow pathways in eyes of humans and four model species provides insight into glaucoma pathogenesis. Proceedings of the National Academy of Sciences. 2020 Apr 27.

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 21, 2020

Last Modified: May 06, 2020

Protocol Integer ID: 35992

Keywords: trabecular meshwork, scRNAseq, single cell, cell dissociation of human trabecular meshwork, fresh trabecular meshwork tissue, human trabecular meshwork, cell dissociation, single cell suspension, cell

Abstract

Fresh trabecular meshwork tissue is dissected from post-mortem human eyes and dissociated into a single cell suspension.

Guidelines

The stock solutions can be prepared up to 6 months in advance and aliquots can be stored at -20°C.



Materials

MATERIALS

- ✕ Sodium bicarbonate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S6014**
- ✕ RNase Zap **Merck MilliporeSigma (Sigma-Aldrich) Catalog #R2020-250ML**
- ✕ Razor blade
- ✕ Papain **Worthington Biochemical Corporation Catalog #LS003126**
- ✕ Ames' Medium **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A1420**
- ✕ Bovine Serum Albumin (Non-acetylated) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #B6917**
- ✕ Bovine Serum Albumin **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A9418**
- ✕ Deoxyribonuclease I from bovine pancreas **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D4527**
- ✕ Trypsin Inhibitor Ovomucoid **Worthington Biochemical Corporation Catalog #LS003087**
- ✕ L-Cysteine hydrochloride **Merck MilliporeSigma (Sigma-Aldrich) Catalog #C1276**
- ✕ 40 µm Falcon® Cell Strainers Sterile Corning **VWR International (Avantor) Catalog #21008-949**
- ✕ Millex-GP Syringe Filter Unit 0.22 µm polyethersulfone 33 mm gamma sterilized **Merck Millipore (EMD Millipore) Catalog #SLGP033RS**
- ✕ 70% EtOH
- ✕ Disposable Petri Dishes
- ✕ Blunt curved Westcott scissors
- ✕ Vannas Microscissors
- ✕ Jewellers Forceps
- ✕ 50 mL Falcon Tubes
- ✕ 15 mL Falcon Tubes
- ✕ 1.7 mL Eppendorf Tubes

STEP MATERIALS

- ✕ L-Cysteine hydrochloride **Merck MilliporeSigma (Sigma-Aldrich) Catalog #C1276**
- ✕ Sodium bicarbonate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S6014**
- ✕ Ames' Medium **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A1420**
- ✕ Deoxyribonuclease I from bovine pancreas **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D4527**
- ✕ Bovine Serum Albumin (Non-acetylated) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #B6917**
- ✕ Bovine Serum Albumin **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A9418**
- ✕ Trypsin Inhibitor Ovomucoid **Worthington Biochemical Corporation Catalog #LS003087**



- ✕ Sodium bicarbonate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S6014**
- ✕ Ames' Medium **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A1420**
- ✕ Papain **Worthington Biochemical Corporation Catalog #LS003126**
- ✕ RNaseZap™ RNase Decontamination Solution **Thermo Fisher Scientific Catalog #AM9780**

Protocol materials

- ✕ Sodium bicarbonate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S6014**
- ✕ Trypsin Inhibitor Ovomucoid **Worthington Biochemical Corporation Catalog #LS003087**
- ✕ Ames' Medium **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A1420**
- ✕ Deoxyribonuclease I from bovine pancreas **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D4527**
- ✕ Disposable Petri Dishes
- ✕ Jewellers Forceps
- ✕ 15 mL Falcon Tubes
- ✕ Razor blade
- ✕ RNase Zap **Merck MilliporeSigma (Sigma-Aldrich) Catalog #R2020-250ML**
- ✕ Millex-GP Syringe Filter Unit 0.22 µm polyethersulfone 33 mm gamma sterilized **Merck Millipore (EMD Millipore) Catalog #SLGP033RS**
- ✕ Papain **Worthington Biochemical Corporation Catalog #LS003126**
- ✕ L-Cysteine hydrochloride **Merck MilliporeSigma (Sigma-Aldrich) Catalog #C1276**
- ✕ 50 mL Falcon Tubes
- ✕ Bovine Serum Albumin **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A9418**
- ✕ RNaseZap™ RNase Decontamination Solution **Thermo Fisher Scientific Catalog #AM9780**
- ✕ 40 µm Falcon® Cell Strainers Sterile Corning **VWR International (Avantor) Catalog #21008-949**
- ✕ Bovine Serum Albumin (Non-acetylated) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #B6917**
- ✕ 70% EtOH
- ✕ Blunt curved Westcott scissors
- ✕ Vannas Microscissors
- ✕ 1.7 mL Eppendorf Tubes



Safety warnings

- ! Universal biosafety precautions should be observed with the handling of all human post-mortem tissue. Instruments should be cleaned and sterilized appropriately after each dissection.

Before start

Fresh Ames' medium should be prepared and oxygenated at the beginning of each dissection to allow for fresh working solutions.

Prepare Stock Solutions

1 Prepare Oxygenated Ames' Medium

- Add 1.9 g sodium bicarbonate to a 1L container (e.g. glass bottle or flask)
- Empty one container of Ames' Medium powder to 1L bottle
- Fill bottle with 1L of ddH₂O, mix well
- Bubble oxygen gas (95% O₂/5% CO₂) through solution

 Sodium bicarbonate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S6014**

 Ames' Medium **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A1420**

2 DNase stock solution

- DNase I 40,000 Units
- oxygenated Ames' solution 3 mL

Filter solution into new 15 mL falcon tube with 0.22 um filter to sterilize

- Store 115 µL aliquots at -20°C.

 Deoxyribonuclease I from bovine pancreas **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D4527**

3 0.4% non-acetylated BSA stock solution

- NA BSA 25 mg
- oxygenated Ames' solution 6.25 mL

Filter solution into new 15 mL falcon tube with 0.22 um filter to sterilize

- Store 100 µL aliquots at -20°C.

 Bovine Serum Albumin (Non-acetylated) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #B6917**

4 Prepare 10X Ovomucoid solution (LO)

- BSA 150 mg
- Ovomucoid 150 mg
- Oxygenated Ames' 10 mL
- dissolve and adjust pH to ~7.4 with a few µL of NaOH (2N)



Filter solution into new falcon tube with 0.22 um filter to sterilize

- Make 1 mL aliquots & store at -20°C.

⊗ Bovine Serum Albumin **Merck MilliporeSigma (Sigma-Aldrich)** Catalog #A9418

⊗ Trypsin Inhibitor Ovomucoid **Worthington Biochemical Corporation** Catalog #LS003087

5 L-cysteine stock solution

- L-cysteine, anhydrous (sigma C1276): 24 mg
- ddH₂O 1 mL

Filter into new 15 mL falcon tube with 0.22 um filter to sterilize

- Store 115 µL aliquots at -20°C.

⊗ L-Cysteine hydrochloride **Merck MilliporeSigma (Sigma-Aldrich)** Catalog #C1276

Prepare Working Solutions and Materials

1h

6 Prepare Fresh Oxygenated Ames' Medium

- Add 1.9 g sodium bicarbonate to a 1L container (e.g. glass bottle or flask)
- Empty one container of Ames' Medium powder to 1L container
- Fill container with 1L of ddH₂O, shake to mix
- Bubble oxygen gas (95% O₂/5% CO₂) through solution

⊗ Sodium bicarbonate **Merck MilliporeSigma (Sigma-Aldrich)** Catalog #S6014

⊗ Ames' Medium **Merck MilliporeSigma (Sigma-Aldrich)** Catalog #A1420

7 Prepare papain solution

- Oxygenated Ames' 10 mL
- Papain 200 Units (ex: for a 44.8 mgP/mL, 26.6 u/mgP batch, use 168 uL)
- DNase I stock solution 100 uL

Filter into new 15 mL falcon tube with 0.22 μ m filter to sterilize

- L-cysteine aliquot (see above, stock solution) 100 μ L

Activate the papain by adding the L-cysteine and incubating at 37 °C for 15 minutes before using for dissociation.

 Papain **Worthington Biochemical Corporation Catalog #LS003126**

8 **Prepare 1X Ovomuroid solution (LO)**

- 10X LO (aliquot) 1 mL
- Oxygenated Ames' 9 mL

9 **Prepare 0.04% non-acetylated BSA**

- Oxygenated Ames' 900 μ L
- 0.4% non-acetylated BSA stock 100 μ L

Tissue collection & dissection

45m

10 **Set up dissection area**

- Clean all instruments and workspace with RNaseZap then 70% EtOH.
- Place all instruments on a Chux Pad or similar sheet to work on

Note: Keep the eye immersed in Ames media or Optisol within a container on ice until dissection; and preferably, the eye should be processed within 6 hours from death.

 RNaseZap™ RNase Decontamination Solution **Thermo Fisher Scientific Catalog #AM9780**

11 **Start Dissection under microscope: Isolate the anterior segment**

- Make a small incision at the pars plana (~4 mm posterior to the limbus) with a razor blade.
- Cut circumferentially with curved scissors to liberate the anterior portion from the rest of the eye.
- Place the anterior segment in oxygenated Ames' medium.

12 Expose the trabecular meshwork

- Free the lens from the ciliary body by carefully cutting the zonules with curved microscissors all the way around the lens
- Gently lift the lens from the rest of the structures with the scissor blades (holding the scissors parallel with the workbench)
- Take care and try to avoid violating or puncturing the lens capsule
- Gently peel the ciliary body and iris from the corneoscleral button with forceps.
- Trim any large chunks of ciliary muscle left behind from the trabecular meshwork with microscissors.
- This is an ideal time to begin activating the papain solution in a 37°C waterbath

13 Isolate the trabecular meshwork

- Identify the TM and insert one prong of the Jewelers forceps into Schlemm's canal
- Grab the TM with the forceps and gently pull the TM out in strips. Some of the ciliary muscle and possibly descemet membrane will come along, and this is okay for the purposes of this protocol
- The back wall of schlemm canal can also be gently scraped so that cells here can be transferred to the papain solution

Dissociate tissue into single cell solution

1h

14 Digest tissue

- Add TM strips to prewarmed, **activated papain solution** in an eppendorf tube
- Incubate for a total of 30 minutes at 37 °C
- Every 7-10 minutes, invert tube to agitate gently

15 Dissociate tissue into a single cell suspension

- When finished incubating, centrifuge digested tissue and papain solution for 30 seconds at 500 rcf
- Carefully remove supernatant to leave behind tissue
- Add 1 mL **1X LO** to the tissue
- Triturate with a 1000 uL pipette
- Filter through a 40 um strainer into a 50 mL falcon tube
- If any large chunks of tissue are left behind, return to eppendorf and add another 1 mL of **1X LO** and repeat.
- Rinse filter with 500 uL - 1000 uL of **1X LO** and collect in the same 50 mL falcon tube

16 Prepare single cell solution for loading onto 10X

- Place the 50ml falcon tube containing the cell/1X LO suspension in centrifuge
- Centrifuge cells suspension at 500 rcf for 5 minutes (preferably at 4 °C)



- Carefully remove supernatant. Note where pellet should be given the falcon tubes' position in the centrifuge; sometimes the pellet is too small to see
- Resuspend the cells in 50 uL - 200 uL **0.04% non-acetylated BSA**, depending on pellet size
- Transfer to a 1.7 mL eppendorf tube
- Count cells with hemocytometer and dilute with **0.04% non-acetylated BSA** to desired concentration based on 10X protocol (we usually aim to achieve a concentration of 1000 cells per microliter and load the 10X for 6000 cells recovered)
- Cell viability assay can be performed at this stage as well prior to loading