

Apr 27, 2022

# Single-Cell Isolation of Human Articular Cartilage

DOI

[dx.doi.org/10.17504/protocols.io.14egn765qv5d/v1](https://dx.doi.org/10.17504/protocols.io.14egn765qv5d/v1)

Hannah Swahn<sup>1</sup>, Martin Lotz<sup>1</sup>

<sup>1</sup>Scripps Research

Human BioMolecular Atlas Program (HuBMAP) Method Development Community  
Tech. support email: [Jeff.spraggins@vanderbilt.edu](mailto:Jeff.spraggins@vanderbilt.edu)



Hannah Swahn

Scripps Research Institute

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn765qv5d/v1>

**Protocol Citation:** Hannah Swahn, Martin Lotz 2022. Single-Cell Isolation of Human Articular Cartilage. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.14egn765qv5d/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**We use this protocol and it's working**

**Created:** March 24, 2022

**Last Modified:** April 27, 2022

**Protocol Integer ID:** 59875

**Keywords:** single cells from healthy knee cartilage, cell isolation of human articular cartilage, single cell rnaseq, isolating single cell, cell isolation, single cell, human articular cartilage, healthy knee cartilage, cell, isolation

## Abstract

This is a protocol that describes the process of isolating single cells from healthy knee cartilage for single cell RNAseq.



## Materials

Dulbecco's Phosphate Buffered Solution [DPBS] (Gibco, 14190-144)

Dulbecco's Modified Eagle Medium [DMEM] (Corning; 10-013-CV)

Bovine Calf Serum [CS] (VWR; 10158-358)

Antibiotic : Antimycotic [Anti-Anti] (GeminiBio; 400-101)

Penicillin-Streptomycin-Glutamine [PSG] (Corning; 30-009-CI)

Collagenase Type II (Worthington; L5004177)

Feather Disposable Scalpel #21 (Electron Microscopy Sciences; 72042-21)

100 um strainer (Fisher Scientific; 22363549)

40 um strainer (Fisher Scientific; 22363547)

50 mL centrifuge tubes (BioPioneer; CNT-50)

Ethylenediamine tetraacetic acid [EDTA] (Fisher Scientific; BP120-500)

DNase I (RNase-free) (Takara; 2270B)

Bovine Serum Albumin [BSA] (Fisher Scientific; 9048-46-8)

## Troubleshooting



- 1 ~1.5 g of articular cartilage from healthy donor knees (grades 0-1) is resected from the medial condyle of the proximal femur. For details regarding the tissue harvesting procedure please see [dx.doi.org/10.17504/protocols.io.14egn7996v5d/v1](https://dx.doi.org/10.17504/protocols.io.14egn7996v5d/v1).
- 2 Cartilage shavings are washed with Room temperature Dulbecco's Phosphate Buffered Solution (DPBS) supplemented with 10% calf serum (CS), 1% Antibiotic : Antimycotic (Anti-Anti) and 1% Penicillin-Streptomycin-Glutamine (PSG).
- 3 Cartilage shavings are minced with a #21 Feather disposable scalpel, and digested in 20mL Dulbecco's Modified Eagle Medium (DMEM) supplemented with 1% Anti-Anti and 2% collagenase type II with 100 rpm shaking at 37 °C for 02:00:00 . 2h
- 4 Cells are gently passed through a 100 µm filter into a 50 mL centrifuge tube followed by gentle passage through a 40 µm filter into a fresh 50 mL centrifuge tube.
- 5 Filtered cells are spun down at 1200 rpm for 00:05:00 at Room temperature 5m
- 6 The supernatant is carefully removed, and the remaining cell pellet is delicately resuspended in 10mL of Room temperature DPBS supplemented with 5% CS and 5 mM Ethylenediamine tetraacetic acid (EDTA).
- 7 Resuspended single cells are treated with DNase I at a concentration of 100 ug/mL for 00:15:00 at Room temperature . 15m
- 8 Single cells are spun down at 1200 rpm for 00:05:00 at Room temperature 5m
- 9 The supernatant is carefully removed, and the remaining cell pellet is delicately resuspended in 10 mL of DPBS supplemented with 0.04% Bovine Serum Albumin (BSA).
- 10 The Invitrogen Countess II FL automated cell counter is used to quantify single cells and determine cell viability. Live cells are determined by trypan blue staining. If >70% cell viability is confirmed, the single cell suspension is diluted to a concentration of  $1 \times 10^6$  cells/mL for single cell RNA-seq library preparation.