

Aug 23, 2022

Single-Cell Isolation of Human Knee Meniscus

 Forked from [Single-Cell Isolation of Human Articular Cartilage](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.n2bvj6k7blk5/v1>

hswahn¹, Martin Lotz¹

¹Scripps Research

Human BioMolecular Atlas Program (HuBMAP) Method Development Community
Tech. support email: Jeff.spraggins@vanderbilt.edu



Merissa Olmer

Scripps Research

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvj6k7blk5/v1>

Protocol Citation: hswahn , Martin Lotz 2022. Single-Cell Isolation of Human Knee Meniscus. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n2bvj6k7blk5/v1>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 10, 2022

Last Modified: August 23, 2022

Protocol Integer ID: 68487

Keywords: single cell isolation, meniscus, single cells from human knee meniscus, cell isolation of human knee meniscus, isolating single cell, human knee meniscus, cell isolation, cell, scrna

Abstract

This is a protocol that describes the process of isolating single cells from human knee meniscus for scRNA-seq.

Protocol materials

⊗ Dulbeccos Phosphate Buffered Solution [DPBS] **Gibco - Thermo Fisher Scientific Catalog #14190-144**

⊗ Bovine Calf Serum [CS] **VWR International (Avantor) Catalog #10158-358**

⊗ Antibiotic-Antimycotic 9100×0 [Anti-Anti] **Thermo Fisher Scientific Catalog #15240062**

⊗ Penicillin-Streptomycin-Glutamine [PSG] **Corning Catalog #30-009-CI**

⊗ DMEM with L-Glutamine 4.5g/L Glucose and Sodium Pyruvate **Fisher Scientific Catalog #MT-10-013-CV**

⊗ Collagenase Type II **Worthington Biochemical Corporation Catalog #L5004177**

⊗ Ethylenediamine tetraacetic acid [EDTA] **Fisher Scientific Catalog #BP120-500**

⊗ Bovine Serum Albumin [BSA] **Fisher Scientific Catalog #9048-46-8**



1 ~1g of tissue from healthy donor knees (Grades 0-1) is collected from the meniscus. For details regarding the tissue harvesting procedure please see dx.doi.org/10.17504/protocols.io.6qpvr614zvmk/v1

2 Meniscal tissue is washed with Room temperature



Dulbeccos Phosphate Buffered Solution [DPBS] **Gibco - Thermo Fisher Scientific Catalog #14190-144**

supplemented with 10%



Bovine Calf Serum [CS] **VWR International (Avantor) Catalog #10158-358** , 1%



Antibiotic-Antimycotic 9100×0 [Anti-Anti] **Thermo Fisher Scientific Catalog #15240062**

and 1% Penicillin-Streptomycin-Glutamine [PSG] **Corning Catalog #30-009-CI** .

3 The tissue is then finely minced with a #21 Feather disposable scalpel, and digested in 20mL



DMEM with L-Glutamine 4.5g/L Glucose and Sodium Pyruvate **Fisher Scientific Catalog #MT-10-013-CV**

supplemented with 1% Anti-Anti and 2%



Collagenase Type II **Worthington Biochemical Corporation Catalog #L5004177**

using 100 rpm shaking at 37 °C for 00:30:00 .

30m

4 Cells are gently passed through a 100 µm filter into a 50 mL centrifuge tube followed by gentle passage through a 40 µm filter into a fresh 50 mL centrifuge tube.

5 Filtered cells are spun down at 1200 rpm for 00:05:00 at Room temperature

5m

6 Carefully remove and collect the collagenase supernatant; the collagenase is reused in subsequent steps.

7 Cells are then resuspended in DMEM supplemented with 10% CS, 1% Anti-Anti and 1% PSG and stored at 37 °C .

8 Repeat the digestion, spin down and filtration steps for 01:00:00 and then 02:00:00 , for a total of 03:30:00 of digestion.

6h 30m



- 9 Combine all collected cells and spin down at 1200 rpm for  00:05:00 at  Room temperature . 5m
- 10 The supernatant is carefully removed, and the remaining cell pellet is delicately resuspended in 10mL of  Room temperature DPBS supplemented with 5% CS and 5 mM  Ethylenediamine tetraacetic acid [EDTA] **Fisher Scientific Catalog #BP120-500** .
- 11 Single cells are spun down at 1200 rpm for  00:05:00 at  Room temperature . 5m
- 12 The supernatant is carefully removed, and the remaining cell pellet is delicately resuspended in 10 mL of DPBS supplemented with 0.04%  Bovine Serum Albumin [BSA] **Fisher Scientific Catalog #9048-46-8** .
- 13 The Invitrogen Countess II FL automated cell counter is used to quantify single cells and determine cell viability. Live cells are determined by trypan blue staining. If >70% cell viability is confirmed, the single cell suspension is diluted to a concentration of 1×10^6 cells/mL for single cell RNA-seq library preparation.