



Feb 08, 2024

## Single-cell RNA-seq

DOI

[dx.doi.org/10.17504/protocols.io.eq2lyjy5mlx9/v1](https://dx.doi.org/10.17504/protocols.io.eq2lyjy5mlx9/v1)

patricia.garcia<sup>1</sup>

<sup>1</sup>The Francis Crick Institute



patricia.garcia

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2lyjy5mlx9/v1>

**Protocol Citation:** patricia.garcia 2024. Single-cell RNA-seq. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.eq2lyjy5mlx9/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** February 08, 2024



**Last Modified:** February 08, 2024

**Protocol Integer ID:** 94895

**Keywords:** single cell, cell rna, seq harvesting, rna, seq, cell

## Abstract

Harvesting and performing single-cell RNA-seq.

## Troubleshooting

## Harvesting cells for single-cell RNA-sequencing

- 1 Samples were trypsinised or scraped from the culture surface and placed in a 15mL conical tube.
- 2 Tubes were centrifuged at 800g in a refrigerated centrifuge for 5 minutes, and the culture media was decanted.
- 3 The pellet was resuspended in 10mL chilled PBS per tube by pipetting, then centrifuged again using the above parameters before decanting the PBS.
- 4 Pellets were resuspended in 1mL PBS, and 100,000 cells were transferred to a 1.5mL falcon tube.
- 5 These were centrifuged at 1000rpm for 3 minutes at 4°C.
- 6 Cells were resuspended in 20µL chilled DPS.
- 7 180µL chilled 100% methanol was added dropwise to the cells while gently vortexing to prevent cell clumps, before fixing the cells on ice for 15 minutes.

## Single-cell RNA-sequencing protocol

- 8 Between 2400 to 4000 cells were loaded for each sample into a separate channel of a Chromium Chip G for use in the 10X Chromium Controller (cat: PN-1000120).
- 9 The cells were partitioned into nanoliter scale Gel Beads in emulsions (GEMs) and lysed using the 10x Genomics Single Cell 3' Chip V3.1 GEM, Library and Gel Bead Kit (cat: PN-1000121).
- 10 cDNA synthesis and library construction were performed as per the manufacturer's instructions.
- 11 The RNA was reversed transcribed and amplified using 12 cycles of PCR.

- 12 Libraries were prepared from 10µL of the cDNA and 13 cycles of amplification. Each library was prepared using Single Index Kit T Set A (cat: PN-1000213) and sequenced on the HiSeq4000 system (Illumina) using 100 bp paired-end run at a depth of 20-50 million reads. Libraries were generated in independent runs for the different samples.

## Single-cell RNA-sequencing data processing

- 13 Reads were aligned to the human reference genome (Ensembl release 93, GRCh38) using Cell Ranger v3.0.2.
- 14 The analysis was carried out using Seurat v3.0 (REF- Butler, 2018 ,Stuart, 2019 #645) following Seurat's standard workflow.
- 15 Cells expressing fewer than 200 genes were excluded from the subsequent analysis. In addition, we excluded cells with more than 3000 detected genes to remove suspected cell doublets or multiplets.
- 16 Due to the nature of the cells used for this experiment, a 10% cut-off was applied for the percentage of mitochondrial genes expressed to filter out likely apoptotic cells.
- 17 Using default parameters of Seurat, data for each sample were log normalised across cells and the 2000 most highly variable genes identified.
- 18 Using the canonical correlation analysis ('CCA') to identify anchors, we integrated the samples using Seurat v 3 (Butler, 2018; Stuart, 2019), followed by regression of the effect of cell cycle and scaling of the data. Dimensional reduction was performed using 50 PCs.
- 19 Clustree v0.4.4 and Seurat's plots functions were used to visualise the expression of astrocytic and neuronal marker genes across different cluster resolutions.
- 20 A clustering resolution of 0.25 was selected.
- 21 The differentially expressed genes between the clusters of interest were identified using Seurat's FindMarkers() and the default 'Wilcox' test.