



Oct 05, 2022

Single-cell RNA-seq for mDA neurons

DOI

dx.doi.org/10.17504/protocols.io.6qpvr4dxpgmk/v1

gurvir.virdi¹

¹UCL Institute of Neurology



gurvir.virdi

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvr4dxpgmk/v1>

Protocol Citation: gurvirdi 2022. Single-cell RNA-seq for mDA neurons. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.6qpvr4dxpgmk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 05, 2022



Last Modified: May 31, 2024

Protocol Integer ID: 70858

Keywords: ASAPCRN, seq for mda neurons harvesting, mda neurons harvesting, derived mda neuron, mda neuron, cell rna, seq on hipsc, rna, seq, cell

Abstract



Harvesting and performing single-cell RNA-seq on hiPSC-derived mDA neurons.

Troubleshooting



Harvesting cells for RNA-seq

5m

- 1 At the desired age of mDA neurons, they are harvested for single-cell RNA-seq:
 - 1.1 mDA neurons are washed 1x in PBS.
 - 1.2 They are incubated with Accutase (Thermo Fisher Scientific) for  00:05:00 at  37 °C .
 - 1.3 mDA neurons are collected as a single cell suspension and diluted 1:3.

5m

Quality and concentration of cells

- 2 The quality and concentration of each single-cell suspension was measured using Trypan blue and the Eve automatic cell counter.
- 3 Each sample presented a concentration between a 1200-1700 cell/ μ l and viability ranged between 55-68%, samples with a viability above 57% were used for sequencing.

single-cell RNA-seq

- 4 Approximately 10000 cells were loaded for each sample into a separate channel of a Chromium Chip G for use in the 10X Chromium Controller (cat: PN-1000120).
- 5 The cells were partitioned into nanoliter scale Gel Beads in emulsions (GEMs) and lysed using the 10x Genomics Single Cell 3' Chip V3.1 GEM, Library and Gel Bead Kit (cat: PN-1000121).
- 6 cDNA synthesis and library construction were performed as per the manufacturer's instructions.
- 7 The RNA was reversed transcribed and amplified using 12 cycles of PCR.



- 8 Libraries were prepared from 10µl of the cDNA and 13 cycles of amplification. Each library was prepared using Single Index Kit T Set A (cat: PN-1000213) and sequenced on the HiSeq4000 system (Illumina) using 100 bp paired-end run at a depth of 65-100 million reads. Libraries were generated in independent runs for the different samples.