



Nov 21, 2025

Single Cell Transcriptional Perturbome in Pluripotent Stem Cell Models: Protocols 1-4

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvw4q29lmk/v1>

elisa.balmas¹

¹University of Torino



elisa.balmas

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvw4q29lmk/v1>

Protocol Citation: elisa.balmas 2025. Single Cell Transcriptional Perturbome in Pluripotent Stem Cell Models: Protocols 1-4. protocols.io <https://dx.doi.org/10.17504/protocols.io.e6nvw4q29lmk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: October 18, 2025

Last Modified: November 21, 2025

Protocol Integer ID: 230162

Keywords: 10X genomics, iPS2-10X-seq, iPS2-CITE-seq, iPS2-MULTI-seq, seq from 10x genomic, conjunction with the appropriate 10x genomics protocol, appropriate 10x genomics protocol, use of cell multiplexing oligo, cell multiplexing oligo, parallel to the standard gene expression library, cell transcriptome, pluripotent stem cell model, 10x genomic, seq relies on microfluidics partitioning, gene expression kit, standard gene expression library, 10x genomics single cell multiome atac, single cell transcriptional perturbome, gene expression library, unique cell, cell doublet, shrna barcode library library, absence of additional multiplexing barcode, cell number, microfluidics chip, single microfluidics channel, microfluidics partitioning, additional multiplexing barcode, shared unique cell, multiplexed sample, manufacturer for barcoded sample pool, treated cell, seq rely, silico matching of uci, barcoded sample pool, based scrna, gene expression

Funders Acknowledgements:

Giovanni Armenise-Harvard Foundation

Grant ID: Career Development Award 2021

ERC starting grant

Grant ID: 101076026

Disclaimer

Contact elisa.balmas@unito.it and or alessandro.bertero@unito.it for further tips and assistance

Abstract

iPS2-10X-seq relies on microfluidics partitioning-based scRNA-seq from 10X Genomics. Despite the higher costs, the ease and reproducibility of this commercial method have led to its widespread adoption. The procedure integrates the manufacturer's protocol to generate an additional NGS library containing UCI-BCs. This is sequenced in parallel to the standard gene expression library to enable in silico matching of UCI-BCs and single-cell transcriptomes through their shared unique cell barcodes. Of note, iPS2-10X-seq provides built-in support for the cost-saving strategy of reasonably overloading the microfluidics chip, since cell doublets are readily identified by the co-expression of two different UCI-BCs and filtered out as part of the standard catchR_10Xcatch pipeline. Thus, it is possible to load the cell number recommended by the manufacturer for barcoded sample pools even in the absence of additional multiplexing barcodes. Other cost-saving strategies include the use of cell multiplexing oligos (CMOs) or barcoded antibodies against housekeeping cell surface antigens to multiplex tet-untreated and tet-treated cells in a single microfluidics channel, enabling robust control versus knockdown experimental design without having to perform two separate reactions.

This Protocol describes the following procedures:

1. Prepare the gene expression library
2. Prepare the shRNA barcode library
3. Library pooling and NGS

This protocol was optimized using the *Chromium Next GEM Single Cell 3' Reagent Kits v3.1* (Dual Index) from 10X Genomics and has been successfully validated with the *next-generation GEM-X Single Cell 3' Reagent Kits v4*, without requiring changes to primer sequences. It must be used in conjunction with the appropriate 10X Genomics protocols: **CG00315 Rev E** (single sample), **CG000388 Rev B** (multiplexed samples), or **CG000419 Rev D** (high-throughput) for v3.1; and **CG000731 Rev A** for v4. Notably, v4 no longer supports cell multiplexing oligo (CMO)-based on sample multiplexing but remains fully compatible with antibody-derived tag (ADT)-based strategies.

A first variation, iPS2-CITE-seq, enables surface protein profiling and ADT-based multiplexing using barcoded TotalSeq-B antibodies from BioLegend. When using antibodies, refer to protocol **CG000149 Rev D** (antibody labeling) in combination with **CG000317 Rev E** for v3.1 or **CG000731 Rev A** for v4.

A second variation, iPS2-multi-seq, adapts the protocol for use with the 10X Genomics Single Cell Multiome ATAC & Gene Expression kit. It should be paired with **CG000365 Rev C** (nuclei isolation) and **CG000338 Rev F** (library construction).

Image Attribution

Figure SP3.1



Materials

- NEBNext High-Fidelity 2X PCR Master Mix
- iPS2-10X-seq_true5 (10 μ M)
- iPS2-10X-seq_inner (10 μ M)
- Preamplified 10X Genomics cDNA (50 ng)
- Nuclease-free water
- Dual Index Kit TT Set A primers
- Qubit dsDNA High Sensitivity kit

Protocol 1: Generation of iPS2-seq plasmids and hiPSCs

- 1 The iPS2-seq pipeline starts with a two-step cloning procedure for a barcoded pool of shRNAs into an AAVS1 locus targeting plasmid to generate all-in-one tet inducible cassette.
This begins with designing single stranded DNA (ssDNA) oligonucleotides each carrying an shRNA with a polythymidine (pT) Pol III terminator sequence, a uniquely matched barcode (BC), and a random sequence serving as a unique clonal identifier (UCI). The shRNA and UCI-BC sequences are separated by a short multicloning site (Sall, Swal, AscI) that is needed at a later stage, and are flanked by short regions of homology to the destination vector to enable assembly.
The oligos are synthesized, pooled, and converted to double stranded DNA (dsDNA) through a single round of isothermal primer extension from the 3' end using a mesophilic DNA polymerase with strong strand displacement activity.
This is essential to overcome the secondary structure of shRNAs in the oligos and minimize representation biases due to differential hairpin strength.
dsDNA fragments are then cloned through Gibson assembly into an AAVS1 targeting plasmid downstream of an inducible H1 Pol III promoter and upstream of a pA signal. The vector contains a puromycin (puro) resistance gene trap and is based on our previously reported plasmid, but was modified to remove repetitive sequences that led to recombination errors in the next cloning step.
The plasmid also includes a novel thymidine kinase (TK) cassette outside of the AAVS1 homology arms. This allows subsequent negative selection of cells carrying random, off-target integrations of the whole plasmid as they become sensitive to fialuridine (FIAU), a thymidine analog preferable to ganciclovir.
The resulting intermediate plasmids are then modified in pool to insert a CAG-OPtettR cDNA cassette between the shRNA and UCI-BC sequence.
This second cloning step reconstitutes the tet-ON system while marking the 3' UTR of the OPtettR with the UCI-BC sequence matched to a given shRNA.
This is performed via restriction digestion and ligation, as the CAG promoter is too GC-rich for efficient and reliable PCR amplification.
An optional strategy built into the system allows depletion of the intermediate plasmid from the final population through restriction digestion with Swal.
The final plasmids are utilized to engineer pools of hPSCs enriched for clones expressing single inducible shRNAs.
Plasmids are transfected alongside vectors encoding for obligate heterodimer zinc finger nucleases (ZFNs) to induce double-strand breaks specifically at the AAVS1 locus and facilitate homology-directed repair (HDR) of the inducible barcoded shRNA cassette. A second AAVS1 targeting plasmid is also co-delivered: this is a filler vector that carries no relevant cargo but has a neomycin (neo) resistance gene trap that allows enrichment

of cells that integrated an shRNA on one AAVS1 allele and the filler sequence on the other allele; these biallelic edited cells are co-selected with puro and neo. Additional negative selection with FIAU kills clones that carry randomly integrated shRNAs.

Collectively, this genome editing strategy is designed to reduce the frequency of double shRNA integrations and maximise the fraction of cells expressing a single shRNA from the AAVS1 locus.

Of note, the first cloning step can result in a low rate of swap between the shRNA and its expected BC based on oligonucleotide design and synthesis.

This issue likely arises from DNA polymerase template switching during dsDNA synthesis and/or Gibson assembly.

The problem can be efficiently countered by next generation sequencing (NGS) analysis of Illumina-compatible libraries spanning the shRNA and UCI-BC region from the intermediate plasmid pool; this enables re-assignment of swapped shRNAs to the observed barcode, leveraging the fact that each plasmid also carries a UCI (also essential for tracking hPSC clones downstream).

2 Plan the screening

Before initiating a screening experiment, it is essential to define three key elements: screening objectives, available budget, and overall feasibility.

The flow chart below summarizes key parameters for optimizing molecular cloning, genome editing, and sequencing efficiency, followed by a worked example: a screen targeting 18 genes with 5 shRNAs per gene plus 10% control constructs.

	Cloning step 1	Cloning step 2	Genome editing	scRNA-seq
Goal:	100 bacterial clones per shRNA	10 bacterial clones per shRNA	5 iPSC clones per shRNA	100 cells per clone
Estimated efficiency:	2,000 bacterial clones per transformation	200 bacterial clones per transformation	50 clones for 1,000,000 iPSCs nucleofected	20,000 cells recovered per reaction 65% of cells with 1 shRNA assigned by catcherR
Example: 100 shRNAs	5 transformations	5 transformations	10 M iPSCs for nucleofection	80,000 target cell recovery* (4 reactions)

Figure SP1.1. Screening design example

Flow chart illustrating how to estimate the number of cells required for sufficient iPS2-seq coverage, accounting for dropout rates from cloning to transfection and single-cell sequencing. The top section defines the screening goal and provides an example estimation based on our protocol; the lower section shows a worked example with 100 shRNAs. scRNA-seq calculations are based the maximum loading capacity per reaction for either the Chromium Next GEM Single Cell 3' HT GEM Kit v3.1 or the Chromium GEM-X Universal 3' v4 Kit.

2.1 Effective screening design should begin with:

Target gene selection — Define the genes to perturb and decide how many shRNAs per gene will be used. We recommend a minimum of 3 shRNAs per gene, and ideally 5, to buffer against off-target effects or poorly performing constructs. Include at least 10% non-targeting control shRNAs.

- 2.2 **Cell editing scale** — Estimate the number of cells to genome-edit in order to achieve robust clonal representation for downstream iPS2-10X-seq analysis. An example calculation is provided in Figure SP1.1, accounting for both bacterial clone yield and the number of hiPSC clones required for sequencing.
- 2.3 **Sequencing depth** — Determine how many cells need to be profiled to capture biological variability and achieve statistical power. Depending on experimental goals and budget, controls may include non-targeting shRNAs (e.g., B2M, SCR), tet-untreated cells, or both. We strongly recommend including tet-untreated cells to control for clonal variability.

3 **Design shRNA oligonucleotides**

ssDNA oligonucleotides for iPS2-seq can be designed using `catcherR_design` as described in Supplemental Protocol 4 - Oligonucleotides design. This function automatically performs the following steps, reducing the chance of manual errors when designing multiple oligos:

1. Identify validated shRNAs or design appropriate sequences. In our lab we use the Broad Institute TRC library, available at the GPP Web Portal. When using this resource take the ~58 bp forward oligos and remove the sequences required for cloning in pLKO vectors (typically the first 4 nucleotides and the last 6 nucleotides).
2. If an shRNA does not start with A or a G, add an initial G (required for efficient transcription by the H1 Pol III promoter).
3. Construct the cloning oligos by adding the following sequences to each shRNA, in order:
 - At the 5', the Gibson homology region 5'-AGTTCCCTATCAGTGATAGAGATCCC-3'
 - At the 3', the polythymidine (pT) Pol III terminator sequence 5'-TTTTTTTTT-3'
 - At the 3', Sall, SwaI, and AclI restriction sites 5'-GTCGACATTTAAATGGCGCGCC-3'
 - At the 3', a random sequence (UCI; i.e., 5'-NNNNNN-3')
 - At the 3', an shRNA-specific barcode (BC; i.e., 5'-CAGTTCCA-3')
 - At the 3', the Gibson homology region 5'-GTAGCTGCGTGATCAGC-3'

The resulting oligos should look like Figure SP1.2

Note: These sequences are passed as arguments to `catcherR_design`. The recommended minimum length for UCIs and BCs is 6 and 8 bp, respectively. Tools such as **BARCOSEL** can be used to generate and select BCs with optimal nucleotide distribution.




Figure SP1.2. Example of iPS2-seq ssDNA oligo

Annotated SMAD2.1 shRNA sequence extracted from TRC ID:TRCN0000010477 and modified by catcheR_design.

1. Modify or replace any cloning oligo that contains BglII (5'-AGATCT-3') or MluI (5'-ACGCGT-3') restriction sites that would interfere with subsequent cloning steps
2. Synthesize oligos as high quality ssDNA. In our lab we use the oPools Oligo Pools service from Integrated DNA Technologies

Pooled cloning step 1: obtain the intermediate plasmid pool

- 4 This section outlines the pooled conversion of barcoded shRNA-containing ssDNA oligos into dsDNA by second-strand synthesis (SSS), and the insertion of such fragments in pAAV-Puro_siKD2.0.
 1. Prepare the oligo pool:
- 4.1 (a) Spin down the lyophilized ssDNA and resuspend it at 50 μ M in nuclease-free water
- 4.2 (b) Dilute the 50 μ M stock 1:10 in nuclease-free water to make a 5 μ M working stock
- 5 2. Convert ssDNA to dsDNA using Bst3.0 DNA polymerase
 - 5.1 (a) Denature and snap cool the ssDNA oligo pool
 - i. Mix the following for 8 reactions in a single tube on ice:

	A	B	C	D
Reagent		one reaction (μ L)	eight reactions (μ L)	Final conc.
ssDNA oligo pool (5 μ M)		5	40	0.5nM
siKDbc_REV Primer (10 μ M)		2.5	20	0.5 μ M
Deoxynucleotide (dNTP) Solution Mix (10mM)		2	16	400 μ M
Nuclease-free water		15.5	124	

A	B	C	D
Total volume	25	200	

Table SP1.1. SSS reaction mix 1

- ii. Distribute 25 μL of the reaction mix into 8 PCR tubes
- iii. Place PCR tubes in a thermocycler and run the following cycling conditions:

A	B	C	D
Steps	Temperature	Time	Cycles
Denaturation	98°C	5 min	1
Hold (rapid ramp down)	4°C	Hold	

Table SP1.2. Cycling conditions for ssDNA denaturation and snap cooling

5.2 (b) Perform SSS

- i. Mix the following for 8 reactions in a single tube on ice:

A	B	C	D
Reagent	One reaction (μL)	Eight reactions (μL)	Final conc.
MgSO ₄ (100mM)	1	8	2mM
Isothermal Amplification Buffer II (10X)	5	40	1X
Bst 3.0 DNA Polymerase (8 U/ μL)	2	16	0.32 U/ μL
Nuclease-free water	17	136	
Reaction from previous step	25	200	
Total volume	25	200	

Table SP1.3. SSS reaction mix 2

- ii. Add 25 μL of this reaction mix to each of the 8 PCR tubes from step 2(a)iii



iii. Place PCR tubes in a thermocycler and run the following cycling conditions:

A	B	C	D
Steps	Temperature	Time	Cycles
Annealing	55°C	5 min	1
Extension	72°C	15 min	1
Inactivation	80°C	5 min	1
Slow ramp down	75°C	-0.1deg/s	1
Hold	75°C	4 min	1
Slow ramp down	70°C	-0.1deg/s	1
Hold	70°C	4 min	1
Slow ramp down	10°C	-0.1deg/s	1
Hold	10°C	Hold	

Table SP1.4. Cycling conditions for SSS

5.3 (c) Remove residual ssDNA

- Add 4 µL of exonuclease I (20 U/µL) to each PCR tube and pipette mix
- Place PCR tubes in a thermocycler and run the following cycling conditions:

A	B	C	D
Steps	Temperature	Time	Cycles
Digestion	37°C	30 min	1
Inactivation	80°C	20 min	1
Slow ramp down	75°C	-0.1deg/s	1
Hold	75°C	4 min	1
Slow ramp down	70°C	-0.1deg/s	1
Hold	70°C	4 min	1
Slow ramp down	10°C	-0.1deg/s	1

A	B	C	D
Hold	10°C	Hold	

Table SP1.5. Cycling conditions for ssDNA removal**5.4** (d) Clean up and concentrate the dsDNA

i. Pool all 8 reactions into a single tube

ii. Column purify using QIAquick PCR Purification Kit, following its manual

iii. Elute in 30 µL elution buffer

Note: the dsDNA shRNA pool can be stored at -20 °C for at least one week

6 Prepare the plasmid backbone**6.1** (a) Prepare a digestion mix and an undigested control mix following the table below:

A	B	C	D
Reagent	Digestion reaction (µL)	Undigested control (µL)	Final conc.
pAAV-Puro_siKD2.0 plasmid	Variable	Variable	5µg
FastDigest Green Buffer (10X)	15	1.5	1X
FastDigest MluI	5	0	
FastDigest BglII	5	0	
FastAP (1 U/µL)	5	0.5	0.03 U/µL
Nuclease-free water	Top up to 150	Top up to 15	
Total volume	150	15	

Table SP1.6. pAAV-Puro_siKD2.0 restriction digestion mix**6.2** (b) Incubate at 37 °C overnight (~16 h), followed by 10 min at 80 °C for heat inactivation**7** Perform electrophoretic separation of the insert (step 2(d)iii) and backbone (step 3b) using 2% and 0.8% (w/v) agarose gels in TBE, respectively. Cast gels including SYBR

Safe DNA gel stain to visualize DNA, preferably using a blue light transilluminator (avoid UV light if possible). Figure SP1.3 shows exemplary electrophoreses

- 7.1 (a) Use 50 bp and 1 kbp DNA ladders for the insert and backbone, respectively
- 7.2 (b) Run 10% of each sample in a separate well and use this to monitor the electrophoretic separation, so as to avoid unnecessary exposure to high-energy light of the gel extracted DNA (particularly important if using UV light by necessity)
- 7.3 (c) Use FastDigest Green Buffer (10X) as loading dye for the insert: commonly used DNA loading dyes containing bromophenol blue and/or xylene cyanol FF can interfere with visualization of the small shRNA pool fragment
- 7.4 (d) If a smear is present in the insert, continue electrophoresis until the 135 bp band is clearly separated
- 7.5 (e) Run 90% of the digested empty plasmid (~140 μ L) in a single well as this will facilitate gel extraction: merge multiple wells using autoclave tape if needed
- 7.6 (f) Load the plasmid undigested control next to 10% of the digestion reaction and continue the electrophoretic separation until the 8715 bp backbone band is completely separated from the undigested control (this can take several hours: use a low voltage and/or change buffer to prevent gel overheating). Incomplete separation can result in contamination of undigested or partially digested plasmid, which can dramatically increase transformation background




Figure SP1.3. Electrophoretic separation and gel extraction during cloning step 1 (A)

Separation of shRNA pool insert (SSS); a PCR-amplified fragment indicates the expected molecular weight. (B) Separation of digested plasmid backbone.

- 8 Use a clean scalpel to excise the DNA bands corresponding to the shRNA pool insert (135 bp) and plasmid backbone (~8.7 kbp) and place them in a microcentrifuge tube
- 9 Proceed with gel extraction using the QIAEX II Gel Extraction Kit following the manual eluting in 20 μ L EB; quantify the DNA concentration using NanoDrop

10 Setup the desired assembly alongside control reactions following the table below:

	A	B	C	D	E
	Reagent	Reaction (μL)	BB only CTRL (μL)	No Assembly CTRL (μL)	Final conc.
	Gel extracted shRNA fragment (135bp)	Variable	0	0	14.34ng (172.4 fmol)
	Gel extracted iPS2-seq backbone (8715bp)	Variable	Variable	Variable	50ng (9.32 fmol)
	NEBuilder HiFi DNA Assembly Master Mix (2X)	10	10	0	1X
	Nuclease-free water	To 20	To 20	To 20	
	Total volume	20	20	20	

Table SP1.7. Pooled cloning step 1 assembly mix

- 10.1 Incubate the reactions at 50 °C for 15 min, then transfer to ice
Note: avoid longer incubation times as they can increase background
- 11 Transform 2 μL of assembly or control reactions in 50 μL NEB 5-alpha competent *E. coli* (high-efficiency), following the manufacturer's instructions
- 12 Plate transformed bacteria
 - 12.1 Pre-warm LB-agar plates containing 100 μg/mL ampicillin to 37 °C
 - 12.2 To facilitate estimation of colony number and cloning efficiency, plate 10% of transformed bacteria on one dish, and the remaining 90% on a second dish
 - 12.3 Incubate the plates at 37 °C overnight

- 13 Estimate the number of bacterial clones by counting the colonies in the plate containing 10% of the transformation of the desired assembly (>200 are to be expected)

Notes:

The number of colonies in the matching backbone-only control and no assembly control plates should be at least 10-fold and 100-fold lower; it is advisable to repeat the procedure should a substantially higher background be observed
Bacterial plates can be stored at 4 °C for up to one week to allow for the subsequent optional quality controls

- 14 **OPTIONAL.** Confirm cloning efficiency by bacterial colony PCR



- 14.1 (a) For each colony to be screened, aliquot in PCR plates a 15 µL PCR mix based on nuclease-free water and containing: iPS2-seq_step1QC_F Primer (200 nM), iPS2seq_seq_QC_R Primer (200 nM), dNTPs (200 µM), MgCl₂ (2 mM), Colorless GoTaq Flexi buffer (1X), and GoTaq DNA polymerase (0.025 U/µL)
- 14.2 (b) Pick an individual bacterial colony from the 10% plate, mix it by pipetting in 10 µL nuclease-free water aliquoted in a PCR plate, and transfer 5 µL of the suspension to the PCR mix for a final volume of 20 µL. The remaining 5 µL can be transferred to a PCR plate containing LB with 100 µg/mL ampicillin and stored at 4 °C for later use. Repeat for all colonies
- 14.3 (c) Place the PCR plate in a thermocycler and run following cycling conditions: (1) 95 °C for 5 min; (2) 95 °C for 30 sec; (3) 60 °C for 30 sec; (4) 72 °C for 1 min; (5) go to step 2 for 34 times; (6) 72 °C for 2 min; (7) hold at 10 °C
- 14.4 (d) Aliquot 5 µL of each colony PCR in a new plate, add 1 µL of 6X DNA loading dye, and perform an electrophoretic run on a 1.5% (w/v) agarose gel in TBE
Note: clones with the correct insert deliver a PCR product of 543 bp.

- 15 **OPTIONAL.** Confirm cloning accuracy by Sanger sequencing

- 15.1 Aliquot 5 µL of each colony PCR in a new plate and add 1.5 µL of a mastermix containing 0.5 µL of exonuclease 1 (20 U/µL) and 1 µL of shrimp alkaline phosphatase (1 U/µL)
- 15.2 Incubate the reaction at 37 °C for 30 min, followed by heat inactivation for 15 min at 80 °C
- 15.3 Add 5 µL of iPS2-seq_QC_R primer (10 µM) to each well, and run Sanger sequencing.
Note: alignment of each electropherogram to the common sequence 5'-TTTTTTTGTGCGACATTTAAATGGCGCGCCNNNNNNNNNNNNNGTAGCTCGCTGATCAGC

-3'

simplifies the extraction of the BC (last 8 Ns); the electropherogram can then be aligned to the expected oligo based on the BC, to verify the correct sequence of the associated shRNA

16 Collect all bacterial colonies from the 90% plate by washing and scraping the surface of LB-agar plates with 5 mL of LB broth with 100 µg/mL ampicillin

17 Inoculate the bacterial suspension in a clean Erlenmeyer flask with 45 mL LB broth containing 100 µg/mL ampicillin, and grow the bacterial cultures at 37 °C at 225 rpm for 16 h

Note: inoculate bacteria in a sterile environment and use a sterile flask to minimize the chances of contamination. Do not grow bacteria for more than 16 h

18 Isolate the intermediate plasmid pool from bacterial cultures using the QIAGEN plasmid midiprep kit following the manual, and quantify the DNA concentration using NanoDrop

Pooled cloning step 2: obtain the final plasmid pool

19 This section outlines the generation of the final plasmid pool by re-inserting the CAG-OPTtetR into the intermediate plasmid pool.

20 Prepare the following restriction digestion reaction of the empty pAAV-Puro_siKD2.0 plasmid with Sall and MluI to extract the CAG-OPTtetR insert:

A	B	C	D
Reagent	Digestion reaction (µL)	Undigested control (µL)	Final conc.
pAAV-Puro_siKD2.0 plasmid	Variable	Variable	5µg
FastDigest Green Buffer (10X)	15	1.5	1X
FastDigest MluI	5	0	
FastDigest Sall	5	0	
Nuclease-free water	Top up to 100	Top up to 10	
Total volume	100	10	

Table SP1.8. pAAV-Puro_siKD2.0 second restriction digestion mix

- 21 Prepare the following restriction digestion reaction of the intermediate plasmid pool from step 15 of the previous section with Sall and Ascl to obtain the backbone:

	A	B	C	D
	Reagent	Digestion reaction (μL)	Undigested control (μL)	Final conc.
	Intermediate iPS2-seq plasmid	Variable	Variable	5μg
	FastDigest Green Buffer (10X)	15	1.5	1X
	FastDigest Sall	5	0	
	FastDigest Sgsl (Ascl)	5	0	
	FastAP (1 U/μL)	5	0.5	0.03 U/μL
	Nuclease-free water	To 150	To 15	
	Total volume	150	15	

Table SP1.9. Intermediate plasmid pool restriction digestion mix

- 22 Incubate at 37 °C overnight (~16 h), followed by 10 min at 80 °C for heat inactivation.
- 23 Perform electrophoretic separation of all fragments using 0.8% (w/v) agarose gel in TBE supplemented with SYBR Safe. Figure SP1.4 shows an exemplary electrophoresis.

Note: follow the same general recommendations described for step 7 of the previous section, paying particular attention to separating the backbone from the undigested plasmid.




Figure SP1.4. Electrophoretic separation and gel extraction during cloning step 2. Separation of the CAG-OPTtetR insert and two digested intermediate plasmid pool backbones.

- 24
- Use a clean scalpel to excise the ~8.8 kbp backbone from the digested intermediate plasmid pool and the ~2.4 kbp CAG-OPTtetR insert from pAAV-Puro_siKD2.0.
- 25
- Proceed with gel extraction using the QIAEX II Gel Extraction Kit following the manual and eluting in 20 μ L EB; quantify the DNA concentrations using NanoDrop.
- 26
- Set up desired ligation alongside a control reaction following the table below:

	A	B	C	D
	Reagent	Reaction (μ L)	BB only CTRL (μ L)	Final conc.
	Gel extracted CAG-OPTtetR fragment (2392bp)	Variable	0	81.59ng (55.38 fmol)
	Gel extracted intermediate iPS2-seq pool (8795bp)	Variable	1	100ng (18.46 fmol)
	Rapid Ligation Buffer (5X)	4	4	1X
	T4 DNA Ligase (5 U/ μ L)	1	1	0.25 U/ μ L
	Nuclease-free water	To 20	To 20	
	Total volume	20	20	

Table SP1.10. Pooled cloning step 2 ligation mix

- 26.1 Incubate the reaction at RT for 1 h.
- 27 Transform 5 μ L of ligation and control reaction in 50 μ L NEB 5-alpha competent *E. coli* (high-efficiency), following the manufacturer's instructions.
- 28 Plate and grow overnight transformed bacteria according to the general recommendations described for step 9 of the previous section, particularly regarding splitting the bacteria on two plates (10% and 90%).
- 29 Estimate the number of bacterial clones by counting the colonies in the plate containing 10% of the transformation of the desired assembly (3e20 are to be expected).
- Notes:**
(a) the number of colonies in the backbone-only control should be at least 5-fold lower.
(b) bacterial plates can be stored at 4 °C for up to one week.
- 30 **OPTIONAL.** Confirm cloning efficiency by bacterial colony PCR
- 30.1 Perform colony PCRs and electrophoresis as described for step 14 of the previous section, except using primers iPS2-seqstep2QCF and iPS2-seqQCR
- Notes:**
i. Clones with the correct insert deliver a PCR product of 312 bp; intermediate plasmids are not amplified, while the parental pAAV-Puro_siKD2.0 plasmid gives a 295 bp band (Figure S1C)
ii. Sanger sequencing can also be performed as described for step 15 of the previous section, to confirm the distribution of shRNA BCs; however, at this final stage, the shRNA and barcode cannot be co-amplified
- 31 Collect all bacterial colonies and inoculate them as described for steps 16-17 of the previous section. Be extra careful to avoid contamination at this stage.
- 32 Isolate the final plasmid pool from bacterial cultures using the QIAGEN plasmid midiprep kit with endotoxin-free buffers, and quantify the DNA concentration using Qubit.
- 33 Perform a diagnostic restriction digestion using EcoRI, MluI, and Swal: the molecular weight of the expected bands is 5669 bp and 5518 bp
- Note:** additional bands at 2002 bp or 4343 bp would indicate contamination with intermediate plasmids or pAAV-Puro_siKD2.0, respectively. These contaminants can be removed by digestion with Swal or BglII and MluI, respectively, followed by re-transformation

OPTIONAL: prepare NGS libraries for quality control of plasmid pools

- 34 This section outlines the optional generation of Illumina-compatible next-generation sequencing (NGS) libraries to determine the composition of intermediate and final plasmid pools (Figure SP1.5A, Figure SP1.6A). Data analysis relies on *catcheRstep1QC* and *catcheRstep2QC*, respectively, as described in Supplemental Protocol 4 - Pooled cloning step 1 QC and Pooled cloning step 2 and hiPSC genome editing QC.

QC n1: Intermediate plasmid pool QC

- 35 **PCR #1:** amplify the shRNA and UCI-BC region from the intermediate plasmid pool while adding the TruSeq read 2 adapter, a diversity index (DI, a random 12 bp sequence that ensures adequate clustering during NGS and can be leveraged to partially demultiplex reads to eliminate duplicates from the subsequent PCR), and partial TruSeq read 1 adapter

- 35.1 (a) Prepare a PCR reaction according to the following calculations:

	A	B	C
	Reagent	Volume (μL)	Final conc.
	Intermediate iPS2-seq plasmid pool	Variable	150 ng
	ThermoPol Reaction Buffer (10X)	2.5	1X
	Deoxynucleotide (dNTP) Solution Mix (10mM)	0.5	200μM
	siKDbc_FW Primer (10μM)	1.5	0.5μM
	siKDbc_REV Primer (10μM)	1.5	0.5μM
	MgSO4 (100mM)	1	4mM
	Deep Vent DNA Polymerase (2 U/μL)	0.25	0.02 U/μL
	Nuclease-free water	To 25	
	Total volume	25	

Table SP1.11. Intermediate plasmid pool QC, PCR #1 mix

- 35.2 (b) Place the PCR tube in a thermocycler and run the following cycling conditions:

A	B	C	D
Steps	Temperature	Time	Cycles
Initial Denaturation	95°C	2 min	1
Denaturation	95°C	30 sec	3
Annealing	54-59°C (increase by 2.5°C each cycle)	15 sec	
Extension	72°C	30 sec	
Denaturation	95°C	30 sec	9
Annealing	63°C	15 sec	
Extension	72°C	30 sec	
Final extension	72°C	5 min	1
Hold (rapid ramp down)	4°C	Hold	

Table SP1.12. Cycling conditions for PCR #1, intermediate plasmid pool QC

Note: a "touch-up" strategy minimizes the annealing of shRNA hairpins during the second part of the reaction: once a sufficient number of molecules have integrated the TruSeq adapters, they anneal to the primers at higher temperatures

- 35.3 (c) Add 25 µL of nuclease-free water to the reaction and perform a double-sided size selection using 0.6X & 1.2X volumes of SPRIselect beads, following the manufacturer's instructions and eluting in 20 µL of elution buffer
- 35.4 (d) Run 2 µL of the undiluted PCR product on a TapeStation High Sensitivity D1000 Screen Tape assay, following the manual. Expect a single peak of ~203 bp. An exemplary result is reported in (Figure SP1.5B)
- 35.5 (e) Quantify the DNA concentration in the range 100-400 bp, and dilute it to 1 ng/µL
- 36 **PCR #2:** prepare dual-indexed Illumina TruSeq sequencing libraries
- 36.1 (a) Prepare a PCR reaction according to the following calculations:

A	B	C
Reagent	Volume (μL)	Final conc.
PCR#1 product (1 ng/ul)	2	2 ng
ThermoPol Reaction Buffer (10X)	2	1X
Deoxynucleotide (dNTP) Solution Mix (10mM)	0.4	200μM
Dual Index primers Kit TT Set A	4,00	
MgSO4 (100mM)	0.8	4mM
Deep Vent DNA Polymerase (2 U/μL)	0.2	0.02 U/μL
Nuclease-free water	10.6	
Total volume	20	

Table SP1.13. Intermediate plasmid pool QC, PCR #2 mix

36.2 (b) Place the PCR tube in a thermocycler and run the following cycling conditions:

A	B	C	D
Steps	Temperature	Time	Cycles
Initial Denaturation	95°C	2 min	1
Denaturation	95°C	30 sec	3
Annealing	45-47°C (increase by 2.5°C each cycle)	15 sec	
Extension	72°C	45 sec	
Denaturation	95°C	30 sec	9
Annealing	63°C	15 sec	
Extension	72°C	30 sec	
Final extension	72°C	5 min	1
Hold (rapid ramp down)	4°C	Hold	

Table SP1.14. Cycling conditions for PCR #2, intermediate plasmid pool QC

- 36.3
- Add 30 μ L of nuclease-free water to the reaction and perform a double-sided size selection using 0.65X & 0.85X volumes of SPRIselect beads, following the manufacturer's instructions and eluting in 20 μ L of elution buffer
- 36.4
- Run 2 μ L of the undiluted PCR product on a TapeStation High Sensitivity D1000 Screen Tape assay, following the manual. Expect a single peak of ~287 bp

37




Figure SP1.5. NGS library generation from intermediate plasmid pools
(A) Library structure and primer sequences used for PCR #1 and PCR #2; read 1 is also indicated.
(B) Exemplary TapeStation High Sensitivity D1000 Screen Tape assays of size-selected DNA from PCRs #1 and #2

QC n2: Final plasmid pool QC

- 38
- PCR #1:** amplify UCI-BCs from the final plasmid pool, adding TruSeq adapters and DIs

39

	A	B	C
	Reagent	Volume (μ L)	Final conc.
	Final iPS2-seq plasmid pool	Variable	150 ng
	Q5 Reaction Buffer (5X)	5	1X
	Deoxynucleotide (dNTP) Solution Mix (10mM)	0.5	200 μ M
	iPS2-seq_step2NGS_Primer (10 μ M)	1.25	0.5 μ M
	iPS2-seq_NGS_R Primer (10 μ M)	1.25	0.5 μ M

A	B	C
Q5 Hot Start High-Fidelity DNA Polymerase (2 U/μL)	0.25	0.02 U/μL
Nuclease-free water	To 25	
Total volume	25	

Table SP1.15. Final plasmid pool QC, PCR #1 mix

40

A	B	C	D
Steps	Temperature	Time	Cycles
Initial Denaturation	98°C	30 sec	1
Denaturation	98°C	10 sec	9 cycles
Annealing	67°C	30 sec	
Extension	72°C	20 sec	
Final extension	72°C	2 min	1
Hold (rapid ramp down)	4°C	Hold	

Table SP1.16. Cycling conditions for PCR #1

41 Purify, quantify, and dilute the PCR product following steps 35.2, 35.3, and 35.4 of the previous section, except for using 0.6X & 1.5X volumes of SPRIselect beads for double-sided size selection. Expect a single peak of ~144 bp (Figure SP1.6B)

42 **PCR #2:** prepare dual indexed Illumina TruSeq sequencing libraries

43

A	B	C
Reagent	Volume (μL)	Final conc.
PCR#1 product (1 ng/ul)	1	1 ng
Dual Index primers Kit TT Set A	4,00	

	A	B	C
	NEBNext high fidelity 2X PCR master mix	10	1X
	Nuclease-free water	5,00	
	Total volume	20	

Table SP1.17. Final plasmid pool QC, PCR #2 mix

44

	A	B	C	D
	Steps	Temperature	Time	Cycles
	Initial Denaturation	98°C	30 sec	1
	Denaturation	98°C	10 sec	8 cycles
	Annealing	58°C	20 sec	
	Extension	72°C	20 sec	
	Final extension	72°C	30 sec	1
	Hold (rapid ramp down)	4°C	Hold	

Table SP1.18. Cycling conditions for PCR #2

45 Purify and quantify the PCR product following steps 36.3 and 36.4 of the previous section, except for using 0.65X & 1.5X volumes of SPRIselect beads for double-sided size selection. Expect a single peak of ~228 bp (Figure SP1.6B)

46 Perform NGS using NextSeq 1000/2000 and 50 cycles reagents, with the following settings: read 1 - 68 cycles; index 1 - 10 cycles; index 2 - 10 cycles (>10,000 reads/shRNA)

47




Figure SP1.6. NGS library generation from final plasmid pools

(A) Library structure and primer sequences used for PCR #1 and PCR #2; read 1 is also indicated.

(B) Exemplary TapeStation High Sensitivity D1000 Screen Tape assays of size-selected DNA from PCRs #1 and #2.

Genome editing: obtain hPSC clonal pools

48 This section outlines the generation of hPSCs genome edited at the *AAVS1* locus with the final iPS2-seq plasmid pool, and the optional genotyping of individual clones as quality control (Figure SP1.1).

49 Preparatory steps:

49.1 Determine the number of nucleofections required to obtain the total target number of hPSC clones, based on the number of shRNAs, the target number of clones per shRNA, and the efficiency of genome editing

Notes:

i. The protocol describes plasmid delivery by nucleofection as this is most efficient, but if the necessary instrumentation is not available, plasmids can instead be transfected as described in the STAR METHODS

ii. While our standard protocol uses ZFNs to target the *AAVS1* locus, alternative programmable nucleases such as TALENs and CRISPR/Cas9 can also be used.

We opted for ZFNs because we previously demonstrated >95% on-target integration efficiency in hPSCs. That said, our targeting vector is fully compatible with other systems, including TALENs and widely used CRISPR/Cas9 reagents (e.g., Addgene #59025, #59026, #129726). Besides plasmid delivery, in side-by-side comparisons, we observed comparable integration efficiencies using 200,000 cells per nucleocuvette strip co-electroporated with the iPS2-seq vector and either 1 µg ZFN plasmids or Cas9 RNPs (1 µg Cas9 protein + 200 ng sgRNA)

- iii. Perform a pilot genome editing experiment to empirically determine the efficiency of genome editing in an hPSC line used for iPS2-seq for the first time
- 49.2 Prepare the following plasmids using the QIAGEN plasmid midiprep kit with endotoxin-free buffers, diluting them to 1 µg/µL (2.5 µg/nucleofection):
- i. **Pool of pAAV-Puro_siKD2.0** with the barcoded shRNAs (step 32 from Pooled cloning step 2: obtain the final plasmid pool)
 - ii. **pAAV-Neo_CAG**
 - iii. **pZFN-AAVS1_ELD**
 - iv. **pZFN-AAVS1_KKR**
- 49.3 Prepare Geltrex-coated 24-well culture dishes (1 well per nucleofection)
- 49.4 Culture hPSCs in Essential 8 at ~40% confluency (1×10^6 /nucleofection)
- 50 16 h before nucleofection, refresh the hPSC culture media to Essential 8 supplemented with 2 µM Thiazovivin. Do not add antibiotics until step 11
- 51 Prepare hPSC Essential 8 supplemented with CEPT [50 nM Chroman 1; 5 µM Emricasan; 0.1% Polyamine Supplement; 0.7 µM trans-ISRIB] (1.5 mL/nucleofection)
- 52 Remove the coating solution from the pre-coated culture dish, add 1 mL of Essential 8 with CEPT per well, and place the dish in the incubator to acclimatize to 37 °C
- 53 Prepare a mix of equal mass of the four plasmids from step 1b (10 µg/transfection) and nucleofection buffers (18 µL P3 Supplement 26 82 µL Nucleofector Solution per reaction)
- 54 Obtain a single cell suspension by treating hPSCs with StemPro Accutase, following the manufacturer's instructions, aliquot 1×10^6 hPSCs/nucleofection in a conical tube, and pellet the cells at 100 g for 5 min
- 55 Remove the supernatant, gently resuspend the pellet with the nucleofection mix from step 5, transfer 110 µL/nucleofection in a Nucleocuvette Vessel (avoiding bubbles), place the cuvette into the 4D-Nucleofector System and run program CA137
- 56 Add 500 µL of Essential 8 with CEPT to each cuvette and use the suction pipette to gently transfer the cell suspension to the 24-well plate from step 4. Place in the incubator

- 57 24 h post nucleofection, refresh media with Essential 8, and prepare Geltrex-coated 100 mm dishes (1 per nucleofection)
- 58 48 h post nucleofection, dissociate cells as small clumps with 0.5 mM EDTA in DPBS, and replat in Essential 8 with 2 μ M Thiazovivin (using one 100 mm dish per well of 24-wp)
- 59 Refresh media daily with Essential 8. Antibiotics can be added from now on
- 60 When hPSCs are ~50% confluent (~5–6 days post-nucleofection), perform dual positive drug selection for 4 days by supplementing Essential 8 with 0.5 μ g/mL puromycin and 25 μ g/mL Geneticin, refreshing media daily; include 2 μ M Thiazovivin for the first 2 days
Note: to optimize the procedure for a specific hPSC line, perform a kill curve with 0.25–2 μ g/mL puromycin and 12.5–100 μ g/mL Geneticin to identify the minimal concentration of each drug that in combination, eliminates all unedited hPSCs within 72 h
- 61 When hPSC colonies are ready for passaging (~1–2 mm in diameter), count the number of colonies (clones), dissociate cells as small clumps with 0.5 mM EDTA in DPBS, and replat in Essential 8 with 2 μ M Thiazovivin at a density of ~1 clones per cm^2
- 62 From 24 h post passaging until 5 days post passaging, refresh media daily with Essential 8 supplemented with 200 nM fialuridine (FIAU) to perform the negative selection
- 63 hPSC clonal pools are now ready for screening experiments and/or cryopreservation
- 63.1 **Notes:**
(a) Minimize the number of passages before the induction of gene knockdown, to minimize the expansion of hPSC clones with a genetic/epigenetic growth advantage
- 63.2 (b) To account for the possible emergence of hPSC clones with neuroectodermal bias (“iPSC-neuro” clones), ensure sufficient clonal representation per perturbation. iPSC2-seq enables robust retrospective control for this source of variability; however, it may still be helpful to screen for expression of early neuroectodermal markers such as ZIC1 or UNC5D. If a given pool contains an overrepresentation of iPSC-neuro clones, consider repeating the targeting or depleting biased clones using, for example, a UNC5D antibody for FACS or MACS
- 63.3 (c) For hPSC growth and differentiation protocols involving media containing FBS and/or BSA, batch-test these reagents for the absence of tetracycline contamination (i.e., by generating a clonal line expressing a validated inducible shRNA, to confirm that shRNA expression is not leaky in the absence of exogenously added tetracycline)

- 63.4 (d) To induce knockdown, add 1 ng/mL tetracycline, refreshing the media daily or every 48 h, when compatible for the specific differentiation protocol

hiPSC pool QC

- 64 This procedure aims to determine clonal diversity and shRNA representation in the genome-edited hiPSC pool by generating Illumina-compatible NGS libraries of UCI-BCs that can be analyzed using *catcher_step2QC* as described in Supplemental Protocol 4 - Pooled cloning step 2 and hiPSC genome editing QC.
- 65 After hPSC passaging, collect the cells from a confluent 6-well plate, and pellet the cells at 100 g for 5 min.
Note: Cell pellets can be frozen at -20 °C if needed.
- 66 Obtain purified genomic DNA (gDNA) of the genome-edited hiPSCs with the Monarch Spin gDNA Extraction Kit.
- 67 Elute samples with 50-100 µL of elution buffer and assess DNA purity and quantity by NanoDrop and Qubit broad range dsDNA.
- 68 Use 500 ng of gDNA per PCR reaction, performing as many reactions as needed to achieve coverage of at least 100-fold the estimated number of clones.
Note: One nanogram of human gDNA contains approximately 152 genome equivalents; therefore, 500 ng corresponds to ~76,000 potential UCI-BC templates at the AAVS1 locus.
- 69 Using the same primers as the Final plasmid pool QC, perform a PCR to amplify UCI-BCs from the genomic DNA, adding TruSeq adapters and DI, according to the adjusted calculations and cycling conditions:

A	B	C
Reagent	Volume (µL)	Final conc.
gDNA iPS2-seq pool	Variable	500 ng
Q5 Reaction Buffer (5X)	10	1X
Deoxynucleotide (dNTP) Solution Mix (10mM)	1	200µM
iPS2-seq_step2NGS_Primer (10 µM)	2.5	0.5µM
iPS2-seq_NGS_R Primer (10 µM)	2.5	0.5µM

	A	B	C
	Q5 Hot Start High-Fidelity DNA Polymerase (2 U/μL)	1,00	0.04 U/μL
	Nuclease-free water	To 50	
	Total volume	50	

Table SP1.19. PS2-seq transfected cells QC, PCR #1 mix

70

	A	B	C	D
	Steps	Temperature	Time	Cycles
	Initial Denaturation	98°C	30 sec	1
	Denaturation	98°C	10 sec	26 cycles
	Annealing	64°C	30 sec	
	Extension	72°C	20 sec	
	Final extension	72°C	2 min	1
	Hold (rapid ramp down)	10°C	Hold	

Table SP1.20. Cycling conditions for PCR #1

- 71 Purify the expected ~144 bp amplification product and purify it using SPRIselect beads with a double side selection using 0.6X & 1.5X bead ratios to remove primer dimers and small fragments.
- 72 Dilute the PCR product to 0.1 ng/μL in EB buffer and perform indexing PCR according to Table SP1.17 and Table SP1.18.
- 73 An indexed library product of ~228 bp is expected. Purify the reaction using a double-sided SPRIselect bead cleanup (0.65X & 1.5X ratios).
- 74 Run 2 μL of the undiluted and 1:10 final library on a TapeStation High Sensitivity D1000 Screen Tape assay, following the manual.
Note: This is the same amplicon described for the final plasmid QC (Figure SP1.6B), but adjusted to work on genomic DNA material.

- 75 Perform NGS using MiSeq or NextSeq 1000/2000 with the following settings: read 1 - 68 cycles; index 1 - 10 cycles; index 2 - 10 cycles (>10,000 reads/clone).

OPTIONAL: hiPSC clone genotyping

- 76 This optional procedure can be used to assess the efficiency of biallelic AAVS1 genome editing in a subset of hPSC clones, or to isolate homozygously edited clones for validation studies. Genotyping leverages various genomic PCRs to determine on-site integration, copy number, and random integrations (Figure SP1.7)




Figure SP1.7. iPS2-seq clone genotyping strategy

Schematic of the genotyping strategies to monitor integration of *pAAV-Puro_siKD2.0* in allele 1 and *pAAV-Neo_CAG* in allele 2. Primer binding sites and directions are shown on the plasmid or integrated cassettes. PCR types and products are numbered as in Table SP1.23.

- 77 After step 13 of Genome editing: obtain hPSC clonal pools, pick and expand individual hPSC clones.
Note: treat clones with 200 nM FIAU to test for random integrations.
- 78 Prepare a culture dish with $2-5 \times 10^5$ cells, and extract gDNA using the Monarch Spin gDNA Extraction Kit, diluting to a concentration of 50-200 ng/ μ L
- 79 Perform PCR reactions according to the following calculations and cycling conditions, employing the primer pairs listed in **Table SP1.23**. Include the following controls: wild type hPSCs; pAAV-Puro_siKD2.0; no template

80

	A	B	C
	Reagent	one reaction (μL)	Final conc.
	Genomic DNA	1	50-200ng
	LongAmp Taq Reaction Buffer (5X)	2	1X
	Deoxynucleotide (dNTP) Solution Mix (10mM)	0,3	300μM
	Forward Primer (10μM)	0,5	0.5μM
	Forward Primer (10μM)	0,5	0.5μM
	DMSO	0,2	2\%
	LongAmp Taq DNA Polymerase (2.5U/μL)	0,4	0.1 U/μL
	Nuclease-free water	To 10	
	Total volume	10	

Table SP1.21. Genotyping mix

81

	A	B	C	D
	Steps	Temperature	Time	Cycles
	Initial Denaturation	94°C	5 min	1
	Denaturation	94°C	15 sec	35 cycles
	Annealing	Tab step SP1.23	30 sec	
	Extension	65°C	Tab step SP1.23	
	Final extension	65°C	5 min	1
	Hold (rapid ramp down)	4°C	Hold	

Table SP1.22. Cycling conditions for genotyping

82 Run half the PCR products on 0.8% (w/v) agarose gels and determine the genotype based on Table SP1.24

83

PCR type NB	PCR type	Primer	Primer name	Primer location	Amplicon wild type ^a	Amplicon target ^b	Amplicon plasmid ^c	Annealing temp. ^d	Extension time ^e
1	AAVS1 locus	FW	Locus_fw	Genomic, 5' to 5' HA	1692 bp	No band	no band	65 °C	1' 30"
		REV	Locus_rev	Genomic, 3' to 3' HA					
2	5' junction (PuroR)	FW	Locus_fw	Genomic, 5' to 5' HA	no band	991 bp	no band	65 °C	1'
		REV	Puro_rev	Puro resistance					
3	3' junction (OPTetR)	FW	OPTetR_fw	OPTetR cDNA	no band	1463 bp	no band	60 °C	1' 30"
		REV	Locus_rev	Genomic, 3' to 3' HA					
4	5' junction (NeoR)	FW	Locus_fw	Genomic, 5' to 5' HA	no band	1032 bp	no band	60 °C	1'
		REV	Neo_rev	Neo resistance					
5	Plasmid (OPTetR)	FW	OPTetR_fw	OPTetR cDNA	no band	no band	1819 bp	60 °C	2'
		REV	Backbone_rev	Backbone, 3' to 3' HA					

Table SP1.23. Genotyping strategies for iPS2-seq clone genotyping

aResult of PCR on wild type AAVS1 allele

bResult of PCR on iPS2-seq-targeted AAVS1 allele

cResult of PCR on pAAV-Puro_siKD2.0 (positive control for off-target plasmid integration)

dVariable parameter in PCR protocol (Table SP1.22)

84

AAVS1 locus	5' junction (PuroR)	3' junction (OPTetR)	5' junction (NeoR)	Plasmid (OPTetR)	AAVS1 genotype	FIAU sensitive?	Number of shRNAs
-	+	+	+	-	Compound heterozygous	No	1
-	+	+	+	+	Compound heterozygous	No	1
+	+	+	-	-	Heterozygous iPS2-seq	No	1
+	+	+	-	+	Heterozygous iPS2-seq	No	1
-	+	+	+	+	Compound heterozygous	Yes	>1
+	+	+	-	+	Heterozygous iPS2-seq	Yes	>1
-	+	+	-	-	Homozygous iPS2-seq	No	2
-	+	+	-	+	Homozygous iPS2-seq	No	2
-	+	+	-	+	Homozygous iPS2-seq	Yes	>2
+	-	-	-	-	Wild type	No	0
+/-	-	+/-	+/-	-	Incorrect iPS2-seq targeting	No	0
+/-	+/-	-	+/-	-	Incorrect iPS2-seq targeting	No	0
+/-	-	+/-	+/-	+	Incorrect iPS2-seq targeting	No	0
+/-	+/-	-	+/-	+	Incorrect iPS2-seq targeting	No	0
+/-	-	+/-	+/-	+	Incorrect iPS2-seq targeting	Yes	1 or more
+/-	+/-	-	+/-	+	Incorrect iPS2-seq targeting	Yes	1 or more

Table SP1.24. Interpretation of iPS2-seq clone genotyping

Protocol 2: iPS2-sci-seq library preparation

85 Decoding single-cell transcriptional perturbomes in iPS2-seq relies on detecting the shRNA-associated UCI-BC within the 3' UTR of the OPTetR transgene (Figure 1A). This

can be in principle achieved with any scRNA-seq method that tags and counts the 3' ends of mRNAs, usually as a result of pA-primed RT.

iPS2-sci-seq is a home-brew adaptation of the 2-level indexing sci-RNA-seq protocol¹⁹, that relies on 2 rounds of transcriptome barcoding through sequential pool-split of defined cell numbers so that the resulting combinatorial barcode in each cell is statistically most likely unique. This approach requires limited custom reagents and supports the analysis of 1,000-10,000 cells per experiment, a scale compatible with the everyday needs of most laboratories. Minor adaptations of this protocol that require additional reagents can enable the 3-level indexing variant of sci-RNA-seq (sci-RNA-seq³), which supports analysis of over ~1 million cells per experiment. Since transcript counts in sci-RNA-seq are notoriously low and zero-inflated (i.e., genes may not be counted at least once in cells that express them), it is key to incorporate a strategy that explicitly enriches for sequences containing the UCI-BC, thereby detecting it with higher sensitivity and specificity. Therefore, in the iPS2-sci-seq protocol, UCI-BCs are enriched at both transcriptome barcoding steps: first during RT, using a set of target-specific primers that share RT barcodes with standard pT primers, well by well, and secondly during PCR, using a set of target-specific primers that share i7 barcodes with standard P7 primers, once again well by well. The resulting multiplexed NGS library pools contain both single-cell transcriptomes and UCI-BCs, which share library structure and carry identical combinatorial indexes, cell by cell.


This Supplemental Protocol describes the following procedures:

1. Reagent setup
2. Nuclei preparation
3. Reverse transcription
4. Nuclei sorting
5. Second-strand synthesis
6. Tagmentation
7. Indexing PCR and NGS

86 Reagent setup: Prepare the RT and PCR oligos using nuclease-free water:

- 86.1 Obtain master stocks of the 96 indexed iPS2-sci-seq_pT primers (250 μ M) and iPS2-sci-seq_tetR primers (100 μ M). Use the RT_indexes indicated in Table SP2.1
- 86.2 Prepare a working stock of multiplexed RT primers by diluting in the same plate iPS2-sci-seq_pT primers (25 μ M) and iPS2-sci-seq_tetR primers (2.5 μ M), matching well positions so that pooled primers share the same RT index
- 86.3 Obtain master stocks of indexed iPS2-sci-seq_P5, iPS2-sci-seq_P7, and iPS2-sci-seq_P7A (all at 100 μ M). Use the i5 and i7 indexes indicated in Table SP2.1.

Note: it is not necessary to obtain 96 indexed oligos of each type, as i5 and i7 can be combined to generate unique combinations.

- 86.4 Prepare a working stock of iPS2-sci-seq_P5 (10 μM).
- 86.5 Prepare a working stock of multiplexed P7 primers by diluting in the same plate iPS2-sci-seqP7 (10 μM) and iPS2-sci-seqP7A (5 μM), matching well positions so that pooled primers share the same i7 index.
- 87 Prepare the following buffers using nuclease-free water, and chill them on ice:
- 87.1 Nuclei buffer (store at 4 °C for up to 1 month, 5 mL/sample): 10 mM Tris-HCl, pH 7.4, 10 mM NaCl, 3 mM MgCl2.
- 87.2 Nuclei lysis buffer (made fresh, 1 mL/sample): nuclei buffer with 0.1-0.3% IGEPAL CA-630, 1% SUPERase In RNase Inhibitor (20 U/μL) and 1% BSA (20 mg/mL).
Note: the optimal concentration of IGEPAL CA-630 for a given cell type should be determined in a pilot experiment. For reference, 0.1% works for hPSCs, but not for hPSC-CMs which are lysed appropriately only at 0.3%.
- 87.3 Nuclei suspension buffer (make fresh, 3 mL/sample): nuclei buffer with 1% SUPERase In RNase Inhibitor (20 U/μL) and 1% BSA (20 mg/mL).
- 87.4 4% PFA (store at 4 °C for up to 1 month, 10 mL/sample). 
- 87.5 DPBS without Ca2+ and Mg2+ (store at 4 °C for up to 1 month, 10 mL/sample). 

88

	Well	RT_index (5'-3')	i7 for primer (5'-3')	i7 for sample sheet (5'-3')	i5 for primer & sample sheet, forward strand (5'-3')	i5 for sample sheet, reverse strand (5'-3')		Well	RT_index (5'-3')	i7 for primer (5'-3')	i7 for sample sheet (5'-3')	i5 for primer & sample sheet, forward strand (5'-3')	i5 for sample sheet, reverse strand (5'-3')



	A 01	TTC TCG CAT G	CCG AAT CCG A	TCG GAT TCG G	CTC CAT CGA G	CTC GAT GGA G		E 01	AGA GAA GGT T	TTC GTT CCA T	ATG GAA CGA A	AGG CGA GAG C	GCT CTC GCC T
	A 02	TCC TAC CAG T	ATA AGC CGG A	TCC GGC TTA T	TTG GTA GTC G	CGA CTA CCA A		E 02	CAT ACT CC GA	TAC CTA ATC A	TGA TTA GGT A	TCA AGA TAG T	ACT ATC TTG A
	A 03	GC GTT GGA GC	CCG GCG GCG A	TCG CCG CCG G	GGC CGT CAA C	GTT GAC GGC C		E 03	GCT AAC TTG C	AAG TAA TAT T	AAT ATT ACT T	TAA TTG ACC T	AGG TCA ATT A
	A 04	GAT CTT ACG C	GGC TTG CCA A	TTG GCA AGC C	CCT AGA CGA G	CTC GTC TAG G		E 04	AAT CCA TCT T	AGC TAA GAA T	ATT CTT AGC T	CAG CCG GCT T	AAG CCG GCT G
	A 05	CTG ATG GTC A	CCG CTA GCT G	CAG CTA GCG G	TCG TTA GAG C	GCT CTA ACG A		E 05	GG CTG AGC TC	GTC GAG GTA T	ATA CCT CGA C	AGA ACC GGA G	CTC CGG TTC T
	A 06	CC GAG AAT CC	CTT ATC CTA C	GTA GGA TAA G	CGT TCT ATC A	TGA TAG AAC G		E 06	CC GAT TCC TG	TTA TTA GTA G	CTA CTA ATA A	GAG ATG CAT G	CAT GCA TCT C
	A 07	GC CG CAA CGA	TGA GCT ACT T	AAG TAG CTC A	CGG AAT CTA A	TTA GAT TCC G		E 07	ACC GC CAA CC	TGC GAA GAT C	GAT CTT CGC A	GAT TAC CGG A	TCC GGT AAT C
	A 08	TGA GTC TGG C	TCA GGA CTT A	TAA GTC CTG A	ATG ACT GAT C	GAT CAG TCA T		E 08	TGG CCT GAA G	AAC TAC GGC T	AGC CGT AGT T	TCG TAA CGG T	ACC GTT ACG A
	A 09	TGC GGA CCT A	CCG CAG CCG C	GCG GCT GCG G	TCA ATA TCG A	TCG ATA TTG A		E 09	AAC CTC ATT C	AAC GGA ACG C	GCG TTC CGT T	TGG CGA CGG A	TCC GTC GCC A
	A 10	ACC TCG TTG A	TGC GCC TGG T	ACC AGG CGC A	GTA GAC CTG G	CCA GGT CTA C		E10	ATA AGG AGC A	GAT GCT ACG A	TCG TAG CAT C	AGT CAT AGC C	GGC TAT GAC T
	A 11	ACG GAG GC GG	AAT CAT ACG G	CCG TAT GAT T	TTA TGA CCA A	TTG GTC ATA A		E11	CGA ACG CC GG	ATC TGC CAA T	ATT GGC AGA T	GTC AAG TCC A	TGG ACT TGA C
	A 12	TAG ATC TAC T	CGC CAA TCA A	TTG ATT GGC G	TTG GTC CGT T	AAC GGA CCA A		E12	GGT ATG CTT G	ATC GTA TCA A	TTG ATA CGA T	ATT CGG AAG T	ACT TCC GAA T



	B 01	AAT TAA GAC T	CAA GGC TTA G	CTA AGC CTT G	GGT ACG TTA A	TTA ACG TAC C		F 01	AAC CTG CGT A	AAC GCC TCT A	TAG AGG CGT T	GTC GGT AGT T	AAC TAC CGA C
	B 02	CCA TTG CGT T	GCG CTC GAC G	CGT CGA GCG C	CAA TGA GTC C	GGA CTC ATT G		F 02	GG CAG ACG CC	ACG GCA ACC A	TGG TTG CCG T	AGG ACG GAC G	CGT CCG TCC T
	B 03	TTA TTC ATT C	TCC AGC AAT A	TAT TGC TGG A	GAT GCA GTT C	GAA CTG CAT C		F 03	TAG CC GTC AT	CAG GCT AAG A	TCT TAG CCT G	CTC CTG GAC C	GGT CCA GGA G
	B 04	ATC TCC GAA C	CAT GAG AAC T	AGT TCT CAT G	CCA TCG TTC C	GGA ACG ATG G		F 04	CCT GGA AGA G	CGC AAT ATC A	TGA TAT TGC G	TAG CCT CGT T	AAC GAG GCT A
	B 05	TTG ACT TCA G	AAC GTA ATC T	AGA TTA CGT T	TTG AGA GAG T	ACT CTC TCA A		F 05	GGA GGT TCT A	TTC GAT AAC C	GGT TAT CGA A	GGT TGA ACG T	ACG TTC AAC C
	B 06	GG CAG GTA TT	ATT CTC CTC T	AGA GGA GAA T	ACT GAG CGA C	GTC GCT CAG T		F 06	CTA GTA GTC T	AAC CTC AAG A	TCT TGA GGT T	AGG TCC TCG T	ACG AGG ACC T
	B 07	AGA GCT ATA A	TCT GCG CGT T	AAC GCG CAG A	TGA GGA ATC A	TGA TTC CTC A		F 07	ATC ATC AAC G	CAG GCG CCA T	ATG GCG CCT G	GGA AGT TAT A	TAT AAC TTC C
	B 08	CTA AGA GAA G	GCT CAT ATG C	GCA TAT GAG C	CCT CCG ACG G	CCG TCG GAG G		F 08	ACG CGA GAT T	AAC TAT TAT A	TAT AAT AGT T	TGG TAA TCC T	AGG ATT ACC A
	B 09	ACT CAA TAG G	AGC GGT AAC G	CGT TAC CGC T	CAT TGA CGC T	AGC GTC AAT G		F 09	GAA GAG GCA T	AAG TTA CCT A	TAG GTA ACT T	AAG CTA GGT T	AAC CTA GCT T
	B 10	CTT GC GC CG C	AAT GAA TAG T	ACT ATT CAT T	TCG TCC TTC G	CGA AGG ACG A		F1 0	GGT ATC CG CC	CGG CAG AGG A	TCC TCT GCC G	TCC GCG GAC T	AGT CCG CGG A
	B 11	AAT CGT AGC G	CCG TAT CTG G	CCA GAT ACG G	TGA TAC TCA A	TTG AGT ATC A		F1 1	AAC TAG GC GC	GCC TCA ATA A	TTA TTG AGG C	TGC GGA TAG T	ACT ATC CGC A



	B 12	GGT ACT GC CT	CCT TAG TCT G	CAG ACT AAG G	TTC TAC CTC A	TGA GGT AGA A		F1 2	TCG CTA AGC A	TTA ACG CCG T	ACG GCG TTA A	TGG CAG CTC G	CGA GCT GCC A
	C 01	TAG AAT TAA C	ACC TAG TTA G	CTA ACT AGG T	TCG TCG GAA C	GTT CCG ACG A		G 01	TAT ATA CTA A	CAT ACG ATG C	GCA TCG TAT G	TGC TAC GGT C	GAC CGT AGC A
	C 02	GC CAT TCT CC	ATA GGA GTA C	GTA CTC CTA T	ATC GAG ATG A	TCA TCT CGA T		G 02	ACT TGC TAG A	AAG CTG ACC T	AGG TCA GCT T	GCG CAA TGA C	GTC ATT GCG C
	C 03	TGC CG GCA GA	CTA CGA CGA G	CTC GTC GTA G	TAG ACT AGT C	GAC TAG TCT A		G 03	AAC CAT TGG A	GAG TCC TTA T	ATA AGG ACT C	CTT AAT CTT G	CAA GAT TAA G
	C 04	TTA CC GAG GC	AGT CGA GTT C	GAA CTC GAC T	GTC GAA GCA G	CTG CTT CGA C		G 04	TCG CG GTT GG	CCT ACG GCA A	TTG CCG TAG G	GGA GTT GCG T	ACG CAA CTC C
	C 05	ATC ATA TTA G	TGG TCC AGT C	GAC TGG ACC A	AGG CGC TAG G	CCT AGC GCC T		G 05	CGT AGT TAC C	AAT ATT CGA A	TTC GAA TAT T	ACT CGT ATC A	TGA TAC GAG T
	C 06	TGG TCA GC CA	ATC TAA GCA A	TTG CTT AGA T	AGA TGC AAC T	AGT TGC ATC T		G 06	TCC AAT CAT C	TTC AAG AAT C	GAT TCT TGA A	GGT AAT AAT G	CAT TAT TAC C
	C 07	ACT ATG CAA T	CGA ATT CGT T	AAC GAA TTC G	AAG CCT ACG A	TCG TAG GCT T		G 07	AAT CGA TAA T	ATG CTC GCA A	TTG CGA GCA T	TCC TTA TAG A	TCT ATA AGG A
	C 08	CGA CG CGA CT	CAG CGA TAG A	TCT ATC GCT G	GTA GGC AAT T	AAT TGC CTA C		G 08	CCA TTA TCT A	GGA GTA AGC C	GGC TTA CTC C	CCG ACT CCA A	TTG GAG TCG G
	C 09	GAT ACG GAA C	GGT CGC TAT G	CAT AGC GAC C	TGC CAG TTG C	GCA ACT GGC A		G 09	TCA ACG TAA G	TTA TCG TAT T	AAT ACG ATA A	GCC AAG CTT G	CAA GCT TGG C
	C 10	TTA TCC GGA T	ATC CGT TAG C	GCT AAC GGA T	CTT AGG TAT C	GAT ACC TAA G		G 10	TCT AAT AGT A	AAG TCT AAT A	TAT TAG ACT T	CAT ATC CTA T	ATA GGA TAT G
	C 11	TAG AGT AAT A	TCG CAA TTA G	CTA ATT GCG A	GAG ACC TAC C	GGT AGG TCT C		G 11	AAC CG CTG GT	CGG CTT ACT A	TAG TAA GCC G	ACC TAC GCC A	TGG CGT AGG T



	C 12	GCA GGT CC GT	GGC TGG CTA G	CTA GCC AGC C	ATT GAC CGA G	CTC GGT CAA T		G 12	GAT CG CTT CT	GAT ATG GTC T	AGA CCA TAT C	GGA ATT CAG T	ACT GAA TTC C
	D 01	TCG GC CTT AC	ACG GTC TTG C	GCA AGA CCG T	GGA GGC GGC G	CGC CGC CTC C		H 01	CTA ACT AGA T	TAG TCG TCC A	TGG ACG ACT A	TGG CGT AGA A	TTC TAC GCC A
	D 02	AGA ACG TCT C	GCT CCA TTC G	CGA ATG GAG C	CCA GTA CTT G	CAA GTA CTG G		H 02	GCT GGA ACT T	TAG CTG CTA C	GTA GCA GCT A	ATT GCG GCC A	TGG CCG CAA T
	D 03	CCA GTT CCA A	ACG ATA AGC G	CGC TTA TCG T	GGT CTC GCC G	CGG CGA GAC C		H 03	AGG TTA GTT C	CTC TTC AAG C	GCT TGA AGA G	TTC AGC TTG G	CCA AGC TGA A
	D 04	GG CGT TAA GG	ACC ATA GCG C	GCG CTA TGG T	GGC GGA GGT C	GAC CTC CGC C		H 04	CAT TCG ACG G	ATG AAC GCG C	GCG CGT TCA T	CCA TCT GGC A	TGC CAG ATG G
	D 05	ACT TAA CCT T	CTC TTA GCG G	CCG CTA AGA G	TAG TTC TAG A	TCT AGA ACT A		H 05	CAT TCA ATC A	GTC GAC GGA A	TTC CGT CGA C	CTT ATA AGT T	AAC TTA TAA G
	D 06	CAA CC GCT AA	TGA TTC AAC T	AGT TGA ATC A	TTG GAG TTA G	CTA ACT CCA A		H 06	CG GAT TAG AA	ACT AAT TGA G	CTC AAT TAG T	GAT TAG ATG A	TCA TCT AAT C
	D 07	GAC CTT GAT A	TAT GGC CGC G	CGC GGC CAT A	AGA TCT TGG T	ACC AAG ATC T		H 07	ATC GG CTA TC	CTT GCA TAA T	ATT ATG CAA G	TAT AGG ATC T	AGA TCC TAT A
	D 08	TCT GAT ACC A	AGA GGT CGC A	TGC GAC CTC T	GTA ATG ATC G	CGA TCA TTA C		H 08	CCT TGA TCG T	TCC TTA CCA A	TTG GTA AGG A	AGC TTA TAG G	CCT ATA AGC T
	D 09	GAA GAT CGA G	AGG AGA TTG A	TCA ATC TCC T	CAG AGA GGT C	GAC CTC TCT G		H 09	ACG AAG TCA A	TGC AGC CTA C	GTA GGC TGC A	GTC TGC AAT C	GAT TGC AGA C
	D 10	AGG AGC GGT A	GGC TAT ATA G	CTA TAT AGC C	TTA ATT AGC C	GGC TAA TTA A		H 10	TTA CCT CGA C	GGA GCT GAG G	CCT CAG CTC C	CGC CTC TTA T	ATA AGA GGC G
	D 11	AAG AAG CTA G	TCG CGT ACT T	AAG TAC GCG A	CTC TAA CTC G	CGA GTT AGA G		H 11	GGA GGA TAG C	GCA GCG GAC T	AGT CCG CTG C	GTT GGA TCT T	AAG ATC CAA C



	D 12	TCC GG CCT CG	AAT AAT AAT G	CAT TAT TAT T	TAC GAT CAT C	GAT GAT CGT A		H 12	GG CTC TCT AT	CAT CGC GCT C	GAG CGC GAT G	GCG ATT GCA G	CTG CAA TCG C

Nuclei preparation

- 89 This is an exemplary protocol for one sample containing $1-5 \times 10^6$ starting cells. Perform all procedures on ice using chilled reagents and using a refrigerated swing bucket centrifuge.
- 90 Obtain a single-cell suspension using a dissociation protocol that preserves >95% cell viability. Eliminate any cell clump using a 40-100 μm cell strainer.
- 91 Extract and fix nuclei:
- 91.1 Centrifuge cells in a 15 mL conical at 150 g for 5 min, discard the supernatant, resuspend the pellet in 1 mL nuclei lysis buffer, and incubate for 5 min.
- 91.2 Centrifuge nuclei at 300 g for 3 min, discard the supernatant, re-suspend the pellet in 100 μL nuclei suspension buffer, add 10 mL 4% PFA, mix and incubate for 15 min.
- 91.3 Centrifuge fixed nuclei at 300 g for 5 min, discard the supernatant, and wash the pellet three times with 1 mL nuclei suspension buffer by repeating this step.
- 91.4 After the final wash, resuspend the pellet in 300 μL nuclei suspension buffer.
- 92 Count nuclei with a hemocytometer with YOYO-1 staining (Figure SP2.1), dilute them to 2,500 nuclei/ μL , prepare single-use aliquots, and either proceed with the next section or snap-freeze nuclei in liquid nitrogen and store them at -80°C for up to one month.




Figure SP2.1. Nuclear extraction QC

Exemplary picture of PFA-fixed nuclei stained with YOYO-1 visualized on a hemocytometer. Nuclei are round, not clumped, and there are minimal cell debris.

Reverse transcription

- 93 Clean the work area and instrumentation to prevent RNase contamination. Unless stated otherwise, keep samples on ice, perform centrifugations at 4 °C, use DNA LoBind 1.5 mL tubes and low-retention tips, and set the thermocycler heated lid at 105 °C to prevent evaporation.
- 93.1 Prepare fresh stop solution (40 mM EDTA and 1 mM spermidine; 600 µL/plate)
- 93.2 Thaw PFA-fixed nuclei in hand or at 37 °C thermal block, mix by flicking, then move to ice
- 93.3 Prepare a mix by combining nuclei suspension (2 µL, 5,000 nuclei, per well) and 10 mM dNTP (0.25 µL per well), and distribute 2.25 µL into each well of a 96-well LoBind plate

- 93.4 Add 1 μ L of multiplexed RT primers to each well (step 1b of Reagent setup). Cap the plate and centrifuge it for 10 seconds at 100 g
- 93.5 Incubate the plate at 55 °C for 5 minutes and immediately place it on ice. After 3-5 minutes on ice, centrifuge the plate for 10 seconds at 100 g
- 93.6 Prepare the following RT mix, distribute 1.75 μ L to each well without mixing, cap the plate, and centrifuge it for 10 seconds at 100 g
- 93.7

	A	B	C
		one well (ul)	one plate
	5X Superscript IV First-strand buffer	1 ul	105 ul
	100mM DTT	0.25 ul	26.25 ul
	Superscript IV reverse transcriptase	0.25 ul	26.25 ul
	RNAseOUT Recombinant Ribonuclease Inhibitor	0.25 ul	26.25 ul

Table SP2.2. iPS2-sci-seq RT mix
Reagent

- 93.8 Incubate the plate at 55 °C for 10 minutes, immediately place it on ice, and add 5 μ L of stop solution into each well to stop the reaction. After 3-5 minutes on ice, centrifuge the plate for 10 seconds at 100 g.
- 93.9 Gently pool all wells into a FACS tube using wide bore tips, passing through a 40 μ m cell strainer, and immediately proceed with the next procedure.

Nuclei sorting

- 94 Prepare collection plate(s) by adding 5 μ L of elution buffer to each well one or more 96-well LoBind plates, cover, and place on ice
- 95 Aliquot 50 μ L of nuclei suspension (8 from Reverse transcription) in a separate tube to be used as unstained control; add 3 μ M DAPI to the remaining volume to stain nuclei

- 96 Setup gating for single nuclei (Fig. SP2.2), and sort 25 nuclei into each wells of collection plate
(assuming 96 indexes in the RT and an accepted rate of ~10% index collisions)
Note: use a cell sorter equipped with a 100 μ m nozzle and minimize sample pressure, chill both sample and collection plate, and mix the sample while sorting
- 97 Centrifuge the plates for 1 minute at 900 g, and either proceed with the next step or store plates at -80 $^{\circ}$ C for up to one week




Figure SP2.2. Nuclei FACS Exemplary gating strategy for single nuclei sorting. Set up the gates based on the negative control (top) and sort the DAPI+ population. DAPI+ single nuclei are back-gated in red; debris are in light blue.

Second strand synthesis

- 99 Prepare reaction master mix, distribute 1 μ L to each well, pipette mix, cap the plate, and centrifuge it for 10 seconds at 100 g

100

	A	B	C
		one well (ul)	one plate
	mRNA Second Strand Synthesis buffer (NEB)	0.6ul	67.5ul
	mRNA Second Strand Synthesis enzyme (NEB)	0.3ul	33.75ul
	mol bio ddH ₂ O	0.1ul	11.25ul

Table SP2.3. iPS2-sci-seq second strand synthesis mix



- 101 Carry out second strand synthesis at 16 °C for 180 minutes, and terminate the reaction by incubating at 75 °C for 20 minutes (lid kept at 80 °C for both steps)
- 102 Centrifuge the plate for 10 seconds at 100 g and either immediately proceed with the next section or store at 4 °C overnight

Tagmentation

- 103 Prepare the tagmentation mix with the Nextera TD buffer (5.75 µL/well) and the TDE1 enzyme (0.01 µL/well), distribute 5.76 µL of mix to each well of the sorting plate, pipette mix, cap the plate, and centrifuge it for 10 seconds at 100 g
- 104 Incubate the plate at 55 °C for 5 min, immediately place it on ice, add 12 µL of Zymo DNA binding buffer to each well, and incubate at room temperature for 5 min.
- 105 Purify tagmented dsDNA using 1.5X volumes of SPRIselect beads (36 µL/well), following the manufacturer's instructions and eluting in 17 µL of elution buffer
- 106 Transfer 16 µL of eluted DNA to a new 96-well LoBind plate and either proceed with the next section or store at 4 °C overnight

Indexing PCR and NGS

- 107 Add 2 µL of 10 µM P5 primer and 2 µL of 10 µM multiplexed P7 primers (step 1e from Reagent setup) to each well, applying a unique combination of i5 and i7
- 108 Add 20 µL NEBNext High-Fidelity 2X PCR Master Mix to each well, pipette mix, cap the plate, and centrifuge it for 10 seconds at 100 g, and run the following PCR:

109

	cycles	temper ature	timing	step
	1	75°C	3 min	
	1	98°C	30 sec	Initial Denaturatio n

	18-22 cycles	98°C	10 sec	Denaturation
	66°C	30 sec	annealing	
	72°C	1 min	extension	
	1	72°C	5 min	final extension

Table SP2.4. iPS2-sci-seq indexing PCR cycling conditions

- 110 Pool all samples of a 96-well plate in a 15 mL Falcon tube and concentrate on 4 columns of the Zymo DNA Clean & Concentrator kit, following the manufacturer's instructions except without performing any washes and eluting in a total volume 25 µL of elution buffer
- 111 Further purify using 0.85X volumes of SPRIselect beads, according to the manufacturer's instructions and eluting in 50 µL of elution buffer
- 112 Run 2 µL of a 1:10 dilution of the NGS library on a TapeStation High Sensitivity 5000 Screen Tape assay, following the manual. Expect a main peak at ~350 bp (Figure SP2.3)
- 113 Perform NGS using NextSeq 1000/2000 and 50 cycles reagents: read 1 - 18 cycles; index 1 & index 2 - 10 cycles each; read 2 - 52 cycles (>25,000 reads/nuclei)
- 114




Figure SP2.3. iPS2-sci-seq NGS library QC

Exemplary TapeStation QC of iPS2-sci-seq NGS library (1:10 dilution).

Protocol 3: iPS2-10X-seq library preparation

115 iPS2-10X-seq relies on microfluidics partitioning-based scRNA-seq from 10X Genomics. Despite the higher costs, the ease and reproducibility of this commercial method have led to its widespread adoption. The procedure integrates the manufacturer's protocol to generate an additional NGS library containing UCI-BCs. This is sequenced in parallel to the standard gene expression library to enable *in silico* matching of UCI-BCs and single-cell transcriptomes through their shared unique cell barcodes. Of note, iPS2-10X-seq provides built in support for the cost-saving strategy of reasonably overloading the microfluidics chip, since cell doublets are readily identified by the co-expression of two different UCI-BCs and filtered out as part of the standard *catcher_10Xcatch* pipeline. Thus, it is possible to load the manufacturer-recommended cell number for barcoded sample pools even without additional multiplexing barcodes. Other cost-saving strategies include using cell multiplexing oligos (CMOs) or barcoded antibodies against housekeeping cell surface antigens to multiplex tet-untreated and tet-treated cells in a single microfluidics channel, enabling robust control over knockdown experimental design without performing two separate reactions. This Supplemental Protocol describes the following procedures:

1. Prepare the gene expression library
2. Prepare the shRNA barcode library
3. Library pooling and NGS

Prepare the gene expression library

- 116 Obtain a single-cell suspension using a dissociation protocol that preserves >95% cell viability. Avoid the use of DNase, and remove any EDTA with at least three washes.
- 116.1 **OPTIONAL:** Label multiple samples using TotalSeq-B hashtags and/or other antibodies with ADTs, following the manufacturer's instructions (CG000149 Rev D).
- 116.2 **OPTIONAL:** Label multiple samples using cell multiplexing oligos (CMOs) in the 3' CellPlex Kit from 10X Genomics, following the manufacturer's instructions (CG000391 Rev B).
- 116.3 **OPTIONAL:** Sort live cells after staining with Fixable Viability Dye eFluor 780.
***Note:** this step increases the fraction of cells that deliver high-quality transcriptomes, and is particularly recommended for organoids that are difficult to dissociate. Avoid the use of other live/dead staining that can intercalate with DNA.*
- 117 Centrifuge cells/cell pools at 150 g for 5 min at 4°C, and resuspend in ice cold DPBS 1% BSA aiming for 1×10^6 cells/mL.
- 118 Manually count cells, and calculate the cell volume needed to load the target cell number per channel of Chromium Next GEM Chip G.
Note: since iPS2-seq cells express shRNA barcodes, chips can be reasonably over-loaded even without any CMO or hashtags (~49,500 cells/reaction at a ~24% doublet rate).
- 119 Follow the 10X Genomics protocols to obtain and quality control the gene expression library (GEX) and, when appropriate, the CMO library and/or the ADT library.
- 120 Save a portion of the pre-amplified full-length cDNA for the next procedure.

Prepare the shRNA barcode library

- 121 Perform a two-step PCR to enrich the UCI-BC region in the 3' UTR of the OPTtetR cDNA.
- 122 Perform the first PCR according to the calculations and cycling conditions below:
 - 122.1

A	B
Reagent	Volume (ul)
NEBNex High-Fidelity 2X PCR Master Mix	10 ul
iPS2-10X-seq_truseq	1 ul
iPS2-10X-seq_inner	1 ul
cDNA 30-60ng/uL	variable
H2O (bring to volume)	to 10ul
Total volume	20 ul

Table SP3.1. iPS2-10X-seq first PCR mix

122.2

A	B	C	D
Steps	Temperature	Time	Cycles
Initial denaturation	98°C	30 sec	1
Denaturation	98°C	10 sec	15
Annealing	60°C	20 sec	
Extension	72°C	20 sec	
Final extension	72°C	30 sec	1

Table SP3.2. Cycling conditions for iPS2-10X-seq first PCR

122.3 **Note:** When performing iPS2-multi-seq, adjust this step by increasing the input material to 200 ng and reducing the cycle amplification to 11 cycles. Adaptations of the current protocol were essential to obtain clean UCI-BC library profiles.

123 Purify the PCR with 1.8X volumes of SPRIselect beads, following the manufacturer's instructions and eluting in 20 µL of elution buffer.

124 Run and quantify 2 µL of undiluted PCR on a TapeStation High Sensitivity 5000 Screen Tape assay, following the manual. Expect a main peak at ~285 bp (Figure SP3.1A).

- 125 Dilute the first PCR to 0.1 ng/μL and use it as a template for the second PCR according to the calculations and cycling conditions below:

125.1

A	B
Reagent	Volume (μL)
NEBNext High-Fidelity 2X PCR Master Mix	10
Dual Index Kit TT Set A primers	4
First PCR (0.1 ng/uL)	variable
Nuclease-free water	to 20ul
Total volume	20

Table SP3.3. iPS2-10X-seq second PCR mix

125.2

A	B	C	D
Step	Temperature	Timing	Cycles
Initial denaturation	98 °C	30 sec	1
Denaturation	98 °C	10 sec	10
annealing	58 °C	20 sec	
extension	72 °C	20 sec	
final extension	72 °C	30 sec	1

Table SP2.7. Cycling conditions for iPS2-10X-seq second PCR

- 126 **Note:** use a different index from the one utilized to prepare the GEX, CMO, and/or library, to allow subsequent demultiplexing of the UCI-BC library.
- 127 Purify and perform a double-sided size selection using 0.7X 26 0.9X volumes of SPRIselect beads, following the manufacturer's instructions and eluting in 20 μL of elution buffer.

- 128 Run 2 μ L of undiluted PCR on a TapeStation High Sensitivity 5000 Screen Tape assay, following the manual. Expect a main peak at $\sim$ 385 bp (Figure SP3.1B).




Figure SP3.1. Exemplary TapeStation QC of iPS2-10X-seq first (A) and second (B) PCR.

Library pooling and NGS

- 130 Quantify the GEX and, when appropriate, CMO and/or ADT libraries (step 3 of Prepare the gene expression library), and the UCI-BC library (step 11 of Prepare the shRNA barcode library) with a Qubit dsDNA High Sensitivity kit.
- 131 Pool GEX, CMO, ADT and UCI-BC libraries at a relative molar ratio of 60:10:10:3.
- 132 Perform NGS using NextSeq 1000/2000 and 100 cycles reagents, with the following settings: read 1 - 28 cycles; index 1 - 10 cycles; index 2 - 10 cycles each; read 2 - 90 cycles (3e30,000 reads/target cells).

Protocol 4: Data processing and analysis with catcherR

- 133 **iPS2-seq design and analysis with catcherR**
This Supplemental Protocol describes the following procedures:
1. Overview
 2. Installation
 3. Oligonucleotides design
 4. Pooled cloning step 1 QC
 5. Pooled cloning step 2 and hiPSC genome editing QC
 6. iPS2-10X-seq perturbation deconvolution
 7. iPS2-sci-seq perturbation deconvolution
 8. Barcode reassignment
 9. Perturbation effect analysis

Overview



- 134 *catcheR* is a comprehensive bioinformatic package for designing and analyzing iPS2-seq experiments. For the complete documentation, see <https://marialuisaratto.github.io/catcheRdocs/>

It comprises the following functions:




Figure SP4.1. Overview of catcherR

Logical relationship and main inputs/outputs for the catcherR analytical pipeline to assign perturbation in iPS2-seq experiments. Functions used for oligo design, plasmid/hiPSC QC, and perturbation effect analyses are not illustrated.

- 134.1 *catcheR_design*, which designs oligonucleotides for Supplemental Protocol 1 - Design shRNA oligonucleotides, facilitating shRNA library cloning.
- 134.2 *catcheR_step1QC*, which analyzes the results of Supplemental Protocol 1 - Intermediate plasmid pool QC, assessing pooled cloning step 1 plasmids for barcode swaps.
- 134.3 *catcheR_step2QC*, which analyzes the results of Supplemental Protocol 1 - Final plasmid pool QC or hiPSC pool QC, assessing pooled cloning step 2 or genome-edited hiPSC pools for shRNA representation.
- 134.4 *catcheR_scicount*, which analyzes 2-level indexing sci-RNA-seq data, facilitating the generation of gene expression matrix for iPS2-sci-seq experiments.
- 134.5 *catcheR_scicatch*, which assigns shRNA perturbations to single nuclei transcriptomes obtained by Supplemental Protocol 2, enabling the primary analysis of iPS2-sci-seq.
- 134.6 *catcheR_10Xcatch*, which assigns shRNA perturbations to single cell transcriptomes obtained by Supplemental Protocol 3, enabling the primary analysis of iPS2-10X-seq.
- 134.7 *catcheR_scicatchQC* and *catcheR_10XcatchQC*, which use the outputs of *catcheR_scicatch* and *catcheR_10Xcatch*, respectively, to fine-tune shRNA assignment thresholds.
- 134.8 *catcheR_filtercatch*, which leverages on the output of *catcheR_scicatchQC* and *catcheR_10XcatchQC* to filter single nuclei/cell transcriptomes expressing a single shRNA.
- 134.9 *catcheR_sortcatch*, which quality controls the cell-by-gene matrix based on the results of *catcheR_step1QC*, reassigning hPSC clones with barcode swaps to the correct shRNA.
- 134.10 *catcheR_scinoatch* and *catcheR_10Xnoatch*, which identify cells expressing no shRNA in iPS2-sci-seq and iPS2-10X-seq experiments, respectively, adding them to the cell-by-gene matrix to be used as additional controls.
- 134.11 *catcheR_load*, which loads gene expression matrices annotated with shRNA perturbations into a Monocle object, preparing the dataset for downstream analysis.
- 134.12 *catcheR_pseudotime*, which analyzes the effects of shRNA perturbations on pseudotime dynamics, highlighting shifts along differentiation trajectories (e.g., Figures 2J, 2K).
- 134.13 *catcheR_modules*, which assesses perturbation-induced changes in gene module expression, such as coordinated activation or repression of functional programs (e.g., Figures 5H–5J).



- 134.14 *catcheR_enrichment*, which quantifies differences in perturbation representation across experimental samples (e.g., Figures 5F, 5G) or cell clusters (e.g., Figure 5E).

Installation

- 135 *catcheR* is available at <https://github.com/alessandro-bertero/catcheR>. The GitHub repository folder "scripts" contains all the bash and R scripts that can be run independently. However, in order to ensure reproducible analyses, we strongly recommend to install *catcheR* package from GitHub, since its functions run all the analysis inside Docker containers.

- 135.1 Install Docker engine, following the instructions at <https://docs.docker.com/engine/install/>

- 135.2 Install *catcheR*
In R ($\geq 3.0.2$) install *devtools*, if not already present

```
install.packages("devtools")  
Install catcheR from GitHub  
install_github("alessandro-bertero/catcheR")  
Install rundocker from GitHub  
install_github("Reproducible-Bioinformatics/rundocker")
```

- 135.3 Load *catcheR* and *rundocker* in your R environment

```
library(catcheR)  
library(rundocker)
```

Oligonucleotides design

- 136 In a new working folder, prepare the following files:
- A comma-separated values (CSV) file with three columns listing: (1) the forward oligos from the TRC shRNA library; (2) the corresponding barcodes (BC); (3) the shRNA names. Below is an example for the shRNA described in Figure SP1.2:



```
CCGGCAAACCTCTCTGCGATTCTCGGAACATCCAGACAGAGACTGTTTTTG, CAGTTCCA, SMAD2
.1
```

- (Optional) - A txt file with a newline-separated list of 5'-3' restriction sites, or other sequences, to be avoided in the shRNAs. By default these are Sall, SwaI and AclI:

```
GTCGAC
ATTTAAAT
GGCGCGCC
```

137 Run `catcherR_design`:

```
catcherR_design(
  group=c("docker", "sudo"),
  folder,
  sequences,
  gibson.five = "AGTTCCTCTACTAGTGATAGAGATCCC",
  gibson.three = "GTAGCTCGGTGATCAGC",
  fixed = "GTCGACATTTAAATGGCGCGCC",
  restriction.sites = NULL)
```

138 *catcherR_design* arguments:

1. group: string with two options: sudo or docker, depending on the user group. For a detailed explanation of Docker user groups, see this page
2. folder: string with the working folder path
3. sequences: string with the CSV file name from step 1a
4. gibson.five: (optional) string with the 5' Gibson homology
5. gibson.three: (optional) string with the 3' Gibson homology
6. fixed: (optional) character string with the multicloning site
7. restriction.sites: (optional) string with the txt file name from Step 25

138.1 Example usage:

```
catcher_design(  
  group = "docker",  
  folder = "path/to/folder",  
  sequences = "filename.csv",  
  restriction.sites = "filename.txt")
```

Pooled cloning step 1 QC

- 139 This step complements Supplemental Protocol 1 - Intermediate plasmid pool QC. In a new working folder, prepare the following files:

- (a) Fastq/fq or fastq.gz files with demultiplexed read 1 from the NGS run
- (b) A CSV file with the shRNA names and their full sequences

```
SMAD2.1, GCAGATCCTCTCTGCTGGATTCTCGGAACATCCAGACAGAGACTG
```

- (c) (Optional) - A txt file with a newline-separated list of clones of interest as "BC_UCI"

```
CAAGACCC_CATCGT
```

- 140 Run `catcher_step1QC`:

```
catcher_step1QC(  
  group=c("docker", "sudo"),  
  folder,  
  fastq.read1,  
  DIs = 100,  
  ratio = 10,  
  plot.threshold = 2000,  
  clones = NULL)
```

- 140.1 *catcher_step1QC* arguments:

- (a) group: string with two options: sudo or docker, depending on the user group (info)
- (b) folder: string with the working folder path
- (c) fastq.read1: string with the read 1 filename from step 1a


- (d) DIs: integer of the minimum number of diversity indexes (DIs, pseudo-unique reads) of the most represented shRNA matching to a given UCI-BC; in combination with "ratio", it selects UCI-BCs for which it is possible to reliably assign an shRNA
- (e) ratio: integer of the minimum ratio between the number of DIs of the most represented and second most represented shRNAs matching to a given UCI-BC
- (f) plot.threshold: integer of the minimum number of DIs per UCI-BC for output 2b
- (g) clones: (optional) a string with the txt file from step 1c

140.2 Example usage:

```
catcher_step1QC (  
  group = "docker",  
  folder = "path/to/folder",  
  fastq.read1 = " filename .fq",  
  clones = " filename .txt ")
```

140.3 *catcher_step1QC* key outputs:

- (a) "reliable_clones_swaps.csv", which lists UCI-BCs with reliable evidence of shRNA-barcode swap, and is used as input for Barcode reassignment
- (b) Bar chart of the number of DIs for each clone above "plot.threshold"
- (c) Bar chart of the number of DIs associated to each BCs, shRNAs, and reliable swaps
- (d) (If "clones" argument is provided) - CSV files and bar charts with the number of DIs for each shRNA matching to each clone of interest (Figure SP4.2)






Figure SP4.2. Example of quantification of DIs for bacterial clones carrying the expected shRNA (A) or with robust evidence of a barcode swap that can be reassigned (B; Figure 1B). These graphs allow visual evaluation of appropriate thresholds for "DI" and "ratio".

Pooled cloning step 2 and hiPSC genome editing QC

- 141 In a new working folder, prepare the following files:
 - (a) Fastq/fq or fastq.gz files with demultiplexed read 1 from the NGS run
 - (b) "rc_barcodes_genes.csv", a CSV file with two columns: (1) the shRNA BCs; (2) the matching shRNA names in the format "GENE.shRNAID"
 - (c) (Optional) - A txt file with clones of interest (step 1c of Pooled cloning step 1 QC)
- 142 Run `catcher_step2QC`:

```
catcheR_step2QC(  
  group=c("docker","sudo"),  
  folder,  
  fastq.read1,  
  DIs = 1000,  
  clones = NULL)
```

142.1 catcheR_step2QC arguments:


- (a) group: string with two options: sudo or docker, depending on the user group
- (b) folder: string with the working folder path
- (c) fastq.read1: string with the read 1 file from step 1c
- (d) DIs: integer of the minimum number of DIs for a given UCI-BC; it selects reliably-measured UCI-BCs
- (e) clones: (optional) a string with the txt file from step 1c

142.2 Example usage:


```
catcheR_step2QC(  
  group = "docker",  
  folder = "path/to/working/folder",  
  fastq.read1 = "filename.fq",  
  clones = "filename.txt")
```

142.3 catcheR_step2QC key outputs:

- (a) Pie charts of DIs per shRNA and per gene targets associated to each clone
- (b) Text file listing all clones above the DI threshold
- (c) Bar chart of the number of DIs for each clone above the DI threshold (Figure SP4.3A)
- (d) Frequency histogram of the percentage of clones above the DI threshold for each shRNA (Figure SP4.3B)



A



B

Figure SP4.3. Example of quantification of DIs for bacterial clones represented by more than 1,000 DIs (A), and of the frequency of bacterial clones for each shRNA. These graphs confirm whether shRNAs are normally distributed in the final plasmid pool.

Obtain iPS2-10X-seq, iPS2-CITE-seq, or iPS2-multi-seq count matrices

- 143 Download and install [cellranger](#) (for iPS2-10X-seq and iPS2-CITE-seq) or [cellranger-arc](#) (for iPS2-multi-seq). Alternatively, Docker containers for [cellranger v7](#), [cellranger v9](#) (recommended for iPS2-CITE-seq), or [cellranger-arc](#) can be pulled from our Docker repository.
 - 143.1 For iPS2-10X-seq and iPS2-CITE-seq, demultiplex Illumina base call (BCL) files using cellranger mkfastq, following the standard 10X Genomics user guide. In the sample sheet CSV, include the index sequences used in *Supplemental Protocol 3* for the GEX and UCI-BC libraries, and, when applicable, CMO and/or ADT libraries.
 - 143.2 For iPS2-multi-seq, use cellranger-arc mkfastq to demultiplex GEX + ATAC dual-index libraries, following the Cell Ranger ARC user guide. Ensure the sample sheet format is compatible with dual-modality sequencing runs, and include the index sequences used for both gene expression and chromatin accessibility libraries.
- 144 Run FastQC to assess the quality of the FASTQ files for each library type.

- 145 Generate cell-by-gene count matrices using cellranger count (for single-sample experiments) or cellranger multi (for multiplexed experiments, including iPS2-CITE-seq). For iPS2-multi-seq, use cellranger-arc count to obtain both gene expression and chromatin accessibility matrices.
- 145.1 Note: In multiplexed experiments (e.g., using CMO or ADT barcodes in iPS2-CITE-seq), count matrices from individual samples can be aggregated using cellranger aggr to produce a unified dataset. This enables joint analysis in a single run of catcheR_10Xcatch by specifying the number of samples via the samples argument.
- 146 Use cellranger mat2csv to convert the sparse matrix output into a dense CSV format for compatibility with downstream tools. For iPS2-multi-seq, run cellranger-arc mat2csv separately for GEX and ATAC outputs, if needed.
After obtaining the count matrix, proceed to *iPS2-10X-seq perturbation deconvolution*.

Obtain iPS2-sci-seq count matrices

- 147 *catcheR_scicount* is a wrapper of the "bbi-sci" pipeline developed by the Brotman Baty Institute for Precision Medicine, which was implemented and dockerized in catcheR to be used on any operating system.

Demultiplex Illumina base calls to fastq files

- (a) Create a "SampleSheet.csv" file with as many sample rows as the PCR wells from Indexing PCR and NGS of Supplemental Protocol 2, where "Sample_ID" is the well identified with format [A-H]01-[12], and index and index2 are the i7 and i5 indexes used in the corresponding row and column (refer to Table SP2.1)
- (b) Run Illumina bcl2fastq following Illumina's manual
- (c) Run fastQC to confirm the quality of the fastq files

- 148 In a new working folder:
- (a) Create the subfolder "fastq", and copy all the demultiplexed fastq.gz files. Ensure all file names begin with the well coordinate (e.g. A01) followed by an underscore
 - (b) Create a tab separated txt file called "sci-RNA-seq-8.RT.oligos" with the association between RT wells and RT barcode sequences (refer to *Table SP2.1*)
 - (c) Create the subfolder "GENOMES", and copy the annotated genome (i.e., **GRCh38**)

- 149 Run catcheR_scicount:

```
catcheR_scicount(  
  group=c("docker", "sudo"),  
  folder,  
  sample.name,  
  UMI.cutoff)
```

149.1 `catcheR_scicount` arguments:

- (a) group: string with two options: sudo or docker, depending on the user group ([info](#))
- (b) folder: string with the working folder path
- (c) sample.name: string with the name of the experiment
- (d) UMI.cutoff: integer of the minimum number of UMI per nucleus needed to consider the single cell transcriptome valid

149.2 Example usage:

```
catcheR_scicount(  
  group = "docker",  
  folder = "path/to/file",  
  sample.name = "experiment",  
  UMI.cutoff = 500)
```

149.3 `catcheR_scicount` outputs, found in the final-output folder: knee-plot of UMI per cell (Figure S1H), statistics files, sparse cell-by-gene count matrix, dense gene matrix "exp_mat.csv" and "exp_mat_no0.csv" (filtered to exclude genes with 0 counts in all cells), also in Rdata format.

After obtaining the count matrix, proceed to *iPS2-sci-seq perturbation deconvolution*.

iPS2-10X-seq perturbation deconvolution

150 shRNA perturbations can be assigned to single cells using *catcheR_10Xcatch*, which identifies NGS reads containing UCI-BCs, matches them to the corresponding transcriptome based on their shared cellular barcodes, applies filters for background noise, and selects cells with robust evidence of a single integration. UCI-BC filtering involves: (1) filtering out UCI-BCs supported by less than a certain number of UMIs in a given transcriptome, to account for PCR and sequencing artifacts; (2) filtering out UCI-BCs that represent less than a certain fraction of all UMIs for UCI-BCs in a given transcriptome, to reduce noise arising from free mRNA; and (3) filtering out UCI-BCs whose UMI fraction in a given transcriptome is not several fold greater than the second most common UCI-BC, to eliminate cells that contain multiple UCI-BCs when the second most common one falls just below one of the first two thresholds. The resulting list of *bona fide* shRNA integrations per cell is then used to filter those expressing a single shRNA. This section outlines a typical analysis workflow for iPS2-10X-seq, iPS2-CITE-seq or iPS2-multi-seq experiments, which follow analogous preprocessing and count matrix structures based on 10X Genomics technology. The pipeline begins with *catcheR_10Xcatch*, which performs a full analysis using automatic thresholding. If needed, these thresholds can be subsequently refined

with *catcherR_10XcatchQC* and applied for re-filtering using *catcherR_filtercatch*. Finally, *catcherR_10Xnocatch* allows incorporation of cells that do not express any shRNA as additional negative controls in the final annotated cell-by-gene matrix.

In a new working folder:

- 150.1 Copy the fastq/fq or fastq.gz files with demultiplexed read 1 and read 2 of the UCI-BC library (step 3 of Obtain iPS2-10X-seq, iPS2-CITE-seq, or iPS2-multi-seq count matrices)
- 150.2 Copy the cell-by-gene count matrix in CSV format (step 6 of Obtain iPS2-10X-seq, iPS2-CITE-seq, or iPS2-multi-seq count matrices)
- 150.3 Create a CSV file named "rc_barcodes_genes.csv" with two columns: (1) the shRNA BCs; (2) the matching shRNA names in the format "GENE.shRNAID"
- 151 Run *catcherR_10Xcatch* to execute the complete analysis with automatic thresholding:



```
catcherR_10Xcatch(  
  group=c("docker","sudo"),  
  folder,  
  fastq.read1, fastq.read2,  
  expression.matrix,  
  reference = "GGGCCGCTCATCTGGGGAGCCG",  
  UCI.length = 6,  
  threads = 2,  
  percentage = 15,  
  mode = "bimodal",  
  ratio = 5,  
  samples = 1,  
  x = 100, y = 400)
```

- 151.1 Example usage:

```
catcherR_10Xcatch(group = "docker",  
  folder = "path/to/folder",  
  fastq.read1 = "R1.fastq.gz", fastq.read2 = "R2.fastq.gz",  
  expression.matrix = "filename.csv",  
  threads = 12)
```

151.2 catcherR_10Xcatch arguments:

- (a) group: string with two options: sudo or docker, depending on the user group ([info](#))
- (b) folder: string with the working folder path
- (c) fastq.read1: string with the filename of read 1 fastq/fq or fastq.gz (step 1a)
- (d) fastq.read2: string with the filename of read 2 fastq/fq or fastq.gz (step 1a)
- (e) expression.matrix: string with the filename of the count matrix file CSV (step 1b)
- (f) reference: string with the reverse complement of the sequence before the shRNA-BC at the start of read 2 (default is the 3' end of the OPTtetR cDNA, optional argument; Figure S2F)
- (g) UCI.length: integer of the UCI length (default is 6, optional argument)
- (h) threads: integer of the threads to be used for parallelization (default is 2, optional argument)
- (i) percentage: integer of the percentage of UMIs supporting a given UCI over the total UMIs supporting all UCIs in a given cell; it is used as threshold to consider the UCI valid (Figure SP4.4D); the recommended default is 15 (optional argument)
- (j) mode: string with two options: bimodal or noise, defining the mode for automatic thresholding the minimum number of UMIs per UCI to consider the UCI valid (optional argument)
 1. "bimodal" (default) sets the threshold at the valley of the bimodal UMIXUCI distribution (Figure SP4.4C)
 2. "noise" sets the threshold at $1.35 \times$ the number of UCIs supported by a single UMI; used if the UMIXUCI distribution is not bimodal
- (k) ratio: the ratio between the number of UMIs supporting the most represented UCI and the number of UMIs supporting the second most represented UCI. The parameter is needed to identify the cell as a single integration cell. The default is 5 (optional argument)
- (l) samples: the number of samples present in the same experiment (cells from different experiments will have the corresponding number after the cell ID in the gene expression matrix, so that multiple reactions of the same experiment can be unified). The default is 1 (optional argument)
- (m) x: an integer indicating the upper limit on the x-axis of the cropped version of the plot "UMIXUCI". Default is 100 (optional argument)
- (n) y: an integer indicating the upper limit on the y axis of the cropped version of the plot "UMIXUCI". Default is 400 (optional argument)

151.3 catcherR_10Xcatch key outputs (will be stored inside the "Result" folder of each sample):

- (a) "log.txt" with the starting number of reads
- (b) "log2.txt" with the number of cells, UMIs, UCIs, and the calculated thresholds of UMI per UCI for both "bimodal" and "noise" (only one is chosen based on the "mode" argument; the other one is provided for a potential reanalysis)
- (c) Bar charts of UMI counts per shRNA and target gene (Figure SP4.4A and Figure SP4.4B)

- (d) Frequency histograms of UCIs supported by a certain number of UMIs (UMIxUCI), to interpret the signal/noise of the experiment and possibly set a custom UMIxUCI threshold for subsequent reanalysis (Figure SP4.4C) **Note:** here and below, copies of the same UCI expressed by different cells are plotted and analyzed separately
- (e) Frequency histogram of UCIs supported by a certain fraction of UMIs over the total number of UMIs supporting all UCIs in a given cell (UMIpercentagexUCI), to further assess signal/noise and possibly adjust the default threshold (Figure SP4.4D)
- (f) Dot plots that combine the UMIpercentagexUCI and UMIxUCI data on the x and y axes, respectively, with each dot representing one or more UCI (quantified by the dot size). Dot colors indicate either the number of valid UCIs (i.e., shRNAs) in the cell containing a given UCI (Figure SP4.4E), or whether said UCI is the only valid one in the relevant cell and is thus assigned to such cell (Figure SP4.4F)
- (g) "log_part3.txt" with how many cells were assigned to a single shRNA or were filtered out due to zero or multiple shRNAs
- (h) "silencing_matrix.csv" is the key output used for the secondary analyses: the cell-by-gene count matrix provided as input, filtered and annotated with the shRNA encoded in each cell. It is also provided in RDS format for easy loading into R. Cell names are modified as follows:

```
cellID_UMIxUCI_BC_GENE_UCI
```


Where:

1. cellID is the original cell name (i.e., the 16 bp of the 10X RT barcode)
2. UMIxUCI is the number of UMI associated with that perturbation
3. BC is the shRNA barcode
4. GENE is the shRNA target
5. UCI is the Unique Clonal Identifier (shared by all cells originating from the same hPSC clone)


Example of a cell name after the analysis:

```
TTCTAACCACAGTCGC_180_CGTGATGC_NKX2 .5 _ACAGTG
```


The annotated matrix can be directly used with the functions implemented in *catcher* (see Perturbation effect analysis) to perform a range of secondary analyses described in the manuscript. Additional scripts are available at our GitHub repository for customized or extended analyses, most of which are now incorporated into *catcher*.




A




B



C



D



E



F

Figure SP4.4. Exemplary key outputs of a *catcherR_10Xcatch* analysis: "barcode_distribution.pdf" (A); "gene_distribution.pdf" (B); "UMIxUCI_400_100.pdf" with the automatic UMIxUCI threshold based on "bimodal" (C); "percentage_of_UMIxUCI_dist.pdf" with the default 15% "percentage" threshold (D); "2D_percentage_of_UMIxUCI_UMI_count_trueorfalse.pdf" (E); and "2D_percentage_of_UMIxUCI_UMI_ValidCells.pdf" (F).

OPTIONAL: fine tune iPS2-10X-seq perturbation assignment

- 152 In the same working folder used to run *catcheR_10Xcatch*, run *catcheR_10XcatchQC*, which uses the txt file outputs to regenerate the quality control plots described in **Figure SP4.4** and suggest new thresholds to enable subsequent filtering.

```
catcheR_10XcatchQC(  
  group=c("docker", "sudo"),  
  folder,  
  reference = "GCGCGCTTCATCTCGGGGAGCCG",  
  mode = "bimodal",  
  sample = 1  
  x = 100, y = 400)
```

- 152.1 Example usage:

```
catcheR_10XcatchQC(  
  group = "docker",  
  folder = "path/to/working/folder",  
  mode = "noise")
```

The function is to be run separately for each sample (specified by argument *sample* with default 1). *catcheR_10XcatchQC* arguments and outputs are the same as *catcheR_10Xcatch*.

- 153 In the same working folder used to run *catcheR_10Xcatch*, run *catcheR_filtercatch*, which uses the output of *catcheR_10Xcatch* to filter cells based on custom thresholds.

```
catcheR_filtercatch(  
  group=c("docker", "sudo"),  
  folder,  
  expression.matrix,  
  UMI.count,  
  percentage = 15,  
  ratio = 5,  
  sample = 1)
```

153.1 Example usage:

```
catcheR_filtercatch(group = "docker",  
  folder = "path/to/working/folder",  
  expression.matrix = "filename.csv",  
  UMI.count = 5, sample = 1)
```

The function is to be run separately for each sample (specified by argument sample with default 1). *catcheR_filtercatch* argument percentage is the minimum percentage of UMI for a given UCI over the total UMIs of UCIs in that cell, to consider the UCI valid.

UMI.count is an integer of the custom UMI threshold. ratio is the minimum ratio between UMIs of top UCI and UMIs of 2nd top UCI. The other arguments and the outputs are the same as *catcheR_10Xcatch*.

OPTIONAL: identify cells expressing no shRNA in iPS2-10X-seq

154 In the same working folder used to run *catcheR_10Xcatch*, run *catcheR_10Xnocatch*:

```
catcheR_10Xnocatch(  
  group=c("docker", "sudo"),  
  folder,  
  expression.matrix,  
  threshold,  
  sample = 1,  
  reference = "TAGCCGTTCACCTCGGGGACGCG")
```

154.1 Example usage:

```
catcheR_10Xnocatch(group = "docker",  
  folder = "path/to/folder",  
  expression.matrix = "filename.csv",  
  threshold = 5)
```

154.2 *catcheR_10Xnocatch* argument threshold is an integer of the minimal UMI count for the "empty" reference (Figure S2F) needed to confidently identify it as having integrated said plasmid. All other arguments are the same as *catcheR_10Xcatch*.



catcheR_10Xnocatch output is an updated "silencing_matrix.csv", also available in RDS format, in which names of empty cells are modified as "cellID__?_empty_NA_empty".

OPTIONAL: merge multiple samples

- 155 In case fine-tuning is applied after catcheR_10Xcatch and multiple samples are present, run catcheR_merge to aggregate the results in a single matrix. This step is done automatically when running catcheR_10Xcatch.

```
catcheR_merge(  
  group=c("docker", "sudo"),  
  folder,  
  samples = 2,  
  empty = T)
```

- 155.1 The argument "empty" determines whether the cells identified by catcheR_nocatch should be added to the matrix.

- 155.2 Example usage:

```
catcheR_merge(group = "docker",  
  folder = "/20tb/ratto/catcheR/test_CM5/",  
  samples = 4)
```

iPS2-sci-seq perturbation deconvolution

- 156 This section describes a typical analysis pipeline similar to the one described in the previous section, but leveraging on the *catcheRscicatch*, *catcheRscicatchQC*, *catcheRfiltercatch*, and *catcheRscinoactch* functions.
In a new working folder:

Create a subfolder called "fastq", and copy all the demultiplexed fastq.gz files (step 2 of Obtain iPS2-sci-seq count matrices). Ensure that all file names begin with the well coordinate (e.g. A01)

- 157 Copy the cell-by-gene expression matrix CSV file obtained with catcheR_scicount function (step 4 of Obtain iPS2-sci-seq count matrices)

- 157.1 Create a CSV file called "rcbarcodesgenes.csv" with two columns: (1) the shRNA BCs; (2) the matching shRNA names in the format "GENE.shRNAID"
- 157.2 Copy the txt file "sci-RNA-seq-8.RT.oligos" used by *catcheR_scicount* (step 3b of Obtain iPS2-sci-seq count matrices)
- 158 Run *catcheR_scicatch* to execute the complete analysis with automatic thresholding; the arguments are the same as for *catcheR_10Xcatch* (step 2 of iPS2-10X-seq perturbation deconvolution), except that no filenames are provided since fastq files must be in the "fastq" subfolder in the working directory:

```
catcheR_scicatch(  
  group=c("docker", "sudo"),  
  folder,  
  expression.matrix,  
  reference = "GGCGGCTTCATCTGGGGGACGCG",  
  UCI.length = 6,  
  threads = 2,  
  percentage = 15,  
  ratio = 5,  
  mode = "bimodal",  
  x = 100, y = 400)
```

- 158.1 Example usage:

```
catcheR_scicatch(group = "docker",  
  folder = "path/to/working/folder",  
  expression.matrix = "filename.csv",  
  threads = 12)
```

- 158.2 *catcheR_scicatch* key outputs are the same as for *catcheR_10Xcatch* (step 2 of iPS2 10X-seq perturbation deconvolution and Figure SP4.4), except that in "silencing_matrix.csv" "cell ID" indicates the PCR well and RT barcode ID (which can be leveraged to identify nuclei from different samples pooled on the same RT plate) Example of cell name after the analysis:

P24_RT_27_7_GCCTGTGT_SCR_ACGGTG

- 158.3 In addition, *catcheR_scicatch* provides two quality control outputs to evaluate experimental biases during sci-RNA-seq library preparation: demux and RT detail about how many cells were identified from each row and column of the PCR and RT plates, respectively

OPTIONAL: fine tune iPS2-sci-seq perturbation assignment

- 159 Run *catcher_scicatch* to execute the complete analysis with automatic thresholding; the arguments are the same as for *catcher_10Xcatch* (step 32.2), except that no filenames are provided since fastq files must be in the "fastq" subfolder in the working directory:

```
catcher_scicatch(  
  group=c("docker", "sudo"),  
  folder,  
  expression.matrix,  
  reference = "GGCGCTTCATCTGGGGGACGCG",  
  UCI.length = 6,  
  threads = 2,  
  percentage = 15,  
  ratio = 5,  
  mode = "bimodal",  
  x = 100, y = 400)
```

159.1 Example usage:

```
catcher_scicatch(group = "docker",  
  folder = "path/to/working/folder",  
  expression.matrix = "filename.csv",  
  threads = 12)
```

- 159.2 *catcher_scicatch* key outputs are the same as for *catcher_10Xcatch* (step 2 of iPS2-10X-seq perturbation deconvolution and Figure SP4.4), except that in "silencing_matrix.csv", "cell ID" indicates the PCR well and RT barcode ID (which can be leveraged to identify nuclei from different samples pooled on the same RT plate).

Example of cell name after the analysis:

```
P24__RT_27_7_GCCTGTGT_SCR_ACCGTC
```

In addition, *catcher_scicatch* provides two quality control outputs to evaluate experimental biases during sci-RNA-seq library preparation: demux and RT detail about how many cells were identified from each row and column of the PCR and RT plates, respectively.

160 Run `catcheR_scicatchQC` following the steps described for `catcheR_10XcatchQC` in *OPTIONAL: fine tune iPS2-10X-seq perturbation assignment*: the functions share the same arguments and outputs, but `catcheR_scicatchQC` used the output of `catcheR_scicatch`.

160.1 Example usage:

```
catcheR_scicatchQC(  
  group = "docker",  
  folder = "path/to/working/folder",  
  mode = "noise")
```

161 Run `catcheR_scinocatch` following the steps described for `catcheR_10Xnocatch` in *OPTIONAL: identify cells expressing no shRNA in iPS2-10X-seq*: the functions share the same arguments and outputs, but `catcheR_scinocatch` uses the output of `catcheR_scicatch`.

161.1 Example usage:

```
catcheR_scinocatch(  
  group = "docker",  
  folder = "path/to/working/folder",  
  expression.matrix = "filename.csv",  
  threshold = 5,  
  reference = "TACGCGTTCATCTGGGGGAG")
```

Barcode reassignment

162 `catcheR_sortcatch` is an optional function that corrects the annotated cell-by-gene count matrix obtained with `catcheR_10Xcatch` or `catcheR_scicatch` reassigning the perturbation of any cell belonging to a HPSC clone with reliable evidence of a shRNA-barcode swap, based on the results of `catcheR_step1QC`.

In a new working folder:

- (a) Copy the annotated count matrix (i.e., "silencing_matrix.csv")
- (b) Copy the CSV file with the list of UCI-BCs with reliable evidence of a shRNA-barcode swap (i.e., "reliable_clones_swaps.csv"; output 2a of Pooled cloning step 1 QC)

163 Run `catcheR_sortcatch`:

```
catcheR_sortcatch(  
  group=c("sudo", "docker"),  
  folder,  
  expression.matrix,  
  swaps)
```

163.1 *catcheR_sortcatch* arguments:

- (a) group: string with two options: sudo or docker, depending on the user group
- (b) folder: string with the working folder path
- (c) expression.matrix: string with the filename of the annotated count matrix CSV
- (d) swaps: a character string with the filename of the txt file (step 1b)

163.2 *catcheR_sortcatch* output is an updated annotated gene expression matrix CSV file called "silencing_matrix_updated.csv", in which "BC" and "GENE" have been modified to reflect the actual shRNA encoded in each cell.

Perturbation effect analysis

164 The second part of *catcheR* provides an exploratory analysis with visual and statistical summaries to highlight perturbation effects.

Annotation

Before proceeding, the gene expression matrix should be annotated with gene symbols using the *scannobyGtf* function from the R package rCASC. As part of quality control, we recommend evaluating the fraction of ribosomal and mitochondrial reads — e.g., using the *mitoRiboUmi* function from the same package — and considering the exclusion of cells with abnormally high proportions, which may indicate poor quality or stress.

Note: After this step, row names of the matrix (the genes) will have the following format:

```
GeneSymbol:EnsemblID
```

Example:

```
ENSG000000000003:TSPAN6
```

Data loading

165 In a new working folder:

- Copy the count matrix annotated with gene names during the previous step (*filtered_annotated_silencing_matrix_complete_all_samples.csv*).
- Copy the file *rc_barcodes_genes.csv* described in step 1a of Pooled cloning step 2 and hiPSC genome editing QC.
- Create a new-line separated file listing the control genes (e.g. SCR, B2M).
- Create a new-line separated file listing the control samples, if any (e.g. 1,3). These are the sample names also used by *aggr* (see CSV file used as input for *aggr*).
- Create a plain text file listing the sample replicate labels, one per line (e.g., batch1, batch1, batch2, batch2). The order of entries must match the order of samples in the input matrix exactly. This file is required for downstream batch-aware analyses.
- If the dataset includes samples from different experiments or processing batches, batch correction is recommended.
- Create a CSV file listing each sample along with its corresponding annotation name, which will be used for display instead of the sample number (this file is mandatory). Example of CSV file to download on GitHub.
- Create a newline-separated file listing the genes of interest whose expression you wish to visualize on the UMAP (optional).

166 Run *catcherR_load*:

```
catcherR_load(  
  group="docker",  
  folder,  
  expression.matrix,  
  control_genes,  
  control_samples = NULL,  
  replicates = NULL,  
  sample_names,  
  resolution = 8e-4,  
  genes = NULL)
```

166.1 Example usage:

```
catcher_load(  
  group="docker",  
  folder="/path/to/working/folder/",  
  expression.matrix =  
    "annotated_silencing_matrix_complete_all_samples.csv",  
  control_genes = "controls.txt",  
  control_samples = "noTET.txt",  
  replicates = "replicates.txt",  
  sample_names = "samples.csv",  
  resolution = 8e-4,  
  genes = "genelist.txt")
```

166.2 The argument `_resolution_` refers to the resolution parameter used by Monocle's `_cluster_cells_` function for clustering.

166.3 `catcher_load` outputs:

- 1. `expression_data.csv` and `cell_metadata.csv`, which can be used to create a Monocle Cell Data Set (CDS) and are also included in the ready-to-load R object `starting_cds.Rdata`
- 2. `UMAP.pdf` plots the dimensionality reduction and `UMAPgeneexpression.pdf` shows the gene expression on the UMAP of the genes provided by the argument "genes"
- 3. `UMAP_clustering.pdf` shows the clustering obtained on the UMAP with the provided resolution
- 4. `processed_cds.RData` is the Cell Data Set after normalization, dimensionality reduction, clustering and calculation of trajectories

166.4 At the end of this step, the data are stored in a format compatible with the `monocle3` package. However, you can switch to other single-cell analysis frameworks such as `Seurat` or `Scanpy` at any point in the workflow. *Note:* For standard iPS2-seq perturbation analysis, always continue using the CDS object generated by *catcher*.

```
library(SeuratWrappers)  
library(Seurat)  
seurat = as.Seurat(cds, assay = NULL)  
scanpy_sce = as.SingleCellExperiment(seurat)
```

166.5 The follow up analysis are: `catcher_pseudotime`, `catcher_modules` and `catcher_enrichment`.

- With the *Pseudotime* function, *catcheR* calculates the cumulative frequency of cells sharing the same gene, shRNA, or clone, compared to each control, along the pseudotime trajectory. *catcheR* can compare the cumulative frequency curves by computing directional Kolmogorov-Smirnov statistics between the target and control
- With the *Genes modules* function, *catcheR* identifies gene modules within the dataset and computes the cumulative frequency of cells with the same gene, shRNA, or clone—relative to each control—based on the expression of each gene module. Also in this case, cumulative frequency curves are compared using a directional Kolmogorov-Smirnov test
- With the *Enrichment / depletion analysis* function, *catcheR* calculates the enrichment of cells with the same gene, shRNA, or clone within clusters, in comparison to each control. *catcheR* can compare each target and control quantity and obtain statistics using the Fisher exact test

Pseudotime

167 Since Monocle pseudotime calculation requires the user to select the starting point interactively, the first step of the pseudotime analysis needs to be done separately within *monocle3*.

After loading the *processed_cds.RData* in R, calculate the pseudotime trajectory and plot it with the R script below:

```
library(monocle3)
cds = order_cells(cds)
pt = as.data.frame(pseudotime(cds))
names(pt) = c("pseudotime")
write.csv(pt, paste0(dir, "/pseudotime.csv"))

plot_cells(cds, color_cells_by = "pseudotime",
label_cellgroups=FALSE, label_leaves=FALSE,
label_branch_points=FALSE, graph_label_size=1.5)
```

168 Run *catcheR_pseudotime*:

```
catcheR_pseudotime( group=c("docker","sudo"), folder, cds,
pseudotime, all = FALSE)
```

168.1 Example usage:

```
catcheR_pseudotime(  
  group="docker",  
  folder="/path/to/working/folder/",  
  cds = "processed_cds.RData",  
  pseudotime = "pseudotime.csv")
```

168.2 Argument "all" is a logical operator indicating whether to perform the Kolmogorov-Smirnov test against the controls together (all = F) or against each control separately (all = T). The default argument is false

168.3 Below is a list of the catcheR_pseudotime outputs:

- (a) The "cumulativefrequency_pseudotime" plots show the number of cells at each given point of the pseudotime for different groups of cells at different comparison levels - gene, shRNA, and clone
- (b) The "ks_statistics" CSV files that include the Kolmogorov-Smirnov test results comparing cumulative the frequency based on the pseudotime between knockdown and different controls
- (c) The Volcano plots of the results of the Kolmogorov-Smirnov test show the significance and fold change between the pseudotime cumulative curves. Different controls and levels are used
- (d) The "correlated_pseudotime_gene_exp.pdf" showing the expression on the UMAP of the genes most correlated with the pseudotime

Gene Modules

169 This function uses Monocle to find gene modules which expression varies in different perturbation groups, i.e., perturbed gene or shRNA or clones.

Run catcheR_modules:

```
catcheR_modules(group=c("docker", "sudo"), folder, cds,  
  resolution=1e-2)
```

169.1 Example usage:

```
catcher_modules(group="docker",
  folder="/30tb/3tb/data/ratto/testing/", cds =
  "processed_cds.RData")
```

- 169.2 The resolution parameter refers to the value used in Monocle's `find_gene_modules` function, which influences the number of gene modules identified and the number of genes included in each module (higher resolution typically results in more, smaller modules)
- 169.3 Below are listed the outputs of `catcher_modules`:
- (a) `gene_modules.csv`, listing the genes present in each module
 - (b) Heatmap plots display the Z-scores of each module, either all modules or the top 10 most variable ones, across perturbation groups
 - (c) The `"modulescells"` folder contains tables with aggregated module expression values for each cell. These can be used as input for `catcherPseudotime` in place of the pseudotime CSV file, allowing the analysis of cumulative frequency based on module expression
- 169.4 Run `catcher_pseudotime` on the CSV files stored in the `"modules_cells"` folder to use the data generated by `catcher_modules` as described in the previous example in the Pseudotime section

Enrichment / depletion

- 170 This function evaluates whether perturbation groups are enriched or depleted in cells (number) or in specific cell subpopulations (clusters).

Run `catcher_enrichment`:

```
catcher_enrichment(
  group=c("docker", "sudo"),
  folder,
  file,
  meta,
  timepoint = "PSC",
  control_gene = "SCR",
  min_cells_cluster = 70,
  min_cells_shRNA = 40)
```

- 170.1 Example usage:

```
catcheR_enrichment(  
  group = "docker",  
  folder = "/3tb/data/ratto/aggr/test/",  
  file = "processedcds.RData",  
  meta = "cell_metadata.csv",  
  timepoint = "PSC",  
  control_gene = "SCR")
```

170.2 The required input is the cds file generated from `catcheR_load`. Timepoint refers to the baseline time point used as a control for statistical analysis. It is required when experiments include multiple time points and aim to assess enrichment or depletion across these time points. `control_gene` specifies the control gene used as a reference for statistical comparisons.

171 Below are the outputs of `catcheR_enrichment`:

- (a) Plots of cells in each perturbation group
- (b) Volcano plot showing enrichment or depletion of cell numbers in perturbation groups compared to the control, based on fold-change ($\log_2$ ratio) and statistical significance. A corresponding bar plot displays the $\log_2$ fold-change for each perturbation
- (c) Barplots showing the distribution of cells from different perturbation groups in the clusters
- (d) Volcano plot showing the results of Fisher's exact test, comparing the distribution of cells from perturbation groups across Monocle-derived clusters. The plot displays $-\log_{10}$ adjusted p-values *versus* $\log_2$ fold-changes relative to the control group
- (e) Table with Fisher's statistics