



Jan 18, 2024

# 🌐 Single-Molecule Antibody Slides For Fluorescence Microscopy

DOI

[dx.doi.org/10.17504/protocols.io.ewov1qbbygr2/v1](https://dx.doi.org/10.17504/protocols.io.ewov1qbbygr2/v1)

Rebecca Andrews<sup>1</sup>

<sup>1</sup>Univeristy of Cambridge



Rebecca Andrews

Univeristy of Cambridge

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov1qbbygr2/v1>

**Protocol Citation:** Rebecca Andrews 2024. Single-Molecule Antibody Slides For Fluorescence Microscopy . **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.ewov1qbbygr2/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** January 18, 2024



**Last Modified:** May 31, 2024

**Protocol Integer ID:** 93745

**Keywords:** ASAPCRN, molecule antibody slides for fluorescence microscopy, molecule antibody slides for fluorescence microscopy, molecule antibody slide, molecule antibody slide, fluorescence microscopy, antibody

## Abstract


This protocol describes how to create single-molecule antibody slides for fluorescence microscopy.

## Troubleshooting



## Coverslip Cleaning

1h 0m 10s



- 1 Coverslips (24 × 50 mm, #1, VWR, Catalogue Number 48404-453) were argon plasma cleaned (Ar plasma cleaner, PDC-002, Harrick Plasma) for  00:30:00 .

30m

## Surface preparation

1h 0m 10s

- 2 A trimmed gasket was placed on top of the slide (CultureWell™ Reusable Gasket, 6mm diameter, Grace Bio-Labs, SKU: 103280).

- 3 Poly-L-Lysine  30 µL (0.01 % w/v PLL, P4707, Sigma-Aldrich) was placed in the wells for  00:30:00 .

30m

- 4 The PLL was removed, the wells washed three times with PBS (pH 7.4, 1x Gibco, Thermo Fisher Scientific, Catalogue number 10010023)


## Depositing Antibodies

1h 0m 10s

- 5 A secondary antibody of choice was added at a final concentration of 0.0002 mg/ml in PBS.

- 6 The antibodies were left in the wells for 5-10 seconds for sufficient surface density before the wells were washed three times with PBS.

10s

- 7 PBS  30 µL was left in the wells for imaging