



Apr 10, 2020

## Single Nucleus Drop-seq (snDrop-seq)

DOI

[dx.doi.org/10.17504/protocols.io.zmvf466](https://dx.doi.org/10.17504/protocols.io.zmvf466)

Song Chen<sup>1</sup>, Blue B. Lake<sup>1</sup>, Sarah Urata<sup>1</sup>, Kun Zhang<sup>1</sup>

<sup>1</sup>University of California, San Diego

Human Cell Atlas Metho...

KPMP



Zhang Lab

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.zmvf466>

External link: <http://genome-tech.ucsd.edu/ZhangLab/>

**Protocol Citation:** Song Chen, Blue B. Lake, Sarah Urata, Kun Zhang 2020. Single Nucleus Drop-seq (snDrop-seq). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.zmvf466>

**Manuscript citation:**

Lake, B.B., Chen, S., Hoshi, M. et al. A single-nucleus RNA-sequencing pipeline to decipher the molecular anatomy and pathophysiology of human kidneys. Nat Commun 10, 2832 (2019). <https://doi.org/10.1038/s41467-019-10861-2>



**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** March 29, 2019

**Last Modified:** April 10, 2020

**Protocol Integer ID:** 21909

**Keywords:** sequencing, single nuclei, DropSeq, seq protocol, seq, original drop, sndrop, mccarroll lab in the department, drop, protocol, mccarroll lab

## Abstract

The protocol presented here is a Drop-Seq protocol modified for single nuclei.

The original Drop-Seq protocol comes from the McCarroll Lab in the Department of Genetics, Harvard Medical School.

<http://mccarrolllab.org/download/905/>

## Materials

### MATERIALS

⊗ DTT Merck MilliporeSigma (Sigma-Aldrich) Catalog #D0632

⊗ 15 mL Falcon tubes

⊗ PCR tubes

⊗ SSC Thermo Fisher Scientific Catalog #AM9765

⊗ 10x PBS Thermo Fisher Scientific Catalog #AM9624

⊗ UltraPure Distilled Water Invitrogen - Thermo Fisher Catalog #10977-015

⊗ TrypLE™ Express Enzyme Thermo Fisher Scientific Catalog #12604013

⊗ TSO: AAGCAGTGGTATCAACGCAGAGTGAATrGrGrG Integrated DNA Technologies, Inc. (IDT)

⊗ SMART PCR primer: AAGCAGTGGTATCAACGCAGAGT Integrated DNA Technologies, Inc. (IDT)

⊗ New-P5-SMART PCR hybrid oligo: AATGATACGGCGACCACCGAGATCTACACGCCT  
GTCCGCGGAAGCAGTGGTATCAACGCAGAGT\* A\*C Integrated DNA Technologies, Inc. (IDT)

⊗ Custom Read 1 primer: GCCTGTCCGCGGAAGCAGTGGTATCAACGCAG AGTAC Integrated DNA Technologies, Inc. (IDT)

⊗ BSA Merck MilliporeSigma (Sigma-Aldrich) Catalog ##A8806

⊗ Barcoded Bead SeqB: 5' -Bead-Linker-TTTTTTTAAGCAGTGGTATCAAC  
GCAGAGTACJJJJJJJJJJJJNNNNNNNN TTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTT-3' Chemgenes

⊗ Inverted Microscope Catalog #Motic AE31

⊗ Syringe Pumps (3) Catalog #Legato 100

⊗ Magnetic Stirrer VP Scientific Catalog #710D2

⊗ Magnetic Mixing Discs VP Scientific Catalog #772DP-N42-5-2

⊗ 3 mL syringes Becton Dickinson (BD) Catalog #BD #309657

⊗ 10 mL syringes Becton Dickinson (BD) Catalog #BD 309695

⊗ Tubing Scientific Commodities Catalog #BB31695PE/2

⊗ Luer lock 26-gage needles Becton Dickinson (BD) Catalog #305111

⊗ PDMS co-flow microfluidic droplet generation device Catalog #CAD file from McCarroll Lab

⊗ 100 micron cell strainers (for beads) VWR International (Avantor) Catalog #21008-950

⊗ 40 micron cell strainers (for cells) VWR International (Avantor) Catalog #21008-949

⊗ Fuchs-Rosenthal hemocytometer INCYTO Catalog #DHC-F01

⊗ Ficoll PM-400 20% in H2O Merck MilliporeSigma (Sigma-Aldrich) Catalog #F5415-50ML

⊗ Nuclease-free H2O Life Technologies Catalog #AM9930



- ✕ Sarkosyl **Merck MilliporeSigma (Sigma-Aldrich) Catalog #L7414**
- ✕ EDTA 0.5 M **Life Technologies Catalog #AM9260G**
- ✕ Tris 2M pH 7.5 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T2944**
- ✕ Perfluorooctanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #370533**
- ✕ 30 um Uberstrainer **pluriSelect Life Science Catalog #43-70030-03**
- ✕ 50 mL centrifuge tubes **VWR International (Avantor) Catalog #734-1876**
- ✕ Droplet Generation Oil **Bio-Rad Laboratories Catalog #186-4006**
- ✕ Mineral Oil **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M5904**
- ✕ Tris 2M pH 8.0 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T3069**
- ✕ 10% SDS **Teknova Catalog #S0288**
- ✕ Tween-20 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P9416**
- ✕ Maxima 5X RT Buffer **Thermo Fisher Scientific**
- ✕ Maxima™ H Minus Reverse Transcriptase **Thermo Fisher Scientific Catalog #EP0753**
- ✕ 10 mM dNTPs **Takara Bio Inc. Catalog #639125**
- ✕ RNase Inhibitor **Lucigen Catalog #30281-1**
- ✕ 10x Exo I Buffer **Thermo Fisher Scientific**
- ✕ Exo I **Thermo Fisher Scientific Catalog #FEREN0582**
- ✕ Kapa HiFi Hotstart Readymix **Kapa Biosystems Catalog #KK2601**
- ✕ Agencourt AMPure XP beads **Beckman Coulter Catalog #A63881**
- ✕ Qubit 4 Fluorometer **Thermo Fisher Scientific Catalog #Q33238**
- ✕ Qubit dsDNA HS Assay Kit **Thermo Fisher Scientific Catalog #Q32854**
- ✕ Nextera XT DNA Library Preparation Kit **Illumina, Inc. Catalog #FC-131-1096**
- ✕ MiSeq v3 (150 cycle) Kit **Illumina, Inc. Catalog #MS-102-3001**

#### STEP MATERIALS

- ✕ Nuclease-free H<sub>2</sub>O **Life Technologies Catalog #AM9930**
- ✕ Sarkosyl **Merck MilliporeSigma (Sigma-Aldrich) Catalog #L7414**
- ✕ EDTA 0.5 M **Life Technologies Catalog #AM9260G**
- ✕ Tris 2M pH 7.5 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T2944**
- ✕ DTT **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D0632**
- ✕ SSC **Thermo Fisher Scientific Catalog #AM9765**
- ✕ Tris 2M pH 8.0 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T3069**



⊗ EDTA 0.5 M **Life Technologies Catalog #AM9260G**

⊗ 10% SDS **Teknova Catalog #S0288**

⊗ EDTA 0.5 M **Life Technologies Catalog #AM9260G**

⊗ Tween-20 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P9416**

⊗ 10x Exo I Buffer **Thermo Fisher Scientific**

⊗ Nuclease-free H<sub>2</sub>O **Life Technologies Catalog #AM9930**

⊗ Exo I **Thermo Fisher Scientific Catalog #FEREN0582**

⊗ Barcoded bead, sequence: TTTTTTTAAGCAGTGGTATCAACGCAGAGTACJJJJJJJJJJJ NNNNNNNNT(30); where J=split-pool oligo; N=random oligo **Chemgenes Catalog #Macosko-2011- 10**

⊗ 100 micron cell strainers (for beads) **VWR International (Avantor) Catalog #21008-950**

⊗ Fuchs-Rosenthal hemocytometer **INCYTO Catalog #DHC-F01**

⊗ Syringe Pumps (3) **Catalog #Legato 100**

⊗ Inverted Microscope **Catalog #Motic AE31**

⊗ Magnetic Stirrer **VP Scientific Catalog #710D2**

⊗ Syringe Pumps (3) **Catalog #Legato 100**

⊗ Inverted Microscope **Catalog #Motic AE31**

⊗ PDMS co-flow microfluidic droplet generation device **Catalog #CAD file from McCarroll Lab**

⊗ Tubing **Scientific Commodities Catalog #BB31695PE/2**

⊗ 15 mL Falcon tubes

⊗ Barcoded bead, sequence: TTTTTTTAAGCAGTGGTATCAACGCAGAGTACJJJJJJJJJJJ NNNNNNNNT(30); where J=split-pool oligo; N=random oligo **Chemgenes Catalog #Macosko-2011- 10**

⊗ Mineral Oil **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M5904**

⊗ Perfluorooctanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #370533**

⊗ Maxima 5X RT Buffer **Thermo Fisher Scientific**

⊗ Nuclease-free H<sub>2</sub>O **Life Technologies Catalog #AM9930**

⊗ SMART PCR primer: AAGCAGTGGTATCAACGCAGAGT **Integrated DNA Technologies, Inc. (IDT)**

⊗ Kapa HiFi Hotstart Readymix **Kapa Biosystems Catalog #KK2601**

⊗ Agencourt AMPure XP beads **Beckman Coulter Catalog #A63881**

⊗ Nextera XT DNA Library Preparation Kit **Illumina, Inc. Catalog #FC-131-1096**

⊗ Nextera XT DNA Library Preparation Kit **Illumina, Inc. Catalog #FC-131-1096**

⊗ Nuclease-free H<sub>2</sub>O **Life Technologies Catalog #AM9930**



- ✕ New-P5-SMART PCR hybrid oligo: AATGATACGGCGACCACCGAGATCTACACGCCT GTCCGCGGAAGCAGTGGTATCAACGCAGAGT\* A\*C Integrated DNA Technologies, Inc. (IDT)
- ✕ Agencourt AMPure XP beads **Beckman Coulter Catalog #A63881**
- ✕ Agencourt AMPure XP beads **Beckman Coulter Catalog #A63881**
- ✕ Ficoll PM-400 20% in H2O **Merck MilliporeSigma (Sigma-Aldrich) Catalog #F5415-50ML**
- ✕ BSA **Merck MilliporeSigma (Sigma-Aldrich) Catalog ##A8806**
- ✕ 10x PBS **Thermo Fisher Scientific Catalog #AM9624**
- ✕ UltraPure Distilled Water **Invitrogen - Thermo Fisher Catalog #10977-015**
- ✕ Nuclease-free H2O **Life Technologies Catalog #AM9930**
- ✕ Maxima 5X RT Buffer **Thermo Fisher Scientific**
- ✕ Ficoll PM-400 20% in H2O **Merck MilliporeSigma (Sigma-Aldrich) Catalog #F5415-50ML**
- ✕ 10 mM dNTPs **Takara Bio Inc. Catalog #639125**
- ✕ RNase Inhibitor **Lucigen Catalog #30281-1**
- ✕ TSO: AAGCAGTGGTATCAACGCAGAGTGAATrGrGrG **Integrated DNA Technologies, Inc. (IDT)**
- ✕ Maxima™ H Minus Reverse Transcriptase **Thermo Fisher Scientific Catalog #EP0753**
- ✕ 50 mL centrifuge tubes **VWR International (Avantor) Catalog #734-1876**
- ✕ 30 um Uberstrainer **pluriSelect Life Science Catalog #43-70030-03**
- ✕ Custom Read 1 primer: GCCTGTCCGCGGAAGCAGTGGTATCAACGCAG AGTAC **Integrated DNA Technologies, Inc. (IDT)**
- ✕ MiSeq v3 (150 cycle) Kit **Illumina, Inc. Catalog #MS-102-3001**
- ✕ Qubit 4 Fluorometer **Thermo Fisher Scientific Catalog #Q33238**
- ✕ Qubit dsDNA HS Assay Kit **Thermo Fisher Scientific Catalog #Q32854**
- ✕ Qubit 4 Fluorometer **Thermo Fisher Scientific Catalog #Q33238**
- ✕ Qubit dsDNA HS Assay Kit **Thermo Fisher Scientific Catalog #Q32854**
- ✕ Tubing **Scientific Commodities Catalog #BB31695PE/2**
- ✕ Luer lock 26-gage needles **Becton Dickinson (BD) Catalog #305111**
- ✕ 10 mL syringes **Becton Dickinson (BD) Catalog #BD 309695**
- ✕ 3 mL syringes **Becton Dickinson (BD) Catalog #BD #309657**
- ✕ Magnetic Mixing Discs **VP Scientific Catalog #772DP-N42-5-2**

## Protocol materials

- ✕ Nuclease-free H2O **Life Technologies Catalog #AM9930**
- ✕ 10x PBS **Thermo Fisher Scientific Catalog #AM9624**
- ✕ MiSeq v3 (150 cycle) Kit **Illumina, Inc. Catalog #MS-102-3001**
- ✕ 10x PBS **Thermo Fisher Scientific Catalog #AM9624**
- ✕ Barcoded bead, sequence: TTTTTTTAAGCAGTGGTATCAACGCAGAGTACJJJJJJJJJJJ NNNNNNNNT(30); where J=split-pool oligo; N=random oligo **Chemgenes Catalog #Macosko-2011- 10**
- ✕ Custom Read 1 primer: GCCTGTCCGCGGAAGCAGTGGTATCAACGCAG AGTAC **Integrated DNA Technologies, Inc. (IDT)**
- ✕ 40 micron cell strainers (for cells) **VWR International (Avantor) Catalog #21008-949**
- ✕ TSO: AAGCAGTGGTATCAACGCAGAGTGAATrGrGrG **Integrated DNA Technologies, Inc. (IDT)**
- ✕ RNase Inhibitor **Lucigen Catalog #30281-1**
- ✕ Qubit dsDNA HS Assay Kit **Thermo Fisher Scientific Catalog #Q32854**
- ✕ Magnetic Mixing Discs **VP Scientific Catalog #772DP-N42-5-2**
- ✕ 50 mL centrifuge tubes **VWR International (Avantor) Catalog #734-1876**
- ✕ Maxima™ H Minus Reverse Transcriptase **Thermo Fisher Scientific Catalog #EP0753**
- ✕ Agencourt AMPure XP beads **Beckman Coulter Catalog #A63881**
- ✕ Fuchs-Rosenthal hemocytometer **INCYTO Catalog #DHC-F01**
- ✕ Qubit dsDNA HS Assay Kit **Thermo Fisher Scientific Catalog #Q32854**
- ✕ 10 mL syringes **Becton Dickinson (BD) Catalog #BD 309695**
- ✕ PDMS co-flow microfluidic droplet generation device **Catalog #CAD file from McCarroll Lab**
- ✕ 30 um Uberstrainer **pluriSelect Life Science Catalog #43-70030-03**
- ✕ PDMS co-flow microfluidic droplet generation device **Catalog #CAD file from McCarroll Lab**
- ✕ Maxima™ H Minus Reverse Transcriptase **Thermo Fisher Scientific Catalog #EP0753**
- ✕ Perfluorooctanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #370533**
- ✕ New-P5-SMART PCR hybrid oligo: AATGATACGGCGACCACCGAGATCTACACGCCT GTCCGCGGAAGCAGTGGTATCAACGCAGAGT\* A\*C **Integrated DNA Technologies, Inc. (IDT)**
- ✕ 10 mM dNTPs **Takara Bio Inc. Catalog #639125**
- ✕ 15 mL Falcon tubes
- ✕ 10 mL syringes **Becton Dickinson (BD) Catalog #BD 309695**
- ✕ UltraPure Distilled Water **Invitrogen - Thermo Fisher Catalog #10977-015**



- ✕ Syringe Pumps (3) **Catalog #Legato 100**
- ✕ Syringe Pumps (3) **Catalog #Legato 100**
- ✕ Qubit 4 Fluorometer **Thermo Fisher Scientific Catalog #Q33238**
- ✕ Agencourt AMPure XP beads **Beckman Coulter Catalog #A63881**
- ✕ Inverted Microscope **Catalog #Motic AE31**
- ✕ Magnetic Stirrer **VP Scientific Catalog #710D2**
- ✕ Ficoll PM-400 20% in H<sub>2</sub>O **Merck MilliporeSigma (Sigma-Aldrich) Catalog #F5415-50ML**
- ✕ Exo I **Thermo Fisher Scientific Catalog #FEREN0582**
- ✕ Nuclease-free H<sub>2</sub>O **Life Technologies Catalog #AM9930**
- ✕ Inverted Microscope **Catalog #Motic AE31**
- ✕ Perfluorooctanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #370533**
- ✕ UltraPure Distilled Water **Invitrogen - Thermo Fisher Catalog #10977-015**
- ✕ 50 mL centrifuge tubes **VWR International (Avantor) Catalog #734-1876**
- ✕ Custom Read 1 primer: GCCTGTCCGCGGAAGCAGTGGTATCAACGCAG AGTAC **Integrated DNA Technologies, Inc. (IDT)**
- ✕ Tubing **Scientific Commodities Catalog #BB31695PE/2**
- ✕ 3 mL syringes **Becton Dickinson (BD) Catalog #BD #309657**
- ✕ Tubing **Scientific Commodities Catalog #BB31695PE/2**
- ✕ Mineral Oil **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M5904**
- ✕ Tris 2M pH 8.0 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T3069**
- ✕ 10% SDS **Teknova Catalog #S0288**
- ✕ EDTA 0.5 M **Life Technologies Catalog #AM9260G**
- ✕ Tris 2M pH 8.0 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T3069**
- ✕ Kapa HiFi Hotstart Readymix **Kapa Biosystems Catalog #KK2601**
- ✕ EDTA 0.5 M **Life Technologies Catalog #AM9260G**
- ✕ EDTA 0.5 M **Life Technologies Catalog #AM9260G**
- ✕ 15 mL Falcon tubes
- ✕ Qubit dsDNA HS Assay Kit **Thermo Fisher Scientific Catalog #Q32854**
- ✕ Sarkosyl **Merck MilliporeSigma (Sigma-Aldrich) Catalog #L7414**
- ✕ Nextera XT DNA Library Preparation Kit **Illumina, Inc. Catalog #FC-131-1096**
- ✕ Ficoll PM-400 20% in H<sub>2</sub>O **Merck MilliporeSigma (Sigma-Aldrich) Catalog #F5415-50ML**



New-P5-SMART PCR hybrid oligo: AATGATACGGCGACCACCGAGATCTACACGCCT  
GTCCGCGGAAGCAGTGGTATCAACGCAGAGT\* A\*C Integrated DNA Technologies, Inc. (IDT)



30 um Uberstrainer **pluriSelect Life Science Catalog #43-70030-03**



Luer lock 26-gage needles **Becton Dickinson (BD) Catalog #305111**



Tris 2M pH 7.5 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T2944**



Maxima 5X RT Buffer **Thermo Fisher Scientific**



Agencourt AMPure XP beads **Beckman Coulter Catalog #A63881**



Barcoded Bead SeqB: 5' -Bead-Linker-TTTTTTTAAGCAGTGGTATCAAC  
GCAGAGTACJJJJJJJJJJNNNNNNNN TTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTT-3' **Chemgenes**



Syringe Pumps (3) **Catalog #Legato 100**



Luer lock 26-gage needles **Becton Dickinson (BD) Catalog #305111**



Tween-20 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P9416**



Magnetic Stirrer **VP Scientific Catalog #710D2**



Qubit 4 Fluorometer **Thermo Fisher Scientific Catalog #Q33238**



Sarkosyl **Merck MilliporeSigma (Sigma-Aldrich) Catalog #L7414**



EDTA 0.5 M **Life Technologies Catalog #AM9260G**



Fuchs-Rosenthal hemocytometer **INCYTO Catalog #DHC-F01**



Inverted Microscope **Catalog #Motic AE31**



BSA **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A8806**



Droplet Generation Oil **Bio-Rad Laboratories Catalog #186-4006**



Nextera XT DNA Library Preparation Kit **Illumina, Inc. Catalog #FC-131-1096**



100 micron cell strainers (for beads) **VWR International (Avantor) Catalog #21008-950**



DTT **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D0632**



SMART PCR primer: AAGCAGTGGTATCAACGCAGAGT **Integrated DNA Technologies, Inc. (IDT)**



Agencourt AMPure XP beads **Beckman Coulter Catalog #A63881**



TSO: AAGCAGTGGTATCAACGCAGAGTGAATrGrGrG **Integrated DNA Technologies, Inc. (IDT)**



100 micron cell strainers (for beads) **VWR International (Avantor) Catalog #21008-950**



10 mM dNTPs **Takara Bio Inc. Catalog #639125**



RNase Inhibitor **Lucigen Catalog #30281-1**



10x Exo I Buffer **Thermo Fisher Scientific**



- ✕ Magnetic Mixing Discs **VP Scientific Catalog #772DP-N42-5-2**
- ✕ Tween-20 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P9416**
- ✕ Nextera XT DNA Library Preparation Kit **Illumina, Inc. Catalog #FC-131-1096**
- ✕ PCR tubes
- ✕ 10% SDS **Teknova Catalog #S0288**
- ✕ DTT **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D0632**
- ✕ Tubing **Scientific Commodities Catalog #BB31695PE/2**
- ✕ SMART PCR primer: AAGCAGTGGTATCAACGCAGAGT **Integrated DNA Technologies, Inc. (IDT)**
- ✕ Nuclease-free H<sub>2</sub>O **Life Technologies Catalog #AM9930**
- ✕ Maxima 5X RT Buffer **Thermo Fisher Scientific**
- ✕ Nuclease-free H<sub>2</sub>O **Life Technologies Catalog #AM9930**
- ✕ SSC **Thermo Fisher Scientific Catalog #AM9765**
- ✕ Qubit 4 Fluorometer **Thermo Fisher Scientific Catalog #Q33238**
- ✕ Exo I **Thermo Fisher Scientific Catalog #FEREN0582**
- ✕ Kapa HiFi Hotstart Readymix **Kapa Biosystems Catalog #KK2601**
- ✕ Maxima 5X RT Buffer **Thermo Fisher Scientific**
- ✕ SSC **Thermo Fisher Scientific Catalog #AM9765**
- ✕ Barcoded bead, sequence: TTTTTTTAAGCAGTGGTATCAACGCAGAGTACJJJJJJJJJJJ NNNNNNNNT(30); where J=split-pool oligo; N=random oligo **Chemgenes Catalog #Macosko-2011- 10**
- ✕ Nuclease-free H<sub>2</sub>O **Life Technologies Catalog #AM9930**
- ✕ Nuclease-free H<sub>2</sub>O **Life Technologies Catalog #AM9930**
- ✕ BSA **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A8806**
- ✕ 3 mL syringes **Becton Dickinson (BD) Catalog #BD #309657**
- ✕ Mineral Oil **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M5904**
- ✕ TrypLE™ Express Enzyme **Thermo Fisher Scientific Catalog #12604013**
- ✕ Tris 2M pH 7.5 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T2944**
- ✕ 10x Exo I Buffer **Thermo Fisher Scientific**
- ✕ Ficoll PM-400 20% in H<sub>2</sub>O **Merck MilliporeSigma (Sigma-Aldrich) Catalog #F5415-50ML**
- ✕ Syringe Pumps (3) **Catalog #Legato 100**
- ✕ PDMS co-flow microfluidic droplet generation device **Catalog #CAD file from McCarroll Lab**
- ✕ Tubing **Scientific Commodities Catalog #BB31695PE/2**



✕ 15 mL Falcon tubes

✕ Inverted Microscope **Catalog #Motic AE31**

✕ Barcoded bead, sequence: TTTTTTTAAGCAGTGGTATCAACGCAGAGTACJJJJJJJJJJJJ NNNNNNNNT(30); where J=split-pool oligo; N=random oligo **Chemgenes Catalog #Macosko-2011- 10**

✕ Ficoll PM-400 20% in H<sub>2</sub>O **Merck MilliporeSigma (Sigma-Aldrich) Catalog #F5415-50ML**

✕ EDTA 0.5 M **Life Technologies Catalog #AM9260G**

✕ Tris 2M pH 7.5 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T2944**

✕ Nuclease-free H<sub>2</sub>O **Life Technologies Catalog #AM9930**

✕ Sarkosyl **Merck MilliporeSigma (Sigma-Aldrich) Catalog #L7414**

✕ DTT **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D0632**

✕ BSA **Merck MilliporeSigma (Sigma-Aldrich) Catalog ##A8806**

✕ SSC **Thermo Fisher Scientific Catalog #AM9765**

✕ EDTA 0.5 M **Life Technologies Catalog #AM9260G**

✕ 10% SDS **Teknova Catalog #S0288**

✕ EDTA 0.5 M **Life Technologies Catalog #AM9260G**

✕ Tween-20 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P9416**

✕ Tris 2M pH 8.0 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T3069**

✕ Syringe Pumps (3) **Catalog #Legato 100**

✕ Inverted Microscope **Catalog #Motic AE31**

✕ Magnetic Stirrer **VP Scientific Catalog #710D2**

✕ Fuchs-Rosenthal hemocytometer **INCYTO Catalog #DHC-F01**

✕ 100 micron cell strainers (for beads) **VWR International (Avantor) Catalog #21008-950**

✕ Barcoded bead, sequence: TTTTTTTAAGCAGTGGTATCAACGCAGAGTACJJJJJJJJJJJJ NNNNNNNNT(30); where J=split-pool oligo; N=random oligo **Chemgenes Catalog #Macosko-2011- 10**

✕ Tubing **Scientific Commodities Catalog #BB31695PE/2**

✕ Luer lock 26-gage needles **Becton Dickinson (BD) Catalog #305111**

✕ 10 mL syringes **Becton Dickinson (BD) Catalog #BD 309695**

✕ 3 mL syringes **Becton Dickinson (BD) Catalog #BD #309657**

✕ Magnetic Mixing Discs **VP Scientific Catalog #772DP-N42-5-2**

✕ 30 um Uberstrainer **pluriSelect Life Science Catalog #43-70030-03**

✕ 50 mL centrifuge tubes **VWR International (Avantor) Catalog #734-1876**



- ✕ Mineral Oil **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M5904**
- ✕ Perfluorooctanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #370533**
- ✕ Maxima 5X RT Buffer **Thermo Fisher Scientific**
- ✕ Ficoll PM-400 20% in H<sub>2</sub>O **Merck MilliporeSigma (Sigma-Aldrich) Catalog #F5415-50ML**
- ✕ 10 mM dNTPs **Takara Bio Inc. Catalog #639125**
- ✕ TSO: AAGCAGTGGTATCAACGCAGAGTGAATrGrGrG **Integrated DNA Technologies, Inc. (IDT)**
- ✕ Maxima™ H Minus Reverse Transcriptase **Thermo Fisher Scientific Catalog #EP0753**
- ✕ Nuclease-free H<sub>2</sub>O **Life Technologies Catalog #AM9930**
- ✕ Maxima 5X RT Buffer **Thermo Fisher Scientific**
- ✕ RNase Inhibitor **Lucigen Catalog #30281-1**
- ✕ Nuclease-free H<sub>2</sub>O **Life Technologies Catalog #AM9930**
- ✕ Exo I **Thermo Fisher Scientific Catalog #FEREN0582**
- ✕ 10x Exo I Buffer **Thermo Fisher Scientific**
- ✕ Nuclease-free H<sub>2</sub>O **Life Technologies Catalog #AM9930**
- ✕ SMART PCR primer: AAGCAGTGGTATCAACGCAGAGT **Integrated DNA Technologies, Inc. (IDT)**
- ✕ Kapa HiFi Hotstart Readymix **Kapa Biosystems Catalog #KK2601**
- ✕ Agencourt AMPure XP beads **Beckman Coulter Catalog #A63881**
- ✕ Qubit 4 Fluorometer **Thermo Fisher Scientific Catalog #Q33238**
- ✕ Qubit dsDNA HS Assay Kit **Thermo Fisher Scientific Catalog #Q32854**
- ✕ Nextera XT DNA Library Preparation Kit **Illumina, Inc. Catalog #FC-131-1096**
- ✕ Nextera XT DNA Library Preparation Kit **Illumina, Inc. Catalog #FC-131-1096**
- ✕ Nuclease-free H<sub>2</sub>O **Life Technologies Catalog #AM9930**
- ✕ New-P5-SMART PCR hybrid oligo: AATGATACGGCGACCACCGAGATCTACACGCCT  
GTCCGCGGAAGCAGTGGTATCAACGCAGAGT\* A\*C **Integrated DNA Technologies, Inc. (IDT)**
- ✕ Agencourt AMPure XP beads **Beckman Coulter Catalog #A63881**
- ✕ Qubit 4 Fluorometer **Thermo Fisher Scientific Catalog #Q33238**
- ✕ Qubit dsDNA HS Assay Kit **Thermo Fisher Scientific Catalog #Q32854**
- ✕ Agencourt AMPure XP beads **Beckman Coulter Catalog #A63881**
- ✕ MiSeq v3 (150 cycle) Kit **Illumina, Inc. Catalog #MS-102-3001**
- ✕ Custom Read 1 primer: GCCTGTCCGCGGAAGCAGTGGTATCAACGCAG AGTAC **Integrated DNA Technologies, Inc. (IDT)**



⌘ 10x PBS Thermo Fisher Scientific Catalog #AM9624

⌘ UltraPure Distilled Water Invitrogen - Thermo Fisher Catalog #10977-015

## Troubleshooting

## Prepare Buffers and Solutions

### 1 **BSA**

Make a 10% stock solution using BSA powder

 BSA Merck MilliporeSigma (Sigma-Aldrich) Catalog #A8806


### 2 **1X PBS**


Make a 1X stock solution using by diluting 10X PBS into nuclease-free distilled water


 10x PBS Thermo Fisher Scientific Catalog #AM9624

 UltraPure Distilled Water Invitrogen - Thermo Fisher Catalog #10977-015

### 3 **Cell Loading Buffer** (makes 1 mL )

Prepare enough for the number of samples (  1 mL per sample ) that will be processed in one experiment

 300  $\mu$ L 20% Ficoll PM-400


 Ficoll PM-400 20% in H<sub>2</sub>O Merck MilliporeSigma (Sigma-Aldrich) Catalog #F5415-50ML

 700  $\mu$ L 1X PBS

### 4 **PBS-BSA:** make this fresh before each experiment

- Cell Loading Buffer
- 0.01% BSA (use the 10% stock)


### 5 **Lysis Buffer** (makes 1.2 mL )

Prepare enough for the number of samples (  1 mL per sample ) that will be processed in one experiment

 960  $\mu$ L H<sub>2</sub>O

 Nuclease-free H<sub>2</sub>O Life Technologies Catalog #AM9930

 12  $\mu$ L 20% Sarkosyl

 Sarkosyl Merck MilliporeSigma (Sigma-Aldrich) Catalog #L7414

 48  $\mu$ L 0.5 M EDTA

 EDTA 0.5 M Life Technologies Catalog #AM9260G

 120  $\mu$ L 2M Tris, pH 7.5



⊗ Tris 2M pH 7.5 Merck MilliporeSigma (Sigma-Aldrich) Catalog #T2944

🧴 60  $\mu$ L 1M DTT

Add this just prior to starting each DropSeq experiment

⊗ DTT Merck MilliporeSigma (Sigma-Aldrich) Catalog #D0632

## 6 **6X SSC**

Make 6X SSC working stock from 20X SSC

⊗ SSC Thermo Fisher Scientific Catalog #AM9765

## 7 **10 mM Tris pH 8.0**

Make 10 mM Tris pH 8.0 working stock from 2M Tris pH 8.0

⊗ Tris 2M pH 8.0 Merck MilliporeSigma (Sigma-Aldrich) Catalog #T3069

## 8 **TE-SDS** (make 🧴 50 mL )

- 10 mM Tris pH 8.0 + 1 mM EDTA
- 0.5% SDS

⊗ EDTA 0.5 M Life Technologies Catalog #AM9260G

⊗ 10% SDS Teknova Catalog #S0288

## 9 **TE-TW** (make 🧴 500 mL )

- 10 mM Tris pH 8.0 + 1 mM EDTA
- 0.01% Tween-20

⊗ EDTA 0.5 M Life Technologies Catalog #AM9260G


⊗ Tween-20 Merck MilliporeSigma (Sigma-Aldrich) Catalog #P9416


## 10 **RT mix** (makes 🧴 200 $\mu$ L )


🧴 75  $\mu$ L H<sub>2</sub>O


⊗ Nuclease-free H<sub>2</sub>O Life Technologies Catalog #AM9930


🧴 40  $\mu$ L Maxima 5x RT Buffer


 Maxima 5X RT Buffer **Thermo Fisher Scientific**

 40 µL 20% Ficoll PM-400


 Ficoll PM-400 20% in H<sub>2</sub>O **Merck MilliporeSigma (Sigma-Aldrich) Catalog #F5415-50ML**


 20 µL 10 mM dNTPs

 10 mM dNTPs **Takara Bio Inc. Catalog #639125**

 5 µL RNase Inhibitor


 RNase Inhibitor **Lucigen Catalog #30281-1**

 10 µL 50 uM Template Switch Oligo

 TSO: AAGCAGTGGTATCAACGCAGAGTGAATrGrGrG **Integrated DNA Technologies, Inc. (IDT)**

 10 µL Maxima H- RTase

**Add after beginning breakage portion of the protocol (Step 27).**

 Maxima™ H Minus Reverse Transcriptase **Thermo Fisher Scientific Catalog #EP0753**

## 11 **Exonuclease Mix** (makes 200 µL )

 20 µL 10x Exo I Buffer

 10x Exo I Buffer **Thermo Fisher Scientific**

 170 µL H<sub>2</sub>O

 Nuclease-free H<sub>2</sub>O **Life Technologies Catalog #AM9930**


 10 µL Exo I

 Exo I **Thermo Fisher Scientific Catalog #FEREN0582**

## Droplet Formation





### 12 Preparing the Barcoded Beads from Dry Resin

Barcoded bead, sequence:


 TTTTTTAAAGCAGTGGTATCAACGCAGAGTACJJJJJJJJJJJ NNNNNNNNT(30);  
where J=split-pool oligo; N=random oligo  
**Chemgenes Catalog #Macosko-2011- 10**



The beads will arrive as a dry resin. To prepare the beads for an experiment from the dry resin, follow these steps:

1. Wash the resin with  30 mL ethanol
2. Wash twice with  30 mL TE-TW
3. Resuspend in  20 mL TE-TW
4. Pass through 100 micron strainer
5. Count the beads using a Fuchs-Rosenthal hemocytometer
6. Store the counted beads at  4 °C

#### Note

To pellet the beads when washing, centrifuge at 1000xg for  00:01:00

#### Note

**Do NOT vortex** the barcoded beads. Vortexing will shear the beads and cause damage. To mix, simply invert tube gently a couple of times until mixed.



100 micron cell strainers (for beads) **VWR International (Avantor) Catalog #21008-950**



Fuchs-Rosenthal hemocytometer **INCYTO Catalog #DHC-F01**

## 13 Arrangement of components



Syringe Pumps (3) **Catalog #Legato 100**




Inverted Microscope **Catalog #Motic AE31**




Magnetic Stirrer **VP Scientific Catalog #710D2**

### Arrangement of components



### Attachment of tubing into the device

Tubing is attached by pressing the tip of the tube into the circular punched holes.



Set up the syringe pumps next to the inverted microscope. Note: arrange the bead pump so that the syringe is pointing downward rather than horizontally. Place the magnet stirrer close to the barrel of the bead syringe.

Source: <http://mccarrolllab.org/download/905/>

## 14 Prepare Pump System


Syringe Pumps (3) **Catalog #Legato 100**


1. Power on the syringe pumps (switches located on the back left of the pump)
2. Set up correct flow rates using the screen on the pump



- Oil: 15,000 ul/hr
- Cells: 4,000 ul/hr
- Beads: 4,000 ul/hr

## 15 Prepare microfluidic cell

 Inverted Microscope **Catalog #Motic AE31**

 PDMS co-flow microfluidic droplet generation device **Catalog #CAD file from McCarroll Lab**


 Tubing **Scientific Commodities Catalog #BB31695PE/2**

 15 mL Falcon tubes

1. Cut and remove protective plastic top layer for one device on the chip
2. Cut the outflow tubing to proper length and connect it to the device so that the collection from the outflow can be safely collected
3. Use clean 15 mL tubes to catch the droplets
4. Use a waste container to catch waste

## 16 Prepare the beads

Barcoded bead, sequence:

 TTTTTTTAAGCAGTGGTATCAACGCAGAGTACJJJJJJJJJJJ NNNNNNNNT(30);  
where J=split-pool oligo; N=random oligo  
**Chemgenes Catalog #Macosko-2011- 10**

1. Take an aliquot of beads (final concentration: 120,000 beads/mL)
2. Spin down in a tabletop centrifuge
3. Remove the TE-TW storage buffer
4. Wash with DTT minus lysis buffer
5. Resuspend in Lysis buffer

## 17 Prepare nuclei according to the protocol "Isolation of single nuclei from solid tissues" steps 1-14.

[dx.doi.org/10.17504/protocols.io.ufketkw](https://dx.doi.org/10.17504/protocols.io.ufketkw)

Dilute nuclei to 100 cells/ul. Use PBS-BSA for final dilution

## 18 Prepare syringes and lines for microfluidic device

Prepare enough syringes and lines for all samples at once (4-6 samples per experiment)

1. Oil syringe: 10 ml
2. Cell/Sample syringes: 3 ml
3. Bead syringes: 3 ml with magnetic mixing disc
4. Measure out tubing for the lines; make sure they are long enough to reach the microfluidic device from where the syringes will sit on the pumps
5. For each line cut one blunt end (for blunt end needle) and one angled end (insert into the device)
6. UV all tubing and needles

 10 mL syringes **Becton Dickinson (BD) Catalog #BD 309695**

 3 mL syringes **Becton Dickinson (BD) Catalog #BD #309657**

 Magnetic Mixing Discs **VP Scientific Catalog #772DP-N42-5-2**

 Tubing **Scientific Commodities Catalog #BB31695PE/2**

 Luer lock 26-gage needles **Becton Dickinson (BD) Catalog #305111**



## 19 Load Syringes

### Note

To load syringe, firmly press the tip of a 1000 ul pipette into the head of the syringe and slowly pull back on the plunger to draw in the solution. Pressing the tip in firmly helps reduce the introduction of bubbles. While holding the syringe in a vertical orientation, gently push out the air bubbles. Affix a 26G needle and tubing.

### Note

Before droplet generation, coat connecting tubing and syringes with 1% BSA to prevent non-specific binding of nuclei to the surface, and then rinse with PBS


1. Flush the tubing with corresponding buffer; keep syringes vertical when flushing
2. Use  1 mL cell buffer to flush syringe/tube for each sample
3. Use  1 mL DTT minus lysis buffer to flush syringe/tube for each bead setup
4. Load the syringes; keep syringes vertical when loading
5. Load cells/nuclei into cell syringes, making sure there are no air bubbles
6. Load beads into bead syringes, making sure there are no air bubbles
7. After beads are loaded in the syringe, keep horizontal to prevent clogging

## 20 System Assembly

1. Place oil syringe in pump and secure; using the screen, hold the free tube over the waste collection, then slowly move the pusher to the plunger until oil drips out
2. Push the angled end of the tubing into the device
3. Place cell/sample syringe in pump and secure; using the screen, hold the free tube over the waste collection, then slowly move the pusher to the plunger until buffer drips out
4. Push the angled end of the tubing into the device
5. Before placing the bead syringe, turn on the magnetic mixer for stirring disc (speed 25-30). Carefully rotate bead syringe by hand to mix up the beads in the lysis buffer.
6. Place bead syringe in pump and secure. Make sure the magnetic disc is circling up and down vertically in the syringe to keep the beads well mixed.
7. Using the screen, hold the free tube over the waste collection, then slowly move the pusher to the plunger until buffer drips out
8. Push the angled end of the tubing into the device

**Attachment of tubing into the device**

Tubing is attached by pressing the tip of the tube into the circular punched holes.



Source: <http://mccarrolllab.org/download/905/>

**21 Run Pumps for Droplet Formation**

1. Start oil pump
2. Start bead pump
3. Start cells/sample pump

**22 Check for Droplet Formation**

1. Have outflow line go to waste
2. Watch the beads move from the syringe to the device
3. Once the beads get to the bead chamber in the device, use a clean microscope slide to catch a drop of output
4. Check slide under a different microscope to see if droplets are uniform in size and shape
5. Once droplets are uniform begin collecting

**23 Collect**

1. Collect uniform droplets in 15 ml Falcon tubes until cells/sample runs out
2. The pump will sound an alarm once the cell/sample runs out; turn off alarm
3. STOP collecting in the 15 ml Falcon tube immediately (place a waste container under the outflow)



## 24 Stop pumps

1. Stop all pumps in any order
2. Place the 15 ml Falcon tube with droplets **on ice**

## 25 When running more than one sample [go to step #20](#)

Repeat Step 20- Step 24 for each sample

## 26 Clean up (**AFTER** Reverse Transcription is set up to incubate)

1. Turn off magnetic mixer
2. Remove syringes from pumps; remove tubing from microfluidic device
3. Power off pumps (switches located on the back left of the pump)
4. Cell/sample outflow tubing goes in biohazardous solid waste
5. All other lines can go in regular waste
6. Cell/sample syringes go in biohazardous solid waste
7. All other syringes can go in regular waste
8. Cell/sample needles go in biohazardous sharps waste
9. All other needles go in regular sharps waste

## STAMP Collection

### 27 **Single cell/nuclei Transcriptomes Attached to MicroParticles (STAMP)**

#### Note

STAMPs are the barcoded beads with mRNA attached; for stability, work **on ice/keep cold** when working with STAMPs/RNA


#### Note

**Add Maxima H- RTase to RT Mix (Step 10)**

Prepare waste and collection tube







- Need one 50 mL Falcon tube for waste (can use same waste from preparing syringes)
- Need one new clean 50 mL Falcon tube for collection
- Need one 30 um Uberstrainer

- Autopipette, 5 mL pipette, and attachment for strainer

 50 mL centrifuge tubes **VWR International (Avantor) Catalog #734-1876**

 30 um Uberstrainer **pluriSelect Life Science Catalog #43-70030-03**

## 28 Breakage


1. Be careful not to disturb the oil/droplet band (direct buffer down the side of the 15 mL Falcon tube and control output so there is smooth flow)
2. Before droplet breakage, add  600 µL mineral oil and incubate in water bath at  72 °C for  00:05:00
3. Let sit on ice for  00:05:00
4. Remove the oil from the bottom of the 15 mL Falcon tube
5. Add  5 mL cold 6X SSC
6. Add  1 mL Perfluorooctanol (PFO) directly to the oil/droplet band to break the droplets; dispense in a circular motion directly over/onto the droplet band
7. Gently roll/rotate the 15 mL Falcon tube on ice to help break the droplets

 Mineral Oil **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M5904**

 Perfluorooctanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #370533**

## 29 Filter/Collection

Set up the 30 um Uberstrainer over the waste tube

1. Place strainer over waste tube (goal: save beads, discard supernatant)
2. Attach autopipette attachment to strainer
3. Angle/tilt tube so that no liquid will go up through the autopipette when using it for a vacuum
4. Use  1 mL cold 6X SSC to moisten the filter


Filter the beads (over waste tube)

1. Pass the aqueous (top layer) phase through the strainer
2. Pass the organic (bottom layer) phase through the strainer
3. Use a couple milliliters cold 6X SSC to rinse the 15 mL Falcon tube to try and get all the beads from the sides; pass the rinse through the strainer




4. Rinse strainer/filter twice with cold 6X SSC

Collect the beads (over new/clean 50 mL Falcon tube)

1. Carefully reverse the strainer over a new/clean 50 mL Falcon tube
2. Wash strainer with  1 mL cold 6X SSC . Repeat four more times (5 washes total)
3. Visually inspect the strainer to make sure all the beads have come off it





## 30 Spin Down


1. Using a  4 °C centrifuge , spin down the 50 mL Falcon tube at 1000xg for

 00:01:00

2. Carefully remove, supernatant, leaving approximately  1 mL cold 6X SSC








## 31 Transfer to 1.5 mL tube

1. Resuspend beads in remaining  1 mL cold 6X SSC
2. Transfer to a new/clean 1.5 mL microfuge tube
3. Spin down (1000xg for  00:01:00 )
4. Wash twice with  1 mL cold 6X SSC
5. Wash with  300 µL 5X RT buffer
6. Remove as much of the 5X RT buffer as possible without taking up any beads

 Maxima 5X RT Buffer **Thermo Fisher Scientific**

## Reverse Transcription

- 32 This step generates cDNA strands on the RNA hybridized to the bead primers. One RT mix is sufficient for processing approximately 90,000 beads.

1. Add  200 µL RT Mix to the beads
2. Incubate at  Room temperature with rotation for  00:30:00
3. Incubate at  42 °C with rotation for  01:30:00
4. Wash beads once with  1 mL TE-SDS
5. Wash beads twice with  1 mL TE-TW



6. If proceeding to Exonuclease I treatment, wash once with

🧴 1 mL 10 mM Tris pH 8.0

#### Note

Stopping point: if stopping, stop at TE-TW wash step. Beads can be stored at

🧊 4 °C in TE-TW

## Exonuclease I Treatment

33 This step chews back the excess bead primers that did not capture an RNA molecule. One Exonuclease Mix is sufficient for processing approximately 90,000 beads.

1. After washing once with 🧴 1 mL 10 mM Tris pH 8.0 , resuspend in

🧴 200 µL Exonuclease Mix

2. Incubate at 🧊 37 °C with rotation for ⌚ 00:45:00

3. Wash beads once with 🧴 1 mL TE-SDS

4. Wash beads twice with 🧴 1 mL TE-TW

5. If proceeding to PCR, wash once with 🧴 1 mL H<sub>2</sub>O

#### Note

Stopping point: if stopping, stop at TE-TW wash step. Beads can be stored at

🧊 4 °C in TE-TW

## PCR


34 This steps utilizes PCR to amplify the cDNA constructed during the Reverse Transcription step.

Bead Count


1. After washing once with 🧴 1 mL H<sub>2</sub>O , spin down to pellet beads

2. Remove supernatant



3. Add  1 mL H<sub>2</sub>O

4. Mix well by pipette to evenly resuspend the beads

5. Quickly remove  20 µL and pipette into a Fuchs-Rosenthal hemocytometer chamber


6. Count all 16 boxes

35 Calculate number of beads to split into PCR tubes

1. Calculate the concentraion (beads/ul) = (#beads counts/16) x 5


2. Apportion 2,000 beads into each PCR tube (want approximately 100 STAMPs per PCR tube)


36 Prepare Master Mix

Master Mix (  50 µL per PCR reaction )

 24.6 µL H<sub>2</sub>O

 Nuclease-free H<sub>2</sub>O **Life Technologies Catalog #AM9930**

 0.4 µL 100 uM SMART PCR Primer

 SMART PCR primer: AAGCAGTGGTATCAACGCAGAGT **Integrated DNA Technologies, Inc. (IDT)**

 25 µL 2x Kapa HiFi Hotstart ReadyMix

 Kapa HiFi Hotstart Readymix **Kapa Biosystems Catalog #KK2601**


37 PCR

1. Spin down tubes


2. Add  50 µL PCR Master Mix to each reaction



3. Mix well by pipette

4. Run the following program on a thermocycler

 95 °C for  00:03:00

**4 cycles of**

 98 °C for  00:00:20

 65 °C for  00:00:45

 72 °C for  00:03:00

**13 cycles of**



🌡️ 98 °C for ⌚ 00:00:20

🌡️ 67 °C for ⌚ 00:00:20

🌡️ 72 °C for ⌚ 00:03:00

**Then**

🌡️ 72 °C for ⌚ 00:05:00

🌡️ 4 °C forever

## Purification and Quantification

### 38 Magnetic Bead Purification

1. Using either KAPA PureBeads or AMPure XP beads, add

🧴 30 µL room temperature magnetic beads to each PCR sample (this is a 0.6x beads to samples ratio based on volume)

2. Purify according to manufacturer's instructions

3. Elute in 🧴 10 µL H<sub>2</sub>O

⊗ Agencourt AMPure XP beads **Beckman Coulter Catalog #A63881**

### 39 Quantification

Can be performed by preferred method (qPCR, Qubit Assay, BioAnalyzer, Tapestation); we used Qubit Assay.

1. Use 🧴 2 µL sample input for this quantification

2. Quantify according to manufacturer's instructions

⊗ Qubit 4 Fluorometer **Thermo Fisher Scientific Catalog #Q33238**

⊗ Qubit dsDNA HS Assay Kit **Thermo Fisher Scientific Catalog #Q32854**



#### Note

The yield for 2000 beads generated from a 50 cell/ul total yield should be approximately 1-100ng.


**QC Cutoff: > 1ng total**

#### 40 Amplified cDNA Gel

1. Confirm average size of cDNA library via gel electrophoresis
2. Average size should be between 1300-2000 bp

### Tagmentation

- 41 This step simultaneously tags and fragments the cDNA library using transposomes













**Figure 2: Nextera Library Preparation Biochemistry**—Nextera chemistry simultaneously fragments and tags DNA in a single step. A simple PCR amplification then appends sequencing adapters and sample indexes to each fragment.

Source: [https://www.illumina.com/content/dam/illumina-marketing/documents/products/datasheets/datasheet\\_nextera\\_dna\\_sample\\_prep.pdf](https://www.illumina.com/content/dam/illumina-marketing/documents/products/datasheets/datasheet_nextera_dna_sample_prep.pdf)

## 42 Nextera Kit

1. Preheat a thermocycler to  $55^{\circ}\text{C}$
2. Working in a cold rack/on ice: for each sample, combine  $600\text{ pg}$  purified cDNA with **H2O** in a total volume of  $5\text{ }\mu\text{L}$

3. To each tube, add  10  $\mu\text{L}$  Nextera TD buffer "Tagment DNA" and  5  $\mu\text{L}$  Amplicon Tagment enzyme "Amplicon Mix" (the total volume of the reaction is now  20  $\mu\text{L}$  )
  4. Mix by pipette approximately 5 times
  5. Spin down
  6. Incubate at  55  $^{\circ}\text{C}$  for  00:05:00
  7. Add  5  $\mu\text{L}$  Neutralization Buffer (the total volume of the reaction is now  25  $\mu\text{L}$  )
  8. Mix by pipette approximately 5 times
  9. Spin down (bubbles are normal)
  10. Incubate at  Room temperature for  00:05:00
-  Nextera XT DNA Library Preparation Kit **Illumina, Inc. Catalog #FC-131-1096**

#### 43 PCR work up of tagmented library


Prepare PCR Master Mix (  24  $\mu\text{L}$  per PCR reaction )


 15  $\mu\text{L}$  Nextera PCR mix

 Nextera XT DNA Library Preparation Kit **Illumina, Inc. Catalog #FC-131-1096**






 8  $\mu\text{L}$  H<sub>2</sub>O

 Nuclease-free H<sub>2</sub>O **Life Technologies Catalog #AM9930**

 1  $\mu\text{L}$  10 uM New-P5-SMART PCR hybrid oligo

New-P5-SMART PCR hybrid oligo: AATGATACGGCGACCACCGAGATCTACACGCCT  
 GTCCGCGGAAGCAGTGGTATCAACGCAGAGT\* A\*C **Integrated DNA Technologies, Inc. (IDT)**

#### 44 PCR

1. Add  24  $\mu\text{L}$  Tagmentation PCR Master Mix to each tube
  2. To each individual tube, add  1  $\mu\text{L}$  10 uM Nextera N7XX oligo The total volume of the reaction is now  50  $\mu\text{L}$
  3. Mix well and run the following program on a thermocycler
-  95  $^{\circ}\text{C}$  for  00:00:30

**12 cycles of**

🔥 95 °C for ⌚ 00:00:10

🔥 55 °C for ⌚ 00:00:30

🔥 72 °C for ⌚ 00:00:30

**Then**

🔥 72 °C for ⌚ 00:05:00

🔥 4 °C forever

## Purification and Quantification

### 45 Magnetic Bead Purification

1. Using either KAPA PureBeads or AMPure XP beads, add

🧴 30 µL room temperature magnetic beads to each PCR sample (this is a 0.6x beads to samples ratio based on volume)

2. Purify according to manufacturer's instructions

3. Elute in 🧴 12 µL to 🧴 15 µL H<sub>2</sub>O

⊗ Agencourt AMPure XP beads **Beckman Coulter Catalog #A63881**

### 46 2nd Magnetic Bead Purification

1. Using either KAPA PureBeads or AMPure XP beads, perform a 0.6x purification

2. Purify according to manufacturer's instructions

3. Elute in 🧴 12 µL to 🧴 15 µL H<sub>2</sub>O

⊗ Agencourt AMPure XP beads **Beckman Coulter Catalog #A63881**

### 47 Quantification

Can be performed by preferred method (qPCR, Qubit Assay, BioAnalyzer, Tapestation); we used Qubit Assay.

1. Use 🧴 1 µL sample input for this quantification

2. Quantify according to manufacturer's instructions

3. Convert ng/ul to nM



 Qubit 4 Fluorometer **Thermo Fisher Scientific Catalog #Q33238**

 Qubit dsDNA HS Assay Kit **Thermo Fisher Scientific Catalog #Q32854**

#### Note

The yield for 2000 beads generated from a 50 cell/ul final cell concentration should be approximately 400-1000 pg/ul.

**QC Cutoff >1ng**

## 48 Sequencing Library Gel

1. Confirm average size of Sequencing library via gel electrophoresis
2. Average size should be between 500-680 bp

#### Note

Smaller-sized libraries will have more polyA reads; larger libraries may have lower sequence cluster density and cluster quality. Although the target size is 500-680 bp, the range can be as broad as 420-700 bp.

## Sequencing

### 49 MiSeq Sequencing - **QC for estimation of library quality and number of nuclei captured**

1. Pool, denature, and dilute to loading concentration according to manufacturer's instructions
2. Sequencing specifications
  - Read 1: 30 bp
  - Read 2: 75~100 bp
  - Read 1 Index: 8 bp
  - Custom Read 1 primer

Hiseq 2500 Sequencing

1. Combine 8-12 snDrop-seq libraries to make a 10  $\mu$ l library pool at 3 nM for denaturation.
2. After final dilution, load a combined library at 12 pM to the sequencer
  - Read 1: 30 bp
  - Read 2: 75~100 bp
  - Read 1 Index: 8 bp
  - Custom Read 1 primer



Custom Read 1 primer: GCCTGTCCGCGGAAGCAGTGGTATCAACGCAG  
AGTAC Integrated DNA Technologies, Inc. (IDT)



MiSeq v3 (150 cycle) Kit Illumina, Inc. Catalog #MS-102-3001

## Sequencing

### 50 snDrop-seq data processing

Mapping, demultiplexing and QC processing:

[https://github.com/chensong611/Dropseq\\_pipeline](https://github.com/chensong611/Dropseq_pipeline)

- Paired-end reads are removed if read 1 had more than four non-T bases in the last ten bases (to remove all non-poly(T)-captured contaminated reads)
- Paired-end reads are removed if read 1 had one or more bases with a poor quality score (<10)
- The right mate of each read pair is trimmed to remove any portion of the SMART adaptor sequence or any large stretches of poly(A) tails (6 consecutive bp or larger)
- The trimmed reads are aligned to the human genome (GENCODE GRCH38) with STAR (e.g. v2.5) with default parameter settings
- Reads that mapped to intronic or exonic regions as per the GENCODE gene annotation are included in gene counts
- Barcode synthesis errors are corrected by inserting N at the last base of the cell barcode for reads in which the first 11 bases of the cell barcode are identical and the last T base of the UMI is the same
- UMI counts for each gene of each nucleus are assigned by collapsing UMI reads that had only 1 edit distance to create a digital expression matrix (genes as rows, cells or nuclei as columns)