

Oct 07, 2025

Single-nucleus RNA (snRNA) sequencing

DOI

<https://dx.doi.org/10.17504/protocols.io.kxygx41ndl8j/v1>

Anita Adami^{1,2}, Raquel Garza^{1,2}

¹Lund University; ²Aligning Science Across Parkinson's (ASAP)



Anita Adami

Lund University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygx41ndl8j/v1>

Protocol Citation: Anita Adami, Raquel Garza 2025. Single-nucleus RNA (snRNA) sequencing. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.kxygx41ndl8j/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 03, 2025



Last Modified: October 07, 2025

Protocol Integer ID: 228932

Keywords: ASAPCRN, nucleus rna, snrna, workflow for nuclei, rna, nuclei, sequencing step, sequencing workflow

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000520

Aligning Science Across Parkinson's

Grant ID: ASAP-024296

Aligning Science Across Parkinson's

Grant ID: ASAP-025170

Abstract

Step-by-step protocol describing the snRNA sequencing workflow for nuclei isolated from post-mortem human brain tissue.



Nuclei isolation from human brain tissue

5h

- 1 Isolate nuclei from human brain tissue following the protocol we have previously published:

Protocol

NAME

Nuclear Extraction from Tissue for FACS

CREATED BY

Raquel Garza

[Preview](#)

Single nuclei preparation for sequencing

30m

- 2 Load the single-nuclei suspension (8500-10000 nuclei/sample) onto the Chromium Next GEM Chip G Single Cell Kit together with the reverse transcription mastermix following the manufacturer's instructions. The protocol for the Chromium GEM single cell 3' kit to generate single-cell gel beads in emulsion can be found on the 10X Genomics website (10X Genomics, PN-1000268). It is possible that new kits and chemistries are released by 10X Genomics. In that case, it can be relevant to use the latest version of the kit and of the protocols provided by the manufacturer.
- 3 Amplify cDNA following the guidelines from 10X Genomics. We used 13 cycles of amplification of the 3' libraries.
- 4 Generate sequencing libraries with unique dual indices (TT set A).

30m

snRNA sequencing

3d

- 5 Pool the libraries from the previous step for sequencing on an Illumina sequencer (we sequenced either on a Novaseq6000 or on a Novaseq X Plus) using a 100-cycle kit and 28-10-10-90 reads.

3d