

Dec 24, 2025

Version 2

Single-step assembly of double guide plasmid (pCas9-Duo) for gene-editing in Plasmodium V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.n2bvjn9d5gk5/v2>

Abhinay Ramaprasad^{1,2}

¹Francis Crick Institute; ²School of Infection and Immunity, University of Glasgow



Abhinay Ramaprasad

University of Glasgow

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvjn9d5gk5/v2>

Protocol Citation: Abhinay Ramaprasad 2025. Single-step assembly of double guide plasmid (pCas9-Duo) for gene-editing in Plasmodium . protocols.io <https://dx.doi.org/10.17504/protocols.io.n2bvjn9d5gk5/v2> Version created by **Abhinay Ramaprasad**

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: December 24, 2025

Last Modified: December 24, 2025

Protocol Integer ID: 235780

Keywords: malaria, Cas9, SHFTiKO, Plasmodium, editing in plasmodium, targeting plasmid, plasmodium falciparum, plasmodium, step assembly of double guide, distinct guide rna, distinct guide rnas in order, mediated modification, successful cas9, pcas9, double guide, editing

Funders Acknowledgements:

Marie Skłodowska Curie Individual Fellowship

Grant ID: 751865

Wellcome Trust

Grant ID: 220318/A/20/Z

European Society of Clinical Microbiology and Infectious Diseases

Grant ID: 2023

MRC Career Development Award

Grant ID: MR/Z504762/1

Abstract

This protocol describes a one-pot GoldenGate assembly of a new dual-guide targeting plasmid (pCas9-Duo) expressing two distinct guide RNAs in order to enhance the chances of a successful Cas9-mediated modification at a target region in *Plasmodium falciparum*. Designed as part of SHFTiKO (frameshift-based trackable inducible knockout) system¹ (based on²).

Attachments



[pCas9-Duo.gb](#)

20KB



Protocol materials

⊗ T4 DNA Ligase Reaction Buffer - 6.0 ml **New England Biolabs Catalog #B0202S**

⊗ T4 Polynucleotide Kinase - 500 units **New England Biolabs Catalog #M0201S**

⊗ Nuclease-free Water

⊗ LB plates with 100 µg/ml ampicillin

⊗ UltraPure ATP 10mM **Promega Catalog #NC2683865**

⊗ BbsI-HF - 300 units **New England Biolabs Catalog #R3539S**

⊗ Cutsmart Buffer

⊗ Hi-T4 DNA Ligase **New England Biolabs Catalog #M2622S**

⊗ NEB 5-alpha Competent E.coli (High Efficiency) - 6×0.2 ml **New England Biolabs Catalog #C2987I**

⊗ XL-10 Gold Ultracompetent cells **Agilent Technologies Catalog #200314**

⊗ M13 Reverse

⊗ M13Forward_reverse



gRNA oligo design

- 1 Add the overhangs "ATTG" and "AAAC" to forward and reverse oligos of gRNA1, and "TTGG" and "TAAA" to gRNA2 respectively. For predicting gRNAs, please refer to this protocol -

Protocol


NAME

Design strategy for DiCre-mediated inducible gene knockout in the malaria parasite

CREATED BY

Abhinay Ramaprasad

Preview



Double guide plasmid containing two gRNA expression cassettes with distinct insertion sites.

Note

When you have lot of targets, you can also use this sheet to create final sequences ready for oligo ordering to avoid inadvertent errors -  [pCas9-Duo-overhangs.xlsx](#) 12.6KB

Anneal gRNA oligos

1h 10m

2 Set up annealing reactions.

10m

 1.0 μ L gRNA.F [M] 100 micromolar (μ M)

 1.0 μ L gRNA.R [M] 100 micromolar (μ M)

 1.0 μ L

 T4 DNA Ligase Reaction Buffer - 6.0 ml **New England Biolabs Catalog #B0202S**



🧪 0.5 µL

🧬 T4 Polynucleotide Kinase - 500 units **New England Biolabs Catalog #M0201S**

🧪 6.5 µL 🧬 Nuclease-free Water

3

Step 1 🌡️ 37 °C ⌚ 00:30:00

1h

Step 2 🌡️ 94 °C ⌚ 00:05:00

Step 3 🌡️ 90 °C to 🌡️ 25 °C RAMP 🌡️ 5 °C per minute

Step 4 🌡️ 4 °C Hold

Single-step GoldenGate assembly

1h 15m

4 Make 1:200 dilution mixture of annealed gRNA1 and gRNA2.

🧪 1 µL gRNA1 annealed reaction

🧪 1 µL gRNA2 annealed reaction

🧪 198 µL 🧬 Nuclease-free Water

5 Set up assembly reaction

🧪 2 µL pCas9 [M] 200 ng/µL

🧪 1 µL 1:200 diluted oligo

🧪 0.5 µL 🧬 BbsI-HF - 300 units **New England Biolabs Catalog #R3539S**

🧪 0.5 µL 🧬 Hi-T4 DNA Ligase **New England Biolabs Catalog #M2622S**

🧪 2 µL 🧬 Cutsmart Buffer

🧪 2 µL 🧬 UltraPure ATP 10mM **Promega Catalog #NC2683865**

🧪 12 µL 🧬 Nuclease-free Water

6

Step 1 🌡️ 37 °C ⌚ 00:05:00

1h

Step 2 🌡️ 25 °C ⌚ 00:05:00

Repeat Step 1-2 for 6 cycles

Step 3 🌡️ 4 °C Hold

Transform, plate and screen colonies

20m 40s



- 7 Add  2 μL of assembly reaction to  15 μL ultracompetent cells (like

 NEB 5-alpha Competent E.coli (High Efficiency) - 6×0.2 ml **New England Biolabs Catalog #C29871**

).

Note

We have observed that the competent cells used in the previous version -

 XL-10 Gold Ultracompetent cells **Agilent Technologies Catalog #200314** do not perform consistently for this reaction.

- 8 Place  On ice for  00:20:00 .

20m

- 9 Heat shock at  42 °C for  00:00:40 and spread transformed cells on

40s

 LB plates with 100 $\mu\text{g}/\text{ml}$ ampicillin .

- 10 Pick colonies, miniprep-isolate plasmids and screen for gRNA1 and gRNA2 insertions by Sanger sequencing using  M13 Reverse (CAGGAAACAGCTATGAC) and  M13Forward_reverse (ACTGGCCGTCGTTTTAC), respectively.

Expected result

Usually we get 2-18 colonies out of which 2 colonies are screened.

Note

If processing a lot of samples, you can use this Shiny app **pCas9-Duo-check** to quickly screen the presence of both gRNA inserts in your colonies. Just ensure that you enter sample and primer names when submitting for Sanger sequencing such that your .ab1 filenames will be a) named uniquely for each colony and b) contain "M13R" or "M13FR" as primer names.

- 11 For some targets, we have observed that whilst gRNA01 insertion is very efficient, gRNA02 repeatedly fails. In these cases, instead of repeating the reaction from scratch or screening more colonies, we have had more success by simply repeating the



GoldenGate reaction (Step 4 onwards) with pCas9-Duo-gRNA01 plasmid obtained in Step 10 as backbone and 1:200-diluted annealed gRNA02 as insert.

Protocol references

Protocol developed as part of

1. Ramaprasad, Abhinay, and Michael J. Blackman. 2024. 'A Scaleable Inducible Knockout System for Studying Essential Gene Function in the Malaria Parasite'. bioRxiv. <https://doi.org/10.1101/2024.01.14.575607>.

Based on

2. Adikusuma, Fatwa, Chandran Pfitzner, and Paul Quinton Thomas. 2017. 'Versatile Single-Step-Assembly CRISPR/Cas9 Vectors for Dual gRNA Expression'. *PLoS ONE* 12 (12): e0187236. <https://doi.org/10.1371/journal.pone.0187236>.