

Feb 12, 2026

Version 2

Single_Stranded_Library_Workflow_Gut_Virome V.2

 Microbiome

DOI

<https://dx.doi.org/10.17504/protocols.io.8epv5r5w5g1b/v2>


Xichuan Zhai¹

¹University of Copenhagen



Xichuan Zhai

University of Copenhagen



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5r5w5g1b/v2>

Protocol Citation: Xichuan Zhai 2026. Single_Stranded_Library_Workflow_Gut_Virome. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.8epv5r5w5g1b/v2> Version created by Xichuan Zhai

Manuscript citation:

<https://doi.org/10.21203/rs.3.rs-4304844/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 12, 2026

Last Modified: February 12, 2026

Protocol Integer ID: 243099

Keywords: Single-stranded library preparation, gut virome, single-stranded-library-workflow-gut-virome protocol, sslr, stranded library, preparation for virome study, virome study

Funders Acknowledgements:

VILLUM FONDEN

Grant ID: 23145

VILLUM FONDEN

Grant ID: 36242

Abstract

Protocol of single-stranded library (SSLR) preparation for virome study

Guidelines

This protocol is designed for virome library preparation with our lab-customed single-stranded library (SSLR). This workflow includes **7 FLEXIBLE** steps, virome isolation and extraction, mock community (positive control) preparation, reverse transcription with the purpose of dsRNA denatures, genome fragmentation, ligation with customized designed adaptors and library preparation for sequencing. Detailed information can be found on the next page.



Materials

SM buffer: lab-prepared according to https://cshprotocols.cshlp.org/content/2006/1/pdb.rec8111.full?text_only=true

PierceTM Universal Nuclease for Cell Lysis: #88701

QIAamp Viral RNA Mini Kit (250): # 52906

SuperScriptTM IV VILOTM Master Mix: #11756050

AMPure XP beads, 60 mL: #A63881

TE, pH 8.0, RNase-free: #AM9849

ET SSB (500 µg/mL): #M2401S

T4 Polynucleotide Kinase (10 U/µL): #EK0031

T4 DNA Ligase (5 U/µL): #EL0012

T4 DNA Ligase Buffer (10X) with 50% PEG-4000: #B69

Forward/reverse adapters: Lab designed and produced by IDT

AccuPrimeTM Taq DNA Polymerase System: #12339016

Nextera XT Index Kit v2 Sets: #FC-131-2001, 2002, 2003, 2004

Qubit 1X dsDNA HS Assay Kit: #Q33231

High Sensitivity D5000 ScreenTape Assay: #5067-5593

Protocol materials

 Fresh 80% Ethanol

Safety warnings

 Not applied

Ethics statement

Not applied

Before start

Fill bucket with crushed ice.



Step1: Virome isolation and extraction

1 Checklist before starting

Note

Where: Fecal lab

Timing: 180 min

Reagents/kits

- SM buffer: lab-prepared according to https://cshprotocols.cshlp.org/content/2006/1/pdb.rec8111.full?text_only=true
- PierceTM Universal Nuclease for Cell Lysis: #88701
- QIAamp Viral RNA Mini Kit (250): # 52906

Equipment

- Centrifuge
- 37°C incubator

Procedure

- 1.1 Make aliquots of  140 µL of enriched virome (procedure can be found from: dx.doi.org/10.17504/protocols.io.b2qagqds).
- 1.2 Add  1 µL of 100-time diluted nuclease (check the stock for the dilution in SM buffer) to each sample and let them incubate at approximately  00:30:00 at  37 °C per sample. Vortex between each sample. 30m
- 1.3 Consider increasing the incubation time if you expect a lot of external DNA.
- 1.4 Immediately after adding  540 µL AVL buffer to inactivate nucleases.
- 1.5 Mix the mixture by pulse vortexing.  
- 1.6 Incubate at  Room temperature for  00:10:00 . 10m

- 1.7 Change gloves.



- 1.8 Briefly centrifuge the mixture with a microcentrifuge.
- 1.9 Add  560 μL of absolute ethanol (96%) to the sample mixture. 
- 1.10 **CRITICAL STEP:** Mix very well by pulse-vortex. 
- 1.11 Briefly centrifuge the mixture with a microcentrifuge.
- 1.12 Add  630 μL of the sample mixture to the spin column.
- 1.13 **CRITICAL STEP:** Do not touch the column rim with the pipet.
- 1.14 Centrifuge at  6,000 x g, Room temperature change the collection tube and repeat steps 1.12 to load all the extracts. 
- 1.15 Add  500 μL AW1 buffer to the spin column the following: 
- 1.16 Centrifuge at  6,000 x g, Room temperature, 00:01:00 , change the collection tube. 

- 1.17 Add  500 μL AW2 buffer to the spin column. 
- 1.18 Centrifuge at  20,000 x g, Room temperature, 00:03:00 , change the collection tube. 

- 1.19 Centrifuge at  20,000 x g, Room temperature, 00:01:00 . Place the spin column in a low-binding RNase-free 1.5 mL tube. 


1.20 Add  30 µL of AVE (elution buffer) to the spin column and incubate

1m

 Room temperature for  00:01:00 .



1.21 **CRITICAL STEP:** Pipet the AVE buffer directly onto the filter membrane without touching it with the pipet tip.



1.22 .

1m

Make sure ethanol 96% has been added to AW1 and AW2 buffer.



Make sure ethanol 96% has been added to AW1 and AW2

bufferCentrifuge at  6,000 x g, Room temperature, 00:01:00 to collect the filtrates.

Note

PAUSE POINT Send for next step or store viral DNA/RNA at -80°C.

CAUTION

- **NOT** interested in **RNA**??? **THEN GO**  .
- Interested in **DNA** and **ssRNA**, **THEN GO**  .
- Interested in DNA, dsRNA and ssRNA, **THEN GO**  .
- Remember to include positive (Mock from extraction) and Negative (H₂O from extraction) controls for each extraction.
- Make sure 96% ethanol has been added to AW1 and AW2 buffers.
- Remember to include a negative control for each extraction.

Step2: Heat treatment

2 Checklist before starting

**Note****Where:** Clean-lab**Timing:** 3 min**Reagents**

No reagents needed

Equipment

- UV-beach
- Thermocycler
- Microcentrifuge

Procedure

2.1 Preheated ThermoCycle machine to 95 °C and sterilized PCR tubes or plates depending on the how many samples are used for reverse transcription.

2.2 Add 16 µL extracted virome in PCR tubes or plates and prepare on ice.

2.3 Briefly centrifuge the mixture with a microcentrifuge.

2.4 Put tubes or plates for heat treatment 95 °C for 00:03:00 .

2.5 **CRITICAL STEP:** Transfer samples to ice immediately after 00:03:00 .

2.6 **PAUSE POINT** Ready for reverse transcription.

Note**CAUTION**

- If you would like to look at the RNA in your samples, it is better to finish extraction, heat treatment, and reverse transcription **on the same day**.
- Pipette, filtered pipette tips, and workbench need to be sterilized by UV lamp for 20 min prior work.
- Remember to include Positive (Mock from extraction Step1) and Negative (SM buffer from extraction) controls for the reaction.

3m



3m



Step3: Reverse transcription (RT)

3 Checklist before starting

Note

Where: Clean-lab

Timing: 25 mim

Reagents/kits

- SuperScript™ IV VILO™ Master Mix: #11756050
- AMPure XP beads, 60 mL: #A63881
- TE, pH 8.0, RNase-free: #AM9849

Equipment

- Thermocycler
- Microcentrifuge
- UV workbench

Procedure

- 3.1 Transfer 16 µL of extracted or heat-treated virome DNA/RNA to a clean PCR plate.
- 3.2 Add  4 µL SuperScript™ IV VILO™ Master Mix to each sample and mix thoroughly, spin them down.
- 3.3 Carefully place the PCR plate into the ThermoCycle machine and select the following program:



	A	B
	Temperature profile	
	25°C	10 min
	50°C	10 min
	85°C	5 min
	4°C	∞



- 3.4 Purified with 1X AMPure XP beads (25 μ L / 25 μ L PCR reaction) and eluted with  20 μ L Tris buffer (10 mM, pH8.0) according to the purification procedures from  [go to step #8](#) .

Note

PAUSE POINT Samples can be stored in the fridge or freezer for downstream preparations.

Step4: Fragmentation

4 Checklist before starting

Note

Where: Basement (using ID card and key for access)

Timing: 15 min

Reagents

No reagent needed

Equipment

- Bioruptor® Pico sonication device (#B01060010)
- Microcentrifuge. Bring the microcentrifuge with you to the basement!
- Autoclaved Bioruptor tubes (volume:  650 μ L)

Procedure

- 4.1 **CRITICAL STEP:** Pre-cooled (4°C) the Biorupter and holder at least  00:30:00 .

30m



- 4.2 Transfer  22 μ L genomic DNA or reverse-transcription products to Bioruptor tubes and briefly spin down to make sure all the liquid is at the bottom of tubes.

- 4.3 Put samples on the ice for at least  00:10:00 .

10m



- 4.4 **CRITICAL STEP:** Set the sonication parameter as following: **15s ON and 90s OFF, using 8 cycles.**
- 4.5 **CRITICAL STEP:** Spin down the samples and carefully load the Bioruptor tubes to sonicator holder and close carefully, make sure there is no liquid at the side of tubes.
- 4.6 **CRITICAL STEP:** Put the tube holder into the sonication chamber and close the lids. 12 samples can be done each time.
- 4.7 Spin down the tubes after shearing.
- 4.8 The fragmented product is ready for ligation.

Note

PAUSE POINT Samples can be stored in fridge or freezer for downstream preparations.

CAUTION

- **Pregnant women** should not stay away from the running machine.
- **DO NOT** turn on the instrument without water.
- **Distilled water** should be used to fill the tank.
- **Always keep 12 tubes** for each run to ensure successful shearing.

Step5: Ligation

5 Checklist before starting

Note

Where: Clean-lab

Timing: 60 min

Reagents

- ET SSB (500 µg/mL): #M2401S
- T4 Polynucleotide Kinase (10 U/µL): #EK0031
- T4 DNA Ligase (5 U/µL): #EL0012
- T4 DNA Ligase Buffer (10X) with 50% PEG-4000: #B69
- Forward/reverse adapters: Lab designed and produced by IDT
- AMPure XP beads, 60 mL: #A63881
- Tris buffer (10 mM, pH 7.5)
- H₂O (autoclaved, UV-sterilized MiliQ)

Equipment

- Thermocycler
- Microcentrifuge
- 96-well PCR plate with seal
- Eppendorf tubes
- 50 mL Falcon tube
- Ice/cooling block

Procedure

5.1 Denature

5.2 According to the number of samples (Include extra 3 for pipetting errors), prepare an ET SSB dilution in a clean tube. In the table below is calculated for a whole plate of 96 wells (Master Mix per 100 samples): diluted ET SSB to 5 ng/µL (500 µg/mL in the stock solution) with Tris and then add same volume of Tris, for example 1 µL ET SSB in 100 µL tris buffer. Mix them, spin them down and place on ice.

	A	B	C
	Reagent	Per 1 sample µL	Per 100 samples µL
	Fragmented DNA	20	-
	ET SSB dilution (10 ng)	2	100
	Total denaturation reaction volume	22	

5.3 Add  20 µL of fragmented DNA to each well.

5.4 Add  2 µL ET SSB solution, pipette several times, and spin down.

5.5 Transfer to ThermoCycle and denature for  00:03:00 at  95 °C .

3m

5.6 Cooling down on ice immediately after heating and set  On ice for at least

5m

 00:05:00 .

5.7 Ligation

5.8 Spin down the heat denatured sample with a microcentrifuge.

5.9 According to the number of samples (Include extra 3 for pipetting errors), prepare a Master Mix containing per sample. In the table below is calculated for a whole plate of 96 wells (Master Mix per 100 samples). Mix them, spin them down and place on ice.

A	B	C
95 °C for 3min, cool down on ice immediately for at least 5 min		
Reagents	Per 1 sample µL	Per 100 samples µL
PEG-4000 (50%)	10	1000
H2O	8	800
T4 ligase buffer (10X final)	5	500
T4 PNK (10,000 units/mL)	1	100
Add the above reagents one by one and then mix for the 1st time		
Forward adapter (1 pmol)	1	100
Reserve adapter (1 pmol)	1	100
Add the above adapters one by one and then mix for the 2nd time		



A	B	C
T4 DNA ligase (400,000 units/mL)	2	200
Add the above reagents one by one and then by finger flicking to mix		
Final volume	50	50000
37 °C for 60 min (set the lid-heating off)		

Note

CRITICAL STEP for prepare master mix:

Make sure the T4 ligase buffer and PEG are dissolve and mix thoroughly;
Add the needed volume of PEG-4000 (sticky, do it slowly), ligase buffer and PNK first and mix thoroughly, then add adapters and mix again, ligase should be added lastly during cooling down the samples and mix thoroughly without vortexing.

5.10 Add  28 µL ligation master mix to the denatured samples, mix thoroughly with pipette, spin briefly.

5.11 Transfer to ThermoCycle for  01:00:00 ligation at  37 °C .

1h

5.12 Purified with 1X AMPure XP beads (25 µL / 25 µL PCR reaction) and eluted with  22 µL MiliQ H₂O according to the purification procedures from [Step7](#).

**Note**

PAUSE POINT Store in fridge for short-time or freezer until beads clean-up.

CAUTION

Ensure that you have booked a ThermoCycle and preheated it to 95°C .

Pipette, filtered pipette tips, and workbench need to be sterilized by UV lamp for

$00:20:00$ prior work.

During preparation, reagents and PCR plate must be on ice/cold block. Reagents are stored at -20°C .

Spin down reagents before using them. Avoid vortexing enzymes.

Remember to include a Positive (Mock from extraction), Negative (H_2O from extraction) &

Blank (H_2O for the ligation) controls for the reaction.

Step6: Post-ligation clean up

6 Checklist before starting

Note

Where: Clean lab

Timing: 30 min

Reagents

- AMPure XP beads, 60 mL: #A63881
- Qubit 1X dsDNA HS Assay Kit: #Q33231
- High Sensitivity D5000
- ScreenTape Assay: #5067-5593
- Freshly prepared 80% ethanol

Equipments/kits

- HulaMixture, Invitrogen™
- Qubit™ 4 Fluorometer
- TapeStation4200 with high-sensitive D5000 Screen Tape (#5067-5592)

Procedure

- 6.1 Place AMPure XP beads into the Hula mixer to resuspend the beads and to equilibrate to

Room temperature for $00:15:00$.

- 6.2 Total volume of ligation product (sample) / well in the PCR plate is $50\ \mu\text{L}$.



- 6.3 Transfer  40 μL (!!! 0.8X) of Beads solution to each PCR product/well and mix with 100 μL pipette tips (10 times up and down), resulting in a  90 μL mixture.
- 6.4 Incubate for  00:05:00 at  Room temperature .
- 6.5 Place the tube or plate with the mixture into the magnetic rack for 2-4 min, until the liquid is clear.
- 6.6 Carefully remove  80 μL the liquid/supernatant with 100 μL pipette tips and discard it, keep  5 μL without disturbing the beads-pellet.
- 6.7 Wash the beads-pellet with  200 μL freshly prepared  Fresh 80% Ethanol by gently dispensing it over the beads with 200 μL pipette tips. Let it rest for  00:00:30 and then remove the liquid.
- 6.8 Repeat washing  [go to step #6.7](#) .
- 6.9 Spin the tubes or plates in a minicentrifuge for  00:00:30 to collect all the residual liquid.
- 6.10 Put the tubes or plates on a magnetic rack and remove the excess with 10 μL pipette tips.
- 6.11 Air dry for approximately  00:00:30 to evaporate the Ethanol.
- 6.12 Remove the tube or plate from the magnetic rack. Add  22 μL nuclease free water and mix with a pipette (>10 times up and down) to resuspend the beads-pellet.
- 6.13 Incubate for  00:02:00 at  Room temperature .
- 6.14 Spin down for  00:00:30 if there is liquid on the wall of the tubes or the well.
- 6.15 Place the tube or plate back on the magnetic rack and wait for  00:03:00 until the liquid clears.



- 6.16 Transfer  22 μL of the liquid/supernatant to a new, labeled Eppendorf tube or a new plate.

Note

PAUSE POINT Store in fridge for short-time or freezer until beads clean-up.

CAUTION

- During the preparation of index PCR, reagents and PCR plate must be on ice/cold block.
- Lead opens while on the cold block.
- Reagents are stored at -20°C .
- Gently mix and spin down reagents before using them. Avoid vortexing enzymes.
- **WARNIN!!!** Change pipette tips and PCR leads during work process!

Step7: Index PCR

7 Checklist before starting

Note

Where: Sequencing lab

Timing: 30 mim

Reagents

- AccuPrimeTM Taq DNA Polymerase System: #12339016
- Nextera XT Index Kit v2 Sets: #FC-131-2001, 2002, 2003, 2004

Equipment

- Thermocycler
- Minicentrifuge

Procedure

- 7.1 Book an available ThermpCycler.
- 7.2 Take the index primer (illumina i5 + i7) plate from the freezer to the fridge before starting.
- 7.3 Remember to add a Blank control for the index PCR.

- 7.4 According to the number of samples (include an extra 3 for pipetting errors), prepare a master mix containing per sample. The table below is calculated for a whole plate of 96 wells (Master Mix per 100 samples):

	A	B	C
		Per 1 sample μL	Per 100 samples μL
	Purified ligated DNA	21.1	
	10× AccuPrime buffer	2.5	250
	AccuPrime DNA polymerase	0.4	40
	Primer P5 (i5)	1	
	Primer P7 (i7)		
	Final volume	25	

- 7.5 **Master Mix:** In a clean EP tube, transfer the needed volume of 10× AccuPrime buffer and AccuPrime DNA polymerase from the stock solution, gently mix, and spin down.

- 7.6 Transfer  2.9 μL Master Mix to each sample.

- 7.7 Transfer  1 μL Primer Mix to each sample.

- 7.8 Carefully mix them and spin them down.

- 7.9 Carefully place the PCR plate into the thermocycle machine and select the following program:



	A	B	C
	Temperature profile		
	95°C	2min	
	95°C	15s	x 20 cycles



	A	B	C
	57°C	30s	
	68°C	30s	
	4°C	∞	

Note

PAUSE POINT Store in fridge for short-time or freezer until beads clean-up.

CAUTION

- During the preparation of index PCR, reagents and PCR plate must be on ice/cold block.
- Lead opens while on the cold block.
- Reagents are stored at -20°C.
- Gently mix and spin down reagents before using them. Avoid vortexing enzymes.
- **WARNIN!!!** Change pipette tips and PCR leads during work process!

Step8: Post-PCR clean up and library quality check

27m

8 Checklist before starting

Note

Where: Sequencing Lab

Timing: 30 min

Reagents

- AMPure XP beads, 60 mL: #A63881
- Qubit 1X dsDNA HS Assay Kit: #Q33231
- High Sensitivity D5000
- ScreenTape Assay: #5067-5593
- Freshly prepared 80% ethanol

Equipments/kits

- HulaMixture, Invitrogen™
- Qubit™ 4 Fluorometer
- TapeStation4200 with high-sensitive D5000 Screen Tape (#5067-5592)

Procedure



- 8.1 Place AMPure XP beads into the Hula mixer to resuspend the beads and to equilibrate to Room temperature for 00:15:00 . 15m
- 8.2 Total volume of PCR product (sample) / well in PCR plate is 25 μL .
- 8.3 Transfer 20 μL (!!! 0.8X) of Beads solution to each PCR product (25 μL) and mix with 100 μL pipette tips (10 times up and down), resulting in a 45 μL mixture.
- 8.4 Incubate for 00:05:00 at Room temperature . 5m
-
- 8.5 Place the tube or plate with the mixture into the magnetic rack for 2-4 min, until the liquid is clear.
- 8.6 Carefully remove 40 μL the liquid/supernatant with 100 μL pipette tips and discard it, keep 5 μL without disturbing the beads-pellet.
- 8.7 Wash the beads-pellet with 175 μL freshly prepared Fresh 80% Ethanol by gently dispensing it over the beads with 200 μL pipette tips. Let it rest for 00:00:30 and then remove the liquid. 30s
- 8.8 Repeat washing [go to step #6.7](#) .
- 8.9 Spin the tubes or plates in a minicentrifuge for 00:00:30 to collect all the residual liquid. 30s
- 8.10 Put the tubes or plates on a magnetic rack and remove the excess with 10 μL pipette tips.
- 8.11 Air dry for approximately 00:00:30 to evaporate the Ethanol. 30s
- 8.12 Remove the tube or plate from the magnetic rack. Add 16 μL nuclease free water and mix with a pipette (>10 times up and down) to resuspend the beads-pellet. 2m
- 8.13 Incubate for 00:02:00 at Room temperature .



8.14 Spin down for  00:00:30 if there is liquid on the wall of tubes or well.

3m 30s

8.15 Place the tube or plate back on the magnetic rack and wait for  00:03:00 until the liquid clears.

8.16 Transfer  15 μL of the liquid/supernatant to a new, labeled Eppendorf tube or a new plate.

8.17 Measure concentration of the cleaned PCR product with Qubit.

8.18 Pool 10 ng per sample into one Eppendorf tube to prepare a pooled library for library quality check with TapeStation following their standard protocols.

Note

PAUSE POINT Send for sequencing or store in the fridge  4 °C for short time or freezer  -80 °C for long time.