


Jul 18, 2020

🌐 Small-scale Expression and Solubility Testing of Proteins in BL21 E. coli

DOI

dx.doi.org/10.17504/protocols.io.261geojkwl47/v1



Vincent A Bielinski¹

¹J. Craig Venter Institute, Synthetic Biology & Bioenergy Group



Vincent A Bielinski

J. Craig Venter Institute, Synthetic Biology & Bioenergy Gro...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261geojkwl47/v1>

Protocol Citation: Vincent A Bielinski 2020. Small-scale Expression and Solubility Testing of Proteins in BL21 E. coli .
protocols.io <https://dx.doi.org/10.17504/protocols.io.261geojkwl47/v1>

Manuscript citation:

Modified from Savitsky et al 2010

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 24, 2020

Last Modified: July 18, 2020

Protocol Integer ID: 36106

Keywords: solubility screen for bacterial overexpression, solubility testing of protein, bacterial overexpression, solubility testing, protein, solubility screen

Abstract

Scaled-down solubility screen for bacterial overexpression constructs.

Attachments



1-s2.0-

S104784771000...

796KB

Guidelines

Use standard asptic techniques and proper PPE. Dispose of waste and used culture medium as appropriate.



Materials

Lysis Buffer

50 mM NaH_2PO_4 (pH 7.5-8.0)

500 mM NaCl

0.1% Triton X-100

10 mM imidazole

1 mg/mL lysozyme

10 uM Beta-mercaptoethanol

1 EDTA-free Protease inhibitor tablet/10mL buffer

Wash Buffer

50 mM NaPO_4 (pH 7.5-8.0)

500 mM NaCl

30 mM imidazole

Elution Buffer

50 mM NaPO_4 (pH 7.5-8.0)


500 mM NaCl

250 mM imidazole

TALON resin (Takara)

Troubleshooting

Safety warnings

 BME can smell noxious and be an irritant, so use a chemical hood when preparing stock for Lysis Buffer.

Before start

Make sure you have fresh buffers and keep samples on ice when noted.

Screening

- 1 Grow and induce BL21 cells in medium of your choice with appropriate temperature and inducer conditions. We often use Terrific Broth and culture cells at 30 °C, with induction carried out at 18 °C overnight (1mM IPTG or 0.5% L-arabinose).
- 2 Harvest 10 mL overnight expression culture by centrifugation for 10 minutes at 3000 x *g*.
- 3 Remove supernatant by decanting and resuspend pellet in 800 uL Lysis Buffer. Transfer cell suspension to a new 1.5 mL Eppendorf tube and store on ice.
- 4 Sonicate as appropriate for complete cell lysis (settings will depend on make and model of your sonicator).
- 5 Clarify lysates by centrifugation at 4 °C for 10 minutes at 12-15000 x *g*.
- 6 While lysates are spinning, equilibrate the Talon resin by washing 3x with 5 volumes of resin bedding.

NOTE: For each pulldown you do, you will need 25 uL of equilibrated Talon resin. This comes supplied as a 50% slurry in ethanol, so calculate how many pulldowns you will be doing and then add 50 uL. If you have 8, estimate you will need enough for 10 pull-downs. To equilibrate, shake the slurry well and remove 500 uL resin slurry with pipette. Transfer slurry to new Eppendorf tube. Pellet resin by spinning in microfuge for 30 seconds at 3000 x *g*, then remove supernatant carefully and add 1 mL of Lysis Buffer. Vortex resin briefly, then pellet again and wash beads 2x more with 1 mL of Lysis Buffer. After the third spin, remove supernatant and add 250 uL of Lysis Buffer to produce 500 uL of equilibrated resin.

- 7 Transfer 750 uL of clarified supernatant to a new Eppendorf tube.
- 8 Add 50 uL equilibrated resin to each Eppendorf tube and incubate for 1 hr at RT with end-over-end mixing.
- 9 After 1 hour, pellet resin with 1 min 3000 x *g* spin at RT, then remove supernatant with pipette.



- 10 Wash resin 3x with 1 mL of Wash Buffer, pelleting resin between each wash and removing supernatant.
- 11 After third wash, remove supernatant and resuspend resin in 50 uL of Elution Buffer.
- 12 Incubate at RT for 5-10 minutes, then quickly spin down beads and transfer eluate to clean Eppendorf tube.
- 13 Repeat elution again, combining eluates in single tube at the end (100 uL eluate).
- 14 Analyze 5-10 uL of eluate on an SDS gel and stain with Coomassie blue.