

Nov 22, 2025

Small-scale test expression of recombinant protein in *E. coli*

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlkydkdg5r/v1>

Ainsley Lederer¹, Prem Shrestha¹

¹University of Pittsburgh



Ainsley Lederer

University of Pittsburgh

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlkydkdg5r/v1>

Protocol Citation: Ainsley Lederer, Prem Shrestha 2025. Small-scale test expression of recombinant protein in *E. coli*. protocols.io <https://dx.doi.org/10.17504/protocols.io.j8nlkydkdg5r/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 04, 2025



Last Modified: November 22, 2025

Protocol Integer ID: 226473

Keywords: protein expression, test expression, protein production, small scale production, scale test expression of recombinant protein, small scale test production of recombinant protein, recombinant protein, small scale test production, scale test expression, procedure

Abstract

This procedure is meant to be used for the small scale test production of recombinant protein in *E. coli* before scaling up to larger volumes.

Materials

Ni-binding buffer (for resuspension):

50 mM Tris-HCl pH 7.5

100 mM NaCl

10% glycerol

5 mM beta-mercaptoethanol



Day 1

3h

- 1 Take LB agar plates with transformed *E. coli*. Select a single colony and add to  4 mL LB media with antibiotic. Place tube in  37 °C shaker.

- 1.1 It is convenient to keep 1000x stocks of antibiotics in the freezer.

	A	B	C	D
	Antibiotic	1000X Concentration	1X Concentration	Solvent
	Ampicillin	100 mg/mL	100 ug/mL	water
	Kanamycin	50 mg/mL	50 ug/mL	water
	Chloramphenicol	34 mg/mL	34 ug/mL	ethanol

- 2 Shake tubes until media becomes cloudy, typically  03:00:00 (2-3 hours) depending on cell strain.

3h

- 3 Take  500 µL of cells and add to  500 µL of 50% glycerol solution (so final of 25% glycerol). Store this at -80°C.

Note

The purpose of creating a glycerol stock at this step is to be able to use this exact colony of cells again in the future if the results are good! Each colony on a plate originated from a single *E. coli*. Yes, theoretically all colonies which grow on the antibiotic plate should act the same, but there are other factors at play that can affect the efficiency of protein production. So, if you plan to scale-up your protein production after the test expression, it is best to use the same stock of cells.

- 4 Take  1 mL of cells and put in a microcentrifuge tube. Centrifuge  14000 rpm, 4°C , 1 minutes . Aspirate the media completely and store the pellet at  -20 °C until day 2. This is your control, **uninduced** sample.



- 5 Add  2.5 μL of 1000X IPTG to the tube. IPTG concentration may vary depending on protein you are trying to express, but  0.5 millimolar (mM) final (0.5M stock) is a good starting point.
- 6 Reduce temperature of shaker to  18 °C and put tubes back in. Allow tubes to shake in incubator overnight.

Day 2

- 7 Take tubes from shaker. Pipette  1 mL of cells into a microcentrifuge tube and centrifuge  14000 rpm, 4°C for 1 minute. Aspirate media completely. This is your **induced** sample.
- 8 Place both uninduced and induced pellet samples on ice. Add  400 μL resuspension buffer. We use Ni-binding buffer (recipe in materials). Pipette to resuspend each pellet.
- 9 Sonicate each cell suspension. For QSonica sonicator, use the smallest tip. The settings are: 1.5 minutes, 03 s on, 08 s off, 30% amplitude.
- 9.1 Each cycle takes about 5 minutes. Make sure that the sonicator tip is submerged in the cell suspension, but not pressed against the side of the tube.
- 10 After all tubes have been sonicated, take  15 μL of each sample and label as "uninduced total" and "induced total". This represents all soluble *and* insoluble protein.
- 11 Spin cells down at  14000 rpm, 4°C for  00:10:00 . Take 15 μL of the supernatant, and label as "uninduced soluble" and "induced soluble". This represents the fraction of protein that is soluble in your buffer.
- 12 Add  5 μL 4X Laemmli sample dye to each tube. Place tubes in  95 °C heat block for 2 minutes. Load samples onto SDS-PAGE gel and analyze by Coomassie staining.

12.1



Note

What would the ideal result be? Ideally, you will see non-specific proteins in the uninduced samples, both total and soluble. In your induced samples, you will hopefully see a thick band at the expected MW for your recombinant protein in *both* the total and soluble samples. If you do not see your target protein, or if you see protein only in the induced total fraction, it indicates a solubility/expression issue that should be addressed before you scale up your protein production!